

SURGERY & ASC GUIDE

2026



Your **One Source** for a **Thriving Practice**



Where Quality, Affordability, and Surgical Confidence *Meet*

Trusted quality and performance
you can count on when it matters most

Ambulatory Surgery Centers are continuing to grow as care shifts to efficient, outpatient settings. At the same time, ASC teams are balancing procedure growth, cost management, and the need to deliver consistent, high-quality patient care. At **Henry Schein Medical**, we understand these challenges and are committed to supporting your success.

This catalog features a focused assortment of medical-surgical supplies and solutions selected specifically for the ASC environment, where reliability, speed, and consistency matter. With decades of experience in non-acute medical distribution, we provide access to trusted products, flexible sourcing options, and dedicated ASC sales support.

Our sales team works closely with centers like yours, providing the products, solutions, and support needed to improve efficiency, manage costs, and maintain clinical excellence. Backed by dependable distribution and responsive service, we help ensure your center has what it needs to operate smoothly.

**Thank you for choosing Henry Schein Medical as your ASC partner.
We look forward to supporting your continued growth.**

INDEX:

AEDs	83	Gloves	13–16	Scales	82
Capital Equipment.....	106	Hypodermics	58–60	Surgical Instruments.....	26–41
Cold Storage.....	61	Infection Prevention	45–48	Surgical Packs.....	7–11
Crash Carts	5–6	Instrument Reprocessing	42–44	Surgical Products	17–25
Cryosurgery	52	IV Administration	56–57	Suture Removal & IV Start Kits.....	12
Diagnostic Equipment.....	80–81	IV Catheters	54–56	Waste Management.....	60
Electrosurgery	49–51	Orthopedics	103–105	Wound Care	93–103
Exam Room Disposables.....	70–75	Patient Monitoring.....	84–85	Wound Closure	86–93
Exam Room Equipment	76–79	Pharmaceuticals	62–66		
Financing	53	Point-of-Care Testing.....	66–69		

HENRY SCHEIN AROUND THE WORLD



OPERATIONS
OR AFFILIATES
IN 33
COUNTRIES
AND
TERRITORIES



SERVING
MORE THAN
1 MILLION
CUSTOMERS



MORE THAN
25,000
TEAM SCHEIN
MEMBERS



AVERAGE OF
142,000
CARTONS
SHIPPED DAILY



51
DISTRIBUTION
AND MANUFACTURING
CENTERS

HISTORICALLY

99%

OF ITEMS ORDERED
ARE SHIPPED ON THE
SAME BUSINESS DAY
THE ORDER IS RECEIVED

MORE THAN
93 YEARS IN
BUSINESS



MORE
THAN
3,200
SUPPLIER
PARTNERS



300,000+
BRANDED AND
HENRY SCHEIN-BRAND
PRODUCTS IN STOCK

Essential Equipment for Your Ambulatory Surgery Center

Henry Schein Medical has partnered with world-class manufacturers to provide a **fully integrated** and **one-stop shop procurement** experience for your ASC's capital equipment needs.



ASC AREA: EQUIPMENT BY LOCATION

A SAMPLE OF OUR EXTENSIVE CAPITAL EQUIPMENT PORTFOLIO

>> ANESTHESIA MACHINES

GE Medical
 Mindray
 DRE Medical Group

>> C-ARMS

GE Medical
 Fuji

>> ECG/EKG

Midmark
 Welch-Allyn

>> ELECTRODES

Henry Schein
 3M
 Cardinal Health
 Nikomed

>> INFUSION PUMPS

B. Braun
 SchureMed

>> IV POLES & STANDS

Brewer
 Dukal
 Medline
 Pedigo

>> MUSCLE STIM ELECTRODES

Ambu
 Cardinal Health
 DJO

>> PAIN TABLES

Brewer
 Midmark
 Oakworks

>> PATIENT MONITORING

Edan Diagnostics
 Mindray
 Fukuda Denshi
 GE Medical
 Phillips
 Raine
 Welch-Allyn

>> STERILIZATION

Midmark
 Tuttnauer

>> ULTRASOUND

GE Medical
 Edan Diagnostics
 B. Braun Philips
 Samsung

>> WAITING ROOM

Furniture
 Bariatric Chairs
 Side Chairs
 Tables
 Benches
 Reception Desks
 Trash Cans
 Furnishings

>> PRE-OP

Medical Stretchers
 Patient Monitors
 Medical Recliners
 Headwalls
 IV Poles
 EKG Machines
 Trash Cans
 Hampers

>> OPERATING ROOM

OR Tables
 Anesthesia Systems
 Patient Monitors
 Lights & Booms
 Anesthesia Carts
 Instrument Tables
 OR Cabinetry
 Electrosurgical Units
 OR Integration
 Video Systems/
 Imaging
 Mayo Stands

>> POST-OP

Medical Stretchers
 Patient Monitors
 Overbed Tables
 Medical Recliners
 Crash Carts
 Back Walls
 Refrigeration
 Nurses Station
 Monitor Mounts
 Electrical, Gas &
 Suction Outlets

>> CENTRAL STERILE

High Vac Sterilizers
 Washer Disinfectors
 Processing Sinks
 Work Tables
 Ultrasonic Cleaners
 Reprocessing Sinks
 Packaging Tables
 Heat Lamps
 Shelving Units
 SCR Wrap Racks

>> CENTRAL SUPPLY

Shelving Units
 Scope Cabinets
 Supply Cabinets
 Case Carts
 Arco Bins
 Barcoding
 Closed Case Carts
 Assorted Shelving
 Central Supply
 Table
 Advanced Cold
 Storage Technology

>> SURGICAL EQUIPMENT

OR Integration
 Insufflators
 Smoke Evacuators
 Camera Systems
 Light Sources
 Grounding Pads
 Endoscopic
 Systems
 Disposables
 PAC Systems
 Electrosurgical Unit
 (ESU)

For more information on Henry Schein Medical's capital equipment options, consult with your dedicated Surgical Field Sales Consultant: ASC@HenrySchein.com

CRASH CART AUTO-REPLENISHMENT

powered by OnTraq

Crash cart management simplified

Tracking expiration dates and reordering crash cart supplies can be a cumbersome task that is subject to compliance risk. To overcome this challenge, the OnTraq Crash Cart Auto-Replenishment service ensures that enrolled emergency medications are tracked and automatically replaced just prior to expiration to maximize shelf life, save money and ensure readiness.

Medication and auto-replenishment service overview

- Provide HealthFirst with a list of medications to be enrolled
- OnTraq notifies you when a medication is expiring and is ready for replenishment
- HealthFirst ships you the expiring medication to replenish your crash cart
- Use the OnTraq online dashboard to view the status of medications at any time



Medications and equipment in low units of measure

- Extensive list of emergency medications available
- HealthFirst is licensed and uniquely qualified to repackage, distribute and track single unit medications
- Medication disposal service included

DSCSA compliant single dose medication distribution



Individual DSCSA Compliant Vial



- QR Code
- RFID
- NDC/Lot Number/Expiration
- No outer packaging to open to access medication



9887 v001 12/2018

OnTraq

POWERING CRASH CART AUTO-REPLENISHMENT

Simplify compliance across facilities with HealthFirst and OnTraq

OnTraq is an online compliance and account management tool. It allows you to view or print each cart's medications, expiration dates, lot numbers, national drug backorder status, and your upcoming medication replenishment shipments. OnTraq is included with the Crash Cart Auto-Replenishment service.

Stay prepared for audits. Tracking and managing compliance across multiple healthcare facilities no longer has to be a time-intensive, manual process. The OnTraq dashboard provides a bird's eye view of the compliance status across all sites within an organization. A single color-coded view displays the status of the crash cart's emergency medications.




Drill down by facility level or compliance area

Deeper insights and details are just a click away. Select a specific site, a specific cart, or compliance area (such as medications status) to see more details.

Easy reporting and task management

Local facility managers can use OnTraq to set up tasks for team members and receive alerts when those tasks are left undone. Summary reports show performance by task, employee and due date.



For more
Surgical Packs
see pages 8–11

Standard & Custom Procedural packs that meet **your needs**

Let our team collaborate with you to ensure that your pack program is optimized for your surgery practice. Custom surgical and procedure packs can help you reduce waste and quicken O.R. turnaround time.



LOOKING FOR MORE SUPPORT AND SAVINGS?

Ask your Henry Schein Surgical Field Sales Consultant to provide a pack assessment.



Breast Procedure Pack



(570-1690) 3/case

Contains: 1 of each: disposal bag, 700-cc emesis basin, 6000-cc ring basin, #15 carbon blade, 32-oz/1000-mL bowl, table back cover, Mayo stand cover, specimen cup with screw lid, chest drape, XL surgical gown, XL surgical gown with towel, waterproof label set, dual-tip skin marker, needle counter, 22 ga x 1½" needle, 25 ga x 1½" Eclipse® needle, 1200-cc pitcher, 6" ruler, Yankauer suction handle, 10-mL Luer-lock syringe, 60-mL irrigating bulb syringe, ¼" x 12" suction tubing & 54" x 54" CSR wrap; 2 of each: #10 carbon blades & 60-cc medicine cups; 4 of each: blue cotton towels, blue low-lint towels & utility drapes; & 10 of each: 4" x 4" 16-ply X-ray gauze sponges & 18" x 18" 4-ply laparoscopic X-ray sponges.



Single Basin Set



(570-1570) 10/case

Contains: 1-6000-cc ring basin & 1-54" x 54" CSR wrap.



Double Basin Set



(570-1571) 10/case

Contains: 1-6000-cc ring basin & 2-54" x 54" CSR wraps.



Major Set-up Procedure Pack



(570-1573) 6/case

Contains: 1 of each: disposal bag, 700-cc emesis basin, 6000-cc ring basin, #15 carbon blade, 32-oz/1000-mL bowl, cautery pencil, cautery-tip cleaner, 60-cc medicine cup, specimen cup with screw lid, scalpel holder, lid for medicine cup, standard-tip skin marker, needle counter, 1200-cc pitcher, Yankauer suction handle, 60-mL irrigating syringe bulb, ¼" x 12" suction tubing, 54" x 54" CSR wrap; 2 of each: #10 carbon blades, flexible light handle cover; 5 of each: 18" x 18" 4-ply X-ray laparotomy sponges; 6 of each: blue cotton towels; & 10 of each: 4" x 4" 16-ply X-ray gauze sponges.

Minor Set-up Procedure Pack



(570-1572) 8/case

Contains: 1 of each: 700-cc emesis basin, 6000-cc ring basin, 32-oz/1000-mL bowl, 60-cc medicine cup; & 54" x 54" CSR wraps.



Basic Eye Procedure Pack



(570-1583) 8/case

Contains: 1 of each: disposal bag, self-closure bag, 32-oz (1000-mL) bowl, back table cover, Mayo stand cover, ophthalmic drape, X-large surgical gown, X-large surgical gown with towel, waterproof label set, standard-tip skin marker, needle counter, 25 ga x 5/8" needle, eye pad, eye shield, 3-mL Luer-lock syringe, 1-cc Luer-slip syringe, 2-oz bulb syringe & instrument wipe; 2-60-cc medicine cups; 4 blue low-lint towels; 10 of each: 3" applicators & 4" x 4" 12-ply gauze sponges; & 20 eye spears.



Head and Neck Procedure Pack



(570-1586) 2/case

Contains: 1 of each: disposal bag, 700-cc emesis basin, 6000-cc ring basin, #10 carbon blade, #15 carbon blade, back table cover, Mayo stand cover, specimen cup with screw lid, half-drape, head drape, X-large surgical gown with towel, waterproof label set, dual-tip skin marker, needle counter, 18 ga x 1 1/2" needle, 25 ga x 1 1/2" needle, 8" x 3" Telfa pad, 1200-cc pitcher, 6" ruler, Yankauer suction handle, 10-mL Luer-lock syringe, 60-mL irrigating bulb syringe, 3/4" x 12" suction tubing & 54" x 54" CSR wrap; 2 of each: 32-oz (1000-mL) bowls, 60-cc medicine cups & X-large surgical gowns; 4 blue low-lint towels; 5-8" x 18" 4-ply X-ray lap sponges, 6 blue cotton towels; & 20-4" x 4" 16-ply X-ray gauze sponges.



ENT Basic Procedure Pack



(570-1585) 6/case

Contains: 1 of each: disposal bag, #15 carbon blade, 16-oz (500-mL) bowl, 32-oz (1000-mL) bowl, back table cover, Mayo stand cover, half-drape, split drape, X-large surgical gown, X-large surgical gown with towel, dual-tip skin marker, needle counter, 27 ga x 1 1/4" aluminum-hub A-bevel needle, 6" ruler, 10-mL Luer-lock control syringe & 3/4" x 12" suction tubing; 3 blue cotton towels; 4-60-cc medicine cups; & 10 of each: blank labels & 4" x 4" 16-ply X-ray gauze sponges.



Podiatry Procedure Pack



(570-1582) 4/case

Contains: 1 of each: 6 1/2" x 12" suture bag, 4" x 108" Esmark bandage, 4" x 5-yd tan cohesive bandage, back table cover, Mayo stand cover, lower-extremity drape, X-large surgical gown, X-large surgical gown with towel, waterproof label set, dual-tip skin marker, needle counter, 18 ga x 1 1/2" needle, 25 ga x 1 1/2" needle, 6" ruler, large impervious stockinette, 10-mL Luer-lock syringe, 30-mL Luer-lock syringe, 60-mL irrigating bulb syringe, 3/4" x 12" suction tubing, 18" x 18" 4-ply X-ray lap sponge, 4" x 4" 16-ply X-ray gauze sponge, 4" x 4" 12-ply gauze sponges; 2 of each: #15 carbon blades, 60-cc medicine cups & blue low-lint towels; 3-6 1/2" ziplock bags; 4 blue cotton towels.



Extremity Procedure Pack

HENRY SCHEIN®

(570-1576) 6/case

Contains: 1 of each: 6½" x 12" suture bag, 4" x 5-yd elastic bandage, 4" x 108" Esmark bandage, 700-cc emesis basin, #10 carbon blade, 32-oz (1000-mL) bowl, cautery pencil with holster, back table cover, Mayo stand cover, 60-cc medicine cup, large drape, lower-extremity drape, impervious split drape, cautery electrode, X-large surgical gown, X-large surgical gown with towel, standard-tip skin marker, needle counter, 23 ga x 1½" needle, 6" ruler, large impervious stockinette, 10-mL Luer-lock syringe, 2-oz bulb syringe, ¼" x 12-ft suction tubing, 4" x 4" 16-ply gauze, 4" x 4" 12-ply gauze; 2 of each: 6½" x 10" ziplock bags, #15 carbon blades, flexible light-handle cover, waterproof label set; 7 of each: 24" x 17" blue cotton towels.



Knee Arthroscopy Procedure Pack

HENRY SCHEIN®

(570-1575) 4/case

Contains: 1 of each: disposal bag, 6" x 5-yd elastic bandage, 6" x 12" Esmark bandage, 4" x 5-yd tan cohesive bandage, #11 carbon blade, 32-oz (1000-mL) bowl, back table cover, Mayo stand cover, 60-cc medicine cup, arthroscopy drape, large drape, impervious split drape, 1" x 8" Xeroform dressing, waterproof label set, standard-tip skin marker, needle counter, 18 ga x 1½" needle, 18 ga x 3½" spinal needle, 6" x 4-yd cast padding, 6" ruler, large impervious stockinette, 10-mL Luer-lock syringe, 30-mL Luer-lock syringe, 9" x 5" x 2" half-deep tray, ¼" x 12" suction tubing & 30" x 30" poly wrap; 2 of each: self-closure bags & breathable X-large surgical gowns; 8 blue cotton towels; & 10 of each: 4" x 4" 16-ply X-ray gauze sponges & 4" x 4" 12-ply gauze sponges.



Upper Extremity Procedure Pack

HENRY SCHEIN®

(570-1580) 4/case

Contains: 1 of each: disposal bag, back table cover, Mayo stand cover, 60-cc medicine cup, lower-extremity drape, three-quarter drape, 1" x 8" Xeroform dressing, X-large surgical gown, X-large surgical gown with towel, standard-tip skin marker, needle counter, 25 ga x 1½" needle, 6" ruler, large impervious stockinette, 10-mL Luer-lock syringe, 60-mL irrigating bulb syringe, 9" x 5" x 2" half-deep tray & ¼" x 12" suction tubing; 2-#15 carbon blades; 6 blue cotton towels; & 10 of each: 4" x 4" 16-ply X-ray gauze sponges & 18" x 18" 4-ply X-ray lap sponges.



Shoulder Arthroscopy Procedure Pack

HENRY SCHEIN®

(570-1578) 6/case

Contains: 1 of each: 6½" x 12" suture bag, 4" x 5-yd tan cohesive bandage, #11 carbon blade, #15 carbon blade, table back cover, Mayo stand cover, large drape, split shoulder drape, impervious split drape, XL surgical gown with towel, standard-tip skin marker, needle counter, 18 ga x 1½" needle, 18 ga x 3½" spinal needle, 6" ruler, medium impervious stockinette, 10-mL Luer-lock syringe, 30-mL Luer-lock syringe, 9" x 5" x 2" half-deep tray, ¼" x 12" suction tubing, U-drape, 4" x 4" 16-ply gauze X-ray sponge; 2 of each: 60-cc medicine cups & XL X-long breathable surgical gowns; 8 blue cotton towels.



Shoulder Beach Chair Procedure Pack



(570-1579) 4/case

Contains: 1 of each: disposal bag, 4" x 5 yd tan cohesive bandage, #11 carbon blade, #15 carbon blade, table back cover, Mayo stand cover, large drape, beach chair shoulder drape, impervious split drape, XL surgical gown with towel, standard-tip skin marker, needle counter, 18 ga x 1½" needle, 18 ga x 3½" spinal needle, 6" ruler, medium impervious stockinette, 10-mL Luer-lock syringe, 30-mL Luer-lock syringe, 9" x 5" x 2" half-deep tray, ¼" x 12" suction tubing; 2 of each: 60-cc medicine cup, XL X-long breathable surgical gown; 8 blue cotton towels; & 10-4" x 4" 16-ply gauze X-ray sponges.

Basic Drape Pack 2




Sterile

Contents include:

- 4 utility drapes with tape, 15" x 26"
- 1 large drape, 60" x 76"
- 1 surgical gown, nonreinforced, set-in sleeve, large, overwrap with towel
- 1 surgical gown, nonreinforced, set-in sleeve, X-large
- 1 back table cover, zone-reinforced, 44" x 90"
- 1 Mayo stand cover, 23" x 54"

(570-1599) 4/case



Basic Drape Pack 3



(570-1600) 6/case

Contains: 1 of each: disposal bag, 4" x 5 yd tan cohesive bandage, #11 carbon blade, #15 carbon blade, table back cover, set-up cover, 60-cc medicine cup, femoral angiography drape, 5" curved Halstead mosquito forcep, sharps holder, waterproof label set, dual-tip skin marker, 18 ga x 1½" needle, 25 ga x 1½" needle, 18 ga x 7 cm needle introducer, 10-mL Luer-lock 3-ring control syringe; 2 each: XL surgical gowns, 20-mL Luer-lock syringes, white paper towels; 4 of each: 10-mL Luer-lock syringes, 5½" forcep towels, blue cotton towels; & 40-4" x 4" 4-ply nonwoven sponges.

Basic Drape Pack 3




Sterile

Contents include:

- 4 utility drapes with tape, 15" x 26"
- 1 half drape, 44" x 60"
- 1 large drape, 60" x 76"
- 1 surgical gown, nonreinforced, set-in sleeve, large, overwrap with towel
- 1 surgical gown, nonreinforced, set-in sleeve, X-large
- 1 back table cover, zone-reinforced, 44" x 90"
- 1 Mayo stand cover, 23" x 54"

(570-1600) 6/case



Cardiac Catheterization Procedure Pack



(570-1588) 6/case

Contains: 1 of each: 36" x 40" band bag, #11 carbon blade, 2500-mL guidewire bowl, 32-oz/1000-mL bowl, table back cover, set-up cover, 60-cc medicine cup, femoral angiography drape, 5" curved Halstead mosquito forcep, sharps holder, waterproof label set, dual-tip skin marker, 18 ga x 1½" needle, 25 ga x 1½" needle, 18 ga x 7 cm needle introducer, 10-mL Luer-lock 3-ring control syringe; 2 each: XL surgical gowns, 20-mL Luer-lock syringes, white paper towels; 4 of each: 10-mL Luer-lock syringes, 5½" forcep towels, blue cotton towels; & 40-4" x 4" 4-ply nonwoven sponges.


EENT Drape Pack



Sterile

- 1 EENT drape, 76" x 124", with 3" x 25" split
- 1 head drape, 50" x 36"
- 1 surgical gown, nonreinforced, set-in sleeve, large
- 1 surgical gown, nonreinforced, set-in sleeve, large, overwrap with towel
- 1 back table cover, zone-reinforced, 44" x 90"
- 1 Mayo stand cover, 23" x 54"
- 1 white hand towel, 15½" x 22"
- 1 disposable bag

(570-1652) 4/case



SUTURE REMOVAL & IV START KITS



Suture Removal Kit with Plastic Handle Instruments



Sterile, Single-Use

Kit features plastic-handled Littauer scissors with plastic thumb forceps. Flatter, thinner-blade scissors provide easier suture removal by allowing blade to slip under suture. Not made with natural rubber latex.

Hard plastic tray includes:

- 1 pair of hook-tipped 4½" Littauer plastic handle scissors
- 1 pair of 4" PosiGrip™ plastic thumb forceps
- 1–3" x 3" gauze sponge

(900-4687) ea



Staple Removal Kit



Used to easily remove surgical staples from patients with minimal discomfort. Not made with natural rubber latex.

Hard plastic tray includes:

- 1 stainless-steel skin staple remover
- 1–3" x 3" gauze sponge

(570-0620) ea



Suture Removal Kit with Iris Scissors



Sterile, Single-Use

Not made with natural rubber latex.

Hard plastic tray includes:

- 1 pair of small-loop, 4½" straight metal iris scissors (sharp/sharp)
- 1 pair of serrated, 4¾" metal forceps
- 1–3" x 3" gauze pad
- 1 large alcohol prep pad

(112-5822) ea



IV Start Kit with Alcohol Prep Pad



- 1 transparent IV dressing
- 1 IV change label
- 2 rolls of 2" x 2" 4-ply gauze
- 1 medium PVP iodophor prep pad
- 1 large alcohol prep pad
- 1 roll of ¾" x 18" transparent tape
- 1 blue ¾" x 18" tourniquet

(570-3407)Ea.



Suture Removal Kit with Littauer Scissors



Sterile, Single-Use

Not made with natural rubber latex.

Hard plastic tray includes:

- 1 pair of hook-tipped 4½" metal Littauer scissors
- 1 pair of 4" metal forceps
- 1–3" x 3" gauze sponge

(570-0619) ea

See pages 85–93 for wound closure products.



GUESS WHAT CAME IN WITH YOUR PATIENTS?

YOU DESERVE THE RIGHT PROTECTION!



FLUIDSHIELD* ASTM LEVEL 3 PROCEDURE MASKS

Fog-Free Foam Strip,
SO SOFT* Lining, Earloops



PURPLE* NITRILE EXAM GLOVES



PROTECTIVE APPAREL

Including Footwear
and Bouffants



AAMI LEVEL 3 SURGICAL GOWNS



QUICK CHECK* STERILIZATION WRAP

In a Wide Range of
Sizes and Weight.



PROCEDURE KITS and CUSTOM TRAYS

Ask your Henry Schein Rep about the full HALYARD* Assortment!

Ansell

HENRY SCHEIN[®]
MEDICAL

**PRECISION YOU CAN FEEL.
DURABILITY YOU CAN COUNT ON.**



GAMMEX[™] PI Hybrid Micro
Hybrid synthetic micro surgical gloves offering the comfort of polyisoprene and strength of neoprene

- **Smooth with micro-textured finish:** Delivers secure grip and handling of instruments
- **High-tactility micro gloves:** Suitable for delicate procedures
- **HYBRID Technology:** Delivers latex-like comfort and high durability

SIZE	ITEM CODE	SIZE	ITEM CODE
5.5	[149-3350]	7.5	[149-3354]
6	[149-3351]	8	[149-3355]
6.5	[149-3352]	8.5	[149-3356]
7	[149-3353]	9	[149-3357]

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.

Ansell, [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending; www.ansell.com/patentmarking. © 2025 Ansell Limited. All Rights Reserved.



CRITERION[™] Latex Surgical

CRITERION

Advanced polymer coating on the inner surface for easy donning when hands are wet, damp, or dry. Microtextured surface and optimized tack level provide excellent grip in all conditions. Supersoft latex film and anatomical design provide exceptional comfort to reduce hand fatigue, even during long procedures. Minimum cuff length of 275 mm for extended coverage and protection. Beaded cuff. 6.3-mil thickness at the palm. White.

Gloves..... 50 pair/box

Specify:

5½	(570-2102)	7½	(570-2106)
6	(570-2103)	8	(570-2107)
6½	(570-2104)	8½	(570-2108)
7	(570-2105)	9	(570-2109)



CRITERION[™] Polychloroprene Surgical

CRITERION

Sterile Polychloroprene Powder-Free Surgeons
Features an anatomically correct design for a more natural fit. Polymer-coating finishes that allow for damp- or dry-hand donning. Textured fingers and palms for enhanced grip and control. Ultralow modulus for enhanced comfort and reduced hand fatigue. Minimum cuff length of 300mm for extended coverage and protection. Beaded cuff to reduce the risk of glove roll-down once worn. Ideal for surgical procedures. 5.91-mil thickness at palm. Sterile. Natural color.

Gloves..... 50 pair/box

Specify:

5½	(112-7079)	7½	(112-7083)
6	(112-7080)	8	(112-7085)
6½	(112-7081)	8½	(112-7086)
7	(112-7082)	9	(112-7088)



CRITERION[™] Polyisoprene Surgical

CRITERION

Sterile Polyisoprene Powder-Free Surgeons
• Powder-free, sterile, and single-use
• Extended and beaded cuff prevents roll-down effect
• Fully microtextured exterior and polymer-coated interior
• Not made with natural rubber latex

Gloves..... 50 pair/box

Specify:

5½	(570-1119)	7½	(570-1123)
6	(570-1120)	8	(570-1124)
6½	(570-1121)	8½	(570-1125)
7	(570-1122)	9	(570-1126)



MICROFLEX® XCEED®
ANSELL

Nitrile Powder-Free Exam
Industry-revolutionizing nitrile formulation and manufacturing process delivers enhanced comfort, strength, and value. Higher-count box. Textured. Packaging is plastic-free and recyclable. 2.8-mil thick. Light Blue.

Gloves..... 250/box
Specify:
X-Small (145-7353) Large (145-7356)
Small (145-7354) X-Large* (145-7357)
Medium (145-7355)

*230/box



MICROFLEX® Neogard™ Touch
73-737 Neoprene
ANSELL

Chloroprene Powder-Free Exam
• Comfortable neoprene (polychloroprene) exam gloves
• Certified for exceptional wet grip
• Ambidextrous
• Beaded cuff
• Internal glove surface is polymer-coated
• Textured fingers
• 3.5-mil fingertip thickness
• 2.8-mil palm thickness
• Green

Gloves..... 200/box
Specify:
X-Small (142-5299) Large (142-5302)
Small (142-5300) X-Large (142-5303)
Medium (142-5301)



MICROFLEX® Ultraform®
ANSELL

Nitrile Powder-Free Exam
Designed with ERGOFORM™ Ergonomic Design Technology, providing musculoskeletal support with 20% less muscle exertion when performing repetitive tasks. Thin design for superior tactile sensitivity. Soft nitrile formulation forms to hands for comfort during extended wear. Protection from Type I latex allergy in health-care workers or patients. Textured fingertips. 9.1" length. 3.9-mil fingertip thickness. Cobalt Blue.

Gloves..... 300/box
Specify:
X-Small (141-1425) Medium (141-1427)
X-Small/Small (142-6336) Medium/Large (142-6338)
Small (141-1426) Large (141-1428)
Small/Medium (142-6337) X-Large* (141-1429)

*250/box



MICRO-TOUCH® NitraFree™
ANSELL

Nitrile Powder-Free Exam
Made without accelerators to protect from allergic reactions. Thin and lightweight to ensure excellent dexterity while working. Tested and approved for use with chemotherapy drugs. Palm thickness 3.1 mil. Finger thickness 5.5 mil. Packaging is plastic-free and recyclable.

Gloves..... 100/box
Specify:
X-Small (685-4447) Large (685-4728)
Small (685-4387) X-Large (685-0140)
Medium (685-7305)



BeeSure® Nitrile
ECOBEE

Nitrile Powder-Free Exam
Ultrasoft nitrile. Optimal tactile sensitivity and contours ergonomically to hands for reduced fatigue. Packaging made from 80% recycled fiber, making it eco-friendly. Finger textured. 3.7-mil thick. Light Blue.

Gloves..... 100/box
Specify:
X-Small (111-8551) Medium (111-8553)
Small (111-8552) Large (111-8554)



INSPIRE™
CRANBERRY

Nitrile Powder-Free Exam
Features InSoft formulation for increased comfort and superior tactile sensitivity. Fingertip textured for enhanced grip in both dry and wet conditions. 2.0-mil thick. Blue.

Gloves..... 300/box
Specify:
X-Small (133-9202) Large (133-9205)
Small (133-9203) X-Large* (133-9206)
Medium (133-9204)

*250/box

NEW

Experience the Softest Criterion™ Chloroprene Exam Glove.



- Chloroprene material provides the comfort and fit of latex but with the synthetic strength of nitrile
- Precision performance provides dexterity and extended-wear comfort
- Puncture and tear resistance for high-risk procedures
- Microtextured fingertips for a secure grip and improved control during procedures
- Interior polymer coating allows for easy donning when hands are wet, damp, or dry
- Beaded cuff helps prevent rollback. Cuff length at 240mm (9.5")
- 3.54 mil/ 0.09 mm thickness at palm
- Ambidextrous
- Conforms to ASTM D3577-19 and ISO 10282-2014

AQUA BLUE COLOR

SIZE	QUANTITY	ITEM CODE
XS	100/box	570-5400
S	100/box	570-5401
M	100/box	570-5402
L	100/box	570-5403
XL	100/box	570-5404



SURGICAL PRODUCTS

HENRY SCHEIN DRAPES

Features:

- Sterile, Single-Use
- Trilayer SMS fabric is soft, drapable and naturally repellent to fluids
- Tri-layer fabric construction provides an excellent barrier to strikethrough
- Excellent antflammability properties provide additional protection for the patient and staff when using lasers and electrocautery devices
- Not made with natural rubber latex



Non-Fenestrated Drapes



Sterile, Polylined, Tissue/Poly/Tissue

- Heat sealed along 2 outer edges to ensure a pinhole-free plastic layer that guarantees maximum protection against bacterial strike-through from below
- Features polyethylene lining between 2 layers of absorbent-strength tissue, virtually eliminating risk of bacterial migration and strike-through
- All sterile fields are available in convenient dispenser boxes of 50 each

White, 18" x 26" 50/box
(900-4686)



Fenestrated Drapes



Sterile, with 3" Round Fenestration
Tissue/Poly/Tissue

- Allows access to the operative site while it protects surrounding area
- 3" round fenestration option is die cut and heat sealed, leaving a smooth, even edge that prevents linting
- Not made with natural rubber latex

White, 18" x 26" 50/box
(900-4685)




Mayo Stand Cover



Sterile

Regular, 23" x 54" 30/pkg
(570-1616)

X-Large, 29½" x 57" 30/pkg
(570-1618)



Split-Sheet Drapes




Sterile

Drape size: 53" x 76", split size: 2½" x 30" with hook-and-loop tube/cord holders, 2" adhesive tape
(570-1647) 26/case

Drape size: 76" x 110", split size: 2½" x 37" with fabric tube/cord holders, 1" adhesive tape
(570-1658) 24/case


Drape size: 76" x 124", split size: 3" x 25" with fabric tube/cord holders and 1" adhesive tape
(570-1648) 22/case



Utility Drapes with Tape



Sterile, 15" x 26" 100/case
(570-1684)



Chest Drapes



Sterile, 100" x 126" 12/case
(570-1626)




Ophthalmic Drapes



With Full Incise, 48" x 51" 30/case
(570-1637)

With Full Incise and Pouch 30/case
(570-1639)



General Drapes



Sterile

Half, 44" x 60" 22/case
(570-1656)

Medium, 44" x 76" 20/case
(570-1617)

Large, 60" x 76" 20/case
(570-1657)

Large, Reverse Fold, 60" x 76" 20/case
(570-1663)

Back Table Covers, Zone Reinforced




Sterile

44" x 75" 20/case
(570-1596)

44" x 90" 20/case
(570-1597)

50" x 90" 20/case
(570-1598)



O.R. Towels



- Prewashed, 100% high-quality cotton
- Highly absorbent and virtually lint-free
- All sterile towels are CSR-wrapped and have a special corner fold that allows for aseptic presentation into a sterile field
- Not made with natural rubber latex
- 17" x 27"
- Blue

Non-Sterile (570-1113)	100/case
Sterile 2s (570-1117)	2/pkg
4s (570-1116)	4/pkg
6s (570-1114)	6/pkg



Laparotomy Sponges



- Prewashed, 100% high-quality cotton with tight weave
- Cross-stitched with super-locked edges to minimize linting and loose threads
- X-ray detectable
- Banded 5s
- Not made with natural rubber latex

18" x 18" (570-1115)	5/pkg
-------------------------------	-------



Surgical Skin Markers



These surgical marking pens provide an accurate and reliable way to mark the incision site prior to surgery.

- Prep-resistant
- Nonabrasive tip
- Large ink reservoir

Non-Sterile, Fine Tip, Ruler, and Up to 9 Labels (900-7018)	100/box
Sterile, Fine Tip, Ruler, and Up to 6 Labels (900-7019)	50/box
Mini (WriteSite Jr.), Dark Blue (900-7020)	100/box

Povidone-Iodine Swabsticks and Prep Pads



Sterile

Topical antiseptic used for skin preparation that helps to reduce bacteria that can cause infection. Saturated with 10% USP povidone-iodine.

Single Pack, 4"L (570-1549)	25 pkg/box
Contains: 1 swabstick per pkg.	
Triple Pack, 4"L (570-1548)	25 pkg/box
Contains: 3 swabsticks per pkg.	
Single Pouches, 1½" x 1½", Folded (570-1543)	1000/case
Contains: 1 per pouch.	



Double-Sided Shave Prep Razors



- Disposable
 - Sharp stainless steel blades for smooth shaving
 - Nonslip grip handle
- (900-7665)

100/box

BZK Towelettes



- Premoistened antiseptic towelettes used for general cleansing for patients and staff
- Cleans and cools simultaneously
- Provides sting-free cleansing of all minor cuts, scrapes, and burns
- Nontoxic and nonirritating
- Alcohol-free
- Unscented

(112-6735)	100/box
------------------	---------




Alcohol Prep Pads



Sterile

- Soft, absorbent, and non-woven pads impregnated with 70% isopropyl alcohol USP
- Primarily used as a topical antiseptic; disinfects prior to administering injections

Medium, 30 mm x 30 mm, Folded (570-2219)	200/box
Large, 70 mm x 62.5 mm (570-1307)	200/box



Surgical Blade Handles



Standard Stainless Steel #3, For Blades #10, #11, #12B, #15 & #15C (100-7520)	ea
German Stainless Steel #3, For Blades #10, #11, #12B, #15 & #15C (104-5814)	ea
#4, For Blades #20, #21, #22 & #23 (100-3973)	ea
#4, For Blades #20, #21, #22 & #23 (104-7296)	ea

*Lightweight, 12 gm (average)

For complete selection of Henry Schein brand surgical blades, see page 28



ChloraPrep™ Patient Preoperative Skin Preparation with Sterile Solution™

BD

- Sterile solution: Patented sterilization process of the antiseptic solution achieves a minimum sterility assurance level (SAL) of 10⁻⁶
- Sterile delivery: The sterile, single-use applicator facilitates aseptic touchless technique by moving the clinician's hand away from the patient's skin, helping to reduce the risk of extrinsic contamination
- Proven efficacy: The FDA-approved dual antiseptic of 2% chlorhexidine gluconate (CHG) + 70% isopropyl alcohol (IPA) kills bacteria more effectively than single agents, providing both immediate and persistent activity
- Preoperative indication: The only CHG and IPA skin antiseptic to offer 2-log reduction for dry sites, and 3-log reduction for wet sites, on its entire portfolio

ChloraPrep™ Applicator

#930480, 1 mL, Clear60/box
 (137-4623)

#930299, FREPP™, 1.5 mL, Clear (137-4627)	20/box
#930400, 3 mL, Clear (137-4626)	25/box
#930415, 3 mL, Hi-Lite Orange™ (137-4625)	25/box
#930700, 10.5 mL, Clear (137-2713)	25/box
#930715, 10.5 mL, Hi-Lite Orange™ (137-4624)	25/box
#930725, 10.5 mL, Scrub Teal™ (137-2711)	100/case
#930800, 26 mL, Clear (137-2716)	25/case
#930815, 26 mL, Hi-Lite Orange™ (137-2717)	25/case
#930825, 26 mL, Scrub Teal™ (137-2715)	25/case



E-Z Scrub® Premium Surgical Brushes

BD

Sterile, 30 Scrub Brushes/Dispenser

This line of products offers a wide range of efficacious antimicrobial formulations specifically designed to enhance performance as they reduce irritation to skin.

- Proprietary polyethylene soft bristles enhance mechanical action while reducing skin irritation
- Soft sponge provides gentle cleaning on irritated areas or when firm mechanical action is not required

With PCMX, Deluxe

(987-0287)30/box

With 4% CHG

(987-6022)30/box

With Iodophor

(780-3379)30/box

ChloraPrep® Patient Preoperative Skin Preparation

BD

Sterile

- 2% chlorhexidine gluconate (CHG) with 70% isopropyl alcohol
- CHG has demonstrated a 50% reduction in the incidence of catheter-related bloodstream infections, compared to povidone iodine
- Rapid activity against gram-positive and gram-negative bacteria
- Nonlatex, nonallergenic tip
- Convenient packaging
- Maintains a 2-year shelf life

1.75-mL Swabstick, Single
 (147-3823)48/box

5.25-mL Swabsticks, Triple
 (147-3826)40/box

Dermablade ACCUTEK

Sterile

This flexible, one-piece shave biopsy instrument has easy-to-pick-up grips that provide stability and control of excision. MicroCoat® surgical blade has a polymer coating to provide sharp, smooth, resistance-free cutting. Individually packaged.

(143-7056)50/box



VenaFlow Elite ENOVIS/AIRCAST

Designed as a prophylaxis for deep-vein thrombosis, the VenaFlow System combines 2 proven technologies— rapid inflation and graduated sequential compression—that work to mimic ambulation and accelerate venous velocity. For the superior emptying of veins, VenaFlow also incorporates asymmetric compression. The breathable and comfortable calf, foot, and thigh cuffs enhance the patient experience and are compatible with one universal pump.

System with Battery

(116-9948)ea

Calf Cuff, 19"

(114-8778) 1 pair

Foam Calf Cuff, 19"

(115-8712) 1 pair

Calf Cuff, X-Large, 22"

(116-1634) 1 pair

Bariatric Calf Cuff, 30"

(116-8326) 1 pair

Foot Cuff, 16"

(116-4740) 1 pair

Johnson & Johnson
MedTech

Do more with less



To see important
safety information,
please scan QR code


SURGICEL™ Powder Absorbable Hemostat stopped bleeding

~50%

faster using

~75%

less product required compared to PerClot® PHS or ARISTA™ AH¹*



* In a comparative quantitative analysis of ARISTA™ AH, PerClot® PHS, and 4DryField PH vs SURGICEL™ Powder in animal models.

1. Stark M, Allen YW, Corrigan B, et al. Comparative analyses of the hemostatic efficacy and surgical device performance of powdered oxidized regenerated cellulose and starch-based powder formulations. 2024. Ethicon, Inc.

Please refer always to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

The third-party trademarks used herein are trademarks of their respective owners.

© Johnson & Johnson and its affiliates 2025. C_US_SRG_BIOS.407276.EN



Go beyond Original

As bleeding needs evolve, go beyond SURGICEL™ Original Absorbable Hemostat with our advanced SURGICEL™ hemostat portfolio.

SURGICEL™
FIBRILLAR™
Absorbable
Hemostat

SURGICEL™
POWDER
Absorbable
Hemostat

SURGICEL™
SNoW™
Absorbable
Hemostat

Johnson & Johnson
MedTech

For complete indications, contraindications, warnings, precautions, and adverse reactions, please reference full package insert.
© Johnson & Johnson and its affiliates 2025. C_US_SRG_BIOS_407276.EN

SURGICEL™ Absorbable Hemostats and SURGIFOAM® Sponges J&J MEDTECH SURGERY

Solutions for a full range of continuous oozing bleeds.

SURGICEL™ #3013SP, Powder Hemostats, 3 g (143-7960).....	5/box
#2081, SNoW™ Hemostats, 1" x 2" (117-9356).....	10/box
#2082, SNoW™ Hemostats, 2" x 4" (117-9357).....	10/box
#2083, SNoW™ Hemostats, 4" x 4" (117-9358).....	10/box
#1961, FIBRILLAR™ Hemostats, 1" x 2" (555-5021).....	10/box
#1962, FIBRILLAR™ Hemostats, 2" x 4" (934-3367).....	10/box



#1963, FIBRILLAR™ Hemostats, 4" x 4" (555-8748).....	10/box
#1946, NU-KNIT™ Hemostats, 6" x 9" (617-7966).....	10/box
#3123SPEA, Endoscopic Applicators (145-7269).....	5/box
SURGIFOAM® #1969, Gelatin Sponges, 1 cm x 1 cm x 1 cm (654-7781).....	24/box
#1972, Gelatin Sponges, 2 cm x 6 cm x 7 mm (654-5663).....	12/box
#1974, Gelatin Sponges, 8 cm x 12.25 cm x 10 mm (654-1617).....	6/box



loban™ 2 Antimicrobial Incise Drapes SOLVENTUM, FORMERLY 3M HEALTH CARE

- Designed with continuous broad-spectrum antimicrobial activity in the drape adhesive
- Clinically shown to help reduce the risk of wound contamination

#6640EZ, 17¼" x 13¾" (777-5899).....	10/box
#6650EZ, 23" x 17" (777-0599).....	10/box
#6651EZ, 33" x 23" (777-2209).....	10/box

DuraPrep™ Surgical Skin Prepping Solution SOLVENTUM, FORMERLY 3M HEALTHCARE

Patient preoperative skin preparation contains iodine povacrylex (0.7% available iodine) and isopropyl alcohol (74% w/w). Includes applicator.
6 mL
(777-4189)ea
26 mL
(777-0115)ea





Bestform Post-Surgical Bra BL INTIMATES

- 2-ply stretch cotton blend (95% cotton, 5% spandex) provides comfort and support
- Unlined cups without wires or padding feel lightweight and breathable
- Front-closure adjustable clasp allows for custom fit
- Wide stretch straps help alleviate shoulder strain and provide support
- White and Nude

White.....ea

Specify:

Size 34	(143-9663)	Size 40	(143-9667)
Size 36	(143-9664)	Size 44	(147-7174)
Size 38	(143-9665)	Size 48	(147-7176)

Nude.....ea

Specify:


Size 34	(148-5603)	Size 42	(148-5607)
Size 36	(148-5604)	Size 44	(148-5608)
Size 38	(148-5605)	Size 48	(148-5610)
Size 40	(148-5606)		

Spinal & Epidural Needles

MYCO MEDICAL

Reli® brand needles are a value alternative to leading brands and offer superior quality in a broad range of sizes.

- Precision-point computerized grinding produces sharp, consistent needles
- Large, clear hubs for enhanced tactile feel and visualization of cerebrospinal fluid (CSF)
- Exact fit between stylet and cannula bevel to minimize coring and trauma
- Color-coded stylet caps for easy identification of needle gauge



Quincke Point Spinal Needles

22g x 3.5" 25/box
(731-1020)
 22g x 5" 25/box
(731-0230)
 22g x 5"–Short Bevel 25/box
(117-9090)
 25g x 3.5" 25/box
(116-8930)
 25g x 5" 25/box
(106-3704)

Tuohy Epidural Needles

18g x 3.5" 25/box
(818-9414)
 20g x 3.5" 25/box
(694-5927)
 22g x 3.5" 25/box
(694-1963)
 20g x3.5" –Fixed Wing 25/box
(116-5010)



Revolutionizing vascular healthcare with affordable, high-quality devices

Delivering innovative vascular medical devices, with a special emphasis on Outpatient and Ambulatory Surgery Centers (OBL/ASCs).

Falexexpand HP PTA Balloon Catheter


High-Pressure with RBP up to 24atm
0.035" OTW, 6F, Dia:3-12mm, BL: 20-200mm, CL: 40cm, 75cm, and 135cm

High-pressure and semi-compliance PTA balloon catheters with soft and tapered tips for better lesion crossing with hydrophilic coating and small profiles for challenging stenosis



Falspeed Semi-Compliance PTA Balloon Catheter

Semi-compliance with Dia: 2.0 to 7.0mm, BL: 20 to 300mm, CL: 150cm (OTW), Tip Crossing: 0.0225", 4F to 6F, NP: 6atm



014" Workhorse Guidewire (300cm)

Distal 40cm high-load Nitinol for fast recovery with Stainless Steel body for predictable tracking.
Shapeable tip design with 0.8g & 1.5g tip weight options; moderate support



018", 032", and 035" Workhorse Guidewire

Nitinol core wire guarantees great tip shape resistance, outstanding kink resistance and 1:1 torque control
Straight, Angled, and Shapeable tips, Standard and Stiff shaft



Balancing Quality & Affordability

Contact your Henry Schein representative or visit www.worldwidevascular.com to learn more

Helping protect your every **move.**



The Bard-Parker® scalpel sheath glides with ease while helping to protect your every pass. For more information or to place an order please contact:

www.henryschein.com/medical

Item #	Description	Size
(780-8145)	Bard-Parker® Stainless Steel Blades	15
(780-3523)	Bard-Parker® Stainless Steel Blades	10
(780-5269)	Bard-Parker® Stainless Steel Blades	11
(780-1424)	Bard-Parker® Carbon Rib-Back™ Blades	15
(780-5031)	Bard-Parker® Carbon Rib-Back™ Blades	10
(987-6084)	Bard-Parker® Carbon Rib-Back™ Blades	11
(987-9248)	Bard-Parker® Protected Disposable Scalpels	15
(987-5136)	Bard-Parker® Protected Disposable Scalpels	11
(987-0484)	Bard-Parker® Protected Disposable Scalpels	10
(780-4962)	Bard-Parker® Disposable Scalpels	15
(780-5593)	Bard-Parker® Disposable Scalpels	11
(780-8028)	Bard-Parker® Disposable Scalpels	10



Bard-Parker® is a registered trademark of Aspen Surgical Products, Inc. Aspen Surgical reserves the right to make changes without notice in design, specifications and models.

©2023 Aspen Surgical Products, Inc. ALL RIGHTS RESERVED.

M-SM-CAT-1083297 rev 2 19-Jan-2023 ENG - US



premierdentals.com/promise

Monsel's Ferric Subsulfate Solution

Pre-Thickened, Ready-to-Use!

Monsel's is a thickened formulation of the popular hemostatic agent that can be applied directly to the surgical site following excision of diseased tissue.

(100-9284) Box of 12 single-application 8ml vials and 12 rayon tip applicators



Also available, for cervical visualization,
(124-9854) Acetic Acid and (100-9046) Lugol's Solution*

*DEA license required.

*For more information on these and other Premier products,
please contact your Henry Schein Sales Consultant.*



SUMMIT MEDICAL™

Our extensive product portfolio of single-use instruments serves ENT, Ophthalmology, Gynecology, and Aesthetics. Speak to your Henry Schein representative for more information.



EAR, NOSE & THROAT



OPHTHALMOLOGY



GYNECOLOGY



AESTHETICS

YOUR TRUSTED SOURCE FOR SINGLE-USE SURGICAL ESSENTIALS.

Scan the QR code to learn more about our products.



TruCath® Indwelling Catheters

The TruCath Foley range includes configurations with various French and balloon sizes and lumens providing suitable options for diverse patient needs, anatomical variations, and medical conditions.

- 100% silicone material provides chemical and thermal stability, low surface tension, and hydrophobicity
- Silicone-elastomer coating helps reduce the potential for encrustation and mucosal irritation
- Balloons are pretested and inspected for integrity and position verification to ensure proper, upright position in the bladder

Item #	Description
(148-6927)	TruCath® 100% Silicone Foley Catheter, 2-Way, 14 Fr, 10 ml Balloon
(148-6928)	TruCath® 100% Silicone Foley Catheter, 2-Way, 16 Fr, 10 ml Balloon
(148-6931)	TruCath® 100% Silicone Foley Catheter, 2-Way, 18 Fr, 10 ml Balloon
(148-4854)	TruCath® Silicone-Elastomer Coated Latex Foley Catheter, 2-Way, 14 Fr, 5 ml Balloon
(148-4852)	TruCath® Silicone-Elastomer Coated Latex Foley Catheter, 2-Way, 16 Fr, 5 ml Balloon
(148-6910)	TruCath® Silicone-Elastomer Coated Latex Foley Catheter, 2-Way, 18 Fr, 5 ml Balloon
(148-8204)	TruCath® Duo 100% Silicone Dual-Balloon Catheter, 14 Fr
(148-8446)	TruCath® Duo 100% Silicone Dual-Balloon Catheter, 16 Fr

Additional Foley catheter configurations available.

For more information, contact your local Henry Schein Medical Representative or visit [henryschein.com/medical](https://www.henryschein.com/medical).

HR HealthCare, the HR HealthCare logo, and TruCath are trademarks or registered trademarks of HR HealthCare. © 2025 HR HealthCare. All rights reserved. LBL-0144 Rev. 0 (12/2025)



SINGLE-USE, ENDLESS OPTIONS

OVER 600 PATTERNS AVAILABLE!



Scan to download our full catalog!

Sklar Econo™ Sterile Surgical Instruments

Please contact your local Henry Schein representative for samples and more product information.

INTEGRA

LIMIT UNCERTAINTY

Surgical Blades for Every Procedure

- Razor-sharp cutting edge
- Individually packaged



Disposable Surgeon Blades – Sterile

HS Item #	Description	Size	Qty/Box
953-3949	Silicone Coated Surgical Blade	#10	100/Bx
953-3950	Silicone Coated Surgical Blade	#11	100/Bx
953-3951	Silicone Coated Surgical Blade	#15	100/Bx
953-1378	Carbon Steel Surgeon Blades	#10	100/Bx
953-1774	Carbon Steel Surgeon Blades	#11	100/Bx
953-1981	Carbon Steel Surgeon Blades	#12	100/Bx
953-2041	Carbon Steel Surgeon Blades	#15	100/Bx
953-1088	Carbon Steel Surgeon Blades	#15c	100/Bx
953-5727	Carbon Steel Surgeon Blades	#20	100/Bx
953-6694	Stainless Steel Surgeon Blades	#10	100/Bx
953-7469	Stainless Steel Surgeon Blades	#11	100/Bx
95-7546	Stainless Steel Surgeon Blades	#12	100/Bx
953-7616	Stainless Steel Surgeon Blades	#15	100/Bx
953-7101	Stainless Steel Surgeon Blades	#15C	100/Bx
953-8126	Stainless Steel Surgeon Blades	#22	100/Bx

Disposable Scalpels & BiopBlade – Sterile

HS Item #	Description	Size	Qty/Box
953-2331	Disposable Scalpels	#10	10
953-1971	Disposable Scalpels	#11	10
953-2844	Disposable Scalpels	#12	10
953-3446	Disposable Scalpels	#15	10
953-1913	Disposable Scalpels	#20	10
953-3959	Disposable Scalpels	#21	10
953-9645	Disposable Scalpels	#22	10
953-6206	Disposable Scalpels	#23	10
953-1831	Disposable Safety Scalpels	#10	10
953-6730	Disposable Safety Scalpels	#11	10
953-7177	Disposable Safety Scalpels	#12	10
953-2142	Disposable Safety Scalpels	#15	10
953-2785	Disposable Safety Scalpels	#15C	10
953-8348	Disposable Safety Scalpels	#22	10
953-7731	BiopBlade		50



Jarit® ASC Solutions Surgical Instrumentation You Can Trust

More than 8,000 instrument patterns, including a robust selection of specialty instrumentation for ASC needs

- Premium grade German crafted
- Fast replenishment and high fill rates
- Available through existing GPO contracts



Scissors

HS Item #	Description	Size	Qty/Box
150-2228	Mayo Dissecting Scissor, Curved Beveled Blade	6¾"	1
150-2229	Mayo Dissecting Scissor, Straight Beveled Blade	6¾"	1
145-6508	Metzenbaum Scissor, Curved	7"	1

Needle Holders

HS Item #	Description	Size	Qty/Box
150-2230	Crile Wood Needle Holder	6"	1
150-2232	Mayo-Hegar Needle Holder	8"	1
150-2233	Mayo-Hegar Needle Holder	6¼"	1

Towel Clamp

HS Item #	Description	Size	Qty/Box
137-4115	Backhaus Towel Clamp	5¼"	1

Knife Handles

HS Item #	Description	Size	Qty/Box
136-4486	Knife Handle	#3	1
133-7343	Knife Handle	#4	1
145-6488	Knife Handle	#7	1
139-1875	Knife Handle	#3L	1

Forceps

HS Item #	Description	Size	Qty/Box
285-0352	Adson Dressing Forcep	4¾"	1
140-3714	Dressing Forcep	5½"	1
145-6686	Dressing Forcep	10"	1
285-0353	Adson Tissue Forcep, 1 x 2 Teeth	4¾"	1
136-5464	Allis Tissue Forcep, 4 x 5 Teeth	6"	1
138-6350	Allis Tissue Forcep, 5 x 6 Teeth	7¼"	1
145-6587	Babcock Tissue Forcep, Straight	6"	1
145-6496	Babcock Tissue Forcep, Straight	8"	1
139-6401	Tissue Forcep, 1 x 2 Teeth	5½"	1
150-2224	Tissue Forcep 1 x 2 Teeth	10"	1
145-6526	Debakay Atraumatic Tissue Forcep, 2.5mm Tips	9½"	1
133-7355	Ferris-Smith Tissue Forcep, 1 x 2 Teeth	7"	1
150-2223	Gemini Artery Clamp, Curved	8¼"	1
145-6465	Halstead Mosquito Forcep, Curved	5"	1
133-7346	Kelly Forcep, Curved	5½"	1
150-2225	Kelly Forcep, Straight	5½"	1
150-2226	Rochester-Carmalt Forcep, Curved	6½"	1
133-7347	Rochester-Ochsner Forcep, Straight, 1 x 2 Teeth	6¼"	1
138-6356	Rochester-Ochsner Forcep, Straight, 1 x 2 Teeth	8"	1
285-0371	Schmidt Hemostatic Forcep, Curved	7½"	1
129-7458	Foerster Sponge Forcep, Straight	9½"	1



For the full line of
Integra Miltex products, visit
www.henryschein.com/medical

HENRY SCHEIN®
BRAND PRODUCTS

DISPOSABLE BIOPSY PUNCH

Precision in Every Punch



Rigid, closed-end translucent guard
protects blade from damage

Seamless, stainless steel finely sharpened blade ensures safe, clean, and easy sample extraction



Ergonomic, ribbed grip handle provides comfort, control and a clear view of the biopsy site

Handles are color-coded
by size for clear identification

Blade size clearly marked
on handle for quick recognition

○ 1.5mm, White.....(570-3435)	● 3mm, Grey.....(570-3438)	● 5mm, Pink.....(570-3441)
● 2mm, Ivory.....(570-3436)	● 3.5mm, Lemon Yellow....(570-3439)	● 6mm, Light Sea Green... (570-3442)
● 2.5mm, Green.....(570-3437)	● 4mm, Light Blue.....(570-3440)	● 8mm, Purple.....(570-3443)

Order today at HENRYSCHEIN.COM/MEDICAL



Surgical Blades

HENRY SCHEIN®

Sterile

- Individual foil packs
- Computer-controlled grinding ensures a consistently sharp edge

Stainless Steel.....100/box

Specify:

#10 (112-6199)	#20 (112-6184)
#11 (112-6190)	#21 (112-6185)
#12 (112-6177)	#22 (112-6195)
#15 (112-6178)	#23 (112-6200)
#15S (112-7294)	#25 (112-6201)

Carbon Steel.....10/pkg

Specify:

#10 (112-6183)	#15S (112-7295)
#11 (112-6175)	#21 (112-6198)
#12 (112-6176)	#22 (112-6197)
#15 (112-6181)	

Carbon Steel, #12B

(112-6180).....10/pkg

#15C

(112-6186).....10/pkg



Disposable Scalpels

HENRY SCHEIN®

- Stainless steel
- Plastic handle with millimeter scale and durable stainless steel blade
- Strong: reduced breakage
- Sharp: clean, true, consistent cut
- Packaging: superior protection and handling
- Individually peel-backed for convenience
- Testing: consistent high quality; ensured quality to international standards

Scalpels.....10/box

Specify:

#10 (112-6196)	#15 (112-6187)
#11 (112-6189)	#15S (112-7296)
#12 (112-6194)	#20 (112-6193)

#12B

(112-6188).....10/box



MAXIMA®

Safety Scalpels

HENRY SCHEIN®

- Sterile
- Single use
- Retractable
- Easily extend or retract the blades while avoiding contact with the sharp blade edges
- Color coded for easy identification of blade size

Safety Scalpels.....10/box

Specify:

Size 10 (570-0243)	Size 12 (570-0244)
Size 11 (570-0245)	Size 15 (570-0246)



Henry Schein Brand Surgical Instruments

Since 1932, Henry Schein has offered an extensive range of products and solutions to optimize your practice and deliver the best care for your patients. Our surgical instruments portfolio is **high quality** and **competitively priced**, available in two grades, and tailored for a broad spectrum of specialties and procedures.

To help you maintain proper infection prevention protocols, we also provide cleaning and care solutions including instrument sterilization and reprocessing. Our ProRepair service includes manufacturer-trained technicians and highly skilled customer service representatives committed to providing quality repair services with a stress-free experience, with complete warranties backed by Henry Schein.

Our available surgical instrument grades include:



- **Maxima® by Henry Schein** – Premium surgical-grade instruments expertly crafted from German Stainless Steel* to conform to rigid specifications and standards that ensure our instruments consistently meet and exceed the high quality and performance expectations of our valued customers. Designed for your general practice needs as well as your more specialized surgical needs, such as **ambulatory surgery centers, plastic surgery, and dermatology practices**.



- **Henry Schein Brand** – High-quality stainless steel instruments, designed for your medical practice/office use. Manufactured in ISO 13485 certified factories, unconditionally guaranteed against defects in quality or craftsmanship.

For our full portfolio of surgical instruments go to www.henryschein.com/medical.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

SCISSORS



Lister Bandage Scissors

HENRY SCHEIN®

Standard Stainless Steel

3½"	(100-9243)	ea
5½"	(100-1087)	ea
7¼"	(100-3334)	ea
8½"	(953-2519)	ea



MAXIMA® Lister Bandage Scissors 8"

HENRY SCHEIN®

German Stainless Steel,
Extra Fine, One Large Finger Ring

(104-5841)	ea
------------	-------	----



MAXIMA® Knowles Bandage Scissors 5½"

HENRY SCHEIN®

German Stainless Steel

Angled

(102-7708)	ea
------------	-------	----

Straight

(104-5504)	ea
------------	-------	----



Universal Scissors 7½"

HENRY SCHEIN®

Standard Stainless Steel

Unusual blade design with 1 serrated edge. Cuts tin, light metal, wire, bandages, and more.

(100-7898)	ea
------------	-------	----



MAXIMA® Bandage and Utility Scissors

HENRY SCHEIN®

German Stainless Steel

6½"	(104-9212)	ea
-----	------------	-------	----

7"	(104-9280)	ea
----	------------	-------	----



MAXIMA® Hercules Heavy-Duty Bandage and Plaster Shears 7½"

HENRY SCHEIN®

German Stainless Steel

Serrated Blade

(104-7448)	ea
------------	-------	----



Wire-Cutting Scissors 4¾"

HENRY SCHEIN®

Standard Stainless Steel

Angled to Side, 1 Serrated Blade

(100-2591)	ea
------------	-------	----



MAXIMA® Lister Bandage Scissors

HENRY SCHEIN®

German Stainless Steel, Extra-Fine

3½"	(104-9717)	ea
-----	------------	-------	----

4½"	(104-8240)	ea
-----	------------	-------	----

5½"	(100-6589)	ea
-----	------------	-------	----



MAXIMA® Esmarch Heavy-Duty Bandage and Cast Shears 8"

HENRY SCHEIN®

German Stainless Steel

(104-7447)	ea
------------	-------	----

SURGICAL INSTRUMENTS

SCISSORS



Gauze Shears 8"

HENRY SCHEIN®

Standard Stainless Steel

(100-4963) ea



**MAXIMA®
Long Stitch Scissors 6¼"**

HENRY SCHEIN®

German Stainless Steel, Extra Long

(100-0109) ea



**MAXIMA®
Littauer Stitch Scissors 5½"**

HENRY SCHEIN®

German Stainless Steel, Curved

(100-2354) ea



**MAXIMA®
Northbent Stitch Scissors 5½"**

HENRY SCHEIN®

German Stainless Steel, Curved

(100-0353) ea



**MAXIMA®
Littauer Jr. Stitch Scissors 4½"**

HENRY SCHEIN®

German Stainless Steel, Slender Pattern, Straight

(100-1420) ea



Littauer Stitch Scissors 5½"

HENRY SCHEIN®

Standard Stainless Steel

(100-5115) ea



**MAXIMA®
Shortbent Stitch Scissors 3½"**

HENRY SCHEIN®

German Stainless Steel, Curved, Delicate

(104-7439) ea



Littauer Jr. Stitch Scissors 4½"

HENRY SCHEIN®

Standard Stainless Steel, Slender Pattern

(100-6306) ea



**MAXIMA®
O'Brien Stitch Scissors 3¾"**

HENRY SCHEIN®

German Stainless Steel, Angled Blade

(104-8241) ea

SCISSORS



Spencer Stitch Scissors 3 1/2"

HENRY SCHEIN®

Standard Stainless Steel

(100-9827) ea



MAXIMA®

Spencer Stitch Scissors 3 1/2"

HENRY SCHEIN®

German Stainless Steel

(100-6024) ea



MAXIMA®

Stevens Tenotomy Scissors

HENRY SCHEIN®

German Stainless Steel

Curved.....ea

Specify:

4 1/2", Blunt Points (104-7440)

4 1/2", Blunt Points (104-7441)

4 1/2", Sharp Points (104-9178)

Straight.....ea

Specify:

4 1/2", Blunt Points (104-7443)

4 1/2", Sharp Points (104-7445)

4 1/2", Sharp Points (104-9639)



MAXIMA®

Gradle Eye Suture Scissors 3 3/4"

HENRY SCHEIN®

German Stainless Steel, Slightly Curved,
Sharp Points

(101-7676) ea



Stevens Tenotomy Scissors

HENRY SCHEIN®

Standard Steel, Blunt Points

Straight, 4 1/2" ea

(100-0973) ea

Curved, 4 1/8" ea

(100-5500) ea



MAXIMA®

Strabismus Scissors 4"

HENRY SCHEIN®

German Stainless Steel

Curved (104-8239) ea

Straight (104-8237) ea

SCISSORS



MAXIMA®
Mayo Dissecting Scissors 5½"
HENRY SCHEIN®
 German Stainless Steel
 Curved
 (100-3105) ea
 Straight
 (100-4018) ea



MAXIMA®
Metzenbaum Scissors 5"
HENRY SCHEIN®
 German Stainless Steel, Curved, Delicate
 (104-7306) ea



MAXIMA®
Operating Room Scissors 4¾"
HENRY SCHEIN®
 German Stainless Steel, Straight
 Blunt Points
 (104-8232) ea
 Sharp Points
 (102-3440) ea
 Sharp/Blunt Points
 (104-8231) ea



MAXIMA®
Mayo Dissecting Scissors 6¾"
HENRY SCHEIN®
 German Stainless Steel
 Curved
 (100-1627) ea
 Straight
 (100-4833) ea



MAXIMA®
Metzenbaum Scissors 8"
HENRY SCHEIN®
 German Stainless Steel, Curved
 (104-8236) ea



MAXIMA®
Operating Room Scissors 6"
HENRY SCHEIN®
 German Stainless Steel, Straight, Sharp/Blunt
 Points
 (104-9283) ea



MAXIMA®
Metzenbaum Scissors 7"
HENRY SCHEIN®
 German Stainless Steel
 Curved
 (100-5142) ea
 Straight
 (101-6892) ea

MAXIMA®
Operating Room Scissors 5½"
HENRY SCHEIN®
 German Stainless Steel
 Curved, Blunt Points
 (104-9957) ea
 Curved, Sharp Points
 (100-3874) ea
 Straight, Blunt Points
 (100-4352) ea
 Straight, Sharp Points
 (100-4549) ea
 Straight, Sharp/Blunt Points
 (100-7561) ea



SURGICAL INSTRUMENTS

SCISSORS

**MAXIMA®
Iris Scissors**
HENRY SCHEIN®

German Stainless Steel

3½", Curved, Delicate Pattern

(102-3439) ea

4", Curved

(101-8111) ea

3½", Straight, Delicate Pattern

(102-5248) ea

4", Straight

(101-6893) ea

4½", Straight

(100-8118) ea



**MAXIMA®
Kelly Angular Scissors 6¼"**
HENRY SCHEIN®

German Stainless Steel, Angled, Sharp Points

Straight, 1 Serrated Blade

(100-0112) ea

Curved, 1 Serrated Blade

(100-3164) ea

Curved, Smooth Blades

(100-7266) ea



Iris Scissors 4½"
HENRY SCHEIN®

Standard Stainless Steel

Curved

(100-8820) ea

Straight

(100-2767) ea

Kelly Scissors 6¼"
HENRY SCHEIN®

Standard Stainless Steel

Angled

(101-4406) ea

Curved, Sharp Points, Smooth Blades

(101-3453) ea

Curved, Serrated Blade

(101-7287) ea



**MAXIMA®
Wagner Plastic Surgery Scissors 4¾"**
HENRY SCHEIN®

German Stainless Steel, Curved, Sharp Points,
1 Serrated Blade

(100-4982) ea

**MAXIMA®
Nail Scissors 4"**
HENRY SCHEIN®

German Stainless Steel, Curved

(104-7810) ea



FORCEPS



MAXIMA®
Adson-Brown Dressing Forceps 4^{3/4}"
HENRY SCHEIN®
 German Stainless Steel, 7 x 7 Side-Grasping Teeth
 (101-6911) ea



MAXIMA®
Micro Adson Dressing Forceps 4^{3/4}"
HENRY SCHEIN®
 German Stainless Steel, Serrated, 0.5-mm Wide Tips
 (102-8861) ea



MAXIMA®
Adson Tissue & Suture Forceps 4^{3/4}"
HENRY SCHEIN®
 German Stainless Steel
 1 x 2 Teeth, Delicate Pattern
 (101-4803) ea
 Cross-Serrated Tips, Delicate Pattern
 (102-3435) ea
 1 x 2 Teeth, with Tying Platform
 (101-6905) ea



Adson-Brown Dressing Forceps 4^{3/4}"
HENRY SCHEIN®
 Standard Stainless Steel, 7 x 7 Side-Grasping Teeth
 (100-5162) ea



MAXIMA®
Dressing Forceps 5^{3/4}"
HENRY SCHEIN®
 German Stainless Steel, Serrated
 (101-1272) ea



MAXIMA®
Hudson Tissue Forceps 4^{3/4}"
HENRY SCHEIN®
 German Stainless Steel, 1 x 2 Teeth
 (104-5547) ea



Adson Dressing Forceps 4^{3/4}"
HENRY SCHEIN®
 Standard Stainless Steel, Delicate Pattern
 (100-5659) ea



MAXIMA®
Adson Tissue & Suture Forceps 4^{3/4}"
HENRY SCHEIN®
 German Stainless Steel, Serrated Tips, Delicate Pattern
 (101-6279) ea



MAXIMA®
Tissue Forceps
HENRY SCHEIN®
 German Stainless Steel, 1 x 2 Teeth
 5^{1/2}", Standard Pattern
 (101-5132) ea
 6^{1/2}", Serrated Handle
 (104-7960) ea

FORCEPS



Adson Tissue Forceps 4 3/4"

HENRY SCHEIN®

Standard Stainless Steel, 1 x 2 Teeth

(100-5891) ea



MAXIMA® Baby Allis Tissue Forceps 5 1/2"

HENRY SCHEIN®

German Stainless Steel, Extra Delicate Jaws,
4 x 5 Teeth

(104-7415) ea



MAXIMA® Swiss Cilia and Suture/Splinter Forceps 3 3/4"

HENRY SCHEIN®

German Stainless Steel, Precision Fitted,
Diagonal Jaws

(102-3436) ea



MAXIMA® Allis Tissue Forceps

HENRY SCHEIN®

German Stainless Steel

6", 4 x 5 Teeth

(101-1848) ea

7", 5 x 6 Teeth

(100-0228) ea



MAXIMA® Tissue Forceps

HENRY SCHEIN®

German Stainless Steel, Standard Pattern,
1 x 2 Teeth

4 1/2"

(104-6618) ea

5"

(104-8243) ea



MAXIMA® Plain Splinter Forceps

HENRY SCHEIN®

German Stainless Steel

3 1/2"

(101-7834) ea

4 1/2"

(101-6895) ea



Allis Tissue Forceps

HENRY SCHEIN®

Standard Stainless Steel

6", 4 x 5 Teeth

(100-6914) ea

6", 5 x 6 Teeth

(100-6095) ea

7 1/4", 5 x 6 Teeth

(100-9430) ea



Hartman-Noyes Alligator Forceps

HENRY SCHEIN®

Standard Steel, Serrated Jaws

3 1/2"

(100-7328) ea

5 1/2"

(100-0886) ea

6 1/2"

(100-6937) ea

SURGICAL INSTRUMENTS

FORCEPS

MAXIMA®
Kelly Forceps 5½"
HENRY SCHEIN®

German Stainless Steel


- Curved (100-8381) ea
- Straight (100-3357) ea



MAXIMA®
Halstead Mosquito Forceps 5"
HENRY SCHEIN®

German Stainless Steel

- Curved (100-2376) ea
- Curved, Extra Delicate (104-5506) ea
- Straight (100-0846) ea
- Straight, Extra Delicate (104-5500) ea



MAXIMA®
Crile Forceps 6¼"
HENRY SCHEIN®

German Stainless Steel, Curved

- (104-7370) ea



MAXIMA®
Hartman Mosquito Forceps 3½"
HENRY SCHEIN®

German Stainless Steel, Curved

- (102-4346) ea



MAXIMA®
Rochester-Pean Hemostat 6¼"
HENRY SCHEIN®

German Stainless Steel, Curved

- (100-6870) ea

MAXIMA®
Hartman Mosquito Forceps 3½"
HENRY SCHEIN®

German Stainless Steel, Straight

- (102-7483) ea



MAXIMA®
Utility Sterilizer Forceps 8"
HENRY SCHEIN®

German Stainless Steel, Straight

- (102-4137) ea



SURGICAL INSTRUMENTS

FORCEPS & TOWEL CLAMPS




**MAXIMA®
Foerster Sponge Forceps 9¾"**

HENRY SCHEIN®

German Stainless Steel, Serrated Jaws

Curved
(104-8859) ea

Straight
(101-1737) ea



**MAXIMA®
Bozeman Sponge Forceps 10¼"**

HENRY SCHEIN®

German Stainless Steel, Curved

(104-7369) ea



**MAXIMA®
Fletcher Sponge Forceps 9½"**

HENRY SCHEIN®

German Stainless Steel, Straight

(104-8249) ea



**MAXIMA®
Backhaus Towel Clamp**

HENRY SCHEIN®

German Stainless Steel

3½"
(101-4934) ea

5¼"
(101-3568) ea



**MAXIMA®
Lorna Nonperforating Towel Clamp 3½"**

HENRY SCHEIN®

German Stainless Steel

(104-8244) ea



**MAXIMA®
Lorna Nonperforating Towel Clamp 5¼"**

HENRY SCHEIN®

German Stainless Steel

(104-8245) ea



Jones Towel Clamp 3½"

HENRY SCHEIN®

Standard Stainless Steel

(100-8003) ea

NEEDLE HOLDERS



Castroviejo Needle Holder 5 1/2"

HENRY SCHEIN®

Standard Stainless Steel, Straight, with Lock
(953-9789) ea



MAXIMA® Derf Needle Holder 4 3/4"

HENRY SCHEIN®

German Stainless Steel and Tungsten Carbide
(104-8597) ea
German Stainless Steel
(104-7419) ea



MAXIMA® Mathieu Needle Holder With Spring Handle 5 1/2"

HENRY SCHEIN®

German Stainless Steel
(104-7424) ea



MAXIMA® Ryder Needle Holder

HENRY SCHEIN®

German Stainless Steel, Tungsten Carbide,
Serrated Jaw
5"
(104-8251) ea
10"
(104-9323) ea



MAXIMA® Webster Needle Holder 5"

HENRY SCHEIN®

German Stainless Steel, Delicate Pattern,
Smooth Jaws
(101-9003) ea
Tungsten Carbide
(104-6385) ea



Webster Needle Holder 5"

HENRY SCHEIN®

Standard Stainless Steel, Delicate Pattern,
Smooth Jaws
(100-6371) ea



MAXIMA® Olsen-Hegar Needle Holder/Scissors Combination 5 3/4"

HENRY SCHEIN®

German Stainless Steel, Tungsten Carbide,
Serrated Jaws
(100-1432) ea




MAXIMA® Crile-Wood Needle Holder

HENRY SCHEIN®

German Stainless Steel
6", Serrated Jaws
(101-7128) ea
6", Tungsten Carbide
(100-7943) ea
6", Tungsten Carbide, Delicate Pattern, Smooth Jaw
(104-5541) ea
7", Tungsten Carbide
(104-8255) ea



NEEDLE HOLDERS & RETRACTORS



Mayo-Hegar Needle Holder

HENRY SCHEIN®

Standard Stainless Steel

- 5" ea
- (100-0959)
- 6" ea
- (100-2570)
- 8" ea
- (100-6577)



MAXIMA® Baumgartner Needle Holder 5½"

HENRY SCHEIN®

German Stainless Steel

- (101-7415) ea
- Tungsten Carbide, Serrated Jaw
- (100-2001) ea



MAXIMA® Halsey Needle Holder 5¼"

HENRY SCHEIN®

German Stainless Steel, Serrated Jaws

- (104-7420) ea



MAXIMA® Mayo-Hegar Needle Holder

HENRY SCHEIN®

German Stainless Steel and Tungsten Carbide

- 5¾" ea
- (104-5517)
- 6" ea
- (100-6116)
- 7¼" ea
- (104-8253)



MAXIMA® Mayo-Hegar Needle Holder

HENRY SCHEIN®

German Stainless Steel

- 5¾", Delicate ea
- (100-8670)
- 6" ea
- (104-5531)
- 7¼" ea
- (104-7548)



MAXIMA® Halsey Needle Holder

HENRY SCHEIN®

Smooth Jaws

- 5", German Stainless Steel and Tungsten Carbide ea
- (104-5514)
- 5¾", German Stainless Steel ea
- (102-7437)



MAXIMA® Weitlaner Retractor 4"

HENRY SCHEIN®

German Stainless Steel, Blunt Points, 2 x 3 Teeth

- (104-9635) ea



MAXIMA® Weitlaner Retractor 5½"

HENRY SCHEIN®

German Stainless Steel, Sharp Points, 3 x 4 Teeth

- (104-7436) ea

SURGICAL INSTRUMENTS

WOMEN'S HEALTH INSTRUMENTS

Disposable Vaginal Specula



Made of hard, clear, nontoxic, high-strength plastic resin.

- Allows examination of vaginal mucosa at a 360° angle without rotation
- Provides a wide range of adjustability (angular and elevation adjustments)
- One-handed manipulation for ease of operation
- Self-holding in various opening positions
- Decreases the risk of contamination



Small (102-5396)	10/pkg
Medium (102-5397)	10/pkg
Large (102-9422)	10/pkg

MAXIMA® Graves Vaginal Specula



German Stainless Steel

Small (104-5540)	ea
Medium (104-5522)	ea
Large (104-5535)	ea




MAXIMA® Pederson Vaginal Speculum



German Stainless Steel, Narrow Blades

Small (104-8971)	ea
Medium (104-5601)	ea
Large (104-5543)	ea
X-Narrow (104-7925)	ea



MAXIMA® Kevorkian-Younger Uterine Biopsy Forceps 9¾"



German Stainless Steel

(104-9098)	ea
------------------	----



PODIATRY INSTRUMENTS

MAXIMA® Nail Splitters



German Stainless Steel, Straight Jaws,
Double Spring

4" (104-5512)	ea
5", Delicate (104-7539)	ea
5" (104-5526)	ea



MAXIMA® Nail Nipper 4¾"



German Stainless Steel, Concave Jaws,
Spring Barrel, Chrome
Smooth barrel-spring action and precise sharp
blades to cut the thickest toenails. Heavy duty.

(102-6341)	ea
------------------	----



Nail Splitters 5"



Standard Stainless Steel, English Anvil, Straight
Jaws, Double Spring

(104-8507)	ea
------------------	----





The Bowie-Dick revolution begins

3M™ Attest™ eBowie-Dick Test System

Automated, digitized sterilization testing is here. Save time, enhance workflow, boost confidence, and get the reassurance that your sterile processing team is operating at its best. The same testing method you've always trusted, now digital.

- **Definitive:** Get pass/fail results without guesswork or subjective interpretation
- **Efficient:** Removes the burden of paperwork and eliminates manual recordkeeping
- **Innovative:** The FDA created a new category for it

Simplify. Standardize. Streamline.

Contact your Henry Schein representative or visit go.solventum.com/ebd to learn more.

3M Health Care is now Solventum

Product size not to scale.

© Solventum 2025. Solventum, the S logo and Attest are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. Other trademarks are the property of their respective owners.



- (148-5846) 3M™ Attest™ eBowie-Dick Auto-reader
- (148-5248) 3M™ Attest™ eBowie-Dick Test Card
- (148-5474) 3M™ Attest™ eBowie-Dick Card Holder

INTEGRA®

Protect Your Instrument Investment.

Integra Sterilization Containers are engineered to safeguard valuable instrument sets while delivering reliable, repeatable sterilization performance.



- 1** Gaskets within retention plates and container perimeter ensure a secure seal of the container, lid and filter.
- 2** Lid and Bottom manufactured from anodized aluminum alloy. Easy to maintain, corrosion resistant, lightweight yet durable, and conducts heat for material drying period.
- 3** Corner design allows for same size container stacking in storage.
- 4** 2 Indicator Card Windows (1 front & 1 back). Indicator card has visible bar for clear indication that steam or gas sterilization has been achieved.
- 5** Tamper Evident Lock: helps to identify a sterility breach.
- 6** Fold-down handles extend and stop at 90° to aid lifting.
- 7** Low profile latch easily and securely closes lid over container bottom.
- 8** 4 ID Tag Windows (2 front & 2 back) ID tags available in 6 colors and can be engraved with up to 18 characters per line.
- 9** Round disposable filters with centering hole to assist in proper alignment.
- 10** Snap on filter plate with a simple locking mechanism that helps provide easy change of disposable filters.

One-Half Size (11¼" x 11¾" x 4⅛")

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4⅛"		953-4325		953-4337
5¼"		953-4326	953-4338	953-4343
6⅞"		953-4330	953-4346	953-4347
7¾"	953-4332	953-4334	953-4350	953-4357
10¼"	953-4335	953-4336	953-4361	953-4362

Three-Quarter Size (18¼" x 11¼" x 4¼")

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4¼"	121-8501	953-4364	953-4369	953-4370
5¼"	953-4365	953-4366	121-9669	953-4371
6"		953-4368	121-9673	953-4372

Full Size (23¼" x 11¼" x 4⅛")

4⅛"	122-9041	953-4376	953-4386	953-4388
5¼"	953-4377	953-4378	953-4389	953-4390
6"		953-4379	121-9668	953-4391
7¾"	953-4380	953-4382	953-4395	953-4396
10¼"	953-4384	953-4385	953-4397	953-4398

Wide-Body (23¼" x 15¾" x 4¾")

4¾"	953-4402	953-4404	953-4408	953-4409
7"	953-4406	953-4407	953-4410	953-4411

For the full line of Integra Miltex products, visit www.henryschein.com/medical.



MAXIZYME® Tabs HENRY SCHEIN®

Enzyme cleaning tablets for ultrasonic cleaning, evacuation system cleaning, and instrument presoak. 2 tablets dissolve easily in 1 gallon of warm tap water. Easy to use and easy to store: requires less storage space.

- Contain protease enzymes that break down hard-to-clean deposits such as blood, bone, and food particles that chemical detergents leave behind
- A single cleaner with 2 uses: ultrasonic and evacuation systems
- Gentle on instruments: nontoxic, neutral pH, biodegradable, and noncorrosive formula that contains no chlorine or harsh alkalines to damage expensive stainless steel instruments
- Refreshing minty scent

(900-0475).....64/box



MAXIZYME® Multi-Enzymatic Detergent HENRY SCHEIN®

Improved multienzymatic detergent that cleans fast in manual, ultrasonic, and automatic cleaning systems. Low-foaming formula rapidly dissolves and carries away organic soils, including blood, protein, tissue, and more. Its corrosion-inhibiting, pH-neutral formula is safe for instruments and equipment, extending the useful life of instruments. Concentrated and economical, each 1-gallon bottle yields 128 gallons of solution. Clean, fresh scent.

1-gal Bottle
(570-3275).....ea



MAXITABS® HENRY SCHEIN®

General Purpose Cleaner
Provide safe and effective cleaning properties in multiple cleaning cycles throughout the day. Does not leave a residue on instruments or inside tank. One tablet makes 1 gallon of solution; 1 box makes 32 gallons.

(900-4088).....32/box



Ultrasonic Cleaning Solutions HENRY SCHEIN®

1-gal Bottles

General Purpose Cleaner

- Nonammoniated concentrate
- Dilute 1:10 for use
- Effectively cleans items such as dried blood, bone, tissue, oxide compounds, and more
- Removes ions from instrument's surface to enhance effectiveness of ultrasonic cleaning unit

(100-3836).....ea

Tartar and Stain Remover ♦

- Used in a beaker in an ultrasonic unit to remove tartar and stains from dentures
- Ready to use, so it's easy and convenient
- Potent stain remover eliminates stains on dentures, crowns, and bridges

(100-4263).....ea

Temporary Cement and Stone Remover ♦

- Ready-to-use formula
- Strong, fast-acting formula removes cement and stone quickly and completely

(100-5168).....ea

Plaster and Stone Remover ♦

(101-3277).....ea

- ♦ Hazardous material



Super Concentrated Ultrasonic Cleaning Solutions HENRY SCHEIN®

- Fast-acting solutions that clean quickly and effectively
- Nontoxic and noncorrosive
- In small, easy-to-store, 16-oz bottles, which are less expensive to ship, handle, and stock
- Metered-dose dispenser top ensures exact dose
- To use: loosen cap ½ turn, squeeze to 1-oz mark, and remove cap and pour into glass beaker
- Available in 3 formulas

16-oz Metered Dose Bottle.....ea

Specify:

General Purpose Cleaner
(Makes 16 gal) Used for very fast cleaning of blood, mucus, and tissue.

(101-4412)

Tartar and Stain Remover (Makes 96 oz) Its citric acid formula is noncorrosive and cleans quickly.

(101-0827)

Enzymatic Cleaning Solution (Makes 16 gal) Effective against blood, tissue, and mucus found on instruments.

(102-6472)



Surgical Milk HENRY SCHEIN®

- Super-concentrated
- Stops corrosion and maintains free-working instruments
- Nontoxic and safe to use
- Permits autoclaving of sharps and carbon steel
- 1 gallon of Surgical Milk diluted with 5 gallons of tap water makes 6 gallons of immersion bath milk

64-fl-oz Bottle

(102-5539).....ea



MAXIWIPE™ Germicidal Wipes



Disinfectant wipe that kills the COVID-19 virus in 1 minute on hard, nonporous surfaces. Provides fast disinfection. Compatible with a broad range of surfaces, including stainless steel, plastic, Formica®, and glass.

- Cleans and disinfects common hard, nonporous surfaces
- Bactericidal, tuberculocidal, and virucidal
- Effective against 16 pathogens, including COVID-1, MRSA, H1N1, TB, HBV, MRSA, VRE, and RSV
- Contains 55% alcohol
- Meets CDC and OSHA regulations

Large—6" x 6¾" 160/canister
(106-6794)

X-Large—7½" x 15" 75/canister
(570-3361)

Packaging label may vary.



MaxiSpray® Plus



Surface Disinfectant Cleaner

- Hospital-level, general-purpose disinfectant/cleaner proven to kill TB in 5 minutes, and HBV, HIV, and HCV in 3 minutes
- EPA-registered as a broad-spectrum disinfectant for both surface and immersion use
- Bactericidal, virucidal, tuberculocidal, and fungicidal

24-oz Trigger-Spray Bottle ea
(101-4114)

1-gal Refill ea
(101-0285)



MAXIWIPE™ Alcohol-Free Germicidal Disposable Wipes



Formulated with care: Alcohol- and fragrance-free wipes for staff and patients with respiratory sensitivities. Efficient, one-step cleaner and disinfectant. Effective bactericidal, tuberculocidal, and virucidal, including SARS-CoV-2 (the virus that causes COVID-19), in 3 minutes. EPA Registration No. 9480-9, now qualified for List Q*. Intermediate level.

Large, 6" x 6¾" 160/cont
(570-3056)

X-Large, 7½" x 15" 75/cont
(570-3057)

*MAXIWIPE™ Alcohol-Free has demonstrated effectiveness against emerging viral pathogens, such as the monkeypox virus, on hard, nonporous surfaces.



BeeSure® Naturals Face Masks Level 3 ECOBEE

Designed with unique, natural elements.

- More than 60% earth-friendly cellulose components
- Quad-Fold® design offers 15% more breathing volume
- 100% non-latex and fiberglass-free filter layers
- Recycled packaging

Earloop Masks 50/box
 Specify:
 Canyon Red **(142-4176)** Cloud White **(142-4178)**
 Glacier Blue **(142-4177)** Forest Green **(142-4180)**



BeeSure® Floral Earloop Masks Level 3 ECOBEE

- Packaging is recyclable
- Quad-Fold® design that offers 15% more breathing volume
- RoundEdge® corners and cushioned earloops for comfort
- Non-latex CoolSilk®, dye-free inner layer for a smooth, nonirritating feeling

Earloop Masks 50/box
 Specify:
 Hibiscus Pink **(142-4158)**
 Plumeria Blue **(142-4159)**
 Orchid Lavender **(142-4163)**
 Fern Green **(142-4165)**
 Daisy Yellow **(142-4164)**



BeeSure® Vibe Face Masks Level 3 ECOBEE

Fluid-resistant mask that features a VA DuoSeal™ nose and chin-adjustable seal, an antifog flap for clearer visibility during procedures, and a Quad-Fold® design for 15% more breathing volume. Lint-free, dye-free inner layer. Recycled packaging

- Fluid resistant
- VA DuoSeal™ nose and chin-adjustable seal
- Antifog flap for clearer visibility during procedures
- Quad-Fold® design for 15% more breathing volume
- Lint-free, dye-free inner layer
- Recycled packaging

Earloop Masks 50/box
 Specify:
 Dazzling Green **(353-0097)**
 Glamorous Pink **(353-0099)**
 Gray **(138-4278)**
 Luminous Blue **(353-0096)**
 White **(138-4277)**

INTRODUCING

CaviCide™ HP

1-Min Hydrogen Peroxide

The power of CaviWipes™ HP is now available in a ready-to-use liquid.

-  1-Minute universal contact time for all listed pathogens
-  Effective disinfection against 70 pathogens including MDRO and **Candida auris**
-  Hydrogen Peroxide with no signal word – no PPE required*
-  1-Step cleaning/disinfection*



Ordering Information

Part No.	Item Description
640-0123	CaviCide HP 24 oz spray bottle
640-0124	CaviCide HP 1 gal refill bottle



Scan the QR code or visit Metrex.com to learn more.

* Follow Product Label for special precleaning and PPE instructions



OF INFECTION PREVENTION

Quality Products. Trusted Solutions. Since 1985.

Metrex™



©2025 Metrex Research, LLC. All trademarks are property of Metrex Research, LLC. All Rights Reserved. MKT-25-1403-HS Rev 0

MAXIWIPE™ ALCOHOL FREE
GERMICIDAL DISPOSABLE WIPE

FORMULATED WITH CARE

Alcohol- and fragrance-free formula for staff and patients with respiratory sensitivities.



LARGE 6" x 6¾" (15.2 x 17.1 cm).....**(570-3056)**
(160/Canister)

EXTRA LARGE 7½" x 15" (19.0 x 38.1 cm)**(570-3057)**
(75/Canister)

Order today at HENRYSCHEIN.COM/MEDICAL

**PROTECTING
WITH PURPOSE**

Clorox Healthcare shares your mission to defeat HAIs.

Studies estimate that up to 60% of Surgical Site Infections are preventable.^{1,2} SSI and other infection-causing pathogens can be introduced into the OR via staff, mobile equipment, surgical instruments, and even the patients themselves. Comprehensive OR cleaning and disinfection has been clinically proven to significantly reduce SSI-causing pathogen transmission.

Scan QR code for product information.

1. Meeks DW, Lally KP, Carrick MM, et al. Compliance with guidelines to prevent surgical site infections: as simple as 1-2-3? Am J Surg. 2011;201(1):76-83. 2. Umscheid CA, Mitchell MD, Doshi JA, Agarwal R, Williams K, Brennan PJ. Estimating the proportion of healthcare-associated infections that are reasonably preventable and the related mortality and costs. Infect Control Hosp Epidemiol. 2011;32(2):101-114.

© 2016 Clorox Professional Products Company NI-66881

Trusted & Trustworthy

It's true. The PURELL® brand is trusted and preferred by healthcare workers.

But more importantly, we work tirelessly to earn that trust.

Our hand hygiene is gentle on hands, even with frequent use. And our powerful surface products can even be used near patients and on electronics.

The Difference Is Clear.



Contact your distributor or GOJO representative today.

Distributed by:



©2023 GOJO Industries, Inc. All rights reserved. | 34993 (05/2023)



142-1402

126-6003

133-4235

137-7898

133-9255

131-1467

The Premium Infection Prevention MedSurg Package



Surface Disinfection



#1 DISINFECTANT WIPES IN HEALTHCARE

Super Sani-Cloth® Wipes Large Canister (138-3942)
 Sani-HP1™ Wipes Large Canister (602-0272)
 Sani-Cloth® Bleach Wipes Large Canister (602-0021)



Skin Prep



PrevaNtic® Swab (602-0035)
 PrevaNtic® Swabstick (602-0036)
 PrevaNtic® Maxi Swabstick (602-0037)



Nasal Decolonization



Profend® Nasal Antiseptic Kit (131-5636)



Hand Hygiene



Sani-Hands® Wipes (610-8876)

Screen Cleaning



Easy Screen® Wipes Canister (125-0248)
 Easy Screen® Wipes Packet (146-6732)

Learn more by contacting your Henry Schein Rep.



*Market leadership claims based on Q3 2025 Clarivate Data.

UPDATE 1225 PDI05259155

Find My PDI Rep



Aspen Surgical® Electrosurgery Solutions

Electrosurgery



SURGI-CENTER PRO™



OR PRO 300™



SMOKE SHARK® II
WITH BOVIE® ORCA®



Contact your Henry Schein Sales Representative at www.henryschein.com/medical to purchase yours today!

- 122-9913 | OR PRO A3350 ELECTROSURGICAL GENERATOR | 1/EA
- 124-9073 | SURGI-CENTER PRO A2350 ELECTROSURGICAL GENERATOR | 1/EA
- 119-3868 | SMOKE SHARK® II SMOKE EVACUATOR | 1/EA
- 139-6511 | BOVIE® ORCA® INTEGRATED SMOKE EVACUATION PENCIL | 10/BOX



©2024 Aspen Surgical Products, Inc., ALL RIGHTS RESERVED M-SM-ADV-1167129, Rev I



Aaron® Smoke Shark™ Smoke Evacuator ASPEN SURGICAL

Designed to be easy to use and is extremely effective at smoke and particle removal at the operative site. While designed for use with all electrocautery generators, it is also excellent for use in controlling laser smoke.

- Easy to operate
- Quiet and effective smoke evacuation
- Long-life, 18-hour filter
- Low cost per procedure
- Excellent particle capture

#SE02
(119-3868) ea



Aaron® Surgi-CenterPro Electrosurgical Generator ASPEN SURGICAL

A 200-watt, high-power operating-room generator. Unit incorporates the power and safety of a premium OR generator at a significant savings to the medical facility. System delivers 12 monopolar and bipolar RF energy modalities to provide a high level of surgical precision.

- Offers pure cut, 4 modes of blend, 3 modes of coagulation, bipolar coagulation, and Bovie FCFS™ technology
- Designed with state-of-the-art safety features
- Backed by a 4-year manufacturer's limited warranty

#A2350
(124-9073) ea

Surgical smoke is toxic. Why is it still in your OR?

Same toxicity as
6 cigarettes
per gram of tissue.¹

20 states have committed to smoke evacuation.
Is your facility ready?

Time to clear the air. Act now with CONMED.
Scan to explore our smoke free portfolio.



150+
hazardous
chemicals.²



CONMED Corporation | 11311 Concept Blvd., Largo, FL 33773 | www.CONMED.com

1 Sahaf, O., Vega-Carrascal, J., Cunningham, F., McGrath, J., Bloomfield, F., (2007). Chemical Composition of smoke produced by high-frequency electrosurgery. Irish Journal of Medical Science, DOI: 10.1007/s11845-007-0068-0 | 2 Pierce JS, Lacey SE, Lippert JF, Lopez R, Franke JE. Laser-generated air contaminants from medical laser applications: a state-of-the-science review of exposure characterization, health effects, and control. J Occup Environ Hyg. 2011;8(7):447-46 | ©2025 CONMED Corporation, 1759289460 11/25

Service & Support



Hall® Titan® Powered Instruments CONMED

Refined for peak performance, the Hall Titan Series is the next generation of powerful, precise, and reliable orthopedic power tools.

- #PRO9300B, Oscillating Saw (142-4717)ea
- #PRO9350B, Primecut+ Oscillating Saw (142-4718) ea
- #PRO9400B, Reciprocating Saw (142-4719) ea
- #PRO9450B, Sternum Saw (142-4720) ea
- #PRO9100B, Single-Trigger Modular Drill (142-4715) ea
- #PRO9200B, Dual-Trigger Modular Drill (142-4716) ea



MicroPower+™ Small Bone Power System CONMED

This system provides a high degree of control and precision.

- Encapsulation and potting protect against moisture
- Rotating sagittal saw head for improved cutting angles
- Indexing oscillating saw collet design for cutting versatility
- Extendable levers with padded tip for more comfort

- #PRO7300SE, Oscillating Saw (142-2763)ea
- #PRO7400SE, Reciprocating Saw (142-2788) ea
- #PRO7200SE, Sagittal Saw (142-2780) ea



High-Temperature (2,200°F), Disposable Cauteries, Cordless



- Lightweight
- Convenient
- Cordless

#AA01, Fine Tip Cautery (112-7177)	ea
#AA03, Loop Tip Cautery (112-7178)	ea
#AA21, Elongated Vasectomy Tip (112-7180)	ea



High-Frequency Desiccator



Designed and engineered to be the most reliable and durable desiccator available today, this unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat. Post-procedure and long-term cleanup are simplified by its nonabsorbing, nonstaining, high-impact plastic case.

- 40-watt high-frequency desiccator
- Adjustable in 1/10-watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100–240-volt AC
- Less to do and less to buy makes this the most cost-effective and user-friendly unit on the market
- This desiccator is a ground reference unit, which means no grounding pad is required
- 4-year warranty

#A942 (570-0583)	ea
------------------------	----

Desiccator Tips



Utilize the standard 3/32" (2.4 mm) stainless steel shaft. All molded electrodes feature safety-grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil.

Sharp Dermal Tips, Non-Sterile (Bulk Packaged) (112-7171)	100/box
Sharp Dermal Tips, Sterile (Packaged Individually) (112-7172)	50/box
Blunt Dermal Tips, Non-Sterile (Bulk Packaged) (112-7173)	100/box
Blunt Dermal Tips, Sterile (Packaged Individually) (112-7175)	50/box
Disposable Handpiece Sheaths, Sterile (112-7176)	25/box
Sharp H-Type Tip, Sterile (Packaged Individually) (570-0888)	50/box
Sharp H-Type Tip, Non-Sterile (Bulk Packaged) (570-0890)	100/box
Blunt H-Type Tip, Sterile (Packaged Individually) (570-0889)	50/box
Sharp H-Type Tip, Non-Sterile (Bulk Packaged) (570-0891)	100/box



Electrosurgical Generator and High-Frequency Desiccator



This electrosurgical generator is a high-frequency desiccator with cutting ability. One half of the unit is the Henry Schein High Frequency Desiccator. The other half is an electrosurgical generator that features 50 watts of cutting power and coag capabilities—plenty of power for most cutting procedures performed in the doctor's office.

- 40-watt high-frequency desiccator
- 2 machines in 1 convenient package
- Adjustable in 1/10-watt increments in the first 10 watts, in monopolar and bipolar
- No grounding pad is required during desiccation mode
- 4-year warranty

#A952 (570-0582)	ea
------------------------	----



Center-Post Colposcopes



- Continuous zoom (3.5X–24X)
- Trinoc tube with built-in C-mount camera adapter
- Dual-diopter adjustment
- Center-post mount
- Quartz halogen, fiber-optic illumination with intensity control and green filter

Binocular Head (570-0148)	ea
Trinocular Head, Camera-Ready Zoom (570-0149)	ea



Portable Cryosurgical H-Type Kits



- 88 mL in each canister
- More fast-charging treatment buds than competing brands
- Safe dispensing system
- Short treatment times
- 3-year shelf life
- Effective treatment for the removal of skin tags, common warts, plantar warts, actinic keratosis, genital warts, molluscum contagiosum, lentigo, and seborrheic keratosis

Single Canister (570-1141)	ea
Contains: 30–5-mm round treatment buds & 20–2-mm arrow treatment buds.	
Dual Canister (570-1137)	ea
Contains: 50–5-mm round treatment buds & 30–2-mm arrow treatment buds.	



Portable Cryosurgical V-Type Kits



- Safe dispensing system
- More fast-charging treatment buds than competing brands
- More treatment cones (6)
- Short treatment times
- 5-year shelf life
- Effective treatment for the removal of skin tags, common warts, plantar warts, actinic keratosis, genital warts, molluscum contagiosum, lentigo and seborrheic keratosis

Single Canister, 162 mL (570-1138)	ea
Contains: 40–5-mm round treatment buds, 20–2-mm arrow treatment buds & 6 treatment cones (3 mm, 5 mm, 7 mm, 9 mm, 12 mm & 16 mm).	
Replacement Canister, 162 mL (570-1139)	ea
Contains: 20 treatment buds & 3 treatment cones.	
Replacement Canister, 236 mL (570-1140)	ea
Contains: 20 treatment buds & 3 treatment cones.	

Equipment Finance Made Easy

Running a successful medical practice requires the right resources and funding. Henry Schein Financial Services offers flexible financing to help you expand, upgrade equipment, and embrace the latest technology, enhancing revenue and efficiency.



Customized Financing Solutions to Expand, Upgrade, and Enhance Your Practice

OUR FINANCING BENEFITS INCLUDE:

- ✓ 100% Financing*
- ✓ Up to 6-month Payment Deferrals
- ✓ Customized Repayment Programs
- ✓ Same-Day Approvals and Documentation
- ✓ Simple, Fast Application Process

PLUS

Take advantage of Section 179 and get a possible tax deduction for up to 100% of your purchase on your 2026 taxes for new equipment, technology, or other qualifying purchases.**

Maximize your practice's potential with **Henry Schein Financial Services**.



Contact us today to schedule a consultation today to **learn more about how we can support your growth and success.**

**SCAN THE QR CODE
TO LEARN MORE**

Disclaimer:

* Rates and programs are subject to change without notice. All transactions are subject to the satisfaction of underwriting guidelines, credit approval by third-party lenders, and documentation requirements, and not all applicants will qualify. Certain other restrictions and additional terms and conditions may apply.

** Neither Henry Schein, Inc. nor Henry Schein Financial Services provides tax advice. Please consult with a qualified professional tax advisor to determine your eligibility for a Section 179 tax deduction. Must purchase and place into service by December 31, 2026.

©2026 Henry Schein, Inc. Neither Henry Schein, Inc. nor Henry Schein Financial Services provides financial advice. Please consult your financial advisor. Neither Henry Schein, Inc. nor Henry Schein Financial Services is a bank and neither represents itself as such, nor conducts banking activities. Henry Schein Financial Services may receive a marketing fee from the vendor for products/services purchased. "World's Most Ethical Companies" and "Ethisphere" names and marks are registered trademarks of Ethisphere LLC.

Introcan Safety[®] 2 IV Catheter with Multi-Access Blood Control

Now enhanced with multi-access blood control to automatically reduce the risk of blood exposure during insertion and every time the hub is accessed*

The Introcan Safety 2 IV Catheter is intended for insertion into a patient's vascular system for short-term use to sample blood, monitor blood pressure or administer fluids and blood intravascularly. The catheters may be used intravascularly with power injectors at a maximum pressure of 325 psi with a Luer lock connection only.

For product availability, please contact your sales representative.




B. Braun Medical Inc. | Bethlehem PA |
bbraunusa.com

* B. Braun Data on File.
Rx only. ©2024 B. Braun Medical Inc., Bethlehem, PA.
All rights reserved. 24-0070_1/24



Scan for eIFU and
more information



Henry Schein® Blood Control Safety IV Catheters Push Button Style, 50/Box


Safety IV Catheter Blood Control (Green), 18Gx1.16"	(570-3712)	Safety IV Catheter Blood Control (Blue), 22Gx1"	(570-3715)
Safety IV Catheter Blood Control (Pink), 20Gx1.16"	(570-3713)	Safety IV Catheter Blood Control (Yellow), 24Gx0.75"	(570-3716)
Safety IV Catheter Blood Control (Pink), 20Gx1"	(570-3714)		

Order today at HENRYSCHEIN.COM/MEDICAL



SAFETY IV CATHETER
From a Name You Can Trust.

Item Number	Description	Quantity
(570-2586)	24G x 0.75"	50/box
(570-2591)	24G x 0.5"	50/box
(570-2587)	22G x 1"	50/box
(570-2588)	20G x 1.16"	50/box
(570-2589)	20G x 1"	50/box
(570-2590)	18G x 1.16"	50/box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



Your Trusted Source for Safety IV Catheters With or Without Blood Control



**TrueSafe
Comfort® BC**



**ClearSafe
Comfort® BC**

For more information visit henryschein.com/medical



I.V. Medication Delivery

From a Name You Can Trust.

Medication Delivery Sets and Connectors

Henry Schein provides market-leading infection control technologies to help minimize bacterial contamination. We carry a wide selection of injection and needle-free medication delivery ports to meet a variety of medication delivery needs.



Features & Benefits:

- Flat top needle-free valve with tight septum seal helps to support effective disinfection
- Metal free for MRI compatibility
- DEHP-free, non-latex and sterile

Item Code	Type	Product Description	Qty/Unit
(570-0881)	Primary I.V. Administration Set	 I.V. Administration Set. 15 drops/mL. 73°. Vented, Slide Clamp, Roller Clamp, Injection Y-site, Rotating Male Luer Lock	Each
(570-0885)	Primary I.V. Administration Set	 I.V. Administration Set. 15 drops/mL. 86°. Vented, Roller Clamp, Needle-free Y-site, Rotating Male Luer Lock	Each
(570-0887)	Primary I.V. Administration Set	 I.V. Administration Set. 15 drops/mL. 105°. Vented, Back Check Valve, Slide Clamp, (2) Needle-free Y-sites, Roller Clamp, Rotating Male Luer Lock	Each
(570-0883)	Secondary I.V. Administration Set	 Secondary I.V. Administration Set. 15 drops/mL. 40°. Vented, Roller Clamp, Hanger, Rotating Male Luer Lock	Each
(570-0884)	Microbore Extension set	 Microbore extension set. 8". Removeable Needle-free Connector, Slide Clamp, Rotating Male Luer Lock	Each
(570-0882)	Needle-free Connector	 Needle-free Connector	Each



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



SafeDAY™ IV Administration and Extension Sets

Cost-effective IV Sets with Needleless Connectors

- Designed to make access clear and easy
- Limited to 24-hour maximum use, ideal for surgical and oncology centers, emergency medicine and physician offices
- Available in convenient lengths and configurations

To learn more, contact your Henry Schein representative.

Rx only. ©2024 B. Braun Medical Inc., Bethlehem, PA. All Rights Reserved. 24-0281_06/24

B | BRAUN
SHARING EXPERTISE



Total ownership. Trusted quality. From design to distribution.

With complete ownership of every step, BD delivers excellence from end-to-end. We are committed to providing a sustainable supply of needles and syringes with safe and reliable performance your clinicians and patients can trust.



Learn more about our syringe portfolio and our commitment to manufacturing quality.

BD, Franklin Lakes, NJ, 07417, U.S.


[bd.com](https://www.bd.com)

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved. BD-139898 (11/24)



HYPODERMIC SAFETY NEEDLES
WITH ISO COLORED PLASTIC HUB

- Two single-handed safety activation methods include hard surface and thumb activation, which may help reduce the risk of accidental needle stick injuries
- Easy-to-use and compatible with all standard luer lock and luer slip-tip syringes
- Hubs and protective shield are color-coded according to ISO standards to facilitate easy and quick identification of needle gauges
- Bevel-up needle allows for low-angle injections
- Protective arm allows for clear visibility of the injection site, ensuring an unobstructed pathway
- Single use



HENRY SCHEIN® HYPODERMIC SAFETY NEEDLES
with ISO Colored Plastic Hub (100/box)

Item Number	Gauge	Color	Needle Length	Bevel Angle	Inner Diameter	Outer Diameter (mm)
570-2709	30	Light Yellow	½"	Long	RW	0.31
570-2706	27	Grey	1 ½"	Long	RW	0.41
570-2707	27	Grey	1 ¼"	Long	RW	0.41
570-2708	27	Grey	⅝"	Long	RW	0.41
570-2704	25	Orange	1 ½"	Long	TW	0.51
570-2703	25	Orange	1"	Long	TW	0.51
570-2705	25	Orange	⅝"	Long	TW	0.51
570-0340	23	Blue	1 ½"	Long	TW	0.63
570-2701	23	Blue	1"	Long	TW	0.63
570-2710	23	Blue	⅝"	Long	TW	0.63
570-2700	22	Black	1 ½"	Long	TW	0.71
570-2699	22	Black	1"	Long	TW	0.71
570-2698	21	Green	1 ½"	Long	TW	0.82
570-2696	21	Green	1"	Long	TW	0.82
570-2697	20	Yellow	1 ½"	Long	TW	0.90
570-2695	20	Yellow	1"	Long	TW	0.90
570-2694	18	Pink	1 ½"	Short	TW	1.28
570-2693	18	Pink	1"	Short	TW	1.28

RW—Regular Wall | **TW**—Thin Wall

Order today at HENRYSCHEIN.COM/MEDICAL

CONVENTIONAL HYPODERMIC NEEDLES
WITH ISO COLORED PLASTIC HUB

- High-quality stainless steel
- Large inside diameter (ID) provides an increased flow rate
- Semi-transparent needle hub is ISO color-coded by size for easy identification and safety
- Needle hub-compatible with all standard Luer-lock and Luer slip-tip syringes
- Sharp needle improves patient comfort
- Wide range of needle sizes for each procedure
- Product conforms to ISO 7864
- Sterile
- Non-pyrogenic



HENRY SCHEIN® CONVENTIONAL HYPODERMIC NEEDLES
with ISO Colored Plastic Hub (100/box)

Item Number	Gauge	Color	Needle Length	Bevel Angle	Inner Diameter	Outer Diameter (mm)
112-6148	30	Light Yellow	½"	Long	RW	0.30
112-6149	30	Light Yellow	1"	Long	RW	0.30
112-7108	27	Grey	½"	Long	RW	0.40
112-7103	27	Grey	1 ¼"	Long	RW	0.40
900-4468	25	Orange	⅝"	Long	TW	0.50
112-7112	25	Orange	⅞"	Long	TW	0.50
900-4470	25	Orange	1"	Long	TW	0.50
900-4472	25	Orange	1 ½"	Long	TW	0.50
900-4473	23	Blue	1"	Long	TW	0.60
112-7107	22	Black	1"	Long	RW	0.70
112-7101	22	Black	1 ½"	Long	RW	0.70
112-7100	21	Green	1"	Long	RW	0.80
112-7111	21	Green	1 ½"	Long	RW	0.80
112-7109	20	Yellow	1"	Long	RW	0.90
112-7110	20	Yellow	1 ½"	Long	RW	0.90
900-4469	18	Pink	1"	Short	TW	1.20
900-4471	18	Pink	1 ½"	Short	TW	1.20

RW—Regular Wall | **TW**—Thin Wall

Order today at HENRYSCHEIN.COM/MEDICAL



SOL-M™ Syringe
SOL-MILLENNIUM

Luer-lock, Without Needle

#P180001PP, Low Dead Space, 1 mL
(138-6993) 100/box



SOL-CARE™ Luer-Lock Syringes
with Safety Needles
SOL-MILLENNIUM

Can be used with any standard Luer-lock or slip-tip syringes. Bevel-up needle allows for low-angle injections. Once activated, the needle remains securely fastened via a double-lock mechanism.

#12558SN, 1 mL, 25 ga x 5/8", Orange
(138-7044) 50/box
#12515SN, 1 mL, 25 ga x 1 1/2", Orange
(139-6166) 50/box
#12310SN, 1 mL, 23 ga x 1", Blue
(139-6165) 50/box
#32558SN, 3 mL, 25 ga x 5/8", Orange
(139-2919) 50/box



SOL-M™ Syringes
with Exchangeable Needles
SOL-MILLENNIUM

- Non-safety
- Regular wall
- Blister pack
- Luer Lock

3-mL Syringe with Needle
#1832558, 25 ga x 5/8" 100/box
#1832510, 25 ga x 1" 100/box
#1832515, 25 ga x 1 1/2" 100/box
#1832310, 23 ga x 1" 100/box
#1832006) 100/box



SOL-GUARD™
Safety Needles
SOL-MILLENNIUM

Designed to help prevent needlestick injury. The robust safety mechanism is intuitive and compatible with all standard Luer-lock and slip-tip syringes.

#SND486, 25 ga x 5/8" 100/box
#SND475, 25 ga x 1" 100/box
#SND463) 100/box



SOL-CARE™ Luer Lock Safety Syringe
with Exchangeable Needles
SOL-MILLENNIUM

- Highly transparent barrel with bold, precise scale for perfect visualization and accurate dosage
- Integrated plunger stopper prevents plunger from being pulled out, avoiding exposure to hazardous fluids

3-mL Syringe
#100078IM, 25 ga x 5/8", Orange
(261-0258) 100/box
#100077IM, 23 ga x 1", Blue
(261-0260) 100/box



Sharps Containers
From a Name You Can Trust.

- For use by professional health-care givers
- Meets all NIOSH and OSHA regulations
- FDA approved
- Puncture-resistant
- Transparent body and/or top allows for easy viewing
- High-impact and shatter-resistant polypropylene



#0310-1500-HS, 1 qt (112-7149) Ea.	#0319-1500-HS, 1-gal Slide Lid (570-0378) Ea.	#0330-150R-HS, 3-gal (570-0212) Ea.
#0310-150T-HS, 1.5 qt (570-0376) Ea.	#0319-150R-HS, 1-gal Rotary Lid (570-0379) Ea.	#0380-150R-HS, 8 gal (570-0380) Ea.
#0322-150R-HS, 2.2 qt (570-0213) Ea.	#0320-1500-HS, 2 gal, Slide Lid (112-7150) Ea.	#0354-150M, M-Style, 5 qt (570-2554) Ea.
#0354-150M-HS, 5 qt (112-7152) Ea.	#0320-150R-HS, 2-gal Rotary Lid (570-0377) Ea.	#0320-1500, Slide Lid, 2 gal (570-2553) Ea.
#0354-150B-HS, 5.4 qt (112-7151) Ea.		



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



General-Purpose Freestanding Undercounter Refrigerator
AMERICAN BIOTECH SUPPLY (ABS)

Intended for general-purpose applications only where precise temperature control and maintenance are not essential and there are no requirements for temperature alarms.

- Swing door with keyed door lock
- 3 shelves (1 fixed/2 adjustable)
- Interior light
- HFC-free refrigerant (R600A), foam insulation, and packaging
- UL/C-UL listed
- Temperature range: 0°C to 10°C
- Capacity: 5 cu ft
- Exterior dimensions: 20¼"W x 22" D x 31½"H
- 1-year parts and labor warranty

#ABT-HC-UCFS-0504W
(125-4697) ea

This unit must have 4" clearance on sides and back for adequate ventilation.



Pharmacy Undercounter Built-In Refrigerator
AMERICAN BIOTECH SUPPLY

- 4.6-cu-ft capacity
- Features a microprocessor temperature controller and display that allows for precise temperature control, verification, and recovery (vs. mechanical-style thermostats)
- Acceptable temperature ranges for vaccine storage are 36°F to 46°F (2°C to 8°C) in the refrigerator and -58°F to 5°F (-50°C to -15°C) in the freezer
- Keyed door lock
- 3 shelves with guardrail on back to prevent product from being pushed against the back wall
- 115V, 60 Hz, ½ hp, 1.74 amps
- 23¾"W x 24½"D x 33¾"H

(125-4626) ea



Premier Pharmacy Dual-Temp Refrigerator & Auto-Defrost Freezer Combination Unit
AMERICAN BIOTECH SUPPLY

- Separate microprocessor temperature controller allows for precise temperature control, verification, and recovery
- Interior fan-driven, forced-draft air-circulation system provides excellent temperature uniformity and recovery after door openings
- Equipped with digital temperature displays for convenient monitoring
- Critical samples and supplies are security backed with an extensive package of audible and visual alarms
- Freezer unit features an automatic defrost operating system
- 3 adjustable shelves for refrigerator, 2 adjustable shelves for freezer
- Shelves have back rails to prevent product from being pushed against back wall
- Capacity: 7 cu ft total (5.2 cu ft for refrigerator, 1.3 cu ft for freezer)
- 53½"H x 23¾"W x 24¼"D
- 1-¾" probe access port for refrigerator, 1-¾" probe access port for freezer
- Keyed door locks

(125-4637)ea



Vaccine Temperature Data Logger
HENRY SCHEIN®

- High accuracy
- 24 million readings
- Maximum of 170 files (.csv)
- Internal (ambient) and external (fridge/freezer) sensors reading display
- Minimum and maximum reading memories
- Low and high alarm limits setting
- Flashing light alarm indicator
- °C or °F scale
- Real-time display
- Waterproof external sensor
- Low battery indicator
- Removable SD memory card or USB interface for data retrieval
- Optional external power

(570-2722)ea

ADRENERGICS

Adrenalin® Injectable PAR STERILE PRODUCTS

1 mg/mL, 1-mL SDV
(248-0883) ea
(890-0257) 25/pkg
1 mg/mL, 30-mL MDV
(129-4847) ea



AUVI-Q® Epinephrine Auto-Injectors KALEO

0.1 mg
(148-1961) 2/pkg
0.15 mg
(148-1959) 2/pkg
0.3 mg
(148-1960) 2/pkg

Compare to Adrenalin® Epinephrine Injectable Ampules BPI LABS

1:1000, 1 mL
(248-0899) ea
(128-9991) 10/box



Epinephrine Auto-Injector TEVA

Authorized generic of EpiPen®.
Adult, 0.3 mg
(137-6239) 2/pkg
Junior, 0.15 mg
(137-6238) 2/pkg

Neffy® Epinephrine Nasal Spray ARS PHARMACEUTICALS

1 mg
(149-3184) 2/pkg
2 mg
(147-7945) 2/pkg



ANESTHETICS

Diprivan® (Propofol) Injectable* FRESENIUS KABI

10 mg/mL
10-mL SDV
(150-0143) 10/box
20-mL SDV
(150-0131) 10/box
50-mL Vial
(110-7244) 20/box



*EDTA as preservative.



Ethyl Chloride Spray GEBAUER

Topical Anesthetic Skin Refrigerant
• Fast: Swab the target area with antiseptic. Spray 4–7 seconds immediately prior to injections, IV start, or venipuncture
• Low cost: pennies per application
Fine-Stream Spray Bottle, 3.9 oz
(138-9019) ea
(138-9016) 4/pkg
Medium-Stream Spray Bottle, 3.9 oz
(138-5777) ea
(138-5778) 4/pkg
(138-5775) 12/case
Mist Spray Can, 3.9 oz
(138-6360) 4/pkg
(138-6330) 12/case



Pain Ease® Spray GEBAUER

Topical Anesthetic Skin Refrigerant
3.9-oz Mist Spray
(791-0634) ea
1-oz Medium-Stream Spray
(132-9467) ea
3.9-oz Medium-Stream Spray
(761-0029) ea
Can be used in facilities that restrict use of flammable components.

Compare to Diprivan® Propofol Injection Emulsion HERITAGE PHARMACEUTICALS

10 mg/mL
20-mL SDV
(142-5726) 20/box
50-mL SDV
(142-5730) 20/box
100-mL SDV
(142-5731) 10/box



Compare to Ultane® Sevoflurane Liquid Inhalation 100% PIRAMAL CRITICAL CARE

250-mL Bottle
(116-9356) ea



Xylocaine® Plain Injectable FRESENIUS KABI

1%, 30-mL SDV MPF
(150-0072) 25/pkg
2%, 5-mL SDV, MPF
(150-0101) 25/pkg
2%, 50-mL MDV
(150-0096) 25/pkg



Xylocaine® With Epinephrine Injectable FRESENIUS KABI

1%, 30-mL SDV
(150-0122) 25/pkg
1%, 50-mL MDV
(150-0073) 25/pkg
2%, 20-mL MDV
(150-0109) 25/pkg
2%, 50-mL MDV
(150-0097) 25/pkg



Zynrelef® Injection HERON THERAPEUTICS

200 mg/6 mg, 10-mL SDV
(148-0933) ea
400 mg/12 mg, 20-mL SDV
(148-0924) ea




BONE RESORPTION INHIBITORS



BILDYOS® (denosumab-nxxp) Injection
ORGANON

1-mL Syringe, 60 mg
(149-5453) ea

Biosimilar to Prolia®



BILPREVDA® (denosumab-nxxp) Injection
ORGANON

1.7-mL SDV, 120 mg
(149-5454) 25/pkg

Biosimilar to Xgeva®



Evenity® Injection
AMGEN

105 mg/1.17 mL
(147-1945) 2/box



Jubbonti® (denosumab-nxxp) Injection
SANDOZ

1-mL Syringe, 60 mg
(149-0249) ea

Biosimilar to Prolia®




Prolia Injectable Prefilled Syringe
AMGEN

60 mg/mL
(147-6037) ea



Xgeva Injection
AMGEN

120 mg/1.7-mL SDV
(149-9424) ea



WYOST® (denosumab-bbdz) Injection
SANDOZ

1.7-mL SDV, 120 mg
(149-0248) ea

Biosimilar to Xgeva®

CNS AGENTS



Byfavo™ for Injection
ACACIA PHARMA

20-mg SDV
(140-0959) 10/box



Botox® Injectable
ABBVIE LABORATORIE

Nonreturnable
100-Unit Vial
(498-2546) ea
200-Unit Vial
(119-9873) ea



Dantrium® IV
PAR STERILE PRODUCTS

20 mg, 70-mL Vial
(115-7589) 6 vials/box



Midazolam HCl Injectables
PFIZER

1 mg/mL, 2-mL SDV
(140-2328) 25/box
1 mg/mL, 5-mL MDV
(248-0885) ea
(140-1490) 10/box
1 mg/mL, 10-mL MDV
(143-6232) 10/box
5 mg/mL, 1-mL MDV
(142-8879) 25/box



Midazolam HCl Injectables–CIV
PFIZER

1 mg/mL, 5-mL SDV
(127-7882) 10/box
5 mg/mL, 1-mL SDV (Nonreturnable)
(248-0822) ea
5 mg/mL, 1-mL SDV
(127-7867) 10/box
5 mg/mL, 2-mL SDV (Nonreturnable)
(248-0863) ea
5 mg/mL, 2-mL SDV
(127-7906) 10/box
5 mg/mL, 5-mL MDV (Nonreturnable)
(248-0884) ea

CORTICOSTEROIDS



Celestone® Soluspan® Injectable Suspension
ORGANON

6 mg/mL, 5-mL MDV
(140-7856) ea



Dexamethasone Sodium Phosphate Injection
HIKMA

4 mg/mL, 1-mL SDV25/pkg
10 mg/mL, 1-mL SDV
(138-6758)25/pkg



Compare to Decadron®
Dexamethasone Sodium Phosphate Injectable
MYLAN INSTITUTIONAL

4 mg/mL, 5-mL Vial
(122-7115) 25/pkg



Dexamethasone Sodium Phosphate Injection
FRESENIUS KABI

4 mg/mL, 30-mL MDV
(147-1174) ea



Compare to Decadron®
Dexamethasone Sodium Phosphate Injectable
FRESENIUS KABI

Preservative-free
10 mg/mL, 1-mL SDV
(406-7616) 25/pkg



Dexamethasone Sodium Phosphate
SOMERSET

10 mg/mL, 10-mL MDV
(138-1136) 10/box



Methylprednisolone Acetate
SANDOZ

40 mg/mL
5-mL MDV
(148-6813) ea



Compare to Depo-Medrol®
Methylprednisolone Acetate Injection
AMNEAL

1-mL SDV
40 mg/mL
(140-9021) ea
(140-9026) 25/box
10-mL MDV
40 mg/mL
(149-0585) ea
80 mg/mL
(140-9031) ea
(140-9048) 25/box



Compare to Kenalog®-40
Triamcinolone Acetonide Injection
AMNEAL

40 mg/mL
1-mL SDV
(142-2167) ea
(142-3667) 25/box
5-mL MDV
(142-3668) ea
10-mL MDV
(142-2227) ea



Zilretta™ Injectable
PACIRA

32-mg SDV
Triamcinolone Acetonide Extended-Release
(142-9357) ea
Series of injections: 1
J Code: J3304

FIBRIN SEALANTS



ARTISS Prefilled Syringe
BAXTER ADVANCED SURGERY

The first and only fibrin sealant custom-designed for tissue adherence in plastic, reconstructive, and burn surgeries. It provides capabilities to adhere autologous skin grafts and skin flaps to wound beds in adults and children.

- 2 mL (147-3491) ea
- 4 mL (147-3494) ea
- 10 mL (147-3495) ea



TISSEEL Injection
BAXTER ADVANCED SURGERY

The only fibrin sealant indicated for both hemostasis and sealing* as well as for use in both heparinized and pediatric patients greater than 1 month of age. It is also effective in patients medicated with antiplatelet drugs.¹

- Needle Kit, 2 mL (147-3489) ea
- Needle Kit, 4 mL (147-3483) ea
- Needle Kit, 10 mL (147-3482) ea
- Prefilled Syringe, 2 mL (147-3492) ea
- Prefilled Syringe, 4 mL (147-3490) ea
- Prefilled Syringe, 10 mL (147-3480) ea

*Sealing indicated to prevent leakage from colonic anastomoses following the reversal of temporary colostomies.


¹ TISSEEL [fibrin sealant] full prescribing information: http://baxterpi.com/pi-pdf/Tisseel_PL.pdf

OPHTHALMIC



Dextenza®
OCULAR THERAPEUTIC
Dexamethasone Ophthalmic Insert

- 0.4 mg (145-5615) ea
- (145-5614) 10/pkg



Susvimo® Implant Injection Kit
GENENTECH

A prescription-only ocular implant for treating wet age-related macular degeneration (AMD), diabetic macular edema (DME), and diabetic retinopathy (DR).

- 10 mg/0.1 mL, 100 mg/mL, 1-mL SDV (149-9996) ea
- Initial Fill Needle, 34 ga x 3/16" (149-9994) ea
- Refill Needle, 34 ga x 1/6" (149-9995) ea
- Implant Insertion Tool (147-6769) ea

OTHER PHARMACEUTICALS



Barhemsys® Injection
ACACIA PHARMA

- 5 mg/2 mL SDV (140-0931) 10/box
- 10 mg/4 mL SDV (140-1706) 10/box



Glycopyrrolate
PIRAMAL CRITICAL CARE

- 0.2 mg/mL, 20-mL MDV (143-6161) 10/box



Ondansetron HCl Injection
APOTEX

- 2 mg/mL, 2-mL SDV (248-0902) ea
- (125-9100) 25/box

Gebauer's Ethyl Chloride®
Instant Topical Anesthetic Spray

First, Best & Always

Gebauer's Ethyl Chloride helps temporarily control needle and minor surgery pain in an instant!

Physicians have trusted Gebauer's Ethyl Chloride topical anesthetic spray for more than 125 years.



Gebauer Company
Ethyl Chloride® | Pain Ease® | Cryo Professional



IMPORTANT RISK & SAFETY INFORMATION

Published clinical trial results support the use in children 3 years of age and older. Ethyl chloride is **FLAMMABLE** and should never be used in the presence of an open flame or electrical cautery equipment.

Use in a well-ventilated area. Intentional misuse by deliberately concentrating or inhaling the contents can be harmful or fatal. Do not spray in eyes. Over application of the product may lead to frostbite and/or altered skin pigmentation. Cutaneous sensitization may occur, but appears to be extremely rare. **CAUTION:** Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

4444 East 153rd Street, Cleveland, OH 44128 To learn more, visit the website at www.Gebauer.com

2351.1

ONESTEP+ DRUG OF ABUSE RAPID URINE TEST CUP

HENRY SCHEIN®
BRAND PRODUCTS

Features & Benefits

- CLIA waived
- Designed to detect a wide range of substances: AMP, BAR, BUP, BZO, COC, MDMA, EDDP, MTD, MET, MOP, OPI, OXY, PCP, PPX, TCA, and THC*
- Simple to use, no need to tilt cup to activate test
- Self-contained cup is ideal for sending preliminary positive specimens to the lab for confirmation
- Lid click feature ensures the cup is sealed properly, maintaining integrity of the test
- Results available in 5 minutes
- Temperature strip helps to monitor specimen tampering
- Adulteration strip offers initial screening and assists in monitoring any tampering of the specimen

*AMP-Amphetamine, BAR-Barbiturate, BUP-Buprenorphine, BZO-Benzodiazepines, COC-Cocaine, MDMA -Ecstasy, EDDP- Methadone Metabolite, MET- Methamphetamine, MTD- Methadone, MOP-Morphine, OPI-Opiates, OXY-Oxycodone, PCP-Phencyclidine, PPX-Propoxyphene, TCA-Tricyclic Antidepressants, THC-Marijuana



Henry Schein® OneStep+ Drug of Abuse Rapid Urine Test Cup 25/Box

Kit Contents: (25) OneStep+ Drug of Abuse Rapid Urine Test Cups, (2) Specimen Validity/Adulteration Strips (if applicable), (30) Disposable Gloves, (25) Security Seal Labels, (1) Instructions For use, (1) Procedure Card

Test Cup, 12 Panel—With Adulteration Strip..... (570-3927)
AMP1000, BAR300, BZO300, COC300, MDMA500, MET1000, MOP300, MTD300, OXY100, PCP25, TCA 1000, THC50 (CR, SG, pH)

Test Cup, 14 Panel..... (570-3928)
AMP1000, BAR300, BUP10, BZO300, COC300, EDDP300, MDMA500, MET1000, MOP300, MTD300, OXY100, PCP25, TCA1000, THC50

Order today at HENRYSCHEIN.COM/MEDICAL



Improve care,
reduce infection risk,
and maintain compliance

With one simple solution



**HemoCue®
Glucose 201 System**

CLIA waived



**Hb 801 Analyzer
HEMOCUE**

The HemoCue® Hb 801 system seamlessly integrates point-of-care testing into your workflow.

- Wireless connectivity through middleware
- Results in under 1 second
- Reagent-free cuvettes, improving open-vial stability
- Large display with intuitive on-screen symbols
- Internal quality control self-test
- USB adapter (included) or 3 AA batteries (not included) power supply

Hb 801 Analyzer
(135-3349) ea

Hb 801 Hemoglobin Microcuvettes
(135-3352) 200/box

Contains: 25 per bottle, 4 bottles per box.



**Glucose 201 Analyzer
HEMOCUE**

The HemoCue Glucose Systems are based on a glucose-dehydrogenase method and consist of a small, dedicated analyzer and a unique, disposable microcuvette. They will produce lab-accurate results for diagnosing or monitoring diabetes, as well as monitoring neonatal blood glucose levels.

- More accurate and 68% higher reimbursement than glucose meters
- FDA-compliant unit when screening and diagnosing
- Analyzer: 0–444 mg/dL (0–24.6 mmol/L)
- Only 5 µL of capillary, venous, or arterial blood is needed
- No inherent limitations related to hematocrit or known drug interferences
- No recalibration needed between cuvette lots
- 3-year warranty

#120706, Analyzer
(104-2290) ea

#110706, Microcuvettes, Glucose
(439-9158) 100/box

HCG

OneStep+ Pregnancy Cassette Tests



- Urine test
- CLIA waived
- 25 mIU/mL detection level
- Sensitivity: 99%
- Room-temperature storage

25-Patient Kit
(900-4073) ea
 Contains: 25 pouched tests per kit & 1 package insert.

100-Patient Kit
(900-4074) ea
 Contains: 100 pouched tests per kit & 1 package insert.
 CPT® Codes: Urine 81025*/Serum 84703*



OneStep+ Pregnancy Combo Cassette Test



- Urine or serum test
- CLIA waived for urine only
- 25 mIU/mL detection level
- Sensitivity: 99%
- Room-temperature storage

25-Patient Kit
(900-4075)25/box
 Contains: 25 pouched tests per kit & 1 package insert.
 CPT® Codes: 81025QW* (Urine); 84703* (Serum)

OSOM® Ultra hCG Combo Test SEKISUI



An easy-to-use lateral-flow immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in serum or urine for the early confirmation of pregnancy. The Urine Control Set is intended for *in vitro* diagnostic use in quality-control testing with the OSOM Ultra hCG Combo Test when a urine specimen is used. The Serum Control Set is intended for *in vitro* diagnostic use in quality-control testing with the OSOM Ultra hCG Combo Test when a serum sample is used.

- Accurate: 100% sensitivity and specificity (urine); 99.5% and 99% (serum)
- Fast: results in 3–5 minutes
- Easy-to-use: single-step procedure with <30 seconds hands-on time
- Objective: detects β-core fragments (data on file)
- QC Inside®: 2 additional tests at no additional charge for external QC testing

#1004, hCG Combo Test
(148-4398) 20/box

#134, Urine Control Set
(327-6476) ea
 Contains: 1–10-mL bottle positive & 1–10-mL bottle negative.

#138, Serum Control Set
(611-1052) ea
 Contains: 1–5-mL bottle positive & 1–5-mL bottle negative.
 CPT® Codes: 81025* (Urine); 84703* (Serum)

TRUE METRIX



True Metrix® Pro Strips



Advanced technologies—the meter, a complex algorithm, chemistry, and electrodes—on the test strip work together as part of the True Metrix® Pro system to produce accurate results. Featuring Triple Sense Technology™, the system provides proven accuracy and confidence in results.

- No-coding technology eliminates the need for coding of meter
- Code is assigned during manufacturing process
- Code is printed onto contacts on end of test strip
- Meter reads code on insertion

#R3099P-450, Strips
(570-0327)50/box

Controls

#R5H01-1, Level 1
(125-4034) ea

#R5H01-2, Level 2
(125-4035) ea

#R5H01-3, Level 3
(125-4036) ea




True Metrix® Pro Meter



- No coding
- Tiny, 0.5-µL sample size
- Stores 500 results
- As fast as 4 seconds
- BG range: 20–600
- Data management
- Download capabilities
- Glucose control detection
- Strip-release button

#RE4099P-43
(570-3418) ea
 Contains: 1 True Metrix® Pro meter, installed 3-volt battery, instructions for use, logbook & compact carrying case.

HBA1C TESTING



A1CNow®+ Test



The handheld monitor provides accurate, in-office A1C results, enabling physicians to make on-the-spot treatment decisions for diabetes patients.

- In-office testing: no waiting for lab results
 - Results in just 5 minutes
 - Hands-on procedure time is less than 1 minute
 - Simple, 3-step procedure
 - CLIA waived
 - No calibration, no daily controls, and no maintenance
 - No refrigeration necessary if used within 3 months
 - No capital equipment required
 - Enables A1C testing in every exam room
- (570-3428)**20/box



DCA™ Vantage Analyzer



Point-of-care immunoassay analyzer for diabetes management. Provides precise and accurate quantitative test results for HbA1c in whole blood during a patient's visit to help improve decision making. With a single urine sample, the analyzer also provides quantitative test results for low levels of albumin, creatinine, and A:C ratio in the physician's office or clinic.

#5075US, Analyzer
(110-0743) ea

#5035C, HbA1c Reagent Kit
(147-8045)10/pkg

#5068A, HbA1c Control Kit
(147-6067) ea
 Contains: 2 vials of each control: high & low.



i-STAT® Handheld System

ABBOTT

Real-time, lab-quality results in minutes—not hours. The i-STAT® Handheld is an advanced, easy-to-use blood analyzer that provides a wide range of tests using just 2 or 3 drops of blood.

- Speed: With components of a full-scale lab in each cartridge, results occur in minutes
- Accuracy: A unique combination of biosensors and calibrants deliver high-quality results
- Reliability: Single-use cartridges reduce many of the problems inherent in multi-use systems

#04J60-20, iSTAT Handheld System (111-3324) ea
 Contains: 1 of each: analyzer, printer kit, downloader/recharger, simulator & rechargeable battery pack.

CLIA Waived Chemistry Tests—Results in 2 Minutes

#03P84-25, Creatine, Crea (111-3394) 25/box

Moderately Complex Chemistry Tests—Results in 2 Minutes

#09P31-26, CHEM 8+: Na, K, Cl, iCa, TCO₂, Glu, BUN/Urea, Crea, Agap*, Hct, Hgb* (*Calculated) (138-1904) 25/box

#03P83-26, Glucose: Glu (146-5249) 25/box

Blood Gas Tests—Results in 2 Minutes

#03P85-51, CG4+: pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Lactate (*Calculated) (146-0531) 25/box

#03P78-26, G3+: pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂* (*Calculated) (139-7158) 25/box

Blood Gas, Electrolytes and Hematology Tests—Results in 2 Minutes

#03P77-25, EG6+: Na, K, pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated) (111-3387) 25/box

#03P76-25, EG7+: Na, K, iCa, pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated) (111-3386) 25/box

#03P88-25, CG8+: Na, K, iCa, Glu, pH, PCO₂, PO₂, TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated) (115-6662) 25/box

Cardiac Marker Tests—Results in 10 Minutes (CK-MB Results in 5 Minutes)

#09P81-25, i-STAT High-Sensitivity Troponin (hsTnl) Test Cartridge, Moderately Complex (148-3789) 25/box

Coagulation Tests—Results in 17 Minutes (PT/INR Results in 5 Minutes)

#03P86-25, Celite ACT (111-3409) 25/box

#03P87-25, Kaolin ACT (111-3412) 25/box

#03P89-50, PT/INR: Prothrombin Time (146-0532) 24/box



Piccolo Xpress® Chemistry Analyzer
ABBOTT

Three steps, 12 minutes. The Piccolo Xpress® chemistry analyzer allows health care practitioners to perform routine multichemistry panels using 0.1 cc—only 3 or 4 drops—of whole blood, serum, or plasma.

- Speed: Add the sample, insert the reagent disc, and in approximately 12 minutes receive a report suitable for the patient file or transferred to a computer, printer, or LIS/EMR
- Ease: A fully automated system requiring no special operating skills
- Exclusivity: Offers the only CLIA-waived comprehensive metabolic panel (CMP) blood test on the market
- Accuracy: Self-calibrates with each run. The onboard continuous intelligent Quality Control (iQCT) monitors the analyzer, reagent reactions, and sample to ensure chemistry and instrument integrity

#07P05-01, Piccolo Xpress® Chemistry Analyzer (106-6368) ea
 Contains: 0.1 cc minipette, minipette tips, power cord, power supply, printer results paper & owners manual.

CLIA-Waived Test Panels

#07P02-08, Comprehensive Metabolic Panel: ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO₂, TP (174-6978) 10/box

#07P02-04, Basic Metabolic Panel: BUN, Ca, Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR* (137-3987) 10/box

#07P02-05, Lipid Panel: CHOL, CHOL/HDL*, nHDLc*, HDL, LDL*, TRIG, VLDL* (* Calculated) (146-0189) 10/box

#07P02-12, Lipid Panel Plus: ALT, AST, CHOL, CHOL/HDL*, nHDLc*, GLU, HDL, LDL*, TRIG, VLDL* (*Calculated) (104-6269) 10/box

#07P02-10, Liver Panel Plus: ALB, ALP, ALT, AMY, AST, GGT, TBIL, TP (104-3998) 10/box

#07P02-01, General Chemistry 6: ALT, AST, BUN, CRE, GGT, GLU, eGFR* (108-5450) 10/box

#07P02-11, General Chemistry 13: ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, GGT, GLU, TBIL, TP, UA, eGFR* (982-0407) 10/box

#07P02-02, Electrolyte Panel: Cl-, K+, Na+, tCO₂ (106-1669) 10/box

#07P02-09, Kidney Check: BUN, CRE, eGFR* (109-7255) 10/box

#07P02-07, Renal Function Panel: ALB, BUN, Ca, Cl-, CRE, GLU, K+, Na+, PHOS, tCO₂, eGFR* (166-0886) 10/box

#07P02-03, MetLyte 8 Panel: BUN, CK, Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR* (114-1284) 10/box

Moderately Complex Test Panels

#07P02-06, Hepatic Function Panel: ALB, ALP, ALT, AST, DBIL, TBIL, TP (195-7347) 10/box

#07P02-13, Basic Metabolic Panel Plus: BUN, Ca, Cl-, CRE, GLU, K+, Na+, tCO₂, Mg, Lactate Dehydrogenase (104-5577) 10/box

#07P02-14, MetLyte Plus CRP: BUN, CK, Cl-, CRE, GLU, K+, Na+, tCO₂, CRP (115-5701) 10/box

#07P02-15, BioChemistry Panel Plus: ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, CRP, GGT, GLU, TP, UA (116-9038) 10/box

#07P02-16, MetLac 12 Panel: ALB, BUN, Ca, Cl-, CRE, GLU, LAC, K+, Mg, Na+, PHOS, tCO₂, eGFR* (631-0050) 10/box



HENRY SCHEIN[®] STRETCHER SHEETS

Features & Benefits

- Easy-to-use, reliable disposable solution that eliminates the need for reusable linens
- Absorbent 2-ply design with a white tissue layer and blue protective poly-backed layer
- Single-use to help reduce cross-contamination
- Flat design for versatile patient care applications
- Smooth shape does not add bulk, providing patient comfort
- Not made with natural rubber latex
- Conveniently packaged with 50 sheets per case
- Measures 40" × 90" to fit most standard stretchers and treatment surfaces



Henry Schein[®] Stretcher Sheets 50 per Case

40" × 90", White/Blue, 2-ply Tissue/Poly(570-5555)

Order today at [HENRYSCHEIN.COM/MEDICAL](https://www.henryschein.com/medical)



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS





Table Paper



- White
- Quiet, absorbent barrier
- Designed and manufactured to provide strength, protection, and value

Smooth, 225 ft
 18" (900-4978) 12 rolls/case
 21" x 225 ft (900-4979) 12 rolls/case



Premium Table Paper Rolls



- Designed and manufactured to provide strength, protection, and value
- White

Crepe, Standard, 125 ft
 18" (100-9865) 12/case
 21" (100-1657) 12/case
 Smooth, Standard, 225 ft
 14" (104-4501) 12/case
 18" (100-9491) 12/case
 21" (100-0177) 12/case
 Smooth, Heavy Duty, 225 ft
 18" (104-6092) 12/case
 21" (104-5273) 12/case



Poly-Perf® Exam Table Paper



- Absorbent top tissue layer with poly back for fluid barrier
- For exam table coverage; can also be used as scale and table liners
- Perforated every 18" or 21"
- 125-ft rolls
- White

18" (570-0728) 9/case
 21" (570-0727) 9/case



Pink Ribbon Table Paper



Smooth, 21" x 225-ft Rolls (900-4276) 12/case



Protection Plus Disposable Underpads



- Made of waterproof, blue poly backing
- Absorbent core of tissue layers
- White, hypoallergenic spunbond facing sheet
- Not made with natural rubber latex

17" x 24" (112-5633) 300/case
 23" x 24" (112-5634) 200/case
 23" x 36" (570-1182) 150/case
 30" x 30" (570-1183) 150/case
 30" x 36" (570-1181) 100/case

EXAM ROOM DISPOSABLES



Everyday Exam Drapes



2-Ply, All-Virgin Tissue
Designed to blanket the patient in lightweight comfort.

White, 36" x 48" (101-3658)	100/case
White, 40" x 48" (100-1797)	100/case
Blue, 40" x 48" (100-3360)	100/case
Mauve, 40" x 48" (100-6428)	100/case
White, 40" x 60" (101-8521)	100/case



Premium Exam Drapes



3-Ply, All-Virgin Tissue
• Soft, opaque, and heavyweight
• Embossed tissue drape with dependable wet strength is ideal for short to moderate exams
• Conforms to the body and offers the ultimate in absorbency

White, 40" x 48" (101-7084)	100/case
--------------------------------------	----------



Pillowcases



Tissue/Poly/Tissue
• Provide linen-like softness for patient comfort and an inner poly barrier for extra protection
• Strong and durable

White, 21" x 30" (104-9206)	100/case
--------------------------------------	----------



Premium Solacel™ Pillowcases



1-Ply Tissue/Poly
• Disposable pillowcase is super soft and strong
• Single-use pillowcase minimizes the risk of cross contamination and infection
• Tissue/poly provides high absorbency and acts as a fluid barrier to protect the pillow

White, 21" x 30" (570-3849)	100/case
--------------------------------------	----------



Pink Ribbon Exam Drapes



2 Ply, 40" x 48"
Provide patient draping for privacy.

(900-4275)	100/case
------------------	----------



Headrest Covers




Tissue/Poly
• Tissue outer surface provides greater comfort for your patients
• Poly lined to help protect your equipment

Small Size 10" x 10"	500/case
-------------------------------	----------

Specify:			
White	(101-7104)	Mauve	(101-3198)
Blue	(101-4595)	Green	(102-8451)
Lavender	(102-5543)		

Regular Size (13" Open Side) 10" x 13"	500/case
---	----------

Specify:			
White	(101-8627)	Mauve	(101-1051)
Blue	(101-8955)	Green	(102-5542)
Lavender	(102-9374)		



Disposable Headrest Covers



Plastic
• Disposable poly barrier that acts as a barrier to protect patient and equipment from cross contamination
• Economical and convenient, space-saver dispenser box

Regular Size 11½" x 10"	250/box
----------------------------------	---------

Specify:			
Clear	(112-4871)	White	(112-4873)

Jumbo Size 14½" x 10"	250/box
--------------------------------	---------

Specify:			
Clear	(112-4872)	White	(112-4874)



Everyday Exam Gowns



3-Ply Tissue, 30" x 42"

- Soft tissue provides patient comfort
- Embossed pattern provides extra strength
- Plastic ties attached

White (100-4471)	50/case
Blue (100-4975)	50/case
Mauve (100-7115)	50/case



Premium Exam Gowns



Tissue/Poly/Tissue, 30" x 42"

- Layer of poly provides strength; layers of tissue absorb moisture
- Opens front or back with attached elastic ties

White (100-2392)	50/case
Blue (100-5728)	50/case
Mauve (100-9652)	50/case



Pink Ribbon Exam Gowns



Tissue/Poly/Tissue, 30" x 42"

- Provides fluid barrier protection and added strength
- Adjustable ties for custom fit

30" x 42" (900-4277)	50/case
-------------------------------	---------



Everyday Exam Capes



3-Ply Tissue, 30" x 21"

- Reliable coverage and strength
- Designed for patient modesty and hygienic protection

White (100-8821)	100/case
Blue (100-2393)	100/case
Mauve (100-2010)	100/case



Premium Exam Capes



Tissue/Poly/Tissue, 30" x 21"

These embossed capes, with a layer of poly film sandwiched between soft absorbent tissue, provide maximum coverage where the exam requires a waist-length garment. The layer of poly provides added strength, making them suitable for most office procedures.

White (101-6908)	100/case
Mauve (101-3029)	100/case
Blue (101-6599)	100/case



Pink Ribbon Exam Capes



Tissue/Poly/Tissue, 30" x 21"

Provides patient modesty for upper-body exams and X-rays.

(900-4278)	100/case
------------------	----------



Disposable Exam Shorts



- Provide protection and security during examination
- Made from SMS material
- Elastic waistband made from latex silk and low-stretch yarn
- Individually wrapped
- Blue

Shorts.....	50/case
-------------	---------

Specify:

Small/Medium (570-2757) 2X-Large (570-2759)

Large/X-Large (570-2758)



Safety Slip Patient Footwear



- Provide slip-resistant, nonskid protection
- Blend of cotton and polyester yarns to provide premium comfort and warmth
- Smooth side, single-tread
- Made to fit left or right foot
- Multiple sizes: all are color coded for easy size identification
- Not made with natural rubber latex

Green, One Size (570-3417)	48 pair/box
-------------------------------------	-------------

Beige, Large (112-6102)	48 pair/box
----------------------------------	-------------

Grey, X-Large (112-6103)	48 pair/box
-----------------------------------	-------------

Mint Green, 2X-Large (112-6105)	48 pair/box
--	-------------



Tray Covers




- Disposable: provide a sanitary means to place instruments on a tray without sterilizing the tray after each patient
- Large variety of sizes: fit all makes and models of trays
- Good quality heavyweight paper: economical
- Not made with natural rubber latex

8½" x 12¼", Ritter B, White
(100-0934) 1000/box

8½" x 12¼", Ritter B
Colors..... 1000/box

Specify:

Aqua	(102-4081)
Beige	(102-8143)
Blue	(101-9336)
Silver-Gray	(101-7620)
Green	(102-7181)
Lavender	(102-4082)
Mauve	(101-4592)
Yellow	(102-4117)



Protex™ Liquid-Proof Tray Covers



- Polycoated, nonseeping paper that prevents cross contamination from paper to paper
- Stops unnecessary waste
- Use poly side down for greater absorbency; use poly side up for quick mixing
- Stops bleed-through
- Economical: only 1 sheet per patient
- Rectangular
- White

8½" x 11"
(100-2254) 500/box

8½" x 12¼"
(100-4202) 500/box

9½" x 12¾"
(100-5098) 500/box

9" x 13½"
(100-5397) 500/box



Plastic Cups



Embossed Drinking Cups with Rolled Lips

- Heavy, sturdy plastic: leakproof and will not crumple in hands
- Releases easily from dispenser: easily accessible
- Rolled lip: easier to drink
- Embossed wall: nonslip
- High-gloss finish and assorted colors: decorative

5 oz..... 1000/case

Specify:

Beige	(570-3734)	Lavender	(570-3728)
Blue	(570-3726)	Peach	(570-3735)
Mauve	(570-3732)	White	(570-3727)
Gray	(570-3731)	Yellow	(570-3733)
Green	(570-3730)		



Paper Cups



Waxed, 5 oz

Pink with Pink Ribbons
(900-4211) 1000/case

White with Pink Hearts
(112-5051) 1000/case



Paper Cups



5 oz, Blue, Polycoated

(949-1498) 1000/case



Graduated Plastic Medicine Cups



Plastic medicine cups with easy-to-read graduations. 1 oz. capacity.

(112-5555) 100/bag



Translucent Plastic Cups



Embossed Drinking Cups with Rolled Lips

- Heavy, sturdy plastic; leakproof and will not crumple in hands
- Releases easily from dispenser; easily accessible
- Rolled lip: easier to drink
- Embossed wall: nonslip

Translucent, 5 oz.
(570-3729) 1000/case



Emesis Bags



- Eliminates exposure to emesis
- Minimizes spillage
- Opening supported by plastic notched ring to provide secure grip and unobstructed access into the bag
- Twist bag and secure notched ring after use
- Holds up to 40 fl oz or 1000 cc/mL
- High-density polyethylene: 1.35-mil thick
- Measures 16.5 cm x 38.1 cm

(570-1178) 24/pkg



Emesis Basins



- Plastic
- Kidney shaped
- 16-oz capacity

Mauve, 8½"
(112-5809) 25/box



Cotton-Tipped Applicators

Applicators are highlighted by a soft cotton tip that is absorbent and will not lint.

Non-Sterile, 3"	(100-6015)	1000/box
Non-Sterile, 6"	(100-9175)	1000/box
Sterile, 6", 100 Packs of 2	(100-9249)	200/box



Wooden Applicator/Mixing Sticks

Plain Wood, 6"	(100-9485)	1000/box
----------------	------------	----------

- International peel symbols: ensure easy opening from proper end of package
- Long shelf life (more than 10 years): economical



Tongue Depressors

- Highly polished and splinter-free
- Smooth
- Designed for oral use or topical application

Non-sterile Junior, 5½"	(100-0758)	500/box
Adult, 6"	(100-2416)	500/box
Sterile Adult, 6"	(100-0238)	100/box



Cotton Balls

Premium quality, absorbent 100% cotton balls designed to effectively prep skin and cleanse wounds. Highest quality, 100% pure-white cotton: soft and absorbent.

Sterile, Medium	(570-3773)	500/box
Non-Sterile, Medium	(570-3771)	4000/case
Non-Sterile, Large	(570-3772)	2000/case



Cotton Bulk Absorbent Rolls, USP

1-lb Bulk Rolls	(570-3769)	ea
Sterile	(570-3770)	ea



Cotton Rolls

Starched, 100% cotton, absorbent nasal roll.

- Starched
- 100% cotton
- Absorbent nasal roll

#2 Medium, Unwrapped, Non-sterile, 1½" x ¾"	(102-5455)	2000/box
#2 Medium, Wrapped, Non-sterile, 1½" x ¾"	(570-3767)	2000/box
#2 Medium, Wrapped, Sterile, 1½" x ¾"	(570-3768)	2000/box



Standard Professional Towels

Specially treated to increase the wet and dry tensile strength of towel. Provides the perfect blend of softness and absorbency. Made in the USA. Not made with natural rubber latex.

3-Ply Paper
13" x 19"..... 500/case

Specify:
White (100-3161) Blue (100-0183)



Dri-Gard® Towels

2-Ply Tissue/1-Ply Poly

- Provide positive moisture barrier and quick absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases the strength of the towel
- Not made with natural rubber latex

13" x 19"..... 500/case

Specify:

Aqua (100-9276)	Lavender (101-0476)
Beige (100-3501)	Peach (100-7384)
Blue (100-7265)	Rose (100-6341)
Gray (100-1313)	White (100-5715)
Green (100-1350)	Yellow (100-9704)



Facial Tissues

Soft, absorbent, value-oriented, white facial tissue: a cost-effective solution for any office. Both patients and employees will notice the softness and quality. 2-ply.

White, 9" x 8"..... 100/box
(900-4214)



Ritter® 225 Barrier-Free® Examination Chair
MIDMARK

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health-care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 225 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

- #225-002, Without Clean Assist Roller System (133-0279)ea
- #225-003, With Clean Assist Roller System (133-0290) ea



463 and 464
Exam/Procedure Chair
MTI

The most feature-packed, ADA compliant* chair on the market. Standard features include ComfortSync™ coordinated footrest, EasyView™ chair-back mounted switches, removable slide-back floating armrests, push-button 333° swivel base and oval articulating headrest.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, recline, and tilt (no tilt on 463)
- 21" seat width
- 18¾"-37" seat height (2" lower than any other ENT chair on the market)
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy, all-steel frame

- #463-002-80, Swivel Base, No Tilt (136-8741) ea
- #464-002-80, Swivel Base, With Tilt (136-8744)ea

*ADA compliant with transfer support accessory for easy patient ingress/egress

Midmark® 626 Barrier-Free® Examination Chair
MIDMARK

Many health care providers are not aware of current standards for proper blood-pressure measurement. If this isn't addressed, it leaves health systems at serious risk. 9.8% of the population is affected by overestimation of high blood-pressure errors. Multiply the cost of treatment for hypertension by the typical patient panel size and the cost of overtreatment could be as much as \$135,000 per each practicing physician. The impact of overmedication can also result in complications such as dizziness, fainting, exhaustion, or even patient falls. The Midmark 626 is the first exam chair of its kind designed to promote proper patient positioning for a more accurate blood-pressure measurement.

- Easily position patients for a better BP measurement
- Powered movement of the back section ensures the patient's back is supported
- A Barrier-Free® low chair-height allows the majority of patients, even female patients with heights in the 3rd percentile, to place their feet flat on the floor

- #626-001, Base Model (133-0301) ea
- #626-002, With Heated Upholstery (133-0312) ea
- #626-003, With IQ Scale (133-0323) ea
- #626-004, With Heated Upholstery and Scale (133-0334) ea



Access High-Low PLUS Power Exam Table
BREWER

Designed to deliver unmatched comfort, convenience, and durability, this table is the ideal choice for modern clinical environments. Engineered to meet the diverse needs of both patients and healthcare professionals with advanced features like ComfortGrip™ patient assist handles and ample storage options. 3-year warranty.

Key features:

- Adjustable height from 18" to 37", fully ADA Access Board compliant
- Extendable table length from 60" to 78" with 18" SafeGlide™ leg extension length
- Mayo Pass-Through Work Surface for multifunctional use
- 6 cubic feet of storage divided between a front drawer (19¾" x 19¼" x 4") and side drawers (22½" x 20¼" x 8")
- 450-lb weight capacity

- #6000PLUS-PR46, Blue Fog (147-3124) ea
- #6000PLUS-PR57, Clamshell (147-3138) ea
- #6000PLUS-PR92, Gunmetal (147-3136) ea



#2135 Revolving Stool Without Back
CLINTON INDUSTRIES

- Nylon casters
- 18"-23" height range
- 14" diameter
- Nylon base
- Pneumatic adjustment
- 4"-thick polyfoam cushioning

(992-0348) ea

Please specify color when ordering: Emerald, Black, Desert Tan, Dove Gray, Mulberry, Burgundy, Royal Blue, or Slate Blue.



425 Air Lift Physician Stool
MIDMARK

Designer Series
Designed for your comfort and support all day, every day. Choose the features you want. Start from the "bottom up" with the frame and seat, then add the back and armrests. Order everything you want at once, or add on later, as you prefer.

(151-4236) ea



Footstool With Hand Rail
CLINTON INDUSTRIES

Chrome-plated steel footstool with 4 nonskid rubber feet.

9"H x 11¼"W x 14¼"L

(766-0390) ea



Midmark® 272 Exam Stool
MIDMARK

Provides optimal comfort, support, and stability. Star-shaped base design allows for closer access to the patient, and stable, easy maneuverability.

- Black, plastic composite, 5-leg caster base and hood
- Pneumatic seat-height adjustment

Exam Stool..... ea

Specify:

#272-001-845, Lunar Gray

(132-2295)

#272-001-814, Stone

(132-2291)



HENRY SCHEIN®
BRAND PRODUCTS

FINGER PULSE OXIMETER

PRECISION AT YOUR FINGERTIPS.

Durable and compact, sized to fit both children and adults.

Oxygen Saturation (SpO₂) Pulse Rate (PR) Pulse Waveform



SpO₂ 70%-100%
Pulse rate measurement range 30-235 bpm

Henry Schein® Finger Pulse Oximeter..... (570-3222)

Each includes: (1) Finger Pulse Oximeter Unit, (1) Lanyard, (1) User Manual, (2) AAA Batteries



Order today at HENRYSCHEIN.COM/MEDICAL

EXAM ROOM EQUIPMENT



Deluxe Sphygmomanometer TECH-MED

Each features a calibrated nylon cuff with a hook-and-loop-type closure. Bulb has deluxe air release valve with built-in filter and extra-large knurled release knob. Two-tube latex bladder. Precision manufactured manometer measures blood pressure from 20 mmHg to 300 mmHg. Lifetime calibration warranty.

Adult, Blue
(100-2396) ea
Contains: 1 of each: sphygmomanometer & instructional manual.
Child, Black
(849-5460) ea.



Mayo Instrument Stand TECH-MED

Ideal under low-clearance units. 13" x 19" removable stainless steel tray. Adjustable from 33" to 52" high. Chrome-plated swivel casters for easy maneuverability.

#4352, Complete
(100-0589) ea



Upholstered Blood Drawing Chair TECH-MED

- High-quality chair with tool-free assembly and padded seat for ultimate patient comfort
- One control knob allows the arms to be raised and lowered effortlessly
- Optional flip-up arm can be adjusted to varying angles with a single knob
- Arms can be raised and lowered in infinite positions and rotate left or right
- Adjustable levelers on all 4 legs allow the chair to sit solidly on an uneven floor
- Vinyl is mildew-resistant, sulfide-stain resistant and static dissipative
- Weight capacity: 500 lb
- Overall dimensions: 38"H x 34"W x 28"D
- Chair dimensions: Seat height to floor: 19", inside arm width: 21", seat depth: 17", height of adjustable arm: 28"-35"

#4382-F
(116-8624) ea



Disposable Stethoscopes DUKAL CORPORATION

Single-Use
Single-head stethoscopes, each with 22" tubing. Not made with natural rubber latex. Individually wrapped.

(138-1784) 10/box



Tech-Med® Mayo Stands and IV Stands



Available in Chrome and Stainless Steel

Mayo Stands

Mayo stands with high-grade stainless steel trays, adjustable heights for versatile use, and bases designed to minimize interference with other equipment. Available with secure locking casters for mobility and weight capacities ranging from 11 lbs to 50 lbs.

HS Item #	Description
(100-0589)	Single Post
(110-6607)	California Base
(122-6910)	Large Tray, California Base, Dual Leg
(100-9541)	Single Post, 5 Caster Base

IV Stands

IV stands feature highly polished chrome or stainless steel poles, adjustable from 47" to 84" with locking clutch collars for secure height positioning. Available with a stable four- or five-caster base for smooth, easy mobility.

HS Item #	Description
(138-7364)	4 Hook, 5 Leg Base, Stainless Steel
(635-4559)	4 Hook, 5 Leg Base, Chrome
(100-5982)	2 Hook, 4 Leg Base, Chrome
(635-5666)	4 Hook, 5 Leg Base, Chrome



Ritter 250-003 LED Exam Light MIDMARK

Features LED technology for longer life, brighter intensity and reduced energy consumption. Provides up to 25,000 hr of lighting—more than 12 times longer than a traditional halogen exam-room light. Light head size: 4½" diam. x 4" long. Height: 36" long, flex length: 22".

- Light source: 12V, 3 LEDs, MR16 type (with glass cover). Illumination @ 18": 1,800 fc. Color temperature: 4,000°K. Life: 25,000 hours
- Electrical: 115V AC, 50–60 Hz, 0.12A. Power cord: 8 ft hospital-grade plug
- Mechanical: Light head size: 4"L x 4½" diameter
- Flex arm: 22"L x 0.585" diameter. Neck material: high carbon steel covered in vinyl
- Base tube: 36"L. Outer tube diameter: 1". Tube material: steel
- Standard base housing: 2¼"H x 13" diameter
- Weight: 12 lb
- Color: Pebble Gray

#250-003, LED Exam Light
(132-3839) ea

Clarus 4™ LED Exam Light MTI

MTI's most economic exam light uses 4 LEDs and is equipped with a high-quality 3" light head. Small light head and flexible neck allow easy movement to many efficient and comfortable positions. Can be used with either a caster base or floor stand. Physical dimensions: head diameter: 3"; caster base height: 73", floor stand height: 71"; caster base diameter: 23¾", floor stand diameter: 15". 2-year limited warranty.

Features:

- 4 white LEDs
- LED life: 40,000 hours
- Color temperature (CCT): 4,300° K
- Light intensity @ 16": 14,100 lux
- Light intensity @ 12": 21,500 lux
- Working distance: 12"
- Spot size at 10": 2"
- Spot size at 12": 3"
- Spot size at 16": 4"
- Spot size at 24": 5"
- CRI: 75
- Control box is dimmable on light head
- Electrical input: 100–240 V, 50–60 Hz, 40W, 12V @ 3.3A
- Power cord: 3-wire grounded cord set, 8' long
- Listings: UL 153, UL 8750, CAN/CS: C22.2, No. 12

Exam Light with Caster Base
(121-8209) ea

Exam Light with Floor Stand
(121-8210) ea



Welch Allyn Green Series™ Exam Light IV BAXTER

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas, the GS Exam Light IV provides you with 3 times the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. And with a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue color rendition during an exam. Utilizing white LEDs, the GS Exam Light IV provides a whiter, brighter light than halogen and consumes less energy. 5-year limited warranty to include the LED light source and power supply only.

- One LED
- Flexible pipe
- Touchless on/off
- Easy to assemble

Light with Table/Wall Mount
(114-2547) ea

Mobile Stand
For great office mobility and space usage. Attaches easily to unit.
(114-4177) ea

Disposable Light Sheaths
(566-0989) 25/box



Welch Allyn Green Series™ 300 General Exam Lights BAXTER

Utilizing a single white LED, the GS 300 provides a whiter, brighter light than halogen and consumes less energy. It's the perfect combination of intense light output and cool operation in a compact design. Features touchless on/off, multiple mounting options for floor space maximization, low cost of ownership due to elimination of bulb replacement and lower energy consumption, and 5-year limited warranty to include the LED light source and power supply.

Light with Mobile Stand
(114-2692) ea

Light with Table/Wall Mount
(114-5632) ea



Clarus 7™ LED Exam Light MTI

MTI's brightest exam light uses 7 LEDs and has a slightly larger 3½" light head, providing maximum illumination and operational ease. Can be used with either a caster base or floor stand. Physical dimensions: head diameter: 3½"; caster base height: 73½"; floor stand height: 71½"; caster base diameter: 23¾", floor stand diameter: 15". 2-year limited warranty.

- LED life: 40,000 hours
- Color temperature (CCT): 4,300° K
- Light intensity @ 16": 33,400 lux
- Light intensity @ 12": 48,500 lux
- Working distance: 16"
- Spot size at 10": 5"
- Spot size at 12": 6"
- Spot size at 16": 8"
- Spot size at 24": 11"
- CRI: 75
- Control box is dimmable on light head
- Electrical input: 100–240 V, 50–60 Hz, 20W, 12V @ 1.67A
- 8-ft power cord: 3-wire grounded cord set
- Listings: UL 153, UL 8750, CAN/CSA: C22.2, No. 12

Light with Caster Base
(121-8211) ea

Light with Floor Stand
(121-8212) ea





Green Series™ 777 Integrated Wall System



All models include a BP aneroid, an ear specula dispenser, and a SureTemp Plus thermometer.

- #777-PM2WAS-HS, With PanOptic Basic LED Ophthalmoscope and MacroView Basic LED Otoscope (570-2666) ea
- #777-PM3WAS-HS, With PanOptic Plus LED Ophthalmoscope and MacroView Plus LED Otoscope for iExaminer (570-2665) ea
- #777-SM2WAS-HS, With Coaxial LED Ophthalmoscope and MacroView Basic LED Otoscope (570-2667) ea
- #777-SM2WAX-HS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser (Does not include SureTemp Plus Thermometer.) (570-2668) ea



The Essentials Series Scope



- Lightweight aluminum chest-piece construction
- Nonchill diaphragm rim for patient comfort
- Proprietary diaphragm enhances acoustic response by 50% over traditional designs
- Adjustable aluminum binaurals with PVC ear tips
- Flexible 22" PVC tubing
- Spare pair of black mushroom ear tips is included
- Inspected and packaged in the USA
- Not made with natural rubber latex
- 1-year warranty

Pediatric, Black (900-4813) ea



Fingertip Pulse Oximeter



- Dual-color LED makes monitoring easy
- Display modes include SpO₂ (range 70%–100%) and pulse rate (range 30–235 bpm) along with pulse bar and plethysmogram
- 2 AAA batteries provide 30 hours of continuous use
- Visual and audible low-battery alarm
- Suitable for both adult and pediatric patients
- Electromagnetic compatibility

(570-3222) ea

Pro Plus Series Clockface Aneroid




Each aneroid features:

- Large, easy-to-read 6" luminescent dial for high visibility in virtually any light condition
- 300mm-Hg, CE marked, no-pin stop movement withstands the heaviest use
- Inflation system includes black nylon cuff, latex inflation bladder, bulb and valve, and 8' length (when extended) of coiled latex tubing
- Contemporary Euro-style ivory casing complements any decor
- Complete with all necessary hardware and assembly instructions
- Available in wall or mobile configurations
- Assembled, inspected, and packaged in the USA from foreign components
- Not made with natural rubber latex
- Wall Aneroid also includes:
 - Wall bracket that swivels through 110° rotation for easy viewing with integral cuff storage
- Mobile Aneroid also features:
 - Height-adjustable Spider-Leg™ caster stand with range of 34" to 48" (from bottom of base to top of unit)
 - Integral ABS cuff-storage compartment
 - Manometer is angled back 15° for easy viewing
 - 18" base diameter

Wall, Adult, Black (112-6069) ea

Mobile, Adult, Black, (112-6070) ea





Welch Allyn SureTemp® Plus 690 Electronic Thermometer
BAXTER

Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes fast, 4–6-second oral temperatures, 10-second pediatric axillary and rectal temperatures, and 15-second adult axillary temperatures.

- Accurate to within ± 0.2°F (0.1°C)
- Interchangeable, removable probe wells to reduce the risk of cross contamination
- Large LCD that can display temperatures in Fahrenheit or Celsius
- Last-temperature recall
- Mode selection and recall buttons
- Four different temperature modes
- Convenient storage for 25 probe covers
- Takes about 6000 readings on 3 AA batteries (included)

Thermometer with Oral Probe, Well, and 4-ft Cord (203-3856).....	ea
Thermometer with Wall Mount, Oral Probe, Well, and 4-ft Cord (178-8374).....	ea
Thermometer with Wall Mount, Oral Probe, and 9-ft Cord (454-2689).....	ea

Probe Covers (566-0460).....	250/sleeve
Wall Holder Mount for Thermometer (206-8985).....	ea
Oral/Axillary Probe, Well, and 9-ft Cord (695-8954).....	ea
Rectal Probe, Well, and 4-ft Cord (741-8045).....	ea



Advantage™ 2200 Fingertip Pulse Oximeter
AMERICAN DIAGNOSTIC CORPORATION

This portable device is perfect for use at home, during sports activities, or in aviation. It offers precise readings with an SpO₂ accuracy of +/-2% and a heart rate accuracy of +/- 5 BPM.

Key features include:

- 2 adjustable display modes
- Bright red LED display
- Compatible with finger thickness of 0.31" to 0.79"
- Includes a carry lanyard and instruction booklet
- Powered by 2 AAA batteries (included)
- 1-year limited warranty

(124-6142).....	ea
-----------------	----

Welch Allyn Braun ThermoScan® PRO 6000 Ear Thermometer
BAXTER

Offers speed and accuracy for all types of patients. Exclusive ExacTemp technology provides an active user feedback system to ensure proper positioning and improved accuracy. Lightweight, simple, fast, and accurate. Preheated probe tip and sensor provide reliable temperature readings. Has automatic probe-cover eject button for quick-and-easy disposal of used probe covers. 3-year limited warranty.



- Innovative PerfecTemp™ technology adjusts for variability in probe placement
- ExacTemp™ technology detects stability of the probe during measurement
- Designed to take quick and accurate readings, to improve patient satisfaction
- 60-second pulse timer assists with manual measurement of pulse rate and respiration
- Electronic security features help prevent theft and loss
- Small or large cradles for probe cover storage that pop up for easy attachment
- Optional charging station stores 200 probe covers, with rechargeable battery
- 6-ft or 9-ft security tether options keep device attached to cradle

#06000-200, Thermometer (124-3558).....	ea
#06000-005, Probe Covers (566-0557).....	20/box



Non-Contact Infrared Thermometer #IR 300

MEDSOURCE INTERNATIONAL

Provides instant temperature readings without worry of physical contact or cross contamination. Measuring distance of 0.4"–2". Temperature displays in <2 seconds. Color digital screen changes color to indicate normal, medium, or high temperature. Automatic power-off in <18 seconds. Memory function records temperature settings for up to 50 groups on 2 AAA batteries. Batteries included. 1-year device warranty.

Model #IR 300 (137-8300).....	ea
-------------------------------	----



Advantage™ 6015N Digital Wrist Blood Pressure Monitor
AMERICAN DIAGNOSTIC CORPORATION

Featuring Smart Logic™ technology and advanced memory capabilities, this monitor measures systolic and diastolic pressure as well as pulse. The Smart Logic™ technology determines the ideal inflation level.

Key features:

- Pressure accuracy: +/- 3 mm Hg
- Pulse range: 40–200 BPM
- Stores the last 99 readings for 2 users
- Wrist cuff fits sizes from 5.3" to 7.7"
- Automatic shut-off feature
- Made without natural rubber latex
- Includes 2 AAA batteries and a compact storage case
- 5-year warranty

(216-1353).....	ea
-----------------	----



seca 777 Physician Scale

SECA

With Eye-Level Display

Combines precision, durability, and safety for high-traffic clinical environments. Featuring a built-in measuring rod with a 4"-90" range, it ensures accurate height and weight readings for a wide range of patients. Designed with an oversized, low-profile, full-metal platform and a tilt-proof safety design, the seca 777 provides stability and reliability without plastic components. Its USB interface enables seamless PC connectivity for efficient data integration. Backed by an 8-year warranty—4 times the industry standard—and compatible with the Panda Bear Bao Bao accessory, this scale is built for long-lasting performance in busy healthcare settings.

- Built-in measuring rod with 4"-90" measuring range
- Oversized, low-profile, full-metal platform (no plastic) with tilt-proof safety design
- 8-year warranty (4 times the industry standard)
- USB interface for PC connectivity
- Compatible with Panda Bear Bao Bao accessory
- Capacity: 550 lb/250 kg
- Graduation: 0.2 lb/100 gm

#7771821004, Scale
(139-2675) ea



Digital Electronic Eye-Height Scale

HENRY SCHEIN®

Compact, eye-height digital scale that features an integrated mechanical height rod, large backlit LCD display, and BMI calculator. Weight capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2½" high. Sta-Bilt stabilizing software ensures an accurate measure, even if the patient is moving. Autohold feature allows the weight to continue to be displayed after the patient has stepped off the scale. Integrated USB interface allows data to be transferred to a computer. Includes 2-year warranty. 22" x 15" x 53"

(112-7154) ea

Eye-Level Digital Physician Scale

HEALTH-O-METER

High-capacity 660 lb/300 kg primary-care scale provides the optimum in quality, reliability, and ease of use.

- Functions: lb/kg conversion, lb/kg lockout, BMI, zero, tare, pre-tare, hold/release, auto zero, auto off, reweigh, recall, and Everlock®
- Low profile, nonskid platform with reinforced base for added stability
- Integrated height rod with 30"-84¼" measuring range
- Wireless connectivity options available

#500KL, Scale
(106-6626) ea

#500HB, Handlebar Accessory
(137-2640) ea



Digital Wheelchair Scale

HENRY SCHEIN®

- Single, stable platform allows scale access with wheelchair
- Easily folded to minimize footprint and moved with casters
- Individual wheelchair weight can be programmed for easy recall
- 32" x 33" x 2½", 1000-lb maximum capacity
- 2-year limited warranty

(570-0597) ea



ZOLL Medical's interactive AEDs guide you through the process of performing high-quality CPR with voice, text, and image prompts. Rigorous, automatic self-testing ensures they are always ready. Each AED has unique features to provide the best solution for your needs.



AED Plus®

- * Real CPR Help® technology
- * Industry leading five-year battery and pad life
- * One-piece pad design

(866-7559) Semi-Automatic Unit
 Package includes: AED, electrodes, batteries, carry case.



Powerheart® G5 AED

- * Intellisense® CPR Feedback (ICPR)
- * Non-polarized pads (can be placed in either position)
- * Dual language functionality

(139-5101) Semi-Automatic Unit
 Package includes: AED, electrodes, battery.



ZOLL R Series®

- * Professional defibrillator/monitor
- * OneStep Complete™ electrodes with integrated three-lead and Real CPR Help® technology

(127-9979) Professional Defibrillator
 Package includes: Monitor/defibrillator, battery, 2 sets of OneStep electrodes.



AED Pro® Automated External Defibrillator
ZOLL

An external defibrillator that provides advanced capabilities for BLS and ALS users with innovative Real CPR Help® technology for rate and depth of chest compressions. See-Thru CPR® capability filters CPR chest compressions, displaying the underlying cardiac rhythm. Unit is unmatched in its ruggedness and durability and is ideal for any environment.

- Available as semiautomatic only or semiautomatic with manual override for ALS professionals
- Built-in internal memory and external USB memory technology allows expanded data storage capacity and maximum flexibility for data management and transfer
- Electrodes and battery compatibility with ZOLL M Series®, E Series®, and R Series® professional defibrillators
- ECG monitoring with 3-lead cable and high-quality display
- Unrivaled battery options (long-life disposable or rechargeable) that meet specific needs

#90110200499991010
(869-6270) ea



Powerheart® G5 AED
ZOLL

Features Intellisense™ CPR, which provides corrective feedback to ensure delivery of high-quality compressions. One-button, dual-language functionality and comprehensive self-tests give rescuers the confidence to act.

- Intellisense™ CPR feedback provides real-time voice and text corrective prompts for CPR quality
- Rescue Ready® conducts regular, automatic self-tests on pads, battery, and internal circuitry
- RescueCoach™ guides users through each critical step of the rescue with user-paced voice/text prompts
- Dual language change from primary to secondary language in the push of a button

#G5S80C-S, Semi-Automatic
(139-5101) ea

#G5A80C-S, Fully Automatic
(139-5102) ea

ZOLL AED Plus®
ZOLL

Bringing Lifesaving Defibrillation to Everyone
 The only AED with CPR feedback and an ECG waveform. By supporting CPR as well as defibrillation, this device can help every unresponsive victim. Features Real CPR Help®, a tool that is able to actually see what user is doing and provides audio and visual prompts to help them do it with confidence and clarity.

- Provides a one-piece electrode for accurate placement
- Uses consumer lithium batteries that are easy to replace, long-lasting, and inexpensive

Semi-Automatic Unit
(866-7559) ea
 Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

Fully Automatic Unit
(119-7910) ea
 Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

CPR-D padz®—Adult Defibrillator Pad
(253-2611) ea

Pedi padz® II—Pediatric Electrode Pad
(203-7054) ea

Type 123 Lithium Batteries
(610-4769) 10/pkg

Soft Carrying Case
(106-1937) ea

3D Wall Sign
(110-5384) ea

AED Plus Simulator
(499-7385) ea



Ambulatory Surgery Center and Clinic



ePM 10M/12M/15M continuous, multi-parameter patient monitors are ideally suited for the ambulatory surgery center where flexibility, efficiency, and value are essential.



The **A4 Advantage** complements clinical workflow with customizable ventilation capabilities and an intuitive touchscreen user-interface. It is well-suited for ambulatory surgery centers and offers outstanding functionality in a comprehensive, cost-effective anesthesia delivery system.

All products offer an industry-leading 3-Year Warranty

To learn more about our entire portfolio of innovative solutions, contact your local Henry Schein Medical representative, or visit www.henryschein.com/medical.



Contact your Henry Schein Rep for more information!

Above the Rest



- (138-9631) X8 Patient Monitor w/ G2
- (138-9695) X8 Patient Monitor w/ Printer
- (138-9709) X8 Patient Monitor w/ Roll Stand
- (138-9710) X8 Patient Monitor w/ Roll Stand and Printer
- (138-9826) X10 Patient Monitor w/ G2
- (138-9825) X10 Patient Monitor w/ Printer
- (138-9827) X10 Patient Monitor w/ Roll Stand
- (138-9828) X10 Patient Monitor w/ Roll Stand and Printer
- (138-8976) X12 Patient Monitor w/ G2
- (138-8900) X12 Patient Monitor w/ Printer
- (138-9992) X12 Patient Monitor w/ Roll Stand
- (138-9993) X12 Patient Monitor w/ Roll Stand and Printer

Edan helps maximize patient care through reliable purchasing and EMR solutions. Edan G2 technology is compatible with generic consumables in order to maximize efficiency across your supply chain. All monitors are out of box ready with accessories.



IQVitals® Zone™ Vital Signs Wireless Monitor MIDMARK

The only wireless vital signs monitor that automatically connects to a computer using zero clicks, importing data directly into the EMR to eliminate manual transcription errors. Using Midmark Zone™ technology, the IQVitals wireless monitor connects to a computer without cables or manual pairing. The Serenity BP™ linear deflation system, designed by Midmark, delivers quick and accurate measurements while ensuring greater comfort for the patient. All functions can be operated at the point of care via the computer or the user-friendly touchscreen. Each IQVitals Zone device includes an Exergen TemporalScanner™ thermometer.



Key features:

- Rechargeable lithium-ion battery
- Non-invasive oscillometric blood pressure method
- Exergen TemporalScanner infrared thermometer
- Masimo or Nellcor SpO₂ measurement options
- 3 mounting options to suit any workflow

- #4-000-0724, Monitor with BP, Temp & Masimo SpO₂ (127-9083) ea
- #4-000-0722, Monitor with BP, Temp & Nellcor SpO₂ (140-5080) ea
- #4-000-0720, Monitor with BP & Temp (141-3660) ea

Connex® Spot Monitor BAXTER

Features an easy-to-use, vivid touch-screen display and provides accurate vital signs measurement, including blood-pressure averaging, spot checking, interval monitoring and custom scoring across patient populations. Device connects wirelessly to an EMR and accurately sends vitals to the patient chart from the point of care to help improve efficiency.



- Designed for the physician-office environment with an easy-to-use, vivid touch-screen display
- More accurate hypertension diagnosis without slowing down your workflow thanks to automated blood-pressure averaging
- Improved capabilities for assessing respiratory conditions such as COPD, asthma, and community-acquired pneumonia (CAP), thanks to integrated professional-grade pulse oximetry
- Increase staff efficiency and eliminate data-entry errors with wireless transfer of vital signs directly into an EMR

- 7100 Base Models with Wired Connectivity
- #71XX-B, SureBP Non-invasive Blood Pressure (566-0464) ea
 - #71XT-B, SureBP, SureTemp (672-0176) ea
 - #71WX-B, SureBP, Nonin SpO₂ (672-0177) ea
 - #71WT-B, SureBP, Nonin SpO₂, SureTemp (566-0466) ea
 - #73WT-B, SurePB, Nonin SpO₂, SureTemp (672-0179) ea
 - #73CT-B, SureBP, Covidien SpO₂, SureTemp (672-0175) ea
 - #73MT-B, SureBP, Masimo SpO₂, SureTemp (566-0467) ea



EXPLORE OUR ROBUST
COLLECTION OF
**PRODUCT
GUIDES**
ONLINE NOW!



www.henryschein.com/productguides




TOPICAL SKIN ADHESIVE — OCTYL

High-Strength 2-Octyl Cyanoacrylate Formulation

Advantages

- Rapid set time
- Flexibility without a plasticizer
- Ergonomic ampule designed for ease of handling and targeted application
- Provides an effective microbial barrier and water-resistant protection
- High tensile strength
- Unique formulation will not solidify in ampule, allowing all .5 mL to be used
- Gravity fed applicator delivery system — no need to squeeze ampule



Additional Information

Polymerization starts when adhesive comes in contact with tissue or fluid.
All 0.5 mL of adhesive can be used on the patient, as adhesive will not solidify in the ampule.

Topical Skin Adhesive will remain attached to the skin surface until natural cellular regeneration occurs (approx. 5 — 7 days). If needed, any standard adhesive remover product may be used to release the adhesive from the skin. Peel off the film, do not pull the skin apart.

Henry Schein® Topical Skin Adhesive — Octyl
0.5mL ampules, Clear, 12/box.....(570-3554)

Order today at HENRYSCHEIN.COM/MEDICAL



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS




TOPICAL SKIN ADHESIVE - BUTYL

Strong Bonding and Reliable Protection

- Sterile, single-use plastic ampules of liquid skin adhesive consisting of n-butyl-2-cyanoacrylate
- Antibacterial barrier protects the wound
- Applicator tips provide precise, controlled application
- Will not clog or dry when opened
- Blue liquid color option available to enable easy control over quantity applied




EASY-TO-USE APPLICATION



Open the package and take out the tip.



Open the ampule by twisting off the ribbed tip.



Apply the tip to opened ampule.



Apply the glue in a thin layer to the wound edges and hold in position for 60 seconds.

Henry Schein® Topical Skin Adhesive

Flexible Ampule with Tip, Blue 0.5 mL, Ea.

(570-0351)

Ampule without Tip, Clear 0.5 mL, 10/box

(112-6106)

Flexible Ampule with Tip, Blue 0.5 mL, 10/box

(112-7158)

Ampule without Tip, Blue 0.5 mL, 10/box

(112-6111)

Order today at HENRYSCHEIN.COM/MEDICAL



Surgical Sutures

From A Name You Can Trust.

Maxima® Sutures by Henry Schein®


RACETRACK FEATURES AND BENEFITS

POLYLITH CARD

Helps reduce moisture content, an advantage with absorbable suture materials that degrade by hydrolysis.

NEEDLE PARK

Two needle parks to accommodate a wide range of needle types, sizes, and double-arm sutures. Needle park also protects needles from damage.



LARGER DIAMETER

Larger overall diameter to ensure less suture memory and minimize the number of suture rotations around track.

PRESS-PUSH NEEDLE LIFT

The tab, when pressed from behind, lifts for easy needle removal.



- Silicone-coated 300 series stainless needles give finer point geometry for superior penetration
- Precision needles are strong enough to resist bending, yet ductile enough to resist breaking
- Racetrack packaging helps reduce material memory as the suture is released from the track
- Large breadth of suture and needle combination for a range of surgeries
- Significant cost savings for your clinic compared to the leading brand

Precision Reverse Cutting Needles			PGA ABSORBABLE SUTURE					
			Material	Color	Length	5-0	4-0	3-0
3/8 circle 13mm	NP-3		Braided	Undyed	18" 45 cm	570-1048 12/box	570-1049 12/box	
3/8 circle 19mm	NPS-2		Braided	Undyed Violet	18" 45 cm		570-1045 12/box 570-1399 12/box	570-1046 12/box
3/8 circle 19mm	NFS-2		Braided	Undyed Violet	30" 75 cm	570-1397 12/box	570-1041 12/box	570-1047 12/box

Precision Reverse Cutting Needles			NYLON NON-ABSORBABLE SUTURE					
			Material	Color	Length	5-0	4-0	3-0
3/8 circle 13mm	NP-3		Monofilament	Black	18" 45 cm	570-1039 12/box	570-1040 12/box	
3/8 circle 19mm	NPS-2		Monofilament	Black	18" 45 cm		570-1042 12/box	570-1044 36/box
3/8 circle 19mm	NFS-2		Monofilament	Black	18" 45 cm	570-1398 12/box	570-1043 12/box	
3/8 circle 24mm	NFS-1		Monofilament	Black	18" 45 cm			570-1050 12/box

Precision Reverse Cutting Needles			POLYPROPYLENE NON-ABSORBABLE SUTURE					
			Material	Color	Length	5-0	4-0	3-0
3/8 circle 13mm	NP-3		Monofilament	Blue	18" 45 cm	570-1591 12/box	570-1593 12/box	
3/8 circle 19mm	NPS-2		Monofilament	Blue	18" 45 cm		570-1592 12/box	570-1594 36/box
3/8 circle 19mm	NFS-2		Monofilament	Blue	18" 45 cm		570-1595 12/box	



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



SUPERIOR POINT

Needles are honed with a unique multiedged point to provide greater ease of penetration.

SUPERIOR ATTACHMENT

Needles are drilled for a seamless needle-suture juncture, providing a smoother transition with less tissue disruption.

SUPERIOR PACKAGING

Unique double-sterile barrier packaging ends costly waste.



NON-ABSORBABLE SUTURES 12/box


MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Silk Black Braided	C-9 24mm 1/2 Circle	18"/45cm				632G* (100-6320)	
	C-43 15mm 1/2 Circle	18"/45cm				(100-7762)	
	C-3 13mm 3/8 Circle	18"/45cm		640G* (101-4148)	641G* (101-2590)		
	C-6 18mm 3/8 Circle	18"/45cm		682G* (101-3207)	683G* (100-6830)	1679H* (100-7772)	
	C-7 24mm 3/8 Circle	18"/45cm				684G* (101-2049)	
	C-22 16mm 3/8 Circle	18"/45cm	(101-2636)				
Nylon-Black Monofilament	C-1 11mm 3/8 Circle	18"/45cm	697G* (900-7482)				
	C-3 13mm 3/8 Circle	18"/45cm	1698G* (900-7473)	698G* (900-7470)	699G* (900-7466)		
	C-6 18mm 3/8 Circle	18"/45cm		661G* (101-7137)	662G* (101-1209)		
		10"/25cm			(900-7464)		
	C-7 24mm 3/8 Circle	18"/45cm				663G* (101-3728)	
	C-8 26mm 3/8 Circle	18"/45cm					664G* (900-7476)
	C-17 12mm 3/8 Circle	10"/25cm	(900-7460)	(900-7459)	(900-7461)		
	C-22 16mm 3/8 Circle	18"/45cm	1665G* (101-6053)				
	PC-31 18mm 3/8 Circle	18"/45cm		1666G* (900-7472)	1667G* (900-7487)		
	DSCM-16 16mm 3/8	18"/45cm		(570-0458)	(570-0459)		
	DSM-19 18mm 3/8 Circle	18"/45cm				(570-0479)	
DSM-24 24mm 3/8 Circle	18"/45cm				(570-0480)		
Polypropylene-Blue Monofilament	C-1 11mm 3/8 Circle	18"/45cm	8697G* (900-7483)				
	C-3 13mm 3/8 Circle	18"/45cm	8695G* (900-7479)	8698G* (900-7484)	8699G* (900-7486)		
	C-6 18mm 3/8 Circle	18"/45cm		8661G* (101-6409)	8683G* (101-7477)		
	C-7 24mm 3/8 Circle	18"/45cm				8684G* (101-7700)	
	C-22 16mm 3/8 Circle	18"/45cm	(101-4378)				
	DS-18 18mm 3/8 Circle	18"/45cm				8665G* (900-7481)	
	PC-31 18mm 3/8 Circle	18"/45cm			8682G* (900-7478)		
	DSM-16 16mm 3/8 Circle	18"/45cm		(570-0471)			
	DSM-19 18mm 3/8 Circle	18"/45cm		(570-0468)		(570-0467)	
	DSM-24 24mm 3/8 Circle	18"/45cm				(570-0470)	
	DSM-26 26mm 3/8 Circle	18"/45cm					(570-0473)
	DSCM-13 13mm 3/8	18"/45cm	(570-0472)				
	DSCM-16 16mm 3/8	18"/45cm		(570-0469)			
Polyester-Green Braided	D-14 17mm 1/2 Circle	18"/45cm		X951H* (100-7773)			

MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Plain Gut	C-9 24mm 1/2 Circle	18"/45cm				612G* (100-6120)	
	C-3 13mm 3/8 Circle	18"/45cm		686G* (101-3192)			
	C-6 18mm 3/8 Circle	18"/45cm			(100-7771)		
	C-6 18mm 3/8 Circle	27"/70cm			H821H* (101-0668)	H822H* (101-2976)	
	DSM-13 13mm 3/8 circle	18"/45cm	1916G* (900-7485)	1915G* (900-7471)			
Chromic Gut	C-3 13mm 3/8 Circle	18"/45cm		687G* (101-8824)			
	C-6 18mm 3/8 Circle	18"/45cm			1637G* (900-7637)	1638H* (900-7636)	
	C-6 18mm 3/8 Circle	27"/70cm			635H* (101-2787)	636H* (101-1582)	
	HR-26 26mm 1/2 Circle	30"/76cm				G122H* (570-0466)	
	DSM-13 13mm 3/8 Circle	18"/45cm			1654G* (570-0475)		
	DSM-19 19mm 3/8 Circle	18"/45cm			1637G* (570-0476)		
Polyglycolic Acid (PGA) Undyed-Braided	C-3 13mm 3/8 Circle	18"/45cm	J492G* (900-7468)	J493G* (900-4417)	J494G* (900-4418)		
	C-6 18mm 3/8 Circle	27"/70cm			J422H* (900-7463)	J423H* (900-7462)	
	PC-30 18mm 1/2 Circle	18"/45cm			J507G* (900-7477)		
	PC-31 18mm 3/8 Circle	18"/45cm			J496G* (900-4419)	J497G* (900-4421)	
	DSM-24 24mm 3/8 Circle	27"/70cm					(570-0482)
	HR-26 26mm 1/2 Circle	27"/70cm				(570-0481)	
Polyglycolic Acid (PGA) Undyed-Braided-Fast Absorbing	DSM-13 13mm 3/8 Circle	18"/45cm		VR493* (570-0477)	VR494* (570-0478)		
Polyglycolic Acid (PGA) Violet	C-6 18mm 3/8 Circle	18"/45cm			J392H* (900-4420)		
	DSM-19 19mm 3/8 Circle	18"/45cm				J393H* (570-0474)	
PGA PCL Undyed-Mono	DSM-13 13mm 3/8 Circle	18"/45cm		Y493G* (900-7467)	Y494G* (900-7465)		
	DSM-19 19mm 3/8 Circle	18"/45cm			Y496G* (900-7488)	Y497G* (900-7474)	
	DSM-13 13mm 3/8 Circle	18"/45cm	Y492G* (570-0465)				
	DSM-19 19mm 3/8 Circle	27"/70cm			(570-0460)	(570-0461)	
	DSM-24 24mm 3/8 Circle	27"/70cm				(570-0463)	
PGA PCL Violet-Mono	DSM-13 13mm 3/8 Circle	18"/45cm		Y463G* (570-0462)	(570-0464)		


*Comparable Ethicon™ material

NEEDLE CHART

REVERSE CUTTING NEEDLES



PRECISION REVERSE CUTTING NEEDLES



PRECISION CONVENTIONAL CUTTING NEEDLES



TAPER CUTTING NEEDLES



HENRY SCHEIN® *Maxima™* PTFE SURGICAL SUTURES

- Maxima PTFE (polytetrafluoroethylene) sutures are an excellent choice for dental bone grafting and implant procedures.
- Proprietary monofilament construction reduces the potential for bacterial migration into the surgical site.
- Ultrasoft, high-density material is designed with a low surface porosity.
- Black premium needles are made from 300-series stainless steel.



ITEM CODE	OSTEOGENICS CYTOPLAST**	NEEDLE	SIZE	LENGTH	RECOMMENDED USAGE	12/BOX
(112-7253)	CS0418	FS-2 (19mm, 3/8 reverse cut)	2-0	18	Larger suture size with a popular needle for procedures where a longer needle is desired	
(112-7254)	CS051819	FS-2 (19mm, 3/8 reverse cut)	3-0	18	A popular suture size and standard sized needle for general dental implant and bone grafting	
(112-7255)	CS0518	FS-3 (16mm, 3/8 reverse cut)	3-0	18	A popular suture size and smaller needle for dental implant and bone grafting	
(112-7256)	CS0618 RC	FS-3 (16mm, 3/8 reverse cut)	4-0	18	Most popular suture size for dental implant and bone grafting	
(112-7257)	CS0618 PERIO	TF (13mm, 1/2 taper point)	4-0	18	Most popular suture size used for soft tissue grafts and delicate tissues requiring taper point needle	
(112-7258)	-	PS-2 (19mm 3/8 premium reverse cut)	4-0	18	Most popular suture size with the most popular sized premium needle	
(112-7259)	CS0618 PREM	P-3 (13mm, 3/8 premium reverse cut)	4-0	18	Most popular suture size used for soft tissue grafts and delicate tissues requiring cutting point needle	
(112-7260)	-	P-3 (13mm, 3/8 premium reverse cut)	5-0	18	Small suture size and premium needle option for delicate tissues to reduce tissue trauma	
(112-7261)	-	P-3 (13mm, 3/8 premium 6-0 reverse cut)	6-0	18	Smaller suture size and premium needle option for delicate tissues and the least amount of tissue trauma	

**Cytoplast is a registered trademark of Osteogenics.

				
FS-2 Reverse Cutting 19mm 3/8 Circle	FS-3 Reverse Cutting 16mm 3/8 Circle	P-3 Premium Reverse Cutting 13mm 3/8 Circle	PS-2 Premium Reverse Cutting 19mm 3/8 Circle	TF Taper Point 13mm 1/2 Circle

Made in USA

Wound Closure Strips

- Adhesive wound closures with an exceptionally porous, spun-bond fiber backing include bilateral nylon filaments for high tensile strength
- To be used in addition to or in place of sutures
- Each size is in a color-coded package for quick and easy size identification
- Medical-grade adhesive properties: Affords excellent adhesion, resistance to moisture, superior longevity, and is gentle to the skin
- Hypoallergenic: Designed for patient comfort and longer-term exposure
- Not made with natural rubber latex

- 1/8" x 3"250/box
Contains: 5 strips per envelope, 50 envelopes per box.
- 1/4" x 3"150/box
Contains: 3 strips per envelope, 50 envelopes per box.
- 3/8" x 4"500/box
Contains: 10 strips per envelope, 50 envelopes per box.



- 1/4" x 1 1/2"300/box
Contains: 6 strips per envelope, 50 envelopes per box.
- 1/2" x 4"300/box
Contains: 6 strips per envelope, 50 envelopes per box.
- 1" x 4"100/box
Contains: 4 strips per envelope, 25 envelopes per box.



Disposable Skin Stapler


- Enables user to work with higher efficiency and assists in avoiding cross infection effectively
- Specially designed for surgical skin suturing
- Novel structure and application security allow quick speed when sewing up during operation, modest histological reaction, and good wound matching and healing
- No pain during removal
- Scarring is minimal

- Skin Stapler, 35 Wide (570-2956)ea
- Skin Stapler, 35 Regular (570-2957)ea
- Skin Staple Remover (900-7011)ea

Johnson & Johnson
MedTech

Every closure counts

Make the most of every closure with
J&J MedTech's advanced wound
closure portfolio.



The third-party trademarks used herein are trademarks of their respective owners.
© Johnson & Johnson and its affiliates 2025. US,SRG,WOUN,xxxxxx



STRATAFIX™ Spiral Knotless Tissue Control Devices

J&J MEDTECH SURGERY

A suturing device that features barbs in a spiral pattern. Available in both bidirectional and unidirectional designs. Barbs allow for tissue approximation without the need to tie surgical knots. Antibacterial technology is available for both designs.

#SXMP1B103, 3-0 Monocryl™ Plus, Undyed, 60 cm, PS-1	(131-9620)	12/box
#SXMP2B412, 3-0 Monocryl™ Plus 12", PS Double Arm, 26 mm	(141-9843)	12/box
#SXPP1B113, 4-0 PDS Plus Violet, 45 cm, PS-2	(131-3098)	12/box
#SXPP2B405, 1 PDS Plus Bidirectional Violet, 14" x 14"	(145-6625)	12/box

STRATAFIX™ Symmetric PDS™ Plus Knotless Tissue Control Device

An antibacterial, absorbable, barbed suture.

#SXPP1A400, 1 Violet, 45 cm, CTX	(128-5753)	12/box
#SXPP1A404, 1 Violet, 45 cm, CT-1	(124-4042)	12/box
#SXPP1A401, 0 Violet, 45 cm, CT-1	(131-5099)	12/box
#SXPP1A445, 1 Violet, 60 cm, CTX	(141-7182)	12/box

ACE® Brand Elastic Bandages

3M

5-yd Rolls

- The elastic has a soft feel, so the elastic bandage wrap is comfortable against your skin and can be washed for reuse
- Intended for reuse, this comfortable wrap conforms to multiple body parts and fits discreetly under clothing to help provide support and customized compression where needed
- Designed with adjustable elastic bandage clips that snap in place and hold securely

2"	(987-0559)	10/box
3"	(987-1590)	10/box
Clip Closure		
#207433, 4"	(847-6319)	10/box
6"	(987-2908)	10/box



DERMABOND™ PRINEO™ Skin Closure System

J&J MEDTECH SURGERY

Combines a sterile, 2-octylcyanoacrylate liquid topical skin adhesive with a self-adhering mesh (not made with natural rubber latex) for closure of surgical incisions and simple lacerations. A noninvasive alternative to skin sutures and staples.

#CLR222US, 22 cm	(122-7672)	2/box
#CLR422US, 42 cm	(136-5543)	2/box
#CLR602US, 60 cm	(116-5081)	2/box

Contains: 2–3.8-mL pens & 2 self-adhering mesh strips.

DERMABOND™ ADVANCED™ and DERMABOND™ Mini Topical Skin Adhesive

DERMABOND™ ADVANCED™ is a unique formulation with an easy-to-use applicator for precise control for wide- and fine-line applications including surgical incisions and lacerations—up to 15 cm long. DERMABOND™ Mini is for smaller incisions and lacerations—up to 4 cm long.

#DNX6, 0.7 mL, DERMABOND ADVANCED	(555-0357)	6/box
#DNX12, 0.7 mL, DERMABOND ADVANCED	(555-0356)	12/box
#DHVM12, 0.36 mL, DERMABOND Mini	(654-0034)	12/box



Nexcare™ Waterproof Bandages 3M

- Offer protection against water, dirt, and germs
- Clear, breathable material
- Ultrathin and comfortable to wear
- Shape designed for better seal around pad
- Easy, one-handed application

2¼" x 1⅞"	(777-2106)	20/box
3½" x 2⅜"	(777-0355)	8/box
Assorted		
(987-0457)		50/box



Ethicon® Plus Antibacterial Sutures J&J MEDTECH SURGERY

A comprehensive portfolio for every need. Available in 3 absorbable polymer sutures: PDS (polydioxanone), Monocryl™ (poliglecaprone 25), and Vicryl™ (polyglactin 910). Choose from a variety of sizes, lengths, and needle types, including multifilament braided, monofilament, and barbed options.

#PDP497G, 3-0 PDS Plus Undyed, 1" x 18", PS-2	(654-0015)	12/box
#PDP148H, 5-0 PDS Plus Violet, 1" 30", RB-2, Double Armed	(128-5177)	36/box
#MCP427H, 3-0 Monocryl™ Plus Undyed, 1" x 27", PS-2	(654-9300)	36/box
#MCP494G, 4-0 Monocryl™ Plus Undyed, 1" x 18", P-3	(654-1377)	12/box
#VCP393H, 3-0 Coated Vicryl™ Plus Undyed, 1" x 18", FS-2	(110-5936)	36/box
#VCP416H, 3-0 Coated Vicryl™ Plus Undyed, 1" x 27", SH	(840-0998)	36/box
#VCP423H, 3-0 Coated Vicryl™ Plus Undyed, 1" x 27", FS-2	(775-2716)	36/box
#VCP422H, 4-0 Coated Vicryl™ Plus Undyed, 1" x 27", FS-2	(175-9721)	36/box

Visit www.henryschein.com/medical to view our full portfolio of products.



Nexcare™ Opticlude™ Orthoptic Eye Patch 3M

Gentle to skin. Ideal for treating amblyopia (lazy eye) and strabismus (crossed eye). Highly breathable backing, comfortable to wear. Hypoallergenic adhesive. Not made with natural rubber latex. Beige.

Regular, 2¼" x 3¼"	(777-1590)	20/box
Junior, 1⅞" x 2½"	(777-0513)	20/box



On the path to new packaging



3M's healthcare business became a standalone company in 2024, but we're still updating the name on all of our products.

The products you rely on are in the process of transitioning to reflect our new company name and packaging. You'll still get the products you love, just with a fresh green look.

© Solventum 2026. All rights reserved.

3M™ Coban™ NL Non-Latex Self-Adherent Wrap SOLVENTUM, FORMERLY 3M HEALTH CARE

5-yd Rolls

- Comfortable: soft, breathable, lightweight, and conformable
- Self-adherent: will not stick to skin, hair, clothing, or drapes
- Residue-free: gentle on fragile/thin skin
- Slip-minimizing cohesion: ideal for active patients
- Available in sterile and nonsterile versions and in a variety of widths and 2 colors

Non-Sterile

#2081T, Tan, 1" Width (142-4831)	72/case
#2081, Blue, 1" Width (142-4830)	72/case
#20815B, Blue, 1½" Width (142-4839)	48/case
#2082, Tan, 2" Width (679-9575)	36/case
#2082B, Blue, 2" Width (142-4847)	36/case
#2083, Tan, 3" Width (395-9243)	24/case
#2084, Tan, 4" Width (313-6816)	18/case
#2086, Tan, 6" Width (608-8687)	12/case

Sterile

#2082S, Tan, 2" Width (777-0573)	36/case
#2083S, Tan, 3" Width (102-5030)	24/case
#2084S, Tan, 4" Width (302-8421)	18/case
#2086S, Tan, 6" Width (777-3803)	12/case

3M™ Transpore™ Surgical Tape SOLVENTUM, FORMERLY 3M HEALTH CARE

- Easy to handle with gloves
- Good adhesion to skin and tubing
- Transparent
- Water-resistant
- Hypoallergenic: not made with natural rubber latex

Standard Rolls, 10 yd

#1527-1, 1" (777-7305)	12/box
#1527-0, ½" (777-9735)	24/box
#1527-2, 2" (777-9422)	6/box
#1527-3, 3" (635-1724)	4/box

Single-Use Rolls, 1½ yd

#1527S-1, 1" (777-2577)	100/box
#1527S-2, 2" (123-8478)	50/box



3M™ Medipore™ H Soft Cloth Surgical Tape SOLVENTUM, FORMERLY 3M HEALTH CARE

- Conformability: bidirectional stretch accommodates swelling, distension, and movement
- Easy to tear at perforations: no scissors necessary
- Not made with natural rubber latex

Standard Rolls, 10 yd

#2861, 1" Width (887-7957)	24/box
#2862, 2" Width (139-3011)	12/box
#2863, 3" Width (777-3265)	12/box
#2864, 4" Width (471-1595)	12/box
#2866, 6" Width (328-0293)	12/box
#2868, 8" Width (110-8807)	6/box

Single-Use Rolls, 2 yd

#2860S-1, 1" Width (137-6301)	72/box
#2860S-2, 2" Width (137-6302)	48/box
#2860S-4, 4" Width (137-6303)	24/box
#2860S-6, 6" Width (147-3865)	16/box



3M™ Durapore™ Medical Tape SOLVENTUM, FORMERLY 3M HEALTH CARE

- Offers strong adhesion: good initial and long-term adhesion to dry skin
- Easy to use: can be torn lengthwise and crosswise to create a variety of sizes as needed
- Adheres to dry skin, tubes, devices, and itself
- Helps reduce risk: individually packaged single-use length rolls can help protect patients by reducing the risk of cross contamination
- Not made with natural rubber latex

Standard Rolls, 10 yd

#1538-0, ½" Width (777-6047)	24/box
#1538-1, 1" Width (777-2892)	12/box
#1538-2, 2" Width (777-0941)	6/box
#1538-3, 3" Width (777-9487)	4/box

Single-Use Rolls, 1½ yd

#1538S-1, 1" Width (777-0731)	100/box
#1538S-2, 2" Width (777-2160)	50/box



3M™ Micropore™ Surgical Tape SOLVENTUM, FORMERLY 3M HEALTH CARE

- Provides good initial and long-term adhesion to skin
- Conformable, breathable paper tape not made with natural rubber latex
- Can be torn by hand
- Provides multi-purpose securement for clinical applications

Standard Rolls, 10 yd, White

#1530-0, ½" (777-5044)	24/box
#1530-1, 1" (777-6080)	12/box
#1530-2, 2" (777-9382)	6/box
#1530-3, 3" (777-9909)	4/box

Single-Use Rolls, 1½ yd, White

#1530S-1, 1" (777-6928)	100/box
#1530S-2, 2" (777-0532)	50/box
Standard Rolls, 10 yd, Tan	
#1533-0, ½" (777-4451)	24/box
#1533-1, 1" (777-5133)	12/box
#1533-2, 2" (777-5407)	6/box



3M™ Micropore™ S Surgical Tape SOLVENTUM, FORMERLY 3M HEALTH CARE

- A gentle tape suitable for use on all patients, including elderly and pediatric patients
- Good initial and long term adhesion to skin
- Tape has a paper-plastic hybrid backing with silicone adhesive and is not made with natural rubber latex
- Individually packaged single-use length rolls (1½ yd) can help reduce cross-contamination risk

1" x 5½ yd (777-0521)	12 rolls/box
2" x 5½ yd (777-0533)	6 rolls/box
1" x 1½ yd (117-3048)	100 rolls/box
2" x 1½ yd (137-5598)	50 rolls/box



3M Tegaderm™ Transparent Film Dressings SOLVENTUM, FORMERLY 3M HEALTH CARE

- Transparent dressing that allows for continuous visibility to the IV insertion site
- Adhesive is gentle on skin
- Dressing flexes with skin for greater patient comfort
- Picture-frame delivery makes accurate placement easy
- Unique shapes and designs for both IV and wound needs
- Hypoallergenic and not made with natural rubber latex
- Provide extra holding power for improved adhesion in challenging conditions

Frame Style

#1622W, 1¾" x 1¾"	(777-9319)	100/box
#1624W, 2¾" x 2¾"	(777-9152)	100/box
#1614, 2¾" x 2¾"	(777-2155)	100/box
#1630, 4" x 4½"	(121-4593)	50/box
#1626W, 4" x 4¾"	(777-9551)	50/box
#1626, 4" x 4¾"	(120-3612)	50/box
#1627, 4" x 10"	(777-2028)	20/box
#1628, 6" x 8"	(777-1556)	10/box
#1629, 8" x 12"	(777-3428)	10/box

Frame Style with Border

#1610, 2" x 2¼"	(129-1117)	100/box
#1614, 2¾" x 2¾"	(777-2155)	100/box
#1633, 2¾" x 3¼"	(777-2176)	100/box
#1635, 3½" x 4¼"	(364-4461)	50/box
#1655, 3½" x 4½"	(600-3007)	50/box
#1616, 4" x 4¾"	(115-2641)	50/box

Tegaderm™ HP (Holding Power)

#9534HP, 2¾" x 2¾"	(926-4993)	100/box
#9536HP, 4" x 4¾"	(777-1776)	50/box
#9546HP, 4" x 4½"	(654-4703)	50/box



3M Tegaderm™ + Pad Film Dressing with Non-Adherent Pad SOLVENTUM, FORMERLY 3M HEALTH CARE

- Provides a barrier to bacteria and viruses
- Enhances patient comfort
- Helps support healing
- Easy to apply
- Convenient for clinicians and practical for patients

#3582, Dressing: 2" x 2¾"; Pad: 1" x 1½"	(777-0359)	50/box
#3584, Dressing: 2¾" x 4"; Pad: 1" x 2¾"	(126-5795)	50/box
#3586, Dressing: 3½" x 4"; Pad: 1¾" x 2¾"	(777-2234)	25/box
#3587, Dressing: 3½" x 4½"; Pad: 1¾" x 2¾"	(113-2888)	25/box
#3589, Dressing: 3½" x 6"; Pad: 1¾" x 4"	(163-8541)	25/box
#3590, Dressing: 3½" x 8"; Pad: 1¾" x 6"	(108-7104)	25/box
#3591, Dressing: 3½" x 10"; Pad: 1¾" x 8"	(777-9136)	25/box
#3593, Dressing: 3½" x 13¾"; Pad: 1¾" x 11¾"	(777-2230)	25/box
#3588, Dressing: 6" x 6"; Pad: 4" x 4"	(272-2407)	25/box



3M Tegaderm™ CHG Chlorhexidine Gluconate IV Securement Dressings SOLVENTUM, FORMERLY 3M HEALTH CARE

- Proven to reduce CRBSIs
- Transparent film and gel pad that allows continuous site visibility to easily assess for early signs of infection
- Integrated design combines infection reduction with site visibility, catheter securement, and consistent application
- Designed to minimize catheter movement and dislodgement
- Provides a waterproof, sterile barrier to external contaminants including liquids, bacteria, and viruses

2¾" x 3¾" (7 cm x 8.5 cm)	(777-0356)	100/case
4" x 4¾" (10 cm x 12 cm)	(777-0357)	100/case
3½" x 4½" (8.5 cm x 11.5 cm)	(777-0301)	25/box
4" x 6½" (10 cm x 15.5 cm)	(114-7418)	100/case

3M Adaptic™ Non-Adhering Dressings SOLVENTUM, FORMERLY 3M HEALTH CARE



- Protects the wound bed and minimizes patient pain during dressing changes
- Enables fluid transfer to a secondary dressing, reducing the risk of exudate pooling and maceration
- Suitable for use under compression
- Can be cut without unraveling or shredding

#2012, 3" x 3"	(140-0916)	50/box
#2015, 3" x 8"	(140-0951)	24/box
#2013, 3" x 8"	(140-0933)	36/box
#2014, 3" x 16"	(140-0937)	36/box
#2019, 5" x 9"	(140-0975)	12/box

3M Fibracol™ Plus Collagen Wound Dressing SOLVENTUM, FORMERLY 3M HEALTH CARE



- Maintains a physiologically moist microenvironment
- Conformable dressing can be cut to the required size to fill irregularly shaped cavities
- Helps protect regenerating tissue
- Minimizes patient pain and trauma during dressing changes
- Absorbent and handles low to highly exuding wounds

#2981, 2" x 2"	(140-0976)	12/box
#2982, 4" x 4¾"	(140-0974)	12/box
#2983, 4" x 8¾"	(140-0979)	6/box
#2984, ¾" x ¾" x 15¾"	(140-1393)	6/box

3M Promogran Prisma™ Collagen Matrix with ORC and Silver SOLVENTUM, FORMERLY 3M HEALTH CARE



- Comprised of 55% collagen, 44% ORC and 1% silver-ORC
- Maintains a physiologically moist microenvironment
- Intended for the management of exuding wounds including: diabetic ulcers, venous ulcers, pressure injuries, traumatic wounds, and dehisced surgical wounds

#MA028, 4.34" x 4.34" Hexagon	(140-0929)	10/box
#MA123, 19.1" x 19.1" Hexagon	(139-8291)	10/box
#MA032, ¾" x ¾" x 12½" Rope	(140-0936)	6/box



Wound Care Solutions

Helping surgeons provide optimal care for over 30 years with effective product solutions.

- › Antiseptics and cleansing
- › Applicators and swabs
- › Bandages and dressings
- › OR towels
- › Surgical drapes
- › Impregnated dressings
- › Medical tapes
- › Small procedure kits
- › Specialty and surgical dressings



Dukal is working with healthcare professionals around the world in manufacturing and supply chain solutions, to mitigate disruptions that impact business and patient health. Applying the expertise of our global workforce, national brand of products, technology investments, data insights, and trusted partnerships to build a better future.

› **Contact your Henry Schein representative for more information**





Flexible Fabric Bandages

DUKAL

- The soft weaved material bends easily and conforms to the body, staying on even when wet
- The highly absorbent pad will not stick to the wound
- Long lasting adhesive gently sticks to skin
- Not made with natural rubber latex
- Sterile

Strip, ¾" x 3"	(127-1285)	100/box
Strip, 1" x 3"	(127-1287)	100/box
Knuckle, 1½" x 3"	(127-1288)	100/box
Fingertip, Large, 1¾" x 3"	(127-1289)	50/box
Patch, 2" x 3"	(127-1290)	100/box
Patch, X-Large, 2" x 4"	(127-1284)	50/box
Assorted Pack	(127-1292)	36/box



Soft Flexible Fabric Bandages

DUKAL

- Consists of a highly absorbent non-stick pad to protect the wound
- Long-lasting adhesive gently sticks to the skin
- Pre-folded wrapper tabs make it simple to open and apply
- Not made with natural rubber latex
- Sterile

Strip, ¾" x 3"	(127-1251)	100/box
Strip, 1" x 3"	(127-1254)	100/box
Spot, ⅞"	(127-1247)	100/box
Oval, ⅞" x 1¼"	(127-1252)	100/box
Toe Shield, 1½" x 2⅝"	(127-1253)	100/box
Knuckle, 1½" x 3"	(127-1246)	100/box
Small Fingertip, 1¾" x 2"	(127-1243)	100/box
Large Fingertip, 1¾" x 3"	(127-1244)	50/box
Patch, 1½" x 2"	(127-1248)	100/box
Patch, 2" x 3"	(127-1357)	50/box
4-Wing, 2⅞" x 3"	(127-1255)	50/box

Stat Strip® Bandages

DUKAL

- Specifically designed for use with gloves
- Unique, easy-opening wrapper allows for aseptic application
- Long lasting adhesive gently sticks to skin
- Not made with natural rubber latex
- Sterile

Flexible Fabric Strips

(127-1279) ¾" x 3"	100/box
(127-1280) 1" x 3"	100/box
Sheer Strips	
(127-1277) ¾" x 3"	100/box
(127-1278) 1" x 3"	100/box

Pediatric Character Adhesive Bandages

DUKAL

- To cover minor cuts, scrapes, and wounds
- Sterile and not made with natural rubber latex
- Long-lasting adhesive
- Nonstick wound pad
- Highly breathable

Bandages, ¾" x 3"	100/box
Specify:	
Strawberry	
Shortcake	(127-1320)
Care Bears	(127-1322)
Transformers	(127-1316)
Candy Land	(127-1319)
Tonka	(127-1318)
My Little Pony	(127-1317)
Scooby Doo	(127-1257)

Designer Adhesive Stat Strip Bandages

DUKAL

- Used to cover and protect minor cuts and wounds
- Sterile and not made with natural rubber latex
- Long-lasting and highly absorbent
- Features fun designs

Strips, ¾" x 3"	100/box
Specify:	
Emoji	(127-5610)
Glitter Stars & Strips	(127-1245)
Blue & Pink Camouflage	(127-1294)
Glitter	(127-1262)
Herbie the Dinosaur	(127-1336)
Farm Friends	(145-8472)
Neon	(127-1338)

Fluff Bandage Rolls

DUKAL

- Made of washed, fluff-dried 100% cotton gauze
- Open, crinkle weave provides quick wicking, maximum absorbency, and greater aeration
- Finished edges reduces linting and loose threads
- Adds bulk and loft that will conform as needed, secure dressings in place, and cushion the wound
- Not made with natural rubber latex

#640, Sterile, 3.4" x 3.1 yd, 6-ply	(115-2712)	ea
#645, Sterile, 4.5" x 4.1 yd, 6-ply	(840-3994)	ea
#650, Sterile, 4.5" x 4.1 yd, 8-ply	(120-7670)	ea
#545, Non-Sterile, 4.5" x 4.1 yd, 6-ply	(378-0740)	ea

Super Hero Stat Strip Bandages

DUKAL

Strip Bandages, ¾" x 3"	100/box
Specify:	
Spider-Man	(127-1324)
Iron Man & Captain America	(127-1325)
Hulk & Thor	(127-1326)
Superman, Wonder Woman & Flash	(127-1311)
Batman, Green Lantern & Aquaman	(127-1312)

Premium Self-Closure Elastic Bandages

DUKAL

- Designed with hook and loop closure at both ends
- Crocheted construction provides a better hook to secure and minimizes pulls on the fabric for repeat use
- Ideal for situations requiring high compression
- Individually wrapped
- Not made with natural rubber latex

#512LF, 2"	(128-4686)	10/box
#513LF, 3"	(681-2732)	10/box
#514LF, 4"	(128-4696)	10/box
#516LF, 6"	(128-4695)	10/box

Cohesive Bandages

DUKAL

- Provides consistent support and maintains compression
- Can easily be torn by hand and is self-adherent; no clips or fasteners are needed
- Lightweight, porous, and conforms well with movement

Non-Latex	
1"	30/box
Specify:	
#8016TLF, Tan	(681-3686)
#8016ASLF, Assorted	(681-3674)
1½"	48/box
Specify:	
#8156TLF, Tan	(146-5621)
#8156ASLF, Assorted	(681-3751)
2"	36/box
Specify:	
#8026TLF, Tan	(681-3699)
#8026ASLF, Assorted	(124-0519)
3"	24/box
Specify:	
#8036TLF, Tan	(681-3712)
#8036ASLF, Assorted	(681-3710)
4"	18/box
Specify:	
#8046TLF, Tan	(681-3743)
#8046ASLF, Assorted	(681-3719)

Type VII Gauze Sponges DUKAL



- USP Type VII gauze, 100% woven cotton
- Highly absorbent and offer exceptional clinical performance
- Used for a number of applications, including general wound care, dressing, prepping, wound packing, and debriding
- Not made with natural rubber latex

Sterile	
2" x 2", 8 Ply (140-8948)	100/box
3" x 3", 12 Ply (140-8915)	100/box
4" x 4", 8 Ply (140-8716)	100/box
4" x 4", 12 Ply (140-8891)	50/box
8" x 4", 12 Ply (140-8970)	50/box
Non-Sterile	
2" x 2", 8 Ply (140-8898)	200/bag
2" x 2", 12 Ply (140-8896)	200/bag
3" x 3", 12 Ply (140-8953)	200/bag
4" x 4", 8 Ply (140-8901)	200/bag
4" x 4", 12 Ply (140-8893)	200/bag
4" x 4", 16 Ply (140-8986)	200/bag

Non-Woven Sponges DUKAL



- Rayon/poly blend construction
- Highly absorbent
- Multipurpose sponges ideal for general wound care and applying ointments
- Virtually nonlinting
- Not made with natural rubber latex

#7224, Sterile, 2" x 2" (120-5675)	50 packs/box Contains: 2 per package, 50 packages per box.
#7334, Sterile, 3" x 3" (120-0238)	40 packs/box Contains: 2 per package, 40 packages per box.
#7444, Sterile, 4" x 4" (120-0685)	50 packs/box Contains: 50 sponges per tray, 2 trays per package.
#7454SP, Sterile, 4" x 4" (600-6705)	800/case Contains: 5 per package, 20 packages per box, 8 boxes per case.
#2102, Non-sterile, 2" x 2" (852-0094)	200/bag Contains: 200 per bag, 20 bags per case
#2103, Non-sterile, 3" x 3" (852-0091)	200/bag Contains: 200 per bag, 20 bags per case
#2104, Non-sterile, 4" x 4" (933-3131)	200/bag

Non-Adherent Pads

DUKAL

- Absorbent pad with nonstick outer facing
- Allows dressing to absorb, without adhering to the wound
- Ease of removal helps minimize trauma to the wound site
- Can be cut to fit any size wound without separating
- Sterile
- Not made with natural rubber latex



Without Adhesive	
#123, 2" x 3" (120-5465)	100/box
#134, 3" x 4" (120-9890)	100/box
#138, 3" x 8" (261-0216)	75/box
With Adhesive	
#124, 2" x 3" (681-3519)	100/box
#125, 3" x 4" (681-3521)	75/box

Caliber™ Non-Woven Island Dressings

DUKAL

- Soft, nonwoven backing with adhesive border is gentle on skin, conformable, breathable, and provides ease of movement
- Absorbent, nonadherent pad collects excess exudates and protects wound from trauma
- Not made with natural rubber latex



#4070, 2" x 3" (933-0082)	50/box
#4075, 4" x 4" (933-0086)	25/box
#4071, 4" x 5" (933-0081)	25/box
#4072, 4" x 8" (933-0083)	25/box
#4073, 4" x 10" (933-0084)	25/box
#4076, 6" x 6" (933-0087)	25/box

Transparent Semi-Permeable Dressings

DUKAL

- Microporous transparent film with high moisture vapor transmission rate (MVTR)
- Ideal for IV sites, minor burns, ulcers, skin tears, and donor sites
- Easy, 2-step application
- Protects against fluid contaminants and bacteria
- Rounded edges prevent roll-up
- Removable patient label for chart
- Not made with natural rubber latex



#4045, 1½" x 2" (600-2488)	100/box
#4067, 2¾" x 2¾" (120-7785)	100/box
#4112, 4" x 4¾" (120-1809)	50/box
#4152, 6" x 8" (108-4765)	10/box

Petrolatum Gauze Non-Adhering Dressings

DUKAL

- Provides soothing non-adherent protection with more convenience and less mess
- Promotes healing by maintaining a moist wound environment
- Clings and conforms to all body contours
- Reduces pain and trauma during dressing changes
- Decreases chances of skin maceration



#204, ½" x 72" (141-9336)	12/box
#210, 1" x 18" (124-9282)	50/box
#206, 3" x 9" (125-2487)	50/box
#207, 3" x 18" (128-6079)	12/box
#208, 3" x 36" (141-8773)	12/box
#209, 6" x 36" (129-4429)	12/box



Xeroform Petrolatum Non-Adhering Dressings

DUKAL

- Provides soothing non-adherent protection with more convenience and less mess
- Promotes healing by maintaining a moist wound environment
- Clings and conforms to all body contours
- Reduces pain and trauma during dressing changes
- Made with Xeroform, a special deodorizing agent containing 3% bismuth tribomphenate

#212, 1" x 8" (156-5251)	50/box
#214, 2" x 2" (128-6077)	25/box
#215, 4" x 4" (128-6078)	25/box
#213, 5" x 9" (125-2489)	50/box
#216, 4" x 3-yd Roll (128-6080)	6/box



All-Gauze Sponges HENRY SCHEIN®

Non-Sterile

- Precision sizes and folding: consistency within functions
- Soft and pliable: nonabrasive touch
- Highly absorbent: soaks up unwanted fluid
- Can be used for general cleaning, debriding, and prepping

8 Ply, 2" x 2" (570-2657).....	200/pkg
12 Ply, 2" x 2" (570-2658).....	200/pkg
8 Ply, 3" x 3" (570-2650).....	200/pkg
12 Ply, 3" x 3" (570-2648).....	200/pkg
8 Ply, 4" x 4" (570-2656).....	200/pkg
12 Ply, 4" x 4" (570-2630).....	200/pkg



Premium All-Gauze Sponges

HENRY SCHEIN®

- Manufactured from 100% USP type-VII cotton gauze, bleached and processed to established standards
- Soft and pliable
- Protects as a wound covering
- Provides maximum absorbency

Non-Sterile	
8 Ply, 2" x 2" (570-2647).....	200/pkg
12 Ply, 2" x 2" (570-2628).....	200/pkg
12 Ply, 3" x 3" (570-2640).....	200/pkg
8 Ply, 4" x 4" (570-2646).....	200/pkg
12 Ply, 4" x 4" (570-2645).....	200/pkg
Sterile	
8 Ply, 2" x 2" (570-2642).....	100/pkg
12 Ply, 3" x 3" (570-2631).....	80/pkg
8 Ply, 4" x 4" (570-2644).....	50/pkg
12 Ply, 4" x 4" (570-2637).....	50/pkg



Compare to Steri-Pad
Gauze Pads

HENRY SCHEIN®

Sterile

- Highly absorbent
- 100% woven cotton, 12-ply
- Can be used for general cleaning, debriding, and prepping
- Sterile 1's in peel-back packages
- Box is perforated at bottom for easier gauze removal

2" x 2" (570-2659).....	100/pkg
3" x 3" (570-2629).....	100/pkg
4" x 4" (570-2649).....	100/pkg

Type VII Gauze

HENRY SCHEIN®

- USP Type VII gauze, 100% woven-cotton
- Not made with natural rubber latex
- Sterile, 4" x 4"

8 Ply (570-2633).....	10/pkg
12 Ply (570-2655).....	10/pkg
16 Ply (570-2638).....	10/pkg
12 Ply, X-ray Detectable (570-2632).....	10/pkg
16 Ply, X-ray Detectable (570-2634).....	10/pkg



Exodontia Sponges

HENRY SCHEIN®

- 100% cotton gauze sponges with maximal absorbency cotton inner liner
- Manufactured under strictest standards: superior and consistent product quality
- Available in 3 sizes: convenient for different uses
- 8-ply
- Not made with natural rubber latex

Non-Sterile-Filled	
2" x 2" (570-2636).....	5000/case
3" x 3" (570-2643).....	4000/case
4" x 4" (570-2651).....	2000/case
Non-Sterile-Filled	
2" x 2" (570-2653).....	200/pkg



Premium Exodontia Sponges

HENRY SCHEIN®

Manufactured from all USP 100% cotton. Soft, absorbent, and pure white. Highest quality. 2" x 2", 8-ply.

Non-Sterile-Filled (570-2635).....	5000/case
Non-Sterile-Filled (570-2652).....	200/pkg



Conforming Stretch Gauze Rolls

HENRY SCHEIN®

4.1-yd Rolls

- Economical knitted bandage in various widths
- Made from high-quality materials: nylon, polyester, and cotton
- Provides light compression as it cushions and protects
- Offers softness, breathability, and alleviation of irritation
- Low lint content reduces contamination to wound site
- Improved packaging is resealable, durable, easy to store, and remains clean
- Not made with natural rubber latex

Non-Sterile

1" (104-2704).....	12 rolls/bag
2" (104-8825).....	12 rolls/bag
3" (104-8997).....	12 rolls/bag
4" (104-2728).....	12 rolls/bag
6" (104-9424).....	6 rolls/bag
Sterile	
2" (104-2739).....	12 rolls/bag
3" (104-9724).....	12 rolls/bag
4" (104-6708).....	12 rolls/bag
6" (104-4293).....	6 rolls/bag



CRITERION™ Non-Woven Sponges



Non-Sterile, 4 Ply

- High-quality, finely constructed, rayon/poly-blend that is soft and gentle to the touch
- Multiuse sponge that can be used for all applications of wound-care treatments
- Thicker and more absorbent than standard sponges
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2"	
(104-2930)	200/pkg
3" x 3"	
(104-2410)	200/pkg
4" x 4"	
(104-2411)	200/pkg



CRITERION™ Non-Woven Sponges



Sterile 2s, 4 Ply

- High-quality, finely constructed, rayon/poly-blend that is soft and gentle to the touch
- Multiuse sponge that can be used for all applications of wound-care treatments
- Thicker and more absorbent than standard sponges
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2"	
(104-7897)	50/pkg
3" x 3"	
(104-4132)	50/pkg
4" x 4"	
(104-8665)	50/pkg



Syngaize 50™



Non-Sterile, 4 Ply

- General-use nonwoven material with a 70/30% blend of rayon and polyester. 47–50% more absorbent than 12-ply gauze
- Absorbs 10 times its weight: greater absorbency reduces the number of sponges required
- Soft and pliable: prevents irritation to mucosa; allows shaping for a better fit
- Sterilizable: lint-free mesh minimizes lint contamination
- Nonraveling: less mess
- Not made with natural rubber latex

2" x 2"	
(100-2125)	200/pkg
3" x 3"	
(101-9668)	200/pkg
4" x 4"	
(101-4583)	200/pkg



Premium Non-Woven Sponges



Non-Sterile, 4 Ply

- General-use, non-woven rayon/polyester-blend material
- Quarter-fold design
- Less linting: minimizes lint contamination
- Soft and pliable: prevents irritation to mucosa
- Greater absorbency: reduces number of sponges required
- Closed weave
- Not made with natural rubber latex

2" x 2"	
(101-9278)	200/pkg
3" x 3"	
(101-8748)	200/pkg
4" x 4"	
(101-4336)	200/pkg



Non-Woven Cotton Sponges



Non-Sterile, 4 Ply

- Highly absorbent 100% cotton
- 4-ply gauze sponges: thick, soft gauze
- Water bleach: no discoloration after sterilization
- See-through hygienic package imprinted with a sterilization indicator: easy identification
- Tear notch: convenient bag opening

2" x 2"	
(100-1090)	100/pkg
3" x 3"	
(100-3835)	100/pkg
4" x 4"	
(100-4420)	100/pkg



Combine ABD Pads



- Highly absorbent multilayer ABD (abdominal) dressing
- Soft, nonwoven outer facing wicks moisture into the cellulose core, which quickly absorbs and disperses fluid
- Moisture-resistant barrier prevents fluid strikethrough
- All edges are sealed to prevent linting and leakage
- Not made with natural rubber latex
- Sterile

5" x 9"	
(570-1470)	25/box
8" x 7½"	
(570-1472)	20/box
8" x 10"	
(570-1471)	20/box



Flexible Fabric Elastic Strips



These strips stretch, conform, and permit free movement. Will not stick to wound. Not made with natural rubber latex.

¾" x 3"	(112-6132)	100/box
1" x 3"	(112-6133)	100/box
Knuckle, ½" x 3"	(112-6144)	100/box



Sheer Plastic Bandage Strips



Sterile
Individually wrapped adhesive bandages for the protection of minor cuts, punctures, and abrasions. Highly absorbent pad will not stick to wound. Air vents let air in and keep moisture out.

¾" x 3"	(112-6138)	100/box
1" x 3"	(112-6142)	100/box
⅞" Spots	(112-6143)	100/box
Professional Tray Pack, ¾" x 3"	(112-6134)	1500/case
Professional Tray Pack, 1" x 3"	(112-6135)	1500/case



Crayon Adhesive Bandages



Sterile
• Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
• Pad will not stick to wounds
• Permits wound to drain
• Not made with natural rubber latex

¾" x 3"	(112-6995)	100/box
---------	------------	---------



Self-Adherent Bandages



Non-Sterile, 5-yd Rolls
• Self-adherent: sticks to itself, not the skin
• Comfortable porous material: pain-free removal that will not damage skin

Latex		
Tan, 1"	(570-0113)	ea
Tan, 2"	(900-4322)	ea
Tan, 3"	(900-4323)	ea
Tan, 4"	(900-4324)	ea
Tan, 6"	(900-4328)	12 rolls/pkg
Multi, 1"	(112-7366)	30 rolls/pkg
Multi, 2"	(900-4331)	36 rolls/pkg
Multi, 3"	(900-4332)	24 rolls/pkg
Non-Latex		
Tan, 1"	(112-7355)	ea
Tan, 2"	(112-7357)	ea
Tan, 3"	(112-7358)	ea
Tan, 4"	(112-7359)	ea
Tan, 6"	(112-7362)	ea
Multicolored, 1"	(112-7363)	30 rolls/pkg
Multicolored, 2"	(112-7364)	30 rolls/pkg
Multicolored, 3"	(112-7365)	24 rolls/pkg



Animal Adhesive Bandages



Sterile
• Exclusive design for Henry Schein
• Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
• Nonstick, highly absorbent pad will not stick to wound
• Air vents let air in and keep moisture out
• Not made with natural rubber latex

Spots, ⅞"	(900-4506)	100/box
¾" x 3"	(900-4507)	100/box



Premium Hypoallergenic Surgical Paper Tape



10-yd Rolls
• Excellent adhesion
• Pinked, nonfraying edges
• Strong woven backing
• Gentle removal
• Not made with natural rubber latex
• Non-sterile

½"	(101-4305)	24/box
1"	(101-9673)	12/box
2"	(101-1332)	6/box



Premium Hypoallergenic Surgical Cloth Tape



10-yd Rolls
• Reliable adhesion
• Highly porous
• Thin, flexible backing
• Gentle removal
• Not made with natural rubber latex
• Non-sterile

½"	(101-3025)	24/box
1"	(101-2323)	12/box
2"	(101-5979)	6/box



Premium Hypoallergenic Transparent Surgical Tape



10-yd Rolls
• Excellent adhesion
• Porous, with flexible backing
• Tears easily in both directions
• Not made with natural rubber latex
• Non-sterile

½"	(101-6655)	24/box
1"	(101-4645)	12/box
2"	(101-9226)	6/box



Elastic Bandages with VELCRO®
HENRY SCHEIN®

Non-Sterile, 4½-yd Rolls

- Elastic bandage with Velcro® brand fasteners
- Provides even compression with strength, comfort, and breathability
- Premium honeycomb-knitted material
- Not made with natural rubber latex
- Velcro® closure

2"		
(900-4681)	10/box
3"		
(900-4682)	10/box
4"		
(900-4683)	10/box
6"		
(900-4684)	10/box



Elastic Bandages with Clips

HENRY SCHEIN®

Non-Sterile, 4½-yd Rolls

- High cotton content to keep patients cool and comfortable
- Durable bandages can be laundered, sterilized, and reused
- Ideal for support and compression applications, such as treatment of varicosities, sprains, strains, and torn or injured ligaments
- Not made with natural rubber latex

2"		
(900-4677)	10 rolls/pkg
3"		
(900-4678)	10 rolls/box
4"		
(900-4679)	10 rolls/box
6"		
(900-4680)	10 rolls/box



Unna's Boot

HENRY SCHEIN®

10-yd Roll

- Cotton gauze formulated with a zinc-oxide emulsion
- Used for treating venous leg ulcers, edema, and dermatitis
- Soothes inflammation and rashes
- Not made with natural rubber latex

3"		
(100-3444)	ea
4"		
(100-3256)	ea
3" with Calamine		
(900-7455)	ea
4" with Calamine		
(900-7456)	ea



Now Medicare Covered Cold Therapy¹

Supported by Clinical Research Demonstrating Reduced Postoperative Opioid Use.²

HS PART NUMBER	DJO PART NUMBER	DESCRIPTION
(106-6950)	11-1421	IceMan Classic® Unit with Cold Pad
(117-3725)	11-0494	IceMan Classic® Unit

¹Traditional Medicare Program coverage is available only through Medicare Hospital Outpatient Prospective Payment System (OPPS) or through the Medicare ASC Payment System.

²Xiao AX, et al. The use of a motorized cold therapy unit in the postoperative period influences pain and narcotic use following arthroscopic rotator cuff repair. J Shoulder Elbow Surg. 2025;[Epub ahead of print]. doi:10.1016/j.jse.2025.08.026

For more information on the NOPAIN Act or Enovis™ Motorized Cold Therapy please contact your Henry Schein Representative



WARNING: This device can be cold enough to cause serious injury. Serious adverse reactions & safety hazards may occur when using this device. Read and understand all warnings and Instructions for Use, including contraindications, before using this device. Additional warnings appear in the Instructions for Use. DO NOT use this device without a prescription from a physician. Your prescription must state how long and how often the device should be used and the length of breaks between uses. Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, legal, or any other type of professional advice. Rather, please consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.



Premium Panel Elastic Binder

ENOVIS

Constructed with panel elastic and soft, cotton-type lining for patient comfort. Universal style, flannel-lined, panel-elastic support with contact closure. Multi-panel design helps prevent panel rolling. Ideal for providing compression and support for post-natal or post-abdominal surgery, abdominal strains, and weakness. 9", 12", or 15" height.

- Waist Measurement: 30"–45", Universal (802-6850)..... ea
- 9", 3-Panel, Waist Measurement: 30"–45", Universal (109-9670).....6/pkg
- Waist Measurement: 45"–62", Universal (213-9079)..... ea
- 9", 3-Panel, Waist Measurement: 45"–62", Universal (109-9671).....6/pkg
- Waist Measurement: 62"–74", X-large (113-5616)..... ea
- Waist Measurement: 72"–84", 2X-large (113-5617)..... ea
- Waist Measurement: 82"–100", 3X-large (153-0722)..... ea
- Waist Measurement: 30"–45", Universal (397-8371).....ea
- 12", 4-Panel, Waist Measurement: 30"–45", Universal (109-9668).....6/pkg
- Waist Measurement: 45"–62", Universal (896-4073).....ea
- Waist Measurement: 45"–62", Universal (109-9669).....6/pkg
- Waist Measurement: 62"–74", X-large (481-0385)..... ea
- Waist Measurement: 72"–84", 2X-large (108-1932)..... ea
- Waist Measurement: 82"–100", 3X-large (102-1525)..... ea
- Waist Measurement: 62"–74", X-large (400-3885)..... ea



Deluxe Abdominal Binders

F. STUBBS

Plush-lined elastic for comfortable feel against the skin. Provide consistent compression during recovery. Fully adjustable with contact closure.

- #1083, 3-Panel 9" X-small, 20"–30" (272-0966).....ea
- #10833, 3-Panel 9" Small, 30"–45" (101-5804).....ea
- #10835, 3-Panel 9" Large, 46"–62" (101-0547)..... ea
- #10843, 4-Panel 12" Small 30"–45" (101-1461).....ea
- #10845, 4-Panel 12" Large, 46"–62" (101-7060).....ea
- #10846, 4-Panel 12" X-large, 62"–75" (272-6577).....ea
- #10847, 4-Panel 12" XX-large, 75"–90" (272-7667).....ea



Blue Vogue Arm Sling

ENOVIS

Generous-size envelope constructed of durable poly-cotton provides patient with comfortable custom fit, yet helps ensure mobility. Fully-padded, 2"-wide foam strap with contact closure and D-rings for easy adjustment and application. Thumb loop helps prevent sling from migrating. Ideal for cast support; injuries to arm, wrist, or hand.

- Small, 7½" x 12½" (380-3055)..... ea
- Medium, 8" x 15" (455-5401).....ea
- Large, 8½" x 17½" (840-2000).....ea
- X-large, 9" x 20" (782-9206).....ea



Deluxe Arm Sling with Pad

ENOVIS

Generously sized envelope to support the arm. Cotton/poly construction with web strap and contact closure. Thumb loop helps prevent sling from migrating. Shoulder pad included.

- Sling..... ea
- Specify:
- #79-84001, XX-small (Ped.), 6½"D x 8½"L (120-5721)
- #79-84002, X-small (Child), 7"D x 11"L (120-6097)
- #79-84003, Small, 7½"D x 12½"L (120-2338)
- #79-84005, Medium, 8"D x 15"L (120-3852)
- #79-84007, Large, 8½"D x 17½"L (120-0628)
- #79-84008, X-large, 9"D x 20"L (120-4568)



Emesis Basin

HENRY SCHEIN®

- Plastic
- Kidney shaped
- 16-oz capacity

Mauve, 8½" (112-5809)..... 25/box



Synthetic Cast Padding

HENRY SCHEIN®

Non-Sterile, 4-yd Rolls

- Made of synthetic fibers
- Provides a safe environment under a hard cast
- Comfortable feel for the patient
- 100% polyester and water-resistant; not waterproof
- Retains consistency under wet or dry conditions
- Conforms to all body contours
- Works with synthetic or plaster casting material
- Not made with natural rubber latex

- 2" (900-4997).....24/bag
- 3" (900-4998)..... 12/bag
- 4" (900-4999)..... 12/bag
- 6" (900-7001).....6/bag




Instant Cold Compresses



Disposable, single-use cold pack that provides immediate cold therapy relief without refrigeration. Ready to activate for quick use in emergencies or first aid situations. Compact, lightweight, and flexible for comfortable application. Ideal for temporary relief of pain and swelling caused by sprains, contusions, minor burns, toothaches, and insect bites. Non-sterile.

- 5" x 5 1/2" (112-6147) 50/case
- 5" x 7" (112-6146) 48/case
- 6" x 9" (15.24 cm x 22.86 cm) (112-6145) 24/case



Reusable Cold or Hot Gel Packs



Versatile gel pack designed for both cold and hot therapy. Simply freeze for cold therapy or microwave for heat. Remains soft and flexible—even when frozen—for added comfort. Reusable and easy to clean for long-lasting value. Provides temporary relief of arthritis, joint pain, muscle strain, and stiffness. Non-toxic and not made with natural rubber latex. Non-sterile.

- 4" x 6 1/2" (10.16 cm x 16.51 cm) (570-2140) ea
- 5 1/2" x 10" (13.34 cm x 25.4 cm) (570-2139) 24/case



Cool Jaw® Soft-Sided Round Gel Packs



Provides uniform cold therapy for an extended period of time. Features one soft-sided surface that shields the skin from direct contact with the frozen pack. Offers a convenient, comfortable, and economical cold therapy application for patients.

- Non-Imprinted
4" ea
- Specify:
- Cool Blue (560-0004) Silver Glitter (560-0079)
 - Pink Ice (560-0005) Black Glitter (560-0073)
 - Grape (560-0054) Purple Glitter (560-0078)
 - Kiwi (560-0055) Clear Glitter (560-0076)
 - White Glitter (560-0081) Teal Glitter (560-0080)
 - Gold Glitter (560-0075)



Aluminum Crutches



- Push-pin adjustment makes it easy to adjust crutch height in 1" increments
- Nonskid, jumbo-size, vinyl-contoured tips provide excellent traction
- Easy wing-nut adjustments for hand grips
- 300-lb weight capacity on aluminum crutch

- Crutches pair
- Specify:
- Youth, Patient Height: 4'6"-5'2", Crutch Height: 37"-46" (570-1466)
 - Adult, Patient Height: 5'2"-5'10", Crutch Height: 45"-53" (570-1464)
 - Tall Adult, Patient Height: 5'10", Crutch Height: 53"-61" (570-1465)



Protection Plus Disposable Underpads



- Made of waterproof, blue poly backing
- Absorbent core of tissue layers
- White, hypoallergenic spunbond facing sheet
- Not made with natural rubber latex

- 17" x 24" (112-5633) 300/case
- 23" x 24" (112-5634) 200/case
- 23" x 36" (570-1182) 150/case
- 30" x 30" (570-1183) 150/case
- 30" x 36" (570-1181) 100/case



Patient Belongings Bags



- Convenient bags for moving and storing patient belongings
- Plastic, opaque white
- Print: Patient Belongings, Name, Room
- 20" x 20" x 4"

- Bags with Handles (570-1179) 250/case
- Bags with Drawstrings (570-1180) 250/case

DawnMist Drawstring Patient Belonging Bags



- Durable for holding patient's clothing and belongings
- Can easily be hung from wheelchair or gurney
- Designated space to write patient's name and room number
- Not made with natural rubber latex

- #DS300, White Bag with Nature Design, 11" x 17" (+4" gusset), 1.5 mil (101-6144) 1000/case
 - Bags, 20" x 20" +4" Gusset 250/case
- Specify:
- #DSPB01, White Bag with Blue Print, 1.3 mil (329-9395)
 - #DSPB01C, Clear Bag with Blue Print, 1.4 mil (840-0751)





HARMONIZED SYSTEMS

CLINICAL CHEMISTRY, IMMUNOASSAY AND INTEGRATED SYSTEMS TO TRANSFORM YOUR LABORATORY



Clinical Chemistry | Immunoassay | Hematology | Transfusion | Molecular | Point of Care | Professional Services



Quality Products, Affordable Prices, and Outstanding Service.

Why choose Henry Schein® Products?

The products you use in your health care facility have to be high quality, whether it's reliable tools you use for a procedure or supplies that help prevent the spread of infection. But because you're running a business, the price of those products is always a factor when you're making a purchasing decision.

That's why we're proud to offer products that support optimal patient care and help protect your bottom line. Additionally, our products are backed by Henry Schein®, a company customers have known and trusted for more than 93 years.




EXCLUSIVE BRANDS
AND PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS



To Order: **1-800-772-4346** 8:00am–8:00pm (et)
henryschein.com/medical

California customers please use website when ordering for Prop 65 information.

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2026 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. †For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.

A FRIENDLY REMINDER TO RECYCLE—
Thank you for helping us care for the environment.

Rely on Us for Non-Acute Supply Chain Expertise

Discover the ultimate convenience of seamless delivery with Henry Schein Medical. Serving global health care providers for more than 90 years, non-acute distribution is our core focus.



Your Single Source for Medical Supplies, Equipment, and Solutions

Get everything you need to deliver exceptional patient care in one order. Our comprehensive portfolio of high-quality products, equipment, and solutions will help you achieve operational and financial excellence.



Medical Equipment



Medical-Surgical Supplies



Pharmaceutical Supplies



Laboratory Supplies



Solutions

Non-Acute Supply Chain Excellence

Designed to support the non-acute market, our 5 mega distribution centers provide the versatility and convenience you demand without additional costs.



99%
Historical Fulfillment Rate



99.91%
Order Accuracy



94%
of Customers Serviced
Next-Day Ground*

- ✓ 5-Day-per-Week Delivery with No Out-of-Cycle Fees
- ✓ No Low Unit of Measure Fees
- ✓ Flexible Stocking Requirements
- ✓ Standard Inside Delivery
- ✓ Low Minimum Order Requirements

Partner with a supplier who prioritizes your non-acute facility. Contact Henry Schein Medical today.

HenrySchein.com/medical | 1-800-772-4346 between the hours of 8:00am–8:00pm ET

***Orders must be received by 5:00pm ET**