

MEDICAL PRODUCT GUIDE

2026-2027



Innovating **Healthcare** for the Future

Your **One Source** for a **Thriving Practice**

Practices are unique and constantly evolving. **Rely on Us** for customized solutions that empower optimal patient & practice outcomes.

www.henryschein.com/medical | 1.800.772.4346 8am-8pm (et)

California customers please use website when ordering for Prop 65 information.

HENRY SCHEIN®

PRODUCTS



Welcome

to our latest Henry Schein® Products Catalog.

Henry Schein® Products is a global product and brand-focused organization with a broad portfolio of dental and medical products. We design, develop, and manufacture our products with attention to performance and reliability.

We leverage our in-depth product expertise, global manufacturing, and supply chain network to offer our customers a comprehensive range of solutions. From everyday clinical supplies to point-of-care products, Henry Schein® Products enable our customers to provide the right care for their patients.

Together we make the world healthier.

Table of Contents:

- Equipment 3-7
- Exam Room Disposables 8-13
- Gloves 14-20
- PPE..... 21-26
- Infection Prevention & Control 27-33
- Hypodermics 34-36
- Blood Collection 37
- IV Products 38
- Waste Management..... 39
- Surgical Instruments 40-71
- Surgical Products 72-81
- Wound Closure 82-86
- Wound Care..... 87-89
- Orthopedics 90-91
- Laboratory 92-99
- Pharmaceuticals 100-101
- Health & Beauty Care 101
- Cold Storage 101
- Anesthetics & Preventives..... 102
- Index 103-105



**Green Series™ 777
Integrated Wall System**

All models include a BP aneroid, an ear specula dispenser, and a SureTemp Plus thermometer.

#777-PM2WAS-HS, With PanOptic Basic LED Ophthalmoscope and MacroView Basic LED Otoscope **(570-2666)**..... ea

#777-PM3WAS-HS, With PanOptic Plus LED Ophthalmoscope and MacroView Plus LED Otoscope for iExaminer **(570-2665)**..... ea

#777-SM2WAS-HS, With Coaxial LED Ophthalmoscope and MacroView Basic LED Otoscope **(570-2667)**.....ea

#777-SM2WAX-HS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser (Does not include SureTemp Plus Thermometer.) **(570-2668)**..... ea



Otoscope Speculum Tips

- Disposable tips carefully manufactured to be free of rough edges or burs
- Designed for a snug fit on Welch Allyn®, Propper®, and Heine® otoscopes
- Will not fit the Welch Allyn® Macroview™

2.5 mm **(112-5527)**.....1000/bag

4.0 mm **(112-5528)**.....1000/bag



Handheld Pulse Oximeter

Can replace the patient monitor and be used as bedside monitoring in hospitals and clinics. Oximeter comes with 1 adult reusable sensor.

- High-resolution TFT screen displays SpO₂ PR, perfusion index (PI), pulse bar, and waveform
- Adjustable audible and visual alarms; battery-low indicator
- Adjustable backlight and voice
- 127 ID setup; 72-hour data storage and review
- MedView software for data analysis
- 3 AA-size alkaline batteries or AC/D adapter (optional)
- Suitable for adult and pediatric use
- Multilanguage menu: English, French, German, Spanish, Italian, Japanese, Russian, and Chinese
- Certification: FDA, CE, CFDA
- Pediatric Sensor sold separately

Pulse Oximeter **(900-8275)**..... ea

Pediatric SpO₂ Sensor Intended for spot checking or continuous monitoring of functional arterial oxygen saturation and pulse rate for pediatric patients. For noninvasive use with oximeter equipment in hospitals, hospital-type facilities, and home environments. **(570-2091)**..... ea



Fingertip Pulse Oximeter

- Dual-color LED makes monitoring easy
- Display modes include SpO₂ (range 70%–100%) and pulse rate (range 30–235 bpm) along with pulse bar and plethysmogram
- 2 AAA batteries provide 30 hours of continuous use
- Visual and audible low-battery alarm
- Suitable for both adult and pediatric patients
- Electromagnetic compatibility

(570-3222).....ea



Digital Ear Thermometer

- For professional or personal use
- Instant reading results in 1 second
- Memory feature recalls last reading
- Fever alarm
- Low battery indicator
- Replaceable lithium battery
- °F or °C measuring mode
- Thermometer includes 15 disposable probe covers

Thermometer **(900-6438)**..... ea

Probe Covers **(900-6439)**.....150/box



PRO Digital Thermometer

- Measures body temperature by oral, axillary, or rectal read
- 60-second response
- 5 probe covers
- More than 200 hours of continuous operation

(570-2055)..... ea



IV Poles & Accessories

- 4 wheels to provide smooth transport and maneuverability
- Chrome-plated steel with coated base for strength, durability, and reduced risk of tipping
- Locking collar for easy height adjustment
- Height range: 120 cm to 216 cm

Chrome IV Pole with 2 Hooks **(570-1484)**..... ea

Chrome IV Pole with 4 Hooks **(570-1486)**..... ea

Mayo Stand **(570-1485)**..... ea

Exam Stool **(570-1487)**.....ea

Overbed Table **(570-1488)**..... ea

HENRY SCHEIN®
BRAND PRODUCTS

DIAGNOSTIC
TOOLS



FINGER PULSE OXIMETER

PRECISION AT YOUR FINGERTIPS.

Durable and compact, sized
to fit both children and adults.



Easy-to-use with easy-to-read results:

- Dual-color OLED displays SpO₂, pulse rate (PR), pulse bar, and waveform
- 6 display modes with level 1-10 adjustable brightness
- Low power consumption with low battery-indicator

Fluid
Resistant
Design



Finger
Skid-Resistance
Design

Shockproof
Occlusion
Design

Battery Cabin Safe
Lock Design



Henry Schein®
Finger Pulse Oximeter
(570-3222)

Order today at [HENRYSCHEIN.COM/MEDICAL](https://www.henryschein.com/medical)



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS



Pro Plus Series Palm Aneroid

This revolutionary palm-style manometer (#804) has an integral bulb and trigger-style air release valve for one-hand convenience.

- Precision-crafted, 300mm-Hg, CE-marked, no-pin stop movement tested to 50,000 cycles
- Movement housed in a zinc safety cage for unsurpassed durability
- Ambidextrous infinite control trigger-release air valve
- Nylon cuff with proprietary size guide marking system, to prevent miscuffing
- Not made with natural latex rubber
- Dip-molded inflation bladder withstands heavy use
- Quick-release Luer-slip connector
- Complete with leatherette carrying case with self-repairing nylon zipper
- Assembled, inspected, and packaged in the USA from foreign components

Adult, Black, Non-Latex (112-6067)	ea
Large Adult (112-6086)	ea
Child (112-6090)	ea
Small Adult (112-6087)	ea
Adult, Cuff and Bladder, 1 Tube (112-6088)	ea
Large Adult (112-6092)	ea
Thigh (112-6091)	ea
Child (112-6093)	ea
Small Adult (112-6094)	ea
Adult (112-6095)	ea
Large Adult (112-6096)	ea
Thigh (112-6097)	ea



Essentials Series Pocket Aneroid

Essentials

- Navy, Latex
- Black enamel, 300mm-Hg, no-pin stop manometer
 - Navy blue nylon cuff with range markings to facilitate selection of correct cuff size, artery label, and gauge holder
 - Latex inflation bladder and bulb
 - Chrome-plated deflation valve
 - Complete with leatherette zippered carrying case and operating instructions
 - Inspected and packaged in the USA from foreign components

Aneroid.....	ea
Specify:	
Child (112-6076)	Large Adult (112-6080)
Adult (112-6077)	



Pro Plus Series Professional Stethoscope

Combination adult-size chest piece is precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Available in 3 striking colors for departmental coding or personal identification.

- Extra-large bell (1 7/8" diameter) for unsurpassed low-frequency response
- Ultrasensitive diaphragm (1 7/8" diameter) for greater amplification and crisper high-frequency response
- Color-coordinated nonchill bell and diaphragm retaining rim for patient comfort
- Stainless steel binaurals with double-leaf internal spring fixed at 15° angle for greater comfort
- Reinforcing yoke molded into flexible 22" PVC "Y" tubing
- Deluxe PVC ear tips (in small and large) for the ultimate in wearing comfort and acoustic seal
- Spare diaphragm included
- Not made with natural rubber latex
- Scope ID tag included
- Lifetime warranty
- Weighs 6.5 oz
- Overall length: 31"

Stethoscope.....	ea
Specify:	
Black (900-4815)	Navy (900-4817)
Burgundy (900-4816)	



EasyOne® Air Spirometer

The all-in-one revolution in lung-function testing. Built for health-care providers large and small, EasyOne Air puts the power of the most advanced spirometry testing tools in the palm of your hand. Revolutionary TrueFlow ultrasound technology provides unprecedented accuracy and reliability. Equipped with ready-to-use EMR connectivity and a user-friendly interface, EasyOne Air makes the best in lung care available anywhere.

- Portable and PC-based solution
- Large, color touch screen
- Wireless connectivity to the PC via Bluetooth
- Instant interpretation based on ATS/ERS guidelines
- Rechargeable battery
- Pediatric incentive screens available with EasyOne Connect
- Calibration-free
- Compatible with EasyOne Connect

#2500-2HS, Spirometer (570-0683)	ea
#2050-1HS, Spirette Mouthpieces (570-0383)	50/case
#2050-5HS, Spirette Mouthpieces (570-0384)	200/case



Nebulizer Compressor System

Quality designed to deliver prescribed medication into aerosol mist for easy inhalation to the bronchial lung passages. Small and lightweight for easy storage and travel mobility. 1.18 kg.

- Operating pressure range: 8–16 psi
- Compressor flow range: 5–8 lpm
- Average nebulization rate: > 0.20 mL/min

(112-7063)	ea
Contains: 1 compressor nebulizer machine, 1 disposable nebulizer, 1 connection cannula or air tube, 1 mouthpiece & 5 filters.	



Nebulizer Kit

Nebulizer Kit, Adult (112-5617)	ea
Contains: 1 of each: mouthpiece, tee adapter, medicine chamber & tubing.	
Nebulizer Kit with Mask, Adult (112-5618)	ea
Contains: 1 of each: adult aerosol mask, medicine chamber & tubing.	
Nebulizer Kit with Mask, Pediatric (112-5619)	ea
Contains: 1 of each: pediatric aerosol mask, medicine chamber & tubing.	



Digital Floor Scale

Compact digital floor scale with a large backlit LCD display. Weighs accurately in both lb and kg, up to 400 lb (200 kg). Features antislip feet and a low profile for added safety. Automatically turns on. Includes 1-year warranty. 14" x 11⁵/₈" x 1"

(112-7153)..... ea



Mechanical Beam Scale with Height Rod

Standing mechanical beam scale with built-in height rod. Weight capacity of 450 lb (200 kg). 1-year warranty. 22³/₄" x 18¹/₄" x 58¹/₄"

(112-7157)..... ea



Digital Remote Indicator Scale

- Convenient 1¹/₂" backlit remote display with 6' long cable connection to base
- Stabilizing software ensures accurate weight even when patient moves
- Auto-Hold feature shows weight even after patient leaves scale
- 15" x 14.5" x 2.3", 500-lb maximum capacity
- 2-year limited warranty

(570-0596)..... ea



Digital Electronic Eye-Height Scale

Compact, eye-height digital scale that features an integrated mechanical height rod, large backlit LCD display, and BMI calculator. Weight capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2¹/₂" high. Sta-Bilt stabilizing software ensures an accurate measure, even if the patient is moving. Autohold feature allows the weight to continue to be displayed after the patient has stepped off the scale. Integrated USB interface allows data to be transferred to a computer. Includes 2-year warranty. 22" x 15" x 53"

(112-7154)..... ea



Digital Wheelchair Scale

- Single, stable platform allows scale access with wheelchair
- Easily folded to minimize footprint and moved with casters
- Individual wheelchair weight can be programmed for easy recall
- 32" x 33" x 2¹/₂", 1000-lb maximum capacity
- 2-year limited warranty

(570-0597)..... ea



Digital Pediatric Table Scale

Small countertop pediatric scale with a curved design to support and protect the infant. Sta-Bilt stabilizing software ensures an accurate measure, even if the infant is moving. Large backlit LCD display measures in lb and oz. Weight capacity of 45 lb. Integrated tape measure to determine the baby's size. Autohold feature allows the weight to be displayed, even after the baby has been taken off the scale. Integrated USB interface allows data to be transferred to a computer. Includes a 2-year warranty. 16³/₄" x 27³/₄" x 4"

(112-7156)..... ea



Digital Waist-High Scale

Compact, waist-high digital scale that features a large, backlit LCD display and shows lb or kg. Capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2¹/₂" high. Wheels for easy movement. Contains BMI calculator and clock. Includes 2-year warranty. 22" x 15" x 36"

(112-7155)..... ea



Digital Hand Rail Scale Height Rod

- Handrails allow patients to hold steady during weighing
- Low-profile design for easy access, platform is only 2³/₄" high
- Contains BMI calculator and clock
- 28" x 23¹/₂" x 2³/₄", 1000-lb maximum capacity
- 2-year limited warranty

(570-0598)..... ea



Center-Post Colposcopes

- Continuous zoom (3.5X–24X)
- Trinoc tube with built-in C-mount camera adapter
- Dual-diopter adjustment
- Center-post mount
- Quartz halogen, fiber-optic illumination with intensity control and green filter

Binocular Head (570-0148)..... ea

Trinocular Head, Camera-Ready Zoom (570-0149)..... ea

Portable Cryosurgical V-Type Kits



- Safe dispensing system
- More fast-charging treatment buds than competing brands
- Short treatment times
- 5-year shelf life
- Effective treatment for the removal of skin tags, common warts, plantar warts, actinic keratosis, genital warts, molluscum contagiosum, lentigo and seborrhic keratosis

Single Canister, 162 mL
Comes with: 40 Round (5 mm) Treatment Buds, 20 Arrow (2 mm) Treatment Buds & 6 Treatment Cones (3 mm, 5 mm, 7 mm, 9 mm, 12 mm & 16 mm) **(570-1138)**..... ea

Replacement Canister, 162 mL
Comes with: 20 Treatment Buds & 3 Treatment Cones **(570-1139)**..... ea

Replacement Canister, 236 mL
Comes with: 20 Treatment Buds & 3 Treatment Cones **(570-1140)**..... ea

Portable Cryosurgical H-Type Kits



- 88 mL in each canister
- More fast-charging treatment buds than competing brands
- Safe dispensing system
- Short treatment times
- 3-year shelf life
- Effective treatment for the removal of skin tags, common warts, plantar warts, actinic keratosis, genital warts, molluscum contagiosum, lentigo, and seborrhic keratosis

Single Canister
Comes with: 30 Round (5 mm) Treatment Buds & 20 Arrow (2 mm) Treatment Buds **(570-1141)**..... ea

Dual Canister
Comes with: 50 Round (5 mm) Treatment Buds & 30 Arrow (2 mm) Treatment Buds **(570-1137)**..... ea

High-Temperature (2,200°F), Disposable Cauteries, Cordless

- Lightweight
- Convenient
- Cordless

#AA01, Fine Tip Cautery **(112-7177)**..... ea
#AA03, Loop Tip Cautery **(112-7178)**..... ea
#AA21, Elongated Vasectomy Tip **(112-7180)**..... ea



High-Frequency Desiccator

Designed and engineered to be the most reliable and durable desiccator available today, this unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat. Post-procedure and long-term cleanup are simplified by its nonabsorbing, nonstaining, high-impact plastic case.

- 40-watt high-frequency desiccator
- Adjustable in 1/10-watt increments in the first 10 watts, in monopolar and bipolar
- More accurate power with Universal Power Control 100-240-volt AC
- Less to do and less to buy makes this the most cost-effective and user-friendly unit on the market
- This desiccator is a ground reference unit, which means no grounding pad is required
- 4-year warranty

#A942 **(570-0583)**..... ea



Electrosurgical Generator and High-Frequency Desiccator

This electrosurgical generator is a high-frequency desiccator with cutting ability. One half of the unit is the Henry Schein High Frequency Desiccator. The other half is an electrosurgical generator that features 50 watts of cutting power and coag capabilities—plenty of power for most cutting procedures performed in the doctor's office.

- 40-watt high-frequency desiccator
- 2 machines in 1 convenient package
- Adjustable in 1/10-watt increments in the first 10 watts, in monopolar and bipolar
- No grounding pad is required during desiccation mode
- 4-year warranty

#A952 **(570-0582)**..... ea



Diagnostic Electrode Resting Tabs

- Conductive adhesive hydrogel for excellent adhesion and clear, reliable tracings
- Silver/silver chloride (Ag/AgCl) sensing system is defibrillation-recoverable
- Not made with natural rubber latex
- DEHP-free, mercury-free

(570-1173)..... 100/pkg



ECG/EKG Paper



- #007976, Z-Fold, Permatrace, Use with Burdick, 216 mm x 150' **(102-8243)**..... 200 sheets/box
- #007981, Z-Fold, Permatrace, Use with Burdick, 216 mm **(102-5507)**..... 300 sheets/box
- #007989, Z-Fold, Permatrace, Use with Burdick, 216 mm **(102-5506)**..... 300 sheets/box
- #M1707A, Z-Fold, Permatrace, Use with Hewlett-Packard M2481A, 216 mm **(102-9711)**..... 200 sheets/box
- #007984, Compare to Burdick#7160237-00, Z-Fold, Permatrace, Use for Atria, Burdick Eclipse Plus, and Quest, 216 mm **(900-4201)**..... 200 sheets/box
- #007868, Z-Fold, Use with Burdick, 216 mm **(102-5508)**..... 200 sheets/box
- #007944, Z-Fold, Use with Burdick or #18416-002, Quinton Instruments Q750B, Q4500, 216 mm x 150' Roll **(100-9529)**..... ea
- #007966, Z-Fold, Use with Burdick, 216 mm **(102-6077)**..... 200 sheets/box
- #007979, Z-Fold, Use with Burdick, 216 mm **(102-7659)**..... 300 sheets/box
- #15022-001, Z-Fold, Use with Quinton Instruments, Q3000, Q4000, 216 mm **(102-5576)**..... 200 sheets/box
- #2157-012A, Z-Fold, Use with Schiller AT10, SP10 or #94010-000, Welch Allyn AT10, SP10, 210 mm **(102-5511)**..... 150 sheets/box
- #2157-017A, Z-Fold, Use with Schiller AT2 or #94002-000, Welch Allyn AT2, 210 mm **(102-8240)**..... 150 sheets/box
- #2157-014A, Z-Fold, Use with Schiller AT1, SP1 or #94001-000, Welch Allyn AT1, SP1, 90 mm **(102-8961)**..... 280 sheets/box
- Use with Marquette (GE) 9402-024, E9001E, DF, 9402-037, 9402-040 **(102-6883)**..... 300 sheets/pad

Henry Schein® Ultrasound Printer Paper

Provides high-quality, video-printed paper with optimal grayscale and resolution. Thermal paper with a heat-sensitive coating that ensures fast, accurate electrocardiographic tracing. Compatible with both Sony® and Mitsubishi® printers.

- Paper..... 5/box
Specify:
Type I (Standard Gloss), 110 mm x 20 mm, 72-µm Thick (±2µm) **(570-5445)**
Type V (High Gloss), 110 mm x 18 mm, 72-µm Thick (±2µm) **(570-5447)**
Type II (High Gloss), 110 mm x 20 mm, 80-µm Thick (±5µm) **(575-5446)**





Everyday Exam Capes

3-Ply Tissue, 30" x 21"

- Durable and soft
- Front or back opening

White (100-8821)	100/case
Blue (100-2393)	100/case
Mauve (100-2010)	100/case



Premium Exam Capes

Tissue/Poly/Tissue, 30" x 21"

These embossed capes, with a layer of poly film sandwiched between soft absorbent tissue, provide maximum coverage where the exam requires a waist-length garment. The layer of poly provides added strength, making them suitable for most office procedures.

White (101-6908)	100/case
Mauve (101-3029)	100/case
Blue (101-6599)	100/case




Pink Ribbon Exam Capes

Tissue/Poly/Tissue, 30" x 21"

Provide patient modesty for upper-body exams and X-rays.

(900-4278)	100/case
-------------------------	----------



Patient Belongings Bags

- Convenient bags for moving and storing patient belongings
- Plastic, opaque white
- Print: Patient Belongings, Name, Room
- 20" x 20" x 4"

Bags with Handles (570-1179)	250/case
Bags with Drawstrings (570-1180)	250/case



Everyday Exam Drapes

2-Ply, All-Virgin Tissue

Designed to blanket the patient in lightweight comfort.

White, 36" x 48" (101-3658)	100/case
White, 40" x 48" (100-1797)	100/case
Blue, 40" x 48" (100-3360)	100/case
Mauve, 40" x 48" (100-6428)	100/case
White, 40" x 60" (101-8521)	100/case



Premium Exam Drapes

3-Ply, All-Virgin Tissue

- Soft, opaque, and heavyweight
- Embossed tissue drape with dependable wet strength is ideal for short to moderate exams
- Conforms to the body and offers the ultimate in absorbency

White, 40" x 48" (101-7084)	100/case
---	----------



Pink Ribbon Exam Drapes

2 Ply, 40" x 48"

Provide patient draping for privacy.

(900-4275)	100/case
-------------------------	----------



Disposable Exam Shorts

- Provide protection and security during examination
- Made from SMS material
- Elastic waistband made from latex silk and low-stretch yarn
- Individually wrapped
- Blue

Shorts.....	50/case
Specify:	
Small/Medium (570-2757)	2X-Large (570-2759)
Large/X-Large (570-2758)	



Everyday Exam Gowns

3-Ply Tissue, 30" x 42"

- Soft tissue provides patient comfort
- Embossed pattern provides extra strength
- Plastic ties attached

White (100-4471)	50/case
Blue (100-4975)	50/case
Mauve (100-7115)	50/case



Premium Exam Gowns

Tissue/Poly/Tissue, 30" x 42"

- Layer of poly provides strength; layers of tissue absorb moisture
- Opens front or back with attached elastic ties

White (100-2392)	50/case
Blue (100-5728)	50/case
Mauve (100-9652)	50/case



Pink Ribbon Exam Gowns

Tissue/Poly/Tissue

Provide patient modesty for full-body exams and X-rays.

30" x 42" (900-4277)	50/case
--------------------------------------	---------



Safety Slip Patient Footwear

- Provide slip-resistant, nonskid protection
- Blend of cotton and polyester yarns to provide premium comfort and warmth
- Smooth side, single-tread
- Made to fit left or right foot
- Multiple sizes: all are color coded for easy size identification
- Not made with natural rubber latex

Green, One Size (570-3417)	48 pair/box
Beige, Large (112-6102)	48 pair/box
Grey, X-Large (112-6103)	48 pair/box
Mint Green, 2X-Large (112-6105)	48 pair/box



Table Paper

- White
- Quiet, absorbent barrier
- Designed and manufactured to provide strength, protection, and value

Smooth, 225 ft
18" 12 rolls/case
(900-4978)
21" x 225 ft
(900-4979) 12 rolls/case



Poly-Perf® Exam Table Paper

- Absorbent top tissue layer with poly back for fluid barrier
- For exam table coverage; can also be used as scale and table liners
- Perforated every 18" or 21"
- 125-ft rolls
- White

18" 9/case
(570-0728)
21" 9/case
(570-0727)



Pink Ribbon Table Paper

Smooth, 21" x 225-ft Rolls
(900-4276) 12/case



Henry Schein® Stretcher Sheets

Reliable, disposable stretcher sheets designed to eliminate the need for reusable liners. Absorbent 2-ply construction features a white tissue layer and a blue poly-backed layer that provides added protection. Flat, smooth design offers comfort without bulk and fits most standard stretchers and treatment surfaces. Single-use design helps reduce cross contamination. Not made with natural rubber latex.

40" x 90", White/Blue, 2-Ply Tissue/Poly
(570-5555) 50 per Case



Premium Table Paper Rolls

- Designed and manufactured to provide strength, protection, and value
- White

Crepe, Standard, 125 ft
18" 12/case
(100-9865)
21" 12/case
(100-1657)

Smooth, Standard, 225 ft
14" 12/case
(104-4501)
18" 12/case
(100-9491)
21" 12/case
(100-0177)

Smooth, Heavy Duty, 225 ft
18" 12/case
(104-6092)
21" 12/case
(104-5273)



Pillowcases

- Tissue/Poly/Tissue
- Provide linen-like softness for patient comfort and an inner poly barrier for extra protection
 - Strong and durable

White, 21" x 30"
(104-9206) 100/case



Premium Solacel® Pillowcases

- 1-Ply Tissue/Poly
- Disposable pillowcase is super soft and strong because it is made from Fabricel® material
 - By bonding supple tissue to polyethylene film with a heat-lamination process, a unique material that is soft, strong, and tear resistant is created
 - Deluxe material prevents soak-through

White, 21" x 30"
(101-2526) 100/case



Headrest Covers

- Tissue/Poly
- Tissue outer surface provides greater comfort for your patients
 - Poly lined to help protect your equipment

Small Size
10" x 10" 500/case
Specify:
White **(101-7104)** Mauve **(101-3198)**
Blue **(101-4595)** Green **(102-8451)**
Lavender **(102-5543)**

Regular Size (13" Open Side)
10" x 13" 500/case
Specify:
White **(101-8627)** Mauve **(101-1051)**
Blue **(101-8955)** Green **(102-5542)**
Lavender **(102-9374)**



Disposable Headrest Covers

- Plastic
- Disposable poly barrier that acts as a barrier to protect patient and equipment from cross contamination
 - Economical and convenient, space-saver dispenser box
 - Made in the USA

Regular Size
11½" x 10" 250/box
Specify:
Clear **(112-4871)** White **(112-4873)**
Jumbo Size
14½" x 10" 250/box
Specify:
Clear **(112-4872)** White **(112-4874)**



Essentials™ Headrest Covers

- Essentials**
- Disposable headrest covers that provide an impervious barrier between the patient and headrest. Available in White and Clear.
- 9½" x 11", Clear Polyethylene
(112-6527) 1000/case



Dri-Gard® Towels

2-Ply Tissue/1-Ply Poly

- Provide positive moisture barrier and quick absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases the strength of the towel
- Made in the USA
- Not made with natural rubber latex

13" x 19" 500/case

Specify:

Aqua	(100-9276)	Lavender	(101-0476)
Beige	(100-3501)	Peach	(100-7384)
Blue	(100-7265)	Rose	(100-6341)
Gray	(100-1313)	White	(100-5715)
Green	(100-1350)	Yellow	(100-9704)



Dri-Gard® Plus Towels

3-Ply Tissue/1-Ply Poly

- Provide a positive moisture barrier with maximum absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases strength of the towel
- Made in the USA
- Not made with natural rubber latex

13" x 19" 500/case

Specify:

Aqua	(100-9851)	Lavender	(101-6810)
Beige	(100-5015)	Peach	(100-7395)
Blue	(100-9825)	Rose	(100-2024)
Gray	(100-9545)	White	(100-7908)
Green	(100-4479)	Yellow	(100-3395)



Dri-Gard® Towels

Teddy Bear Design, 2-Ply Tissue/1-Ply Poly

- Add color to the operator and bring smiles to younger patients
- Made in the USA
- Not made with natural rubber latex

13" x 19" 500/case
(100-9712)



Standard Professional Towels

Unique properties of wet or dry tensile strength, combined with facial softness and super-absorbency.

3-Ply Paper

13" x 19" 500/case

Specify:

White	(100-3161)	Blue	(100-0183)
-------	-------------------	------	-------------------



Contour Bibs

Plastic-Backed, 3-Ply Paper/1-Ply Poly

Provide 30% more coverage than standard dental bibs. 3-ply tissue offers extra absorption, while the poly backing provides durable moisture barrier. Perfect for facial exam procedures.

18" x 22" 400/case

Specify:

Blue, Ultimate	(102-8707)
Green, Ultimate	(100-6332)
Lavender, Ultimate	(103-2185)
White, Ultimate	(100-8260)



Essentials Patient Bibs

Essentials

- Provide a protective barrier from fluids and quick absorbency
- 2-ply poly pattern protects patients' clothing and equipment from moisture
- Available in a cost-effective package of 500 bibs and 10 designer colors
- Not made with natural rubber latex

(2-Ply Tissue/1-Ply Poly)

13" x 18" 500 per Case

Specify:

Aqua	(112-6748)	Lavender	(112-6742)
Beige	(112-6739)	Peach	(112-6741)
Blue	(112-6749)	Silver Gray	(112-6743)
Dusty Rose	(112-6747)	White	(112-6744)
Green	(112-6746)	Yellow	(112-6740)

13" x 19" 800/case

Specify:

Aqua Blue	(112-6521)	White	(112-6526)
Lavender	(112-6524)		



Tray Covers

Manufactured from a heavyweight paper stock that impedes the flow of moisture to minimize liquid seep-through. Size B tray covers are available in 10 designer colors that can be color coordinated in the operator with Henry Schein brand matching patient towels, plastic cups, and headrest covers. Made in the USA. Not made with natural rubber latex.

- Disposable: provide a sanitary means to place instruments on a tray without sterilizing the tray after each patient
- Large variety of sizes: fit all makes and models of trays
- Good quality heavyweight paper: economical

8½" x 11", White
(100-5912) 1000/box

8½" x 12¼", Ritter B, Blue
(101-9336) 1000/box

8½" x 12¼", Ritter B, Mauve
(101-4592) 1000/box

8½" x 12¼", Ritter B, Silver-Gray
(101-7620) 1000/box

8½" x 12¼", Ritter B, White
(100-0934) 1000/box

8½" x 12¼", Ritter B
Colors 1000/box

Specify:

Aqua	(102-4081)
Beige	(102-8143)
Green	(102-7181)
Lavender	(102-4082)
Yellow	(102-4117)

9½" x 12¾", Chayes, Weber "A", White
(100-2977) 1000/box

9" x 13½", Certified, SS White, White
(100-8128) 1000/box

9½" x 13½", Midwest, ADEC E, White
(101-1309) 1000/box

10¼" x 15¾", SS White D, White
(100-9489) 1000/box

11" x 17¼", Weber C, White
(100-2550) 1000/box

11" x 17¼", Weber C, Blue
(101-4675) 1000/box

5" x 8", Mini F, White
(101-0839) 1000/box



Protex™ Liquid-Proof Tray Covers

- Polycoated, nonseeping paper that prevents cross contamination from paper to paper
- Stops unnecessary waste
- Use poly side down for greater absorbency; use poly side up for quick mixing
- Stops bleed-through
- Economical: only 1 sheet per patient
- Rectangular
- White

8½" x 11"	(100-2254)	500/box
8½" x 12¼"	(100-4202)	500/box
9½" x 12¾"	(100-5098)	500/box
9" x 13½"	(100-5397)	500/box



Graduated Plastic Medicine Cups

Plastic medicine cups with easy-to-read graduations. 1 oz. capacity.

(112-5555)	100/bag
------------	---------



Essentials™ Disposable Plastic Cups

Essentials

Embossed plastic cup provides enhanced grip and cup strength. Rolled lip allows easy dispensing and drinking.

5 oz.	1000 per Case
Specify:	
Blue (112-6792)	White (112-6791)
Lavender (112-6793)	



Paper Cups

5 oz. Blue, Polycoated	1000/case
(949-1498)	



Translucent Plastic Cups

Embossed drinking cups with rolled lips that are made from heavy, sturdy plastic to prevent crumpling and leaking. Each cup releases effortlessly from dispenser. Embossed wall eliminates slippage, and rolled lip makes drinking easier.

Translucent, 5 oz.	1000/case
(570-3729)	



Paper Cups

Waxed, 5 oz

Pink with Pink Ribbons	1000/case
(900-4211)	



White with Pink Hearts	1000/case
(112-5051)	

HENRY SCHEIN®
BRAND PRODUCTS

DISPOSABLE PLASTIC CUPS



Fits standard dispensers.

Special embossed grip provides greater strength and enhances the holding and dispensing of cups. Unique double lips allow for consistent, single-cup dispensing. Color coordinates with our Henry Schein brand of matching patient towels and tray covers. Not made with natural rubber latex.

- Heavy, sturdy plastic: leakproof; will not crumple in hands
- Releases easily from dispenser: easily accessible
- Embossed wall: nonslip
- High-gloss finish and assorted colors: decorative

5oz Cup - 1000 per Case

Beige	(570-3734)	Lavender	(570-3728)
Blue	(570-3726)	Mauve	(570-3732)
Clear	(570-3729)	Peach	(570-3735)
Gray	(570-3731)	White	(570-3727)
Green	(570-3730)	Yellow	(570-3733)

Order today at HENRYSCHEIN.COM/MEDICAL

Cotton Balls

- 100% pure-white cotton
- Soft and absorbent



Sterile, Medium
(570-3773) 500/box

Non-sterile, Medium
(570-3771) 4000/case

Non-sterile, Large
(570-3772) 2000/case

Cotton Bulk Absorbent Rolls, USP

1-lb Bulk Rolls

Sterile
(570-3769) ea

Non-sterile
(570-3770) ea



Tongue Depressors

- Highly polished
- Smooth
- Designed for oral use or topical application



Non-sterile
 Junior, 5½"
(100-0758) 500/box

Adult, 6"
(100-2416) 500/box

Sterile
 Adult, 6"
(100-0238) 100/box

Cotton Rolls

Starched, 100% cotton, absorbent nasal roll.

#2 Medium, Unwrapped, Non-sterile, 1½" x ¾"
(102-5455) 2000/box

#2 Medium, Wrapped, Non-sterile, 1½" x ¾"
(570-3767) 2000/box

#2 Medium, Wrapped, Sterile, 1½" x ¾"
(570-3768) 2000/box



Essentials™ Cotton Rolls
✓Essentials

Soft, pliable, and easy to cut with excellent absorbency. Available in a cost-effective package of 2500.

Medium #2-1½" x ¾"
 Non-sterile
(112-6917) 2500/bag



Alcohol Prep Pads

Sterile

- Soft, absorbent, and non-woven pads impregnated with 70% isopropyl alcohol USP
- Primarily used as a topical antiseptic; disinfects prior to administering injections

Medium, 30 mm x 30 mm, Folded
(570-2219) 200/box

Large, 70 mm x 62.5 mm
(570-1307) 200/box



Premium Cotton Rolls

- Made from 100% cotton for greater absorbency
- More absorbent than the leading brands
- Longer fibers reduce linting
- Longer bleaching process produces a whiter cotton roll
- Can easily be cut into smaller pieces with scissors for effective use in tight oral-cavity areas
- Nonlinting, non-sterile, and pliable
- Maintains shape when wet
- Will not adhere to mucous membranes

Non-sterile
(896-9286) 2000/box

Sterile
(104-8192) 2000/box



Cotton-Tipped Applicators

Applicators are highlighted by a soft cotton tip that is absorbent and will not lint.

Non-sterile, 3"
(100-6015) 1000/box

Non-sterile, 6"
(100-9175) 1000/box

Sterile, 6", 100 Packs of 2
(100-9249) 200/box



Wooden Applicator/Mixing Sticks

Plain Wood, 6"

- International peel symbols: ensure easy opening from proper end of package
- Long shelf life (more than 10 years): economical

(100-9485) 1000/box



Emesis Basins

- Plastic
- Kidney shaped
- 16-oz capacity

Mauve, 8½"
(112-5809) 25/box





C-Fold Towels

- 1-ply, strong, and absorbent
- Individually folded towels
- Whitened without chemicals containing chlorine bleach
- Fits a wide range of universal dispensers
- 150 sheets/pack

White, 12.87" x 10.5"2400/case
(570-0704).....2400/case
 Contains: 150 sheets per pack, 16 packs per case.



Multifold Towels

- White
- 1-ply, strong and absorbent
- Interlocking technology
- Whitened without chemicals containing chlorine bleach
- Fits wide range of universal dispensers
- 250 sheets/pack

White, 9½" x 9¼"4000/case
(570-0703).....4000/case
 Contains: 250 sheets per pack, 16 packs per case.



Disposable Airlaid Washcloths

- Airlaid is made from synthetic staple fibers and natural cellulose pulp
- Offers absorbency, strength, and softness
- Ideal for incontinence cleanups, bathing, and other patient applications

White, 10" x 10"1000/case
(570-2403).....1000/case
 Contains: 50 washcloths per bag, 20 bags per case.



Facial Tissues

- 2-Ply
- Soft, absorbent, value-oriented, white facial tissue.

White, 9" x 8"100/box
(900-4214).....100/box



Protection Plus Disposable Underpads

- Made of waterproof, blue poly backing
- Absorbent core of tissue layers
- White, hypoallergenic spunbond facing sheet
- Not made with natural rubber latex

17" x 24"300/case
(112-5633).....300/case
 23" x 24"200/case
(112-5634).....200/case
 23" x 36"150/case
(570-1182).....150/case
 30" x 30"150/case
(570-1183).....150/case
 30" x 36"100/case
(570-1181).....100/case



HENRY SCHEIN® STRETCHER SHEETS

- Easy-to-use, reliable disposable solution that eliminates the need for reusable linens
- Absorbent 2-ply design with a white tissue layer and blue protective poly-backed layer
- Single-use to help reduce cross-contamination
- Not made with natural rubber latex



Henry Schein® Stretcher Sheets(570-5555)
 50 per Case, 40" x 90", White/Blue, 2-ply Tissue/Poly



Order today at HENRYSCHEIN.COM/MEDICAL

FROM HERE



TO HERE

WEAR IT WITH CONFIDENCE.

CRITERION™

Criterion offers a wide range of quality gloves for your everyday needs. From nitrile to surgical, we've got you covered.



by HENRY SCHEIN®



ORDER NOW AT HENRYSCHEIN.COM/CRITERION

CRITERION™

Gloves Comparison Chart

GLOVE NAME	TYPE	COLOR	STERILE STATUS	NOT MADE WITH NATURAL RUBBER LATEX	SURFACE TEXTURE	THICKNESS AT PALM (MM/MIL)	ADDITIONAL FEATURES	ITEM CODES
NITRILE								
CRITERION™ Nitrile Exam N100 Gloves	Exam	Blue	Non-Sterile	Yes	Finger Textured	0.08/3.15	Chemotherapy tested to ASTM D6978 standard.	100/Box: XS (111-8533), S (111-8535), M (111-8536), L (111-8537) 90/Box: XL (111-8539)
CRITERION™ Nitrile Exam N200 Gloves	Exam	Periwinkle Blue	Non-Sterile	Yes	Finger Textured	0.06/2.36	Chemotherapy tested to ASTM D6978 standard.	200/Box: XS (900-7437), S (900-7438), M (900-7439), L (900-7440) 180/Box: XL (900-7441)
	Exam	Sky Blue	Non-Sterile	Yes	Finger Textured	0.06/2.36	Chemotherapy tested to ASTM D6978 standard.	200/Box: XS (570-2097), S (570-2098), M (570-2099), L (570-2100) 180/Box: XL (570-2101)
	Exam	Black	Non-Sterile	Yes	Finger Textured	0.06/2.36	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard.	200/Box: XS (570-3487), S (570-3488), M (570-3489), L (570-3490) 180/Box: XL (570-3491)
CRITERION™ Nitrile Exam N300 Gloves	Exam	Ice Blue	Non-Sterile	Yes	Finger Textured	0.05/1.97	Chemotherapy tested to ASTM D6978 standard. Eco-friendly Packaging.	300/Box: XS (112-6829), S (112-6830), M (112-6831), L (112-6832) 250/Box: XL (112-6833)
	Exam	Electric Blue	Non-Sterile	Yes	Finger Textured	0.05/1.97	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard. Eco-friendly Packaging.	300/Box: XS (570-3540), S (570-3541), M (570-3542), L (570-3543) 250/Box: XL (570-3544)
	Exam	Black	Non-Sterile	Yes	Finger Textured	0.05/1.97	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard. Eco-friendly Packaging.	300/Box: XS (570-3482), S (570-3483), M (570-3484), L (570-3485) 250/Box: XL (570-3486)

by HENRY SCHEIN®

For additional information visit: HENRYSCHEIN.COM/CRITERION

GLOVE NAME	TYPE	COLOR	STERILE STATUS	NOT MADE WITH NATURAL RUBBER LATEX	SURFACE TEXTURE	THICKNESS AT PALM (MM/MIL)	ADDITIONAL FEATURES	ITEM CODES
NITRILE (CONT.)								
CRITERION™ Nitrile Exam Colloidal Oatmeal Gloves	Exam	White	Non-Sterile	Yes	Finger Textured	0.05/1.97	Colloidal Oatmeal Coated.	<u>200/Box:</u> XS (112-6220), S (112-6221), M (112-6222), L (112-6223), XL (112-6224)
CRITERION™ Nitrile Exam Low Derma Gloves	Exam	Slate Blue	Non-Sterile	Yes	Finger Textured	0.05/1.97	Fentanyl tested. Chemotherapy tested to ASTM D6978 standard. Low Dermatitis Potential.	<u>200/Box:</u> XS (570-0630), S (570-0631), M (570-0632), L (570-0633) <u>180/Box:</u> XL (570-0634)
CRITERION™ Nitrile Exam Aloe + Vitamin E Gloves	Exam	Magenta	Non-Sterile	Yes	Finger Textured	0.05/1.97	Aloe Vera + Vitamin E Coated. Polymer Coated.	<u>200/Box:</u> XS (112-7143), S (112-7144), M (112-7146), L (112-7147)
	Exam	Lavender	Non-Sterile	Yes	Finger Textured	0.05/1.97	Aloe Vera + Vitamin E Coated. Polymer Coated.	<u>200/Box:</u> XS (570-3545), S (570-3546), M (570-3547), L (570-3548), XL (570-3549)
CRITERION™ Nitrile Mint Scented	Exam	Blue	Non-Sterile	Yes	Finger Textured	0.08/3.15	Mint Scented.	<u>100/Box:</u> XS (570-3535), S (570-3536), M (570-3537), L (570-3538), XL (570-3539)
POLYCHLOROPRENE								
CRITERION™ Chloroprene Exam Ultrathin Gloves	Exam	Lime	Non-Sterile	Yes	Finger Textured	0.05/1.97	Polymer Coated.	<u>100/Box:</u> XS (112-7306), S (112-7307), M (112-7308), L (112-7309), XL (112-7310)
CRITERION™ Chloroprene Exam Gloves	Exam	Green	Non-Sterile	Yes	Finger Textured	0.08/3.15	Polymer Coated.	<u>100/Box:</u> XS (112-6857), S (112-6858), M (112-6859), L (112-6860) <u>90/Box:</u> XL (112-6861)
	Exam	Pink	Non-Sterile	Yes	Finger Textured	0.08/3.15	Polymer Coated.	<u>100/Box:</u> XS (570-3532), S (570-3533), M (570-3534), L (570-3572) <u>90/Box:</u> XL (570-3573)
CRITERION™ Chloroprene Soft Exam Gloves	Exam	Aqua	Non-Sterile	Yes	Micro-textured fingertips	0.09/3.54	Polymer Coating, Criterion's softest chloroprene exam glove	<u>100/Box:</u> XS (570-5400), S (570-5401), M (570-5402), L (570-5403), XL (570-5404)
CRITERION™ Polychloroprene Surgical Gloves	Surgical	White	Sterile	Yes	Palm to Fingertip Textured	0.15/5.91	Polymer Coated. Extended Cuff.	<u>50 Pair/Box:</u> 5.5 (112-7079), 6 (112-7080), 6.5 (112-7081), 7 (112-7082), 7.5 (112-7083), 8 (112-7085), 8.5 (112-7086), 9 (112-7088)
CRITERION™ Polyisoprene Surgical Gloves	Surgical	White	Sterile	Yes	Palm to Fingertip Textured	0.21/8.27	Polymer and Polyurethane Coated. Extended Cuff.	<u>50 Pair/Box:</u> 5.5 (570-1119), 6 (570-1120), 6.5 (570-1121), 7 (570-1122), 7.5 (570-1123), 8 (570-1124), 8.5 (570-1125), 9 (570-1126)
LATEX								
CRITERION™ Latex Aloe & Vitamin E Gloves	Exam	White	Non-Sterile	No	Finger Textured	0.08/3.15	Aloe Vera + Vitamin E Coated. Polymer Coated. Chlorinated.	<u>100/Box:</u> XS (570-3448), S (570-3449), M (570-3450), L (570-3451) <u>90/Box:</u> XL (570-3452)
CRITERION™ Latex Exam Polymer-Coated Gloves	Exam	White	Non-Sterile	No	Finger Textured	0.08/3.15	Polymer Coated. Chlorinated.	<u>100/Box:</u> XS (570-3412), S (570-3413), M (570-3414), L (570-3415) <u>90/Box:</u> XL (570-3416)
CRITERION™ Latex Exam Chlorinated Gloves	Exam	White	Non-Sterile	No	Palm to Fingertip Textured	0.11/4.33	Chlorinated.	<u>100/Box:</u> XS (102-5354), S (102-8296), M (102-5355), L (102-7763), XL (102-5357)
CRITERION™ Latex Extra Grip Gloves	Exam	White	Non-Sterile	No	Palm to Fingertip Textured	0.10/4.33	Polymer Coated. Non Chlorinated.	<u>100/Box:</u> XS (570-3527), S (570-3528), M (570-3529), L (570-3530) <u>90/Box:</u> XL (570-3531)
CRITERION™ Latex Surgical Gloves	Surgical	White	Sterile	No	Smooth	0.14/5.51	Polymer Coated. Extended Cuff. Non Chlorinated.	<u>50 Pair/Box:</u> 5.5 (570-4495), 6 (570-4496), 6.5 (570-4497), 7 (570-4498), 7.5 (570-4499), 8 (570-4500), 8.5 (570-4501), 9 (570-4502)
VINYL								
CRITERION™ Vinyl Exam Gloves	Exam	Clear	Non-Sterile	Yes	Smooth	0.08/3.15	Smooth. Accelerator Free.	<u>100/Box:</u> XS (570-0087), S (570-0088), M (570-0089), L (570-0090) <u>90/Box:</u> XL (570-0091)



**CRITERION™
Nitrile N100**

Nitrile Powder-Free Exam
Durable nitrile construction for reliable protection and tactility. Chlorinated for reduced friction and easy donning. Soft and stretchy, this nonchlorinated glove is the ideal alternative to latex. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 3.15-mil thickness at the palm. Blue.

Gloves.....100/box
Specify:
X-Small **(111-8533)** Large **(111-8537)**
Small **(111-8535)** X-Large* **(111-8539)**
Medium **(111-8536)**

*90/box



CRITERION™ Nitrile N200

Nitrile Powder-Free Exam
Slightly thinner than the N100, the durable nitrile construction provides reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.36-mil thickness at the palm.

Periwinkle Blue200/box
Specify:
X-Small **(900-7437)** Large **(900-7440)**
Small **(900-7438)** X-Large* **(900-7441)**
Medium **(900-7439)**

Sky Blue.....200/box
Specify:
X-Small **(570-2097)** Large **(570-2100)**
Small **(570-2098)** X-Large* **(570-2101)**
Medium **(570-2099)**

*180/box



**CRITERION™
Nitrile N200**

Nitrile Powder-Free Exam
Slightly thinner than the N100, this glove's durable nitrile construction provides reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy and fentanyl tested to the ASTM D6978 standard. Beaded cuff. 2.7-mil thickness at the palm. Black.

Gloves.....200/box
Specify:
X-Small **(570-3487)** Large **(570-3490)**
Small **(570-3488)** X-Large* **(570-3491)**
Medium **(570-3489)**

*180/box

**CRITERION™
Nitrile N300**

Nitrile Powder-Free Exam
Ultrathin formulation for excellent tactile sensitivity and dexterity. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy and fentanyl tested to the ASTM D6978 standard. Ecofriendly packaging—with 300 gloves in a box—reduces waste and saves space. Beaded cuff. 2.0-mil thickness at the palm. Black.

Gloves.....250/box
Specify:
X-Small **(570-3482)** Large **(570-3485)**
Small **(570-3483)** X-large* **(570-3486)**
Medium **(570-3484)**

*250/box



CRITERION™ Nitrile N300

Nitrile Powder-Free Exam, Non-Sterile
Ultrathin formulation for excellent tactile sensitivity and dexterity. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy tested to the ASTM D6978 standard. Ecofriendly packaging. With 300 gloves in a box, reduces waste and saves space. Beaded cuff. 2.0-mil thickness at the palm.

Ice Blue250/box
Specify:
X-small **(112-6829)** Large **(112-6832)**
Small **(112-6830)** X-large* **(112-6833)**
Medium **(112-6831)**

Electric Blue.....250/box
Specify:
X-small **(570-3540)** Large **(570-3543)**
Small **(570-3541)** X-large* **(570-3544)**
Medium **(570-3542)**

*250/box

**CRITERION™ Nitrile
Colloidal Oatmeal**

Nitrile Powder-Free Exam
Made from 100% all-natural colloidal oatmeal, clinically proven to act as a skin protectant and repair damaged skin. Provides continuous protection from prolonged glove use and hand washing. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Gluten-free. Beaded cuff. 2.0-mil thickness at the palm. White.

Gloves.....200/box
Specify:
X-small **(112-6220)** Large **(112-6223)**
Small **(112-6221)** X-large **(112-6224)**
Medium **(112-6222)**



**CRITERION™
Nitrile Low Derma**

Nitrile Powder-Free Exam
Made from a soft, durable nitrile formulation without the use of chemical accelerators or sulfur, preventing Type-IV allergic reactions. Reduces the potential for dermatitis for glove wearers, yet still provides softness and comfort and maintains superior protection. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.0-mil thickness at the palm. Slate Blue.

Gloves.....200/box
Specify:
X-small **(570-0630)** Large **(570-0633)**
Small **(570-0631)** X-large* **(570-0634)**
Medium **(570-0632)**

*180/box



CRITERION™ Nitrile Aloe + Vitamin E

Nitrile Powder-Free Exam
Extra-soft stretch with an enhanced inner coating of aloe and vitamin E to soothe, moisturize, and rejuvenate hands. Fast and easy donning. Comfortable fit and feel to provide a superior range of motion. Textured fingers promote sure handling of instruments and materials. Beaded cuff. 1.97-mil thickness at the palm.

Magenta.....200/box
Specify:
X-small **(112-7143)** Medium **(112-7146)**
Small **(112-7144)** Large **(112-7147)**

Lavender.....200/box
Specify:
X-small **(570-3545)** Medium **(570-3547)**
Small **(570-3546)** Large **(570-3548)**

**CRITERION™
Mint Scented**

Nitrile Powder-Free Exam
Mint-scented glove that is chlorinated for easy donning. Non-sterile and powder-free. Textured fingertips provide increased grip of instruments and materials. Beaded cuff. 3.15-mil/0.08mm thickness at the palm. Blue.

Gloves.....100/box
Specify:
X-small **(570-3535)** Large **(570-3538)**
Small **(570-3536)** X-large **(570-3539)**
Medium **(570-3537)**



Essentials

VALUE-PRICED PRODUCTS FOR EVERYDAY USE.



Essentials™ Nitrile 200

Nitrile Powder-Free Exam

Great fit, feel, tactile sensitivity, strength, and value in a latex-free glove. Textured fingertips. 2.4-mil thickness. Periwinkle.

Gloves.....200/Box

Specify:

XS	(112-6765)	L	(112-6768)
S	(112-6766)	XL**	(112-6769)
M	(112-6767)		

**180/Box



Essentials™ 300 Nitrile

Nitrile Powder-Free Exam

Soft and strong nitrile glove that offers increased dexterity while providing maximum protection. Easy Release Technology allows for single-glove dispensing with fewer wrinkles. Textured fingertips deliver improved tactile sensitivity. Chemotherapy tested: appropriate for handling cytotoxic chemicals. 300-count box saves space and reduces waste. Not made with natural rubber latex. 2.0-mil thickness. Indigo.

Gloves.....300/box

Specify:

XS	(570-3631)	L	(570-3635)
S	(570-3632)	XL	(570-3634)
X-Large (250/box)	(570-3633)		

CRITERION

COMFORT IN PROTECTION™



CRITERION™ Chloroprene Ultrathin

Powder-Free Chloroprene Exam, Non-Sterile Ultrathin chloroprene glove that offers the best of both worlds: the comfort and fit of latex, but with the synthetic strength of nitrile. Excellent dexterity for precision procedures and extended-wear comfort. Latex-like feel for enhanced softness and better fit. Strong puncture-resistance and tear-resistance for high-risk procedures. Textured fingertips for a superior grip. Proprietary polymer coating for "glide-on" donning. 2.0-mil thickness at the palm. Beaded cuff. Lime.

Gloves.....100/box

Specify:

X-small	(112-7306)	Large	(112-7309)
Small	(112-7307)	X-large	(112-7310)
Medium	(112-7308)		

CRITERION

COMFORT IN PROTECTION™



CRITERION™ Chloroprene

CRITERION

Chloroprene Powder-Free Exam

Synthetic gloves that offer the comfort and fit of latex, but with the strength of nitrile. Plus, they provide enhanced dexterity for precise procedures and extended-wear comfort. Latex-like feel: the softest CRITERION chloroprene glove. Strong puncture- and tear-resistance for high-risk procedures. Textured fingertips for grip. Polymer-coated for easy, glide-on donning, even when hands are wet. 3.54-mil thick at the palm. Beaded cuff is 240 mm. Not made with natural rubber latex. Aqua Blue.

Gloves.....100/box

Specify:

X-Small	(570-5400)	Large	(570-5403)
Small	(570-5401)	X-Large	(570-5404)
Medium	(570-5402)		



CRITERION™ Chloroprene

Chloroprene Powder-Free Exam

Offers the best of both worlds: the comfort and fit of latex, but with the synthetic strength of nitrile. Excellent dexterity for precision procedures and extended-wear comfort. Latex-like feel for enhanced softness, extended wear comfort, and better fit. Strong puncture-resistance and tear-resistance for high-risk procedures. Textured fingertips for superior grip. Proprietary polymer coating for "glide-on" donning. 3.15-mil thickness at the palm. Beaded cuff.

Green.....100/box

Specify:

X-small	(112-6857)	Large	(112-6860)
Small	(112-6858)	X-large*	(112-6861)
Medium	(112-6859)		

Pink.....100/box

Specify:

X-small	(570-3532)	Large	(570-3572)
Small	(570-3533)	X-large*	(570-3573)
Medium	(570-3534)		

*90/box



CRITERION™ Polychloroprene Surgical

Sterile Polychloroprene Powder-Free Surgeons

Features an anatomically correct design for a more natural fit. Polymer-coating finishes that allow for damp- or dry-hand donning. Textured fingers and palms for enhanced grip and control. Ultralow modulus for enhanced comfort and reduced hand fatigue. Minimum cuff length of 300mm for extended coverage and protection. Beaded cuff to reduce the risk of glove roll-down once worn. Ideal for surgical procedures. 5.91-mil thickness at palm. Sterile. Natural color.

Gloves.....50 pair/box

Specify:

5½	(112-7079)	7½	(112-7083)
6	(112-7080)	8	(112-7085)
6½	(112-7081)	8½	(112-7086)
7	(112-7082)	9	(112-7088)



CRITERION™ Polyisoprene Surgical

Sterile Polyisoprene Powder-Free Surgeons

Made from synthetic polyisoprene, these gloves feature skin-friendly accelerators to reduce allergy risk while providing the fit and feel of natural rubber latex. Anatomical shape minimizes hand fatigue, and the microrough texture enhances grip and control. Patented multilayer coating ensures easy donning, with an extended, beaded cuff and 8.27-mil thickness at the palm. Sterile. Individually wrapped. White.

Gloves.....50 pair/box

Specify:

5½	(570-1119)	7½	(570-1123)
6	(570-1120)	8	(570-1124)
6½	(570-1121)	8½	(570-1125)
7	(570-1122)	9	(570-1126)

CRITERION™
COMFORT IN PROTECTION™



CRITERION™ Latex Aloe + Vitamin E

Latex Powder-Free Aloe Exam
Enhanced inner coating of aloe and vitamin E to soothe and moisturize. Fast and easy donning. Comfortable fit and feel provide superior range of motion. Textured finish for an enhanced grip. Beaded cuff. 3.15-mil thickness at the palm. White.

Gloves.....100/Box

Specify:
X-small **(570-3448)** Large **(570-3451)**
Small **(570-3449)** X-large* **(570-3452)**
Medium **(570-3450)**

*90/box



CRITERION™ Latex Chlorinated

Latex Powder-Free Exam
Offers a comfortable fit and feel for a superior range of motion and a textured finish for an enhanced grip. Chlorinated for reduced friction and lower protein levels. Beaded cuff. 4.33-mil thickness at the palm. White.

Gloves.....100/box

Specify:
X-small **(102-5354)** Large **(102-7763)**
Small **(102-8296)** X-large* **(102-5357)**
Medium **(102-5355)**

*90/box



CRITERION™ Latex Polymer Coated

Latex Powder-Free Exam, Non-Sterile
Finished with a smooth interior polymer coating for fast and easy donning. Comfortable fit and feel for a superior range of motion. Textured finish for an enhanced grip. Beaded cuff. 3.15-mil thickness at the palm. White.

Gloves.....100/box

Specify:
X-small **(570-3412)** Large **(570-3415)**
Small **(570-3413)** X-large* **(570-3416)**
Medium **(570-3414)**

*90/box



CRITERION™ Latex Exam Extra Grip

Latex Powder-Free Exam
Fully textured powder-free latex examination gloves for enhanced grip. Polymer-coated for easy donning. Low protein content minimizes allergies. Non-sterile and ambidextrous for convenience. Beaded cuff. 3.94-mil/0.10-mm thick at the palm. White.

Gloves.....100/box

Specify:
X-small **(570-3527)** Large **(570-3530)**
Small **(570-3528)** X-large **(570-3531)**
Medium **(570-3529)**



CRITERION™ Latex Surgical

Sterile Latex Powder-Free Surgical
Features microtextured fingertips for a secure grip and improved control. Made with soft, elastic latex with an anatomical, curved finger. Provides comfort and minimizes hand fatigue, even during long procedures. Polymer-coated inner surface allows for easy donning when hands are wet, damp, or dry. Beaded cuff prevents rollback. Made with natural rubber latex. 5.51-mil thickness at palm. Minimum cuff length of 280 mm (11") for extended coverage and protection. Hand specific in individually packaged pairs. Natural color.

Gloves.....50 pair/box

Specify:
5½ **(570-4495)** 7½ **(570-4499)**
6 **(570-4496)** 8 **(570-4500)**
6½ **(570-4497)** 8½ **(570-4501)**
7 **(570-4498)** 9 **(570-4502)**



CRITERION™ Vinyl

Powder-Free Vinyl Exam
Non-latex formulation, making it an ideal solution for those at risk of latex sensitivity. Accelerator-free and nonchlorinated. Smooth finish with a supple touch for maximum tactile sensitivity. 3.15-mil thickness at the palm. Beaded cuff. Clear.

Gloves.....100/box

Specify:
X-small **(570-0087)** Large **(570-0090)**
Small **(570-0088)** X-large* **(570-0091)**
Medium **(570-0089)**

*90/box

Glove-Sizing Assistant

X-SMALL/ Sterile: 5–5.5


SMALL/ Sterile: 6–6.5

MEDIUM/ Sterile: 7–7.5

LARGE/ Sterile: 8–8.5

X-LARGE/ Sterile: 9

Place the widest part of your hand over the size grid as shown





CRITERION® Wire Glove Box Dispensers

Formerly known as Maxi-Dispensz®

Durable steel-wire frame with white-powdered finish can be mounted on wall for easy dispensing of gloves. Selection of 5 different choices meets any office needs. Hardware is included. Designed for standard 9½" x 5" x 3" glove boxes.

- A. Dual Box Dispenser**
Holds 1 box horizontally or 2 boxes vertically.
10¾"W x 4"H x 3-1/6"D
(112-5471)..... ea
 - B. Double Horizontal Box Dispenser**
Holds 2 boxes horizontally.
11"W x 11½"H x 4"D
(112-5473)..... ea
 - C. Double Vertical Box Dispenser**
Holds 2 boxes vertically.
11"W x 11½"H x 4"D
(112-5472)..... ea
 - D. Triple Horizontal Box Dispenser**
Holds 3 boxes horizontally.
11"W x 14"H x 4"D
(112-5475)..... ea
 - E. Triple Vertical Box Dispenser**
Holds 3 boxes vertically.
16"W x 8"H x 4"D
(112-5474)..... ea
- (Designed for standard 9½" x 5" x 3" glove boxes.)



CRITERION® Acrylic Glove Box Dispensers

Formerly known as Maxi-Dispensz®

Clear acrylic glove box holder. Available in 4 sizes to meet the needs of any office or facility. Accommodates most brands and sizes. Option to wall mount with predrilled holes (screws NOT included) or place on countertop.

- A. Single Box Holder**
10 15/16"W x 4¾"H x 3½"D
(102-6467)..... ea
- B. Double Box Holder**
12"W x 10"H x 3¾"D
(102-6468)..... ea
- C. Triple Box Holder**
11"W x 14½"H x 3½"D
(102-8344)..... ea
- D. Quadruple Box Holder**
11"W x 19¼"H x 3½"D
(102-7201)..... ea



CRITERION™ by HENRY SCHEIN®



EXPERIENCE THE SOFTEST

Chloroprene Exam Glove Criterion™ has to offer.

- Precision performance provides dexterity and extended-wear comfort
- Microtextured fingertips for a secure grip and improved control during procedures



AQUA BLUE COLOR

SIZE	QUANTITY	ITEM CODE
● XS	100/box	570-5400
● S	100/box	570-5401
● M	100/box	570-5402
● L	100/box	570-5403
● XL	100/box	570-5404

Order today at HENRYSCHEIN.COM/CRITERION

✓ Essentials

VALUE-PRICED GLOVES FOR EVERYDAY USE.

Providing reliable protection and user comfort.



Essentials Nitrile Powder-Free Exam Gloves

Color: Periwinkle Blue

SIZE	QUANTITY	ITEM CODE
XS	200/Box	(112-6765)
S	200/Box	(112-6766)
M	200/Box	(112-6767)
L	200/Box	(112-6768)
XL	180/Box	(112-6769)

Color: Dark Blue

SIZE	QUANTITY	ITEM CODE
XS	300/Box	(112-7321)
S	300/Box	(112-7322)
M	300/Box	(112-7323)
L	300/Box	(112-7324)
XL	250/Box	(112-7325)


Contact your local Henry Schein Representative or visit HENRYSCHEIN.COM/ESSENTIALS to learn about our complete assortment of products.

FROM HERE



TO HERE

WEAR IT WITH
CONFIDENCE.



CRITERION™

COMFORT IN PROTECTION™

A range of smart PPE. Strong and reliable. Easy to use and comfortable to wear. Focus more on your patients and not on your PPE. **Count on it. All the way.**

ORDER NOW AT [HENRYSCHEIN.COM/CRITERION](https://www.henryschein.com/criterion)

by  HENRY SCHEIN®



Not all masks are created equal.

Know the key facts about ASTM-rated mask levels so you can ensure you have the right protection.

The ASTM (American Society for Testing and Materials) updated its medical face mask standards in response to the Covid-19 pandemic. The key changes were aimed at improving the performance criteria of medical masks, ensuring healthcare professionals are able to select the right level of protection.

To learn more about these changes, please visit www.astm.org.

ASTM Mask Level

CHARACTERISTIC	LEVEL 1	LEVEL 2	LEVEL 3
Bacterial Filtration Efficiency (BFE) (%)	≥ 95	≥ 98	≥ 98
Sub-Micron Particulate Filtration Efficiency (PFE) (%)	≥ 80	≥ 85	≥ 85
(ΔP) (mm H ₂ O/cm ²)	< 4.0	< 5.0	< 6.0
Resistance to Penetration by Synthetic Blood, (mmHg)	80	120	160
Flame Spread	Class 1	Class 1	Class 1

Pediatric ASTM Mask Level

CHARACTERISTIC	RESULTS
Bacterial Filtration Efficiency (BFE) (%)	≥ 95
Sub-Micron Particulate Filtration Efficiency (PFE) (%)	≥ 80
(ΔP) (mm H ₂ O/cm ²)	< 5.0
Flame Spread	Class 1

Note: The Pediatric/Child's face mask was not tested to the fluid barrier parameters of ASTM F1862 (Test Method for Resistance of Medical Face Masks to Penetration by Synthetic Blood).



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS

CRITERION™
COMFORT IN PROTECTION™

CRITERION™ Earloop Masks Level 1



Light, comfortable masks for extended wear, these 3-ply procedural face masks offer ASTM Level 1 protection, safeguarding against low levels of moisture, spray and splatter. Each mask features a soft, nonirritating inner and outer layer.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks.....50/box

Specify:
Black (570-3272) White (104-6611)
Blue (104-3809) Teal (570-3280)
Lavender (104-8600) Yellow (104-2849)
Pink (104-3730)

Blue—Made in the USA (570-2443).....50/box

CRITERION™ Pink Ribbon Earloop Masks Level 1



Same great features as the Criterion™ Earloop Masks Level 1, but with a pink-ribbon design to support breast cancer awareness.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

White with Pink Ribbons (570-2829).....50/box

CRITERION™ Earloop Masks Level 2



Same great features as the Criterion™ Earloop Masks Level 1, but with ASTM Level 2 protection.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for secure fit
- Not made with natural rubber latex
- Made in the USA

Masks.....50/box

Specify:
Blue (104-7799) Pink (104-9200)
Lavender (104-8071)

CRITERION™ Pediatric Masks



Designed to fit smaller faces, these masks offer optimal respiratory protection for children ages 4-12. Each mask features a soft, non-woven material to provide enhanced breathability and comfort that children need.

- Fluid resistant
- Flexible nosepiece
- Pleated design for enhanced breathability
- Not made with natural rubber latex

Blue (570-2268).....50/box

CRITERION™ Earloop Masks Level 3



Same great features as the Criterion (TM) Earloop Masks Level 1, but with the highest barrier of protection. These non-sterile, disposable face masks are suitable for longer-duration procedures and surgeries where fluid, splash and splatter protection are essential.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit.
- Not made with natural rubber latex

Masks.....50/box

Specify:
Black (570-3273) Teal (570-2042)
Blue (570-2044) White (570-2043)
Lavender (570-2047) Yellow (570-2046)
Pink (570-2045)

Blue—Made in the USA (570-2448).....50/box

CRITERION™ Fitted Earloop Masks



Comfortable, nonirritating, 3-ply procedural face masks that offer ASTM Level 1 and Level 3 protection. Each mask features an aluminum nosepiece and chinpiece that provides a secure seal around any face shape.

- Fluid resistant
- Earloops and a flexible nosepiece and chinpiece for secure fit
- Not made with natural rubber latex
- Made in the USA

Level 1 Masks.....50/box

Specify:
Pink (900-7452) Blue (900-7451)

Level 3 Masks.....50/box

Specify:
Blue (570-3278) Pink (570-3279)

CRITERION™ Extra Breathing Room Earloop Masks



With a contoured design that enhances airflow and breathability, these 3-ply procedural face masks offer ASTM Level 1 and Level 3 protection. Each mask features a soft, nonirritating inner and outer layer.

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Level 1 Masks.....50/box

Specify:
Blue (570-0722) Lavender (570-0724)
White (570-0723)

Level 3 Masks.....50/box

Specify:
Blue (570-3289) Lavender (570-3291)
White (570-3290)

CRITERION™ Anti-Fog Earloop Masks Level 3



Lightweight and comfortable, these 3-ply procedural face masks offer ASTM Level 3 protection. Each mask features an integrated antifog strip to prevent eyewear from fogging and a soft, nonirritating inner and outer layer.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks.....50/box

Specify:
Blue (570-3281) White (570-3283)
Lavender (570-3282) Pink (570-3284)

CRITERION™ Surgical Tie-On Masks



3-ply surgical masks with cloth-like, ample-length ties for a secure and comfortable fit. Each mask features a soft, nonirritating inner and outer layer.

- Fluid resistant
- Flexible nosepiece
- Pleated design for enhanced breathability
- Not made with natural rubber latex

Blue Level 1 (570-3115).....50/box

Level 3 (570-3294).....50/box

CRITERION™ Sensitive Skin Earloop Masks



Lightweight, comfortable, 3-ply procedural face masks that offer ASTM Level 1 and Level 3 protection. Each mask features a hypoallergenic inner layer that will not lint, tear, or shred for users with sensitive skin.

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Hypoallergenic and free of chemicals, inks, and dyes
- Not made with natural rubber latex
- Made in the USA

Level 1 Masks.....50/box

Specify:
Blue (104-8072) Pink (104-8073)

Level 3 Masks.....50/box

Specify:
Blue (570-3276) Pink (570-3277)

Essentials™ Earloop Procedure Masks
Essentials



Light and comfortable, 3-ply procedural face mask. Meets ASTM 2100 filtration standard for primary, low-barrier performance characteristics. ASTM Level 1.

Masks.....50/box
 Specify:
 Blue (570-3113) Pink (570-3114)

Disposable Face Shields



CRITERION
 Non-Latex
 Optically clear, see-through, full-length visor provides increased eye protection against splash, splatter and aerosols. Features a 1 1/2" thick, 9" long foam headband that keeps the face shield off the face to provide maximum comfort and protection. Ultrasonically welded, 13" elastic headband for added strength and increased wearer comfort. Special fog-free coating on inside and outside of visor. Not made with natural rubber latex. Made in the USA.
 Standard, 13"W x 7 1/2"L
 (101-2254)..... 24/box

CRITERION™ Earloop Masks with Eye Shields Level 3



CRITERION
 4-ply procedural face masks that offer ASTM Level 3 protection. Each mask features an eye shield treated on both sides with antifog material, offering maximum barrier protection.
 • Fluid resistant
 • Pleated design for enhanced breathability
 • Earloops and a flexible nosepiece for a secure fit
 • Not made with natural rubber latex
 Teal
 (570-3116)..... 25/box

Protective Eyewear

Lightweight, wraparound, unisex protective eyewear. 10-base curve for complete coverage. Excellent optical quality. Ideal for patient and staff use. Gray lens and temples.
 (570-0435).....ea



Maxi-Gard® Protective Eyewear

Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 9-base curve allows for great wraparound protection and coverage. Stop hinge prevents temple tip from scratching lens. Protective case included. Metal.
 811 Series, Metal
 (900-4437).....ea



Maxi-Gard® Protective Eyewear

Stylish, wraparound protective eyewear for greater comfort and protection. Clear, shatter-resistant polycarbonate lens. 10-base curve for complete coverage. Includes a scratch-resistant, antifog coating and foam inserts for added comfort. Meets ANSI Z87+ standards.
 Eyewear..... ea
 Specify:
 Black Frame (570-0434) Pink Frame (570-0431)
 Blue Frame (570-0433) Purple Frame (570-0432)



Maxi-Gard® Protective Eyewear

Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 10-base curve allows for complete coverage. Sure-grip rubber temple insert provides extra comfort.
 • For eye protection in operator or lab
 • Lightweight and comfortable
 • Unisex
 806 Series, Plastic..... ea
 Specify:
 Black (900-4436) Pink (900-4779)
 Blue (900-4488) Red (900-4781)



Safety Glasses

Clear plastic lens with side shields for eye protection in the operator or lab. Unisex.
 • Clear plastic lenses
 • Side shields for eye protection in the operator or lab
 • Unisex
 #5148002
 (100-7382)..... ea



Essentials™ Protective Eyewear
Essentials

Wraparound, unisex protective eyewear. Clear, lightweight lenses and temples with a 10-base curve for complete coverage. Excellent optical quality. Ideal for patient and staff use.
 (570-0436)..... 6/box



CRITERION™ Bouffant Caps
CRITERION

Full-sized cut provides total hair coverage and allows the elastic to be comfortable. Prevents hair loss from contaminating sensitive area or site. Elastic that is not made with natural rubber latex minimizes allergic reaction. Available in White and Blue.
 White
 21" (102-4151)..... 100/box
 24" (102-6614).....100/box
 Blue
 21" (104-6032).....100/box
 24" (104-6058).....100/box



CRITERION™ Shoe Covers
CRITERION

Good tensile strength, made from non-woven, spunbonded polypropylene. Stretchy elastic openings for easy on and off. Nonconductive with a nonskid design. Single-use and non-sterile. Not made from natural rubber latex. Blue.
 One Size (Fits up to Size 11)
 (104-9452)..... 50 pair/box
 (104-6350).....150 pair/case
 X-Large (Fits Sizes 12-14)
 (570-3834)..... 50 pair/box





CRITERION

CRITERIONSM SMS Jackets
Formerly Maxi-Gard

Compare to ValuMax
Extra-Safe Jackets and Lab Coats

Premium protective jacket available in 9 fashionable, attractive colors. Stylish hip-length jacket is made from high-quality SMS fabric that meets OSHA fluid-resistance requirements for protection. Lightweight, breathable material for maximum wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure. Not made with natural rubber latex.

Sizes S–X-Large..... 10/pkg.
5 Pkg.@

Specify:

Small		Medium		Large		X-Large	
Black	(570-2288)	Black	(570-2298)	Black	(570-2327)	Black	(570-2318)
Blueberry	(570-2289)	Blueberry	(570-2299)	Blueberry	(570-2328)	Blueberry	(570-2319)
Ceil Blue	(570-2290)	Ceil Blue	(570-2300)	Ceil Blue	(570-2329)	Ceil Blue	(570-2320)
Light Pink	(570-2291)	Light Pink	(570-2301)	Light Pink	(570-2297)	Light Pink	(570-2321)
Purple	(570-2292)	Purple	(570-2302)	Purple	(570-2278)	Purple	(570-2322)
Raspberry	(570-2293)	Raspberry	(570-2303)	Raspberry	(570-2330)	Raspberry	(570-2323)
Sky Blue	(570-2294)	Sky Blue	(570-2304)	Sky Blue	(570-2331)	Sky Blue	(570-2324)
Teal	(570-2295)	Teal	(570-2305)	Teal	(570-2332)	Teal	(570-2325)
White	(570-2296)	White	(570-2306)	White	(570-2333)	White	(570-2326)

Size 2X-Large–Ceil Blue
(570-2307)..... 10/pkg.
5 Pkg.@

Size 3X-Large–Ceil Blue
(570-2308)..... 10/pkg.
5 Pkg.@

CRITERION

CRITERIONSM SMS Lab Coats
Formerly Maxi-Gard

Stylish knee-length coat is manufactured from lightweight, breathable SMS fabric with optimal tensile strength, making them suitable for everyday use. Detailed finishing includes a 3-pocket design and snap-front closure for easy removal. Available in 9 attractive colors. Not made from natural rubber latex.

Sizes S–X-Large..... 10/pkg.
5 Pkg.@

Specify:

Small		Medium		Large		X-Large	
Black	(570-2334)	Black	(570-2346)	Black	(570-2279)	Black	(570-2309)
Blueberry	(570-2335)	Blueberry	(570-2347)	Blueberry	(570-2280)	Blueberry	(570-2310)
Ceil Blue	(570-2336)	Ceil Blue	(570-2348)	Ceil Blue	(570-2281)	Ceil Blue	(570-2311)
Light Pink	(570-2337)	Light Pink	(570-2349)	Light Pink	(570-2282)	Light Pink	(570-2312)
Purple	(570-2338)	Purple	(570-2350)	Purple	(570-2283)	Purple	(570-2313)
Raspberry	(570-2339)	Raspberry	(570-2351)	Raspberry	(570-2284)	Sky Blue	(570-2315)
Sky Blue	(570-2340)	Sky Blue	(570-2352)	Sky Blue	(570-2285)	Raspberry	(570-2314)
Teal	(570-2341)	Teal	(570-2353)	Teal	(570-2286)	Teal	(570-2316)
White	(570-2342)	White	(570-2354)	White	(570-2287)	White	(570-2317)

Size 2X-Large–Ceil Blue
(570-2343)..... 10/pkg.
5 Pkg.@

Size 3X-Large–Ceil Blue
(570-2344)..... 10/pkg.
5 Pkg.@





Compare to ValuMax Dual-Fabric Surgical Cover Gowns

CRITERION™ Cover Gowns



Knee-length cover gowns with complete front coverage. Dual-fabric design offers protection as it maximizes comfort. Sleeves and front of gown are made from quality SMS fabric that meets OSHA fluid-resistance requirements for protection; back of gown is made of a lightweight, breathable material for wearer comfort. Open-back design includes neck and waist ties for individual adjustments. Knitted cuffs. Not made with natural rubber latex. Blue.

- Regular (112-6245).....10/pkg
- X-Large (570-3461).....10/pkg



CRITERION™ Polycoated Polypropylene Aprons



Disposable, full-length coverage with a fabric-like feel for medical and dental settings. Lightweight, breathable, and fluid-resistant polycoated polypropylene material. Adjustable fit with neck strings and waist ties.

- White (570-3462).....25/pkg



Isolation Gowns

Single-use, trilayered isolation gowns made with SMS material. Each features antistatic properties, elastic cuffs, and waist ties. AAMI levels 1–3. Blue and Yellow.

- AAMI–L1.....10/pkg
 - Specify:
 - Medium–Blue (570-2357)
 - Medium–Yellow (570-2358)
 - Large–Blue (570-2359)
 - Large–Yellow (570-2360)
 - X-Large–Blue (570-2361)
 - X-Large–Yellow (570-2362)
- AAMI–L2.....10/pkg
 - Specify:
 - Medium–Blue (570-2375)
 - Medium–Yellow (570-2376)
 - Large–Blue (570-2377)
 - Large–Yellow (570-2378)
 - X-Large–Blue (570-2363)
 - X-Large–Yellow (570-2364)
- AAMI–L3.....10/pkg
 - Specify:
 - Medium–Blue (570-2371)
 - Medium–Yellow (570-2372)
 - Large–Blue (570-2373)
 - Large–Yellow (570-2374)
 - X-Large–Blue (570-2355)
 - X-Large–Yellow (570-2356)



CRITERION™ Polypropylene Cover Gowns



Full-length, lightweight, and protective polypropylene-coated fabric that provides coverage during medical and dental procedures to minimize fluid exposure. Adjustable fit with a full-cut design, open neck, open back, and waist ties. Sewn neck binders and knit cuffs for enhanced protection.

- Gowns.....10/pkg
- Specify:
- Small/Medium, Pink (570-3466)
- Medium/Large, Blue (570-3465)
- Large, Pink (570-3464)



Maxi-Gard™ CPE Isolation Gowns, AAMI Level 2

Made of single-layer CPE (chlorinated polyethylene) film, not affected by liquid. Provides economical, comfortable, and reliable protection with full-length maximum coverage and fixed-ring cuff band.

- Flame and chemical resistant
 - Knee-length with long sleeves
 - Cuff bands with thumb hooks
 - Pullover design with waist ties
 - One size fits most
 - Not made with natural rubber latex
 - Blue
- (570-2122).....10/bag

BRAND PRODUCTS

FOR COMPLETE OFFERING OF RAPID TESTING ITEMS

— SEE PAGES 96-97 —



Hand Sanitizer

Waterless, alcohol-based gel hand sanitizer. Ideal for use where soap and water are not readily available. Rinse-free: no water required. Kills 99.9% of most common disease-causing germs spread by hands. Soothing emollients and a pleasant, fresh, clean scent. Leaves hands soft and clean. 66.5% ethyl alcohol. Made in the USA.

8-oz Pump Bottle
(900-4443)..... ea



16-oz Pump Bottle
(900-4442)..... ea



Maxiclens™ Instant Hand Sanitizer

- Kills 99.99% of germs in 15 seconds without the use of soap and water
- Contains aloe vera to add moisture with every use
- Henry Schein's hand sanitizer helps meet OSHA, APIC, and CDC hand-washing recommendations
- Citrus scent

800-mL Bag for Manual Dispenser
(112-5988)..... ea

800-mL Bag for Automatic Dispenser
(112-5991)..... ea

Automatic Dispenser, Soap
(112-5985)..... ea

Manual Dispenser
(112-5984)..... ea

Automatic Dispenser, Hand Sanitizer
(112-5986)..... ea



Hand Lotion

Eliminates rough skin and dryness that comes from overexposure to soaps and sanitizers. Contains cucumber and aloe emollients that help to leave hands feeling luxuriously soft.

- Luxurious hand cream
- Soothing to skin
- Contains emollients to leave hands feeling soft
- Pleasant Cucumber and Aloe fragrance

8-oz Pump Bottle
(900-4446)..... ea



Foaming Antibacterial Hand Soap

- Provides a rich, dense lather
- Rinses easily using minimal water
- Leaves hands feeling soft and smooth
- Pleasant Spring Water Clean fragrance
- Contains 0.6% chloroxylenol (PCMX) for antimicrobial action, and aloe vera to add moisture with each use

18.6-oz Pump Bottle
(900-4441)..... ea



Quick Wipes

- Soft, strong towelettes: allow easy and gentle cleansing of hands and face
- Effective and gentle formulation that contains aloe: formulated for fast drying and gentleness
- Convenient packaging: towelettes easily dispensed through orifice; each towelette is saturated with cleaning liquid

Pop-Up Container, 9" x 5½"
(104-4176)..... 150/cont



Antibacterial Liquid Hand Soap

- Provides a luxurious feel, leaving hands smooth and soft
- Pleasant Cucumber and Aloe fragrance
- Contains 0.6% chloroxylenol (PCMX) for antimicrobial action, and aloe vera to add moisture with each use

8-oz Pump Bottle
(900-4439)..... ea

16-oz Pump Bottle
(900-4438)..... ea

1-gal Refill Bottle
(900-4440)..... ea

Maxiclens™ Antiseptic Skin Cleanser

4.0% chlorhexidine gluconate formulation. Fast acting in just 15 seconds, while inhibiting bacterial regrowth for up to 6 hours. Skin emollients protect against irritation, keeping hands soft and soothed. Residue-free, nonsticky formulation with no added dyes. All-in-one solution is tough on bacteria, yet gentle on the skin. Pleasant floral scent.



1-qt Bottle
(570-3268)..... ea



1-gal Refill
(570-3269)..... ea

MAXIWIPE™ GERMICIDAL DISPOSABLE WIPES

Choosing the Right Wipe for Your Practice

MaxiWipe by Henry Schein® offers broad spectrum, one-step intermediate-level disinfecting wipes in alcohol and alcohol-free formulations. Select the best fit for your healthcare setting, staff, and patient needs.



**MaxiWipe™
Germicidal Wipes**

Ideal for daily use in fast-paced environments that require short contact times and broad coverage of microorganisms.

(106-6794) Large – 6" x 6¾" (15.2 x 17.1 cm)
(570-3361) XL – 7½" x 15" (19.1 x 38.1 cm)



**MaxiWipe™ Alcohol Free
Germicidal Wipes**

Alcohol- and fragrance-free formula ideal for use around patients and staff with respiratory sensitivities.

(570-3056) Large – 6" x 6¾" (15.2 x 17.1 cm)
(570-3057) XL – 7½" x 15" (19.0 x 38.1 cm)

 <p>KEY ACTIVE INGREDIENT</p>	<p>Quaternary ammonium and isopropyl alcohol (IPA)</p>		<p>Quaternary ammonium</p>
 <p>CONTACT TIME</p>	<p>1-MINUTE CONTACT TIME Effective against SARS-CoV-2 and Tuberculosis</p>	<p>2-MINUTE OVERALL CONTACT TIME Effective Bactericidal and Virucidal against Hepatitis, HIV-1, Influenza A (H1N1) and more</p>	<p>3-MINUTE OVERALL CONTACT TIME Effective Bactericidal, Tuberculocidal, and Virucidal, including SARS-CoV-2 (virus that causes COVID-19), and Influenza B</p>
 <p>SURFACE & EQUIPMENT COMPATIBILITY</p>	<p>Compatible with a broad range of surfaces including stainless steel, plastic, glass, and Formica®</p>		<p>Compatible with equipment sensitive to alcohol and a variety of hard, nonporous surfaces. Low-streak formula is great on glass, mirrors, dark countertops, and stainless steel</p>
 <p>FORMULATION</p>	<p>Free of fragrance, acetone, bleach, phenol, ammonia, and acid</p>		<p>Free of alcohol, fragrance, acetone, bleach, phenol, ammonia, and acid</p>
 <p>EPA LIST COMPLIANCE</p>	<p>EPA List N, Registration No. 9480-4</p>		<p>EPA List N and now also qualified for List Q*, Registration No. 9480-9 <small>*MaxiWipe™ Alcohol Free has demonstrated effectiveness against emerging viral pathogens such as the Monkeypox virus on hard, nonporous surfaces.</small></p>

Formica® is registered trademarks of The Diller Corporation.



EXCLUSIVE BRANDS & PRODUCTS



QUALITY YOU CAN RELY ON



SMART SAVINGS



MaxiSpray® Plus

Surface Disinfectant Cleaner

- Hospital-level, general-purpose disinfectant/cleaner proven to kill TB in 5 minutes, and HBV, HIV, and HCV in 3 minutes
- EPA-registered as a broad-spectrum disinfectant for both surface and immersion use
- Bactericidal, virucidal, tuberculocidal, and fungicidal

24-oz Trigger-Spray Bottle
(101-4114).....ea

1-gal Refill
(101-0285).....ea



MAXIWIPE™ Germicidal Wipes

Disinfectant wipe that kills the COVID-19 virus in 1 minute on hard, nonporous surfaces. Effective against 16 pathogens, including COVID-1, MRSA, and H1N1. Provides fast disinfection. Compatible with a broad range of surfaces, including stainless steel, plastic, Formica, and glass.

- Cleans and disinfects common hard, nonporous surfaces
- Bactericidal, tuberculocidal, and virucidal
- Effective against 13 microorganisms including TB, HBV, MRSA, VRE, and RSV
- Contains 55% alcohol
- Meets CDC and OSHA regulations

Large—6" x 6¾"
(106-6794).....160/canister

X-Large—7½" x 15"
(570-3361).....75/canister

Packaging label may vary.



MAXIWIPE™ Alcohol-Free Germicidal Disposable Wipes

Formulated with care: Alcohol- and fragrance-free wipes for staff and patients with respiratory sensitivities. Efficient, one-step cleaner and disinfectant. Effective bactericidal, tuberculocidal, and virucidal, including SARS-CoV-2 (the virus that causes COVID-19), in 3 minutes. EPA Registration No. 9480-9, now qualified for List Q*. Intermediate level.

Large, 6" x 6¾"
(570-3056).....160/cont

X-Large, 7½" x 15"
(570-3057).....75/cont

*MAXIWIPE™ Alcohol-Free has demonstrated effectiveness against emerging viral pathogens, such as the monkeypox virus, on hard, nonporous surfaces.



MAXICIDE® 3% Glutaraldehyde

This advanced 3% glutaraldehyde high-level disinfection solution is twice as fast as leading competitors, taking only 25 minutes to process. It provides effective sterilization in 10 hours and can be reused for up to 28 days. The moderate pH is compatible with heat-sensitive medical devices and contains a rust inhibitor to prevent corrosion. Once the powdered activator is added, the solution turns blue and has a mild lemon scent. Makes 1 gallon with powdered activator bottle.

1-gal Bottle
(570-3274).....ea



Surgical Milk

Super-concentrated formula that stops corrosion and maintains free-working instruments. Nontoxic and safe to use. Permits autoclaving of sharps and carbon steel. One gallon of Surgical Milk diluted with 5 gallons of tap water makes 6 gallons of immersion bath milk.

1-qt Bottle
(115-9589).....ea

64-fl-oz Bottle
(102-5539).....ea



MAXICIDE® Activated

Preactivated 2.65% Glutaraldehyde Solution
Can be used as a sterilant and high-level disinfectant for up to 30 days. Tuberculocidal in 45 minutes. Ready to use; has a fresh, pleasant scent.

1-gal Bottle
(104-4681).....ea



MAXICIDE® OPA 28 High-Level Disinfectant

Nongluteraldehyde 0.575% orthophthaldehyde (OPA) solution. Ready-to-use high-level disinfectant that achieves disinfection in just 10 minutes at room temperature (20°C) and offers a reuse period of up to 28 days. Significantly reduces the risk of occupational hazards including eye, skin, and respiratory tract irritations and sensitization. Test strips should be used to verify the minimum recommended concentration (MRC) prior to each use.

Solution
(112-6976).....ea

Test Strips
(570-3271).....50/bottle



Self-Seal Sterilization Pouches with Dual Indicators

Superior quality pouch with dual indicators that comply with the recommendations of the CDC for effective sterilization. Preprinted internal and external indicators change color when processed with steam and EtO. Made with medical-grade paper and an easy-to-open thumb notch. Lead-free and not made with natural rubber latex.

2¼" x 4" (112-4859)	200/box
2¾" x 9" (112-4856)	200/box
3½" x 5¼" (112-4857)	200/box
3½" x 9" (112-4853)	200/box
3½" x 9" (112-4855)	500/box
3¾" x 12" (112-4866)	200/box
4¼" x 11" (112-4863)	200/box
5¼" x 10" (112-4854)	200/box
5¼" x 12" (112-4861)	200/box
7½" x 13" (112-4858)	200/box
8" x 13" (112-4864)	100/box
9" x 15" (112-4860)	100/box
10½" x 17" (112-4862)	100/box
13" x 20" (112-4865)	100/box

*The use of weekly biological indicators is recommended.



Dual-Purpose Nylon Film Pouches

All-nylon, self-sealing film pouches with indicators. Can be used in steam or dry heat.

- All-nylon
- Self-sealing film pouches with indicators
- Can be used in steam or dry heat

2" x 9½" (101-3363)	100/box
3" x 10" (101-3373)	100/box
4" x 10" (101-3383)	100/box
6" x 11" (101-3393)	100/box



MAXIMA® Multiparameter Self-Seal Sterilization Pouches

- Confirms that all 3 parameters for sterilization have been met: time, temperature, and presence of steam
- Internal and external multiparameter indicators (Type 4, Class 1) are printed directly on the pouch
- Medical-grade paper with triple seals; chevron safety seal at base for added protection; thumb notch for easy opening
- Not made with natural rubber latex
- Made in the USA

2¼" x 4" (900-7838)	200/box
2¾" x 9" (900-4602)	200/box
3½" x 5¼" (900-4604)	200/box
3½" x 9" (900-4601)	200/box
3½" x 9" (570-3718)	500/box
5¼" x 6½" (570-3719)	200/box
5¼" x 10" (900-4603)	200/box
7½" x 13" (900-4782)	200/box
12" x 18" (900-4783)	100/box



Autoclave Bags

Built-in Indicator/Paper, 2½" Wide, 1½" Opened Sides

- Wet-strength, sterilizable paper permits sterilant penetration and acts as a barrier after sterilization
- Chemical indicator on each bag verifies exposure to sterilization process
- Indicator ink will have noticeable color change when exposed to sterilization conditions
- Safety fold design to enhance effectiveness of sterile field
- Cost-effective alternative to paper/poly pouch; biodegradable

6" (100-0166)	1000/box
8" (100-0529)	1000/box
10" (100-1686)	1000/box



Essentials™ Self-Seal Sterilization Pouches with Dual Indicators
Essentials

Self-sealing sterilization pouches for sterilization with autoclaves and ethylene oxide (EtO) chambers. Dual (internal and external) process steam indicators will change color when a sterilization process has begun (steam/EtO). To confirm sterilization, use a biological monitor. Leak-resistant. Maintains sterilization for up to 12 months (once sterile and items remain sealed in the pouch). Made with medical-grade paper and an easy-to-open thumb notch. When used properly, it assists with CDC infection control guidelines compliance.

- Lead-free and made without natural rubber latex
- Suitable for sterilization with dental autoclaves and ethylene-oxide (EtO) chambers
- Perforated fold for easy and accurate folding of adhesive strip to form an airtight seal
- Meets ANSI/AAMI ST79, ISO 11607, ISO 13485, and ISO 11140-1 standards

2.25" x 4" (570-4241)	200 per Box
2.75" x 9" (570-4238)	200 per Box
3.5" x 5.25" (570-4240)	200 per Box
3.5" x 9" (570-4236)	200 per Box
(570-4244)	500 per Box
4.25" x 11" (570-4239)	200 per Box
5.25" x 10" (570-4237)	200 per Box
7.5" x 13" (570-4243)	200 per Box
9" x 15" (570-4242)	100 per Box



Self-Seal Autoclave Bag

Self-seal sterilizer bag with autoclave and chemical vapor sterilization indicator ink. Heavy-duty surgical-grade kraft paper. Dual indicators for steam or chemical vapor. Self-sealing feature ensures high-integrity seal.

- Indicator ink shows noticeable color change when exposed to sterilization conditions
- Wet-strength, sterilizable-grade, white paper permits sterilant penetration and acts as a poststerilization barrier
- Chemical indicator on each bag verifies exposure to sterilization process
- Convenient, user-friendly, self-seal closure: no heat sealer necessary, no release liner to peel off to get to the adhesive
- Safety-fold design to enhance effectiveness of sterile field
- Cost-effective alternative to paper/poly pouch; biodegradable

2½" x 1½" x 10½" (101-5915)	1000/box
---	----------



Dual-Purpose Nylon Sterilization Tubing

- 100-ft Rolls
- Transparent: contents are easy to identify.
- Sturdy: resists tearing and puncturing. Dual purpose: can be used in steam or dry heat.
- 100 feet per roll.
- Continuous Clear, 2" (101-2224)..... ea
- Continuous Clear, 3" (101-2225)..... ea
- Continuous Clear, 4" (101-2226)..... ea
- Continuous Clear, 6" (101-2227)..... ea
- Continuous with Indicator, 2" (101-1111)..... ea
- Continuous with Indicator, 3" (101-1110)..... ea
- Continuous with Indicator, 4" (101-2222)..... ea
- Continuous with Indicator, 6" (101-2223)..... ea
- Perforated with Indicator, 2" x 9½" (101-3323)..... ea
- Perforated with Indicator, 3" x 10" (101-3333)..... ea
- Perforated with Indicator, 4" x 10" (101-3343)..... ea
- Perforated with Indicator, 6" x 11" (101-3353)..... ea



MaxiTest® Steam Sterilization Integrator

- Class 5 integrating indicator provides users with immediate readout results with no incubation
- Monitors both gravity displacement and prevacuum steam sterilizers and is equivalent in performance to a biological integrator
- Precise, easy-to-read format provides permanent, physical proof of sterilization conditions
- (570-0310).....100/bag
- (570-0311).....1000/bag



Dry Heat Sterilizer Indicator Labels

- For dry heat sterilizers. Indicator ink changes from green to black, indicating sterilization process has been completed.
- (102-7982)..... 1000/box




Lead-Free Steam Indicator Tape

- Resilient and pressure-sensitive
- Lead-free
- Adhesive tape designed for use on disposable wraps
- The backing is a natural-colored saturated crepe paper that provides secure wraps throughout sterilization expansion
- Indicator lines will change from beige color to dark brown/black when proper levels of moisture and temperature have been achieved
- ISO 11140-1:2014 Type 1
- ¾" Roll (570-1510)..... ea
- 1" Roll (570-1509)..... ea



Indicator Strips

- For steam and chemical sterilizers
- Indicator ink changes from white to black, indicating sterilization process has been completed
- (102-6332).....250/box



MaxiTest® Chemical Vapor Indicator Strips

- Makes a distinct color change when exposed to the chemical vapor sterilization process
- An internal chemical indicator strip should be used within each package or tray to be sterilized
- 4" strip can be written on and kept as permanent record
- (570-0382)..... 250/box



Instrument Cassettes

- 5-Instrument Cassette
- 7.1" x 2.8" x 1.3"..... ea
- Specify:
- Blue (112-4963) Red (112-5297)
- Green (112-5301) Yellow (112-5442)
- Purple (112-5336)
- Replacement Rail..... ea
- Specify:
- Blue (112-5455) Red (112-5456)
- Green (112-5457) Yellow (112-5459)
- Purple (112-5458)



- 10-Instrument Cassette
- 7.1" x 5.6" x 1.3"..... ea
- Specify:
- Blue (112-4964) Red (112-5299)
- Green (112-5302) Yellow (112-5443)
- Purple (112-5439)
- Replacement Rail..... ea
- Specify:
- Blue (112-5460) Red (112-5461)
- Green (112-5462) Yellow (112-5463)
- Purple (112-5464)



- 16-Instrument Cassette
- 11.3" x 7.1" x 1.3"..... ea
- Specify:
- Blue (112-4965) Red (112-5300)
- Green (112-5303) Yellow (112-5444)
- Purple (112-5441)
- Replacement Rail..... ea
- Specify:
- Blue (112-5465) Red (112-5466)
- Green (112-5468) Yellow (112-5470)
- Purple (112-5469)



MaxiTest® Value Test Pack

- Used to test steam, dry heat, chemical vapor, and EtO processes
- Mail-in service with results reported quarterly
- Immediate and confidential notification of test failures
- Test results available electronically
- Postage not included

(101-4715)..... ea
 Contains: 52 test packets, each containing 1 test strip & 1 control strip.

Not for export.



MaxiTest®

Mail-in Service

- Easy to use: place MaxiTest strips within an average load; process according to normal procedure; when cycle is complete, remove load from sterilizer and remove test strips; forward to bacteriological laboratory for culturing
- Suitable for use with all sterilizers: dual-organism spore strips contain both *Bacillus stearothermophilus* and *Bacillus subtilis* so they can be used to monitor all sterilizers, including unsaturated chemical vapor
- 2 test strips and 1 control strip: provide increased safety to patients and staff; 2 test strips placed in separate areas of the sterilizer provide greater accuracy of test results
- Biological Sterilization Monitoring Service: MaxiTest includes test sets, culturing service, laboratory reports, all postage and handling, and any device concerning sterilization procedures; the lab cultures are completed in 7 days; however, if any test strip cultures positive (spore growth indicating sterilization failure), the office will be telephoned immediately; if the test strips do not culture positive, the office will receive a laboratory report in approximately 2 weeks, indicating sterilization was achieved

12-Test Package
(101-5830)..... ea
 Contains: 12 packets, each containing 2 test strips & 1 control strip.

48-Test Package
(101-1253)..... ea
 Contains: 48 packets, each containing 2 test strips & 1 control strip.

Not for export.



MaxiTest® In-Office Biological Monitoring System

- In-office system for monitoring steam sterilizers
- Easy to use and interpret
- Results in 24–48 hours
- Simply process vial in a normal sterilizer cycle, then activate vial, incubate, and record results
- If any spores survive the process, the vial medium will turn yellow
- MaxiTest Incubator is 24 volts and holds up to 14 indicators

Starter Kit
(102-6076)..... ea
 Contains: 1 of each: box of 25 indicators, incubator, crusher & record keeper.

Indicator Refill
(102-5804)..... 25/box
(112-6848)..... 100/box



MaxiTest® Steam BI Test Pack

- This test pack, with the MaxiTest® Class 5 Integrator, is indicated for use in routine monitoring and sterilizer qualification testing, and allows users to immediately release processed loads
- Per ANSI/AAMI standards, steam sterilizers should be monitored at least weekly, preferably every day the sterilizer is run

(570-0312)..... 25/case



MaxiTest® Bowie-Dick Test Pack

- Disposable pack used to evaluate the efficacy of air removal and steam penetration in prevacuum steam sterilizers
- Includes frequency of use, instructions for use, and storage/shelf-life information preprinted on the box for ease of use and staff training
- Lead-free and easy to interpret

(570-0309)..... 30/case



MaxiTest® Sterilizer Test Pack

- This test pack, with the MaxiTest® Class 5 Integrator, provides immediate readout results
- Equivalent in performance to an AAMI BI test pack
- No incubation required
- For use in all steam sterilizer cycles
- For the highest level of sterility assurance, use in every wrapped, nonimplant sterilizer load

(570-0313)..... 25/case



Compare to Omni® Cleaner XL and Midmark SpeedClean® Autoclave Cleaner

- Mild, slightly alkaline liquid that acts as a powerful, concentrated cleaner to remove grease, grime, and accumulation of scale deposits
- Prolongs the life of the autoclave with routine maintenance
- Mix 4 oz. of cleaner into 64 oz. of distilled water
- Excellent for use with steam

32 oz
(112-5839)..... ea



Super Concentrated Ultrasonic Cleaning Solutions

Fast-acting solutions that clean quickly and effectively. Nontoxic and noncorrosive. In small, easy-to-store, 16-oz bottles, which are less expensive to ship, handle, and stock. Metered-dose dispenser top ensures exact dose. To use: loosen cap 1/2 turn, squeeze to 1-oz mark, and remove cap and pour into glass beaker. Available in 3 formulas.

16-oz Metered Dose Bottles..... ea
Specify:

General Purpose Cleaner
(Makes 16 gal) Used for very fast cleaning of blood, mucus, and tissue. (101-4412)

Tartar and Stain Remover (Makes 96 oz)
Its citric acid formula is noncorrosive and cleans quickly. (101-0827)

Enzymatic Cleaning Solution (Makes 16 gal)
Effective against blood, tissue, and mucus found on instruments. (102-6472)



Ultrasonic Cleaning Solutions

General Purpose Cleaner

- Nonammoniated concentrate
- Dilute 1:10 for use
- Effectively cleans items such as dried blood, bone, tissue, oxide compounds, and more
- Removes ions from instrument's surface to enhance effectiveness of ultrasonic cleaning unit

1-gal Bottle
(100-3836)..... ea



Unit-Dose Ultrasonic General-Purpose Cleaner

- Single-dose, superconcentrated powder for use in ultrasonic cleaning units, auxiliary pans, and beakers for cleaning instruments directly in the ultrasonic cleaning tank
- Eliminates lifting and storing heavy gallon containers
- 1 packet of powder makes 1 gallon of solution (101-0377).....24/box



MAXIZYME® Multi-Enzymatic Detergent

Improved multienzymatic detergent that cleans fast in manual, ultrasonic, and automatic cleaning systems. Low-foaming formula rapidly dissolves and carries away organic soils, including blood, protein, tissue, and more. Its corrosion-inhibiting, pH-neutral formula is safe for instruments and equipment, extending the useful life of instruments. Concentrated and economical, each 1-gallon bottle yields 128 gallons of solution. Clean, fresh scent.

1-gal Bottle
(570-3275)..... ea



Ultra Clean®

General-purpose, nonammoniated ultrasonic cleaner. Removes dried blood, tissue, soil, plaster, stone, pumice, and other light contaminants. Nonflammable and nonfoaming. Supplied as powder in container. Concentrated solution: dilute 1 part to 9 parts of water.

1-gal Bottle
(102-2255)..... ea



Ultra Clean-L

General-purpose, nonammoniated, and super concentrated. Easier to dilute than powders. Extra cleaning power and pleasant scent. Add water to fill container; dilute 1:10 for working-strength solution.

1-gal Bottle
(100-3185)..... ea



MAXITABS®

General Purpose Cleaner
Provide safe and effective cleaning properties in multiple cleaning cycles throughout the day. Does not leave a residue on instruments or inside tank. One tablet makes 1 gallon of solution; 1 box makes 32 gallons.

(900-4088).....32/box



MAXIZYME® Tabs

Enzyme Cleaning Tablets

- Contain protease enzymes that break down hard-to-clean deposits such as blood, bone, and food particles that chemical detergents leave behind
- A single cleaner with 2 uses: ultrasonic and evacuation systems
- Gentle on instruments: nontoxic, neutral pH, biodegradable, and noncorrosive formula that contains no chlorine or harsh alkalines to damage expensive stainless steel instruments
- Refreshing minty scent
- Easy to store

(900-0475).....64/box



Safety Needle and Syringe Combinations

Features and benefits:

- 3 single-handed safety activation methods: hard-surface, thumb, or finger
- Double-locking mechanism with an audible click
- Needle hub-compatible with all standard Luer-lock and Luer slip-tip syringes
- ISO color coded: easy and quick identification of the needle gauge
- Bevel-up needle: allows for low-angle injections

Normal Wall

- 1 cc, 27 ga x ½", Gray (900-7648).....100/box
- 1 cc, 25 ga x 1", Orange (900-7663).....100/box
- 3 cc, 25 ga x 1", Orange (900-7644).....100/box
- 3 cc, 25 ga x ⅝", Orange (900-7646).....100/box
- 3 cc, 23 ga x 1", Blue (900-7643).....100/box
- 3 cc, 21 ga x 1½", Green (900-7640).....100/box

Normal Wall

- 1 cc, 27 ga x ½", Gray, for Allergy (900-7659).....25/tray
- 1 cc, 26 ga x ⅜" IDB, Tan, for Allergy (900-7660).....25/tray



HYPODERMIC SAFETY NEEDLES
WITH ISO COLORED PLASTIC HUB

- Two single-handed safety activation methods include hard surface and thumb activation, which may help reduce the risk of accidental needle stick injuries
- Easy-to-use and compatible with all standard luer lock and luer slip-tip syringes
- Hubs and protective shield are color-coded according to ISO standards to facilitate easy and quick identification of needle gauges
- Bevel-up needle allows for low-angle injections
- Protective arm allows for clear visibility of the injection site, ensuring an unobstructed pathway
- Single use



HENRY SCHEIN® HYPODERMIC SAFETY NEEDLES
with ISO Colored Plastic Hub (100/box)

Item Number	Gauge	Color	Needle Length	Bevel Angle	Inner Diameter	Outer Diameter (mm)
570-2709	30	Light Yellow	½"	Long	RW	0.31
570-2706	27	Grey	1 ½"	Long	RW	0.41
570-2707	27	Grey	1 ¼"	Long	RW	0.41
570-2708	27	Grey	⅝"	Long	RW	0.41
570-2704	25	Orange	1 ½"	Long	TW	0.51
570-2703	25	Orange	1"	Long	TW	0.51
570-2705	25	Orange	⅝"	Long	TW	0.51
570-0340	23	Blue	1 ½"	Long	TW	0.63
570-2701	23	Blue	1"	Long	TW	0.63
570-2710	23	Blue	⅝"	Long	TW	0.63
570-2700	22	Black	1 ½"	Long	TW	0.71
570-2699	22	Black	1"	Long	TW	0.71
570-2698	21	Green	1 ½"	Long	TW	0.82
570-2696	21	Green	1"	Long	TW	0.82
570-2697	20	Yellow	1 ½"	Long	TW	0.90
570-2695	20	Yellow	1"	Long	TW	0.90
570-2694	18	Pink	1 ½"	Short	TW	1.28
570-2693	18	Pink	1"	Short	TW	1.28

RW—Regular Wall | **TW**—Thin Wall

Order today at HENRYSCHEIN.COM/MEDICAL



CONVENTIONAL HYPODERMIC NEEDLES
WITH ISO COLORED PLASTIC HUB

- High-quality stainless steel
- Large inside diameter (ID) provides an increased flow rate
- Semi-transparent needle hub is ISO color-coded by size for easy identification and safety
- Needle hub-compatible with all standard Luer-lock and Luer slip-tip syringes
- Sharp needle improves patient comfort
- Wide range of needle sizes for each procedure
- Product conforms to ISO 7864
- Sterile
- Non-pyrogenic



HENRY SCHEIN® CONVENTIONAL HYPODERMIC NEEDLES
with ISO Colored Plastic Hub (100/box)

Item Number	Gauge	Color	Needle Length	Bevel Angle	Inner Diameter	Outer Diameter (mm)
112-6148	30	Light Yellow	½"	Long	RW	0.30
112-6149	30	Light Yellow	1"	Long	RW	0.30
112-7108	27	Grey	½"	Long	RW	0.40
112-7103	27	Grey	1¼"	Long	RW	0.40
900-4468	25	Orange	⅝"	Long	TW	0.50
112-7112	25	Orange	⅞"	Long	TW	0.50
900-4470	25	Orange	1"	Long	TW	0.50
900-4472	25	Orange	1½"	Long	TW	0.50
900-4473	23	Blue	1"	Long	TW	0.60
112-7107	22	Black	1"	Long	RW	0.70
112-7101	22	Black	1½"	Long	RW	0.70
112-7100	21	Green	1"	Long	RW	0.80
112-7111	21	Green	1½"	Long	RW	0.80
112-7109	20	Yellow	1"	Long	RW	0.90
112-7110	20	Yellow	1½"	Long	RW	0.90
900-4469	18	Pink	1"	Short	TW	1.20
900-4471	18	Pink	1½"	Short	TW	1.20

RW—Regular Wall | **TW**—Thin Wall

Order today at HENRYSCH EIN.COM/MEDICAL



Safety Insulin Syringes with Needles

- Transparent locking safety shield protects the health-care worker during transportation and after use of the device. Shield can be locked with either a right or left twist. An audible click is heard to confirm the final locked position
- Revolutionary orbit cover eliminates the need to recap the needle after draw-up of medication
- Low dead space reduces medication waste
- Clear and bold gradations are easy to read

½ cc, 29 ga x ½" (900-7661).....100/box
 1 cc, 29 ga x ½" (900-7662).....100/box



Syringes without Needles

Sterile Luer-lock and slip-tip syringes. Low drag coefficient and low penetration force. Clear, visible scale marking for accuracy and dosage control. Compatible with all Luer-style and slip-tip needles and hubs.

Slip-tip Syringe
 1 cc (900-4480).....100/box
Luer-lock Syringe
 3 cc (900-4475).....100/box
 5 cc (900-4477).....100/box
 10 cc (900-4476).....100/box
 20 cc (112-6151).....50/box
 30 cc (112-6153).....50/box
 60 cc (112-6152).....50/box



Blunt Fill Needles

Used for draw-up of medication, reconstitution of medication, and accessing IV bags.

- 45° bevel to help reduce accidental needlesticks
- 18 ga for ease in draw-up of medications
- 5-micron filter provides effective filtration of glass particles
- Hub is color coded for identification

18 ga x 1½" (112-7097).....100/box
 18 ga x 1½" with 5-Micron Filter (112-7098).....100/box



Luer Lock Syringes with Exchangeable Needle Tray

Features and benefits:

- 50 syringes per tray reduces set-up time
- Clear and bold graduations are easy to read
- ISO color-coded for easy and quick identification of the needle gauge
- Ultra-sharp needle makes it easy to access vial

3 cc, 23 ga x 1", Blue (900-7844)	50/tray
3 cc, 22 ga x 1", Black (900-7843)	50/tray
3 cc, 21 ga x 1", Green (900-7842)	50/tray



Tuberculin Syringes with Attached Needles

- 1 cc
- Slip tip
- Clear, bold graduations are easy to read

27 ga x 1/2" (570-2120)	100/box
26 ga x 3/8" (570-2118)	100/box
25 ga x 3/8" (570-2119)	100/box



Conventional Syringes with Needles

- Sterile, disposable Luer-lock syringes with needles
- Low drag coefficient and penetration force
- Sharp lancet point needle
- Excellent for vaccine delivery
- Thin wall

3-cc Syringes.....	100/box
Specify:	
25 ga x 5/8" (900-4458)	21 ga x 1" (900-4462)
25 ga x 1" (900-4461)	20 ga x 1" (112-7104)
23 ga x 1" (900-4460)	20 ga x 1 1/2" (112-7115)
22 ga x 1" (112-7102)	18 ga x 1 1/2" (112-7113)
22 ga x 1 1/2" (112-7105)	
5-cc Syringes	
21 ga x 1" (112-7099)	100/box
10-cc Syringes.....	100/box
Specify:	
21 ga x 1" (900-4466)	20 ga x 1" (900-4465)

MAXIWIPE™ ALCOHOL FREE
GERMICIDAL DISPOSABLE WIPE

FOR OUR COMPLETE OFFERING OF DISINFECTANTS SEE PAGE 29

FORMULATED WITH CARE

Alcohol- and fragrance-free formula for staff and patients with respiratory sensitivities.





Primary IV Administration Set

- Flat-top, needle-free valve with tight septum seal helps to support effective disinfection
- Metal-free for MRI compatibility
- Not made with natural rubber latex
- DEHP-free
- Sterile
- Single use

15 drops/mL, 73", vented, slide clamp, roller clamp, injection Y-site, rotating male Luer lock.

(570-0881).....ea



15 drops/mL, 86", vented, roller clamp, needle-free Y-site, rotating male Luer lock.

(570-0885).....ea



15 drops/mL, 105", vented, back check valve, slide clamp, 2 needle-free Y-sites, roller clamp, rotating male Luer lock.

(570-0887).....ea



Secondary IV Administration Set

- Flat-top, needle-free valve with tight septum seal helps to support effective disinfection
- Metal-free for MRI compatibility
- Not made with natural rubber latex
- DEHP-free
- Sterile
- Single use

15 drops/mL, 40", vented, roller clamp, hanger, rotating male Luer lock

(570-0883).....ea



Safety IV Catheters

- Rapid-Flash notched needles in 20 ga to 24 ga give immediate visualization of successful vein placement
- Tri-beveled back-cut needle features industry-leading sharpness, resulting in less trauma during insertion
- Proprietary PUR catheter material softens upon insertion, providing a better dwelling experience
- Box-shape grip allows better hold of the device, preventing rolling
- Raised activation button helps reduce accidental activation

24 ga x 1/2"	(570-2591)	50/box
24 ga x 3/4"	(570-2586)	50/box
22 ga x 1"	(570-2587)	50/box
20 ga x 1.16"	(570-2588)	50/box
20 ga x 1"	(570-2589)	50/box
18 ga x 1.16"	(570-2590)	50/box

Blood Control Safety IV Catheters

Ergonomically designed with contoured catheter housing and a textured hub that creates a nonslip grip for easy handling. These catheters feature a triple-facet, bevel needle with a fully encapsulating needle shield that encases the needle upon full retraction, providing safety without compromising functionality.

- Triple facet, back cut, beveled needle provides ease of insertion with increased patient comfort at insertion site
- Gently tapered catheter tip provides smooth transition from needle to catheter during insertion
- Ridged catheter hub offers superior grip and stability during connection of IV administration sets
- Blood control with active check valve technology inhibits blood flow when valve is activated
- Push-off tab assists with one-handed catheter advancement
- Push-button style needle retraction safety mechanism retracts needle fully and safely into the catheter body
- Contoured catheter body is ergonomically designed for one-handed technique
- Flashback chamber allows for clear visual confirmation of successful venipuncture
- Audible click lock safety mechanism: audible confirmation ensures the needle is safely encapsulated in the catheter body
- VenaPUR® catheter material is designed to remain kink-free and soften within the body for increased patient comfort

24 ga x 3/4"	(570-3716)	50/box
22 ga x 1"	(570-3715)	50/box
20 ga x 1"	(570-3714)	50/box
20 ga x 1.16"	(570-3713)	50/box
18 ga x 1.16"	(570-3712)	50/box



Safety Winged Blood Collection Needles

- Flexible wings secure device to patient's skin
- Thin-walled needle minimizes discomfort
- 3 ways to activate safety mechanism
- DEHP-free
- Sterile
- Single use

3/4" x 12" Needles.....	50/box
Specify:	
25 ga (570-0735)	21 ga (570-0733)
23 ga (570-0734)	

Microbore Extension Set

- Flat-top, needle-free valve with tight septum seal helps to support effective disinfection
- Metal-free for MRI compatibility
- Not made with natural rubber latex
- DEHP-free
- Sterile
- Single use

8", removable Luer-lock connector, slide clamp, rotating male Luer lock

(570-0884).....ea



Needle-Free Connector

- Flat-top, needle-free valve with tight septum seal helps to support effective disinfection
- Metal-free for MRI compatibility
- Not made with natural rubber latex
- DEHP-free
- Sterile
- Single use

(570-0882).....ea



Pressure-Activated Safety Lancets

- Application: Glucose and many other *in vitro* tests
- User safety: Needle is shielded before and after use
- Convenience: Preloaded, pressure-activated design
- Comfort: Ultrasharp and beveled needle or blade
- Efficacy: Various depths to ensure accurate sample


Safety Lancets.....	100/box
Specify:	
28 ga x 1.8", Purple	(900-7991)
21 ga x 1.8", Light Blue	(900-7992)
21 ga x 2.2", Orange	(900-7980)
21 ga x 2.8", Rose	(900-7982)
18 ga x 1.8", Green	(900-7993)





BRAND PRODUCTS

BLOOD CONTROL SAFETY IV CATHETERS PUSH BUTTON STYLE



ADDITIONAL FEATURES:

- Ground stainless-steel needle
- Straight hub
- Single-use
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant
- Radiopaque

Henry Schein® Blood Control Safety IV Catheters Push Button Style, 50/Box

	Safety IV Catheter Blood Control (Green), 18Gx1.16"	(570-3712)
	Safety IV Catheter Blood Control (Pink), 20Gx1.16"	(570-3713)
	Safety IV Catheter Blood Control (Pink), 20Gx1.0"	(570-3714)
	Safety IV Catheter Blood Control (Blue), 22Gx1.0"	(570-3715)
	Safety IV Catheter Blood Control (Yellow), 24Gx0.75"	(570-3716)



Sharps Containers

Shatter-resistant, puncture-resistant, and transparent lid allows for easy viewing. Ideal for safe disposal of needles, needles with syringes, vacutainers, lancets, and blades (e.g., razors and scalpels). Counterbalanced lid is designed to prevent infection/cross contamination by limiting access to container contents and facilitating one-handed operation. Constructed of puncture-resistant polypropylene. Safe to autoclave or incinerate and meets or exceeds FDA standards, including OSHA and EPA requirements for local biohazard waste-collection services

#0310-1500-HS, 1 qt (112-7149)	ea
#0310-150T-HS, 1.5 qt (570-0376)	ea
#0354-150B-HS, 5.4 qt (112-7151)	ea
#0354-150M-HS, 5 qt (112-7152)	ea
#0319-1500-HS, 1 gal, Slide Lid (570-0378)	ea
#0320-1500-HS, 2 gal, Slide Lid (112-7150)	ea
#0330-150R-HS, 3 gal (570-0212)	ea
#0380-150R-HS, 8 gal (570-0380)	ea
#0354-150M, M-Style, 5 qt (570-2554)	ea
#0320-1500, Slide Lid, 2 gal (570-2553)	ea



Infectious Waste Collection Bags

Constructed of high-grade polyethylene for safe collection, handling, and disposal of infectious wastes destined for incineration or landfill.

- Red and labeled with a biohazard symbol
- Twist ties provided

Small, 10 gal, 15" x 9" x 24" (100-9305)	100/box
Medium, 18 gal, 15" x 9" x 31" (100-8878)	100/box
Large, 33 gal, 23" x 8" x 41" (100-3416)	100/box



Sharps Containers with Rotary Lid

Rotating lid allows for maximum storage and safe disposal of larger sharps. High-impact polypropylene material. Fracture- and puncture-resistant. Lid locks for final disposal.

#0322-150R-HS, 2.2 qt (570-0213)	ea
#0319-150R-HS, 1 gal, Rotary Lid (570-0379)	ea
#0320-150R-HS, 2 gal, Rotary Lid (570-0377)	ea
#0330-150R-HS, 3 gal (570-0212)	ea
#0380-150R-HS, 8 gal (570-0380)	ea



Emesis Bags

- Eliminates exposure to emesis
- Minimizes spillage
- Opening supported by plastic notched ring to provide secure grip and unobstructed access into the bag
- Twist bag and secure notched ring after use
- Holds up to 40 fl oz or 1000 cc/mL
- High-density polyethylene: 1.35-mil thick
- Measures 16.5 cm x 38.1 cm

(570-1178)	24/pkg
-------------------------	--------

Trash Liners

These liners are a strong, durable, and affordable option, and fit many different shapes and sizes of trash cans.

- Available in Clear, Black, or White
- 11 different colors, sizes, and thicknesses
- Made in the USA



Clear

7–10 gal, 6 mic (570-1340)	1000/case
15–20 gal, 8 mic (570-1346)	1000/case
20–30 gal, 0.39 mil (570-1336)	250/case
20–30 gal, 10 mic (570-1341)	250/case
33 gal, 13 mic (570-1342)	250/case
40–45 gal, 12 mic (570-1343)	250/case



Black

33 gal, 0.51 mil (570-1337)	250/case
40–45 gal, 0.51 mil (570-1338)	250/case
60 gal, 1.4 mil (570-1339)	100/case



White

13 gal, 0.8 mil (570-1344)	100/case
20–30 gal, 0.85 mil (570-1345)	250/case



Henry Schein Brand Surgical Instruments

Since 1932, Henry Schein has offered an extensive range of products and solutions to optimize your practice and deliver the best care for your patients. Our surgical instruments portfolio is **high quality** and **competitively priced**, available in two grades, and tailored for a broad spectrum of specialties and procedures.

To help you maintain proper infection prevention protocols, we also provide cleaning and care solutions including instrument sterilization and reprocessing. Our ProRepair service includes manufacturer-trained technicians and highly skilled customer service representatives committed to providing quality repair services with a stress-free experience, with complete warranties backed by Henry Schein.

Our available surgical instrument grades include:




- **Maxima[†] by Henry Schein** – Premium surgical-grade instruments expertly crafted from German Stainless Steel* to conform to rigid specifications and standards that ensure our instruments consistently meet and exceed the high quality and performance expectations of our valued customers. Designed for your general practice needs as well as your more specialized surgical needs, such as **ambulatory surgery centers, plastic surgery and dermatology practices**.



- **Henry Schein Brand** – High-quality stainless steel instruments, designed for your medical practice/office use. Manufactured in ISO 13485 certified factories, unconditionally guaranteed against defects in quality or craftsmanship.

*Excludes products marked by †



Surgical Blades

- Sterile
- Individual foil packs
 - Computer-controlled grinding ensures a consistently sharp edge
 - Stainless steel

Stainless Steel.....100/box

Specify:

#10	(112-6199)	#20	(112-6184)
#11	(112-6190)	#21	(112-6185)
#12	(112-6177)	#22	(112-6195)
#15	(112-6178)	#23	(112-6200)
#15S	(112-7294)	#25	(112-6201)

Carbon Steel.....10/pkg

Specify:

#10	(112-6183)	#20	(112-6182)
#11	(112-6175)	#21	(112-6198)
#12	(112-6176)	#22	(112-6197)
#15	(112-6181)	#23	(112-6192)
#15S	(112-7295)		

Carbon Steel.....10/pkg

Specify:

#12B	(112-6180)	#15C	(112-6186)
------	-------------------	------	-------------------



Surgical Blade Handles

Standard Stainless Steel
 #3, For Blades #10, #11, #12B, #15 & #15C
(100-7520)..... ea

German Stainless Steel
 #3, For Blades #10, #11, #12B, #15 & #15C
(104-5814)..... ea

#4, For Blades #20, #21, #22 & #23
(100-3973)..... ea

#4, For Blades #20, #21, #22 & #23
(104-7296)..... ea

Air Collection* #3
(900-8098)..... ea

Air Collection* #4
(900-8099)..... ea



Disposable Biopsy Punch

- Seamless, stainless steel, sharpened blade for safe, clean, and easy sample extraction.
- Rigid, closed-end translucent guard protects blade from damage
 - Ergonomic, ribbed-grip handle provides comfort, control, and a clear view of the biopsy site
 - Handle is color coded for clear identification, and blade size is clearly marked
 - Free from phthalates
 - Sterile
 - Single use
 - Not made from natural rubber latex

Biopsy Punch..... ea

Specify:

1.5 mm, White	(570-3435)	4 mm, Blue	(570-3440)
2 mm, Ivory	(570-3436)	5 mm, Pink	(570-3441)
2.5 mm, Green	(570-3437)	6 mm, Sea Green	(570-3442)
3 mm, Gray	(570-3438)	8 mm, Purple	(570-3443)
3.5 mm, Yellow	(570-3439)		



Disposable Scalpels

- Stainless steel
- Plastic handle with millimeter scale and durable stainless steel blade
- Strong: reduced breakage
- Sharp: clean, true, consistent cut
- Packaging: superior protection and handling
- Individually peel-backed for convenience
- Testing: consistent high quality; assured quality to international standards

Scalpels.....10/box

Specify:

#10	(112-6196)	#15	(112-6187)
#11	(112-6189)	#15S	(112-7296)
#12	(112-6194)	#20	(112-6193)

#12B
(112-6188).....10/box

MAXIMA® Safety Scalpels



- Sterile
- Single use
- Retractable
- Easily extend or retract the blades while avoiding contact with the sharp blade edges
- Color coded for easy identification of blade size

Safety Scalpels.....10/box

Specify:

Size 10	(570-0243)	Size 12	(570-0244)
Size 11	(570-0245)	Size 15	(570-0246)



Instrument ID Tape Kit

Autoclavable, self-adhering strips for instrument identification. Available in 8 colors.

- Complete Set
(570-1979).....ea
 Contains: 1–3-ft. roll of each color: orange, blue, neon pink, green, white, neon purple, red & yellow.
- Refill Roll–10 ft.....ea
 Specify:
- | | | | |
|-----------|-------------------|-------------|-------------------|
| Orange | (570-1976) | White | (570-1971) |
| Blue | (570-1972) | Neon Purple | (570-1977) |
| Neon Pink | (570-1978) | Red | (570-1974) |
| Green | (570-1973) | Yellow | (570-1975) |



Color Code Instrument Rings

Brightly colored silicone rings. Autoclavable up to 392°F and can be reused. Standard size fits regular hand instruments. Large size has a larger diameter hole for use with the larger ergonomic handle instruments.

- Standard Size
 Assorted Kit
(570-1711).....ea
 Contains: 8 compartment cases with 600 standard-size code rings in 8 assorted colors.

- Standard Size
 Single Color..... 50/pkg
 Specify:
- | | | | |
|-------|-------------------|--------|-------------------|
| Black | (570-1707) | Purple | (570-1706) |
| Blue | (570-1702) | Red | (570-1701) |
| Green | (570-1704) | White | (570-1705) |
| Mauve | (570-1708) | Yellow | (570-1703) |

- Large Size
 Assorted Kit
(570-1720).....ea
 Contains: 8 compartment cases with 320 large-size code rings in 8 assorted colors.

- Large Size
 Single Color..... 50/pkg
 Specify:
- | | | | |
|-------|-------------------|--------|-------------------|
| Black | (570-1718) | Purple | (570-1717) |
| Blue | (570-1713) | Red | (570-1712) |
| Green | (570-1715) | White | (570-1716) |
| Mauve | (570-1719) | Yellow | (570-1714) |



Instrument Tray

Grooved instrument trays designed for maximum convenience and efficiency. Autoclavable plastic. Size: 4" x 8".

- Light White.....ea
 Specify:
- | | | | |
|-----|-------------------|----|-------------------|
| 16 | (570-1568) | 18 | (570-1555) |
| 16A | (570-1553) | 19 | (570-1556) |
| 17 | (570-1554) | 20 | (570-1557) |



Disposable Tray Inserts

Premolded disposable tray inserts that offer an organized, efficient, and effective method for eliminating cross contamination.

- Standard Size
 Fits Instrument Trays 11¼" x 7¼" x ¾"
(900-0937)..... 50/pkg
- Mini Size
 Fits Instrument Trays 5¾"W x 7¼"D x ¾"H
(900-7598).....100/pkg

CRITERION™ by HENRY SCHEIN®

COMFORT AND CONTROL:
 Designed for Surgical Procedures

CRITERION™ LATEX SURGICAL GLOVES

Natural color, 50 pair/box, Individually Packaged Pairs

SIZE	ITEM CODE	SIZE	ITEM CODE
5½"	570-4495	7½"	570-4499
6"	570-4496	8"	570-4500
6½"	570-4497	8½"	570-4501
7"	570-4498	9"	570-4502



Order today at HENRYSCHEIN.COM/CRITERION



Lister Bandage Scissors

Standard Stainless Steel

- 3½" (100-9243)..... ea
- 5½" (100-1087)..... ea
- 7¼" (100-3334)..... ea
- 8½" (953-2519)..... ea



Universal Scissors 7½"

Stainless Steel

Unusual blade design with 1 serrated edge. Cuts tin, light metal, wire, bandages, and more. Ideal for the lab, office, home, or first aid kit.

- (100-7898)..... ea



Wire-Cutting Scissors 4¾"

Standard Stainless Steel

Angled to Side, 1 Serrated Blade

- (100-2591)..... ea



MAXIMA® Knowles Bandage Scissors 5½"

German Stainless Steel

- Angled (102-7708)..... ea
- Straight (104-5504)..... ea



MAXIMA® Bandage and Utility Scissors

German Stainless Steel

- 6½" (104-9212)..... ea
- 7" (104-9280)..... ea



MAXIMA® Hercules Heavy-Duty Bandage and Plaster Shears 7½"

German Stainless Steel, Serrated Blade

- (104-7448)..... ea



MAXIMA® Lister Bandage Scissors

German Stainless Steel, Extra-Fine

- 3½" (104-9717)..... ea



- 4½" (104-8240)..... ea



- 5½" (100-6589)..... ea



- 7¼" (104-7364)..... ea



8", One Large Finger Ring

- (104-5841)..... ea



MAXIMA® Esmarch Heavy-Duty Bandage and Cast Shears 8"


German Stainless Steel

- (104-7447)..... ea



Northbent Stitch Scissors 5½"

German Stainless Steel, Curved
(100-0353)..... ea



MAXIMA® Gauze Shears 8"

Standard Stainless Steel
(100-4963)..... ea



MAXIMA® Shortbent Stitch Scissors 3½"

German Stainless Steel, Curved, Delicate
(104-7439)..... ea



Long Stitch Scissors 6¼"

German Stainless Steel, Extra Long
(100-0109)..... ea



MAXIMA® Littauer Jr. Stitch Scissors 4½"

German Stainless Steel, Slender Pattern
 Straight
(100-1420)..... ea



MAXIMA® Littauer Stitch Scissors 5½"

German Stainless Steel, Curved
(100-2354)..... ea



Littauer Jr. Stitch Scissors 4½"

Standard Stainless Steel, Slender Pattern
(100-6306)..... ea



Littauer Stitch Scissors 5½"

Standard Stainless Steel
(100-5115)..... ea



Spencer Stitch Scissors 3 1/2"

Standard Stainless Steel

(100-9827).....ea



MAXIMA® O'Brien Stitch Scissors 3 3/4"

German Stainless Steel, Angled Blade

(104-8241).....ea



MAXIMA® Stevens Tenotomy Scissors

German Stainless Steel

Curved.....ea

Specify:

4 1/2", Blunt Points **(104-7440)**

4 1/8", Blunt Points **(104-7441)**

4 1/8", Sharp Points **(104-9178)**

Scissors, Straight.....ea

Specify:

4 1/2", Blunt Points **(104-7443)**

4 1/2", Sharp Points **(104-7445)**

4 1/8", Sharp Points **(104-9639)**



MAXIMA® Spencer Stitch Scissors 3 1/2"

German Stainless Steel

(100-6024).....ea



MAXIMA® Gradle Eye Suture Scissors 3 3/4"

German Stainless Steel, Slightly Curved, Sharp Points

(101-7676).....ea



Stevens Tenotomy Scissors

Standard Steel, Blunt Points

Straight, 4 1/2" **(100-0973)**.....ea

Curved, 4 1/8" **(100-5500)**.....ea



FOR COMPLETE OFFERING OF INSTRUMENT CASSETTES SEE PAGE 31

←

HENRY SCHEIN®
BRAND PRODUCTS



MAXIMA® Dean Scissors 7"

German Stainless Steel, Sharp Points, 1 Serrated Blade
(100-8553)..... ea



MAXIIMA® Mayo Dissecting Scissors 6 3/4"

German Stainless Steel
 Straight
(100-4833)..... ea
 Curved
(100-1627)..... ea



Dean Scissors 7"

Standard Stainless Steel, Sharp Points, 1 Serrated Blade
(101-2187)..... ea



MAXIMA® Metzenbaum Scissors 5"

German Stainless Steel, Curved, Delicate
(104-7306)..... ea



MAXIIMA® Mayo Dissecting Scissors 5 1/2"

German Stainless Steel
 Straight
(100-4018)..... ea
 Curved
(100-3105)..... ea



Metzenbaum Scissors 5 1/2"

German Stainless Steel
 Curved
(100-8335)..... ea
 Straight
(101-8823)..... ea



MAXIMA® Metzenbaum Scissors 7"

German Stainless Steel

Curved
(100-5142)..... ea

Straight
(101-6892)..... ea



Metzenbaum Scissors 8"

German Stainless Steel, Curved

(104-8236)..... ea



MAXIMA® Strabismus Scissors 4"

German Stainless Steel

Curved
(104-8239)..... ea

Straight
(104-8237)..... ea



MAXIMA® Operating Room Scissors 5½"

HENRY SCHEIN®

German Stainless Steel

Curved, Blunt Points
(104-9957)..... ea



Curved, Sharp Points
(100-3874)..... ea



Straight, Blunt Points
(100-4352)..... ea



Straight, Sharp Points
(100-4549)..... ea



Straight, Sharp/Blunt Points
(100-7561)..... ea



MAXIMA® Operating Room Scissors 4 3/4"

German Stainless Steel, Straight

- Blunt Points
(104-8232)..... ea
- Sharp Points
(102-3440)..... ea
- Sharp/Blunt Points
(104-8231)..... ea



Goldman-Fox Scissors 5"

Standard Stainless Steel

- (101-2805)**..... ea



MAXIMA® Goldman-Fox Scissors 5"

German Stainless Steel

- (100-9837)**..... ea



MAXIMA® Operating Room Scissors 6"

German Stainless Steel, Straight, Sharp/Blunt Points

- (104-9283)**..... ea



Operating Room Scissors 6 1/2"

German Stainless Steel, Straight, Sharp Points

- (100-4104)**..... ea



Operating Room Scissors 6 1/2"

German Stainless Steel, Straight, Sharp/Blunt Points

- (100-4502)**..... ea



Operating Room Scissors, 5 1/2"

Curved, Sharp Points

- (100-7833)**..... ea



Operating Room Scissors 5 1/2"

Straight, Blunt Points

- (100-6022)**..... ea



Operating Room Scissors 5 1/2"

Curved, Sharp/Blunt Points

- (100-3875)**..... ea



Operating Room Scissors 5 1/2"

Straight, Sharp/Blunt Points

- (100-2372)**..... ea



MAXIMA® Iris Scissors

German Stainless Steel

- 3½", Curved, Delicate Pattern
(102-3439).....ea
- 4", Curved
(101-8111).....ea
- 3½", Straight, Delicate Pattern
(102-5248).....ea
- 4", Straight
(101-6893).....ea
- 4½", Straight
(100-8118).....ea



Iris Scissors 4½"

Standard Steel

- Curved
(100-8820).....ea
- Straight
(100-2767).....ea



MAXIMA® Wagner Plastic Surgery Scissors 4¾"

German Stainless Steel, Curved, Sharp Points, 1 Serrated Blade

- (100-4982)**.....ea



Kelly Angular Scissors 6¼"

German Stainless Steel, Sharp Points, 1 Serrated Blade

- Straight
(100-0112).....ea



- Curved
(100-3164).....ea



Kelly Scissors 6¼"

Standard Stainless Steel


- Angled
(101-4406).....ea
- Curved, Sharp Points, Smooth Blades
(101-3453).....ea
- Curved, Serrated Blade
(101-7287).....ea




MAXIMA® Nail Scissors 4"

German Stainless Steel, Curved

- (104-7810)**.....ea



Ruskin Bone-Cutting Forceps 6"
 Standard Stainless Steel, Standard Blades, Double Action
(104-6986)..... ea



Adson-Brown Dressing Forceps 4 3/4"
 Standard Steel, 7 x 7 Side-Grasping Teeth
(100-5162)..... ea



Adson Dressing Forceps 4 3/4"
 Standard Steel, Delicate Pattern
(100-5659)..... ea



MAXIMA® Adson-Brown Dressing Forceps 4 3/4"
 German Stainless Steel, 7 x 7 Side-Grasping Teeth
(101-6911)..... ea



MAXIMA® Adson Tissue & Suture Forceps 4 3/4"
 German Stainless Steel, Serrated Tips, Delicate Pattern
(101-6279)..... ea



MAXIMA® Micro Adson Dressing Forceps 4 3/4"
 German Stainless Steel, Serrated, 0.5-mm Wide Tips
(102-8861)..... ea



Dressing Forceps 5 3/4"
 German Stainless Steel, Serrated
(101-1272)..... ea



MAXIMA® Hudson Tissue Forceps 4 3/4"
 German Stainless Steel, 1 x 2 Teeth
(104-5547)..... ea



Thumb Forceps
 Standard Stainless Steel
 5 1/2" Standard Pattern
(100-5152)..... ea



MAXIMA® Adson Tissue & Suture Forceps 4 3/4"
 German Stainless Steel
 1 x 2 Teeth, Delicate Pattern
(101-4803)..... ea
 Cross-Serrated Tips, Delicate Pattern
(102-3435)..... ea
 1 x 2 Teeth, with Tying Platform
(101-6905)..... ea



6" Fenestrated Handle
(100-6424)..... ea



Adson Tissue Forceps 4 3/4"
 Standard Stainless Steel, 1 x 2 Teeth
(100-5891)..... ea



MAXIMA® Tissue Forceps
 German Stainless Steel, 1 x 2 Teeth
 5 1/2", Standard Pattern
(101-5132)..... ea
 6 1/2", Serrated Handle
(104-7960)..... ea



Allis Tissue Forceps

Standard Stainless Steel

- 6", 4 x 5 Teeth
(100-6914)..... ea
- 6", 5 x 6 Teeth
(100-6095)..... ea
- 7¼", 5 x 6 Teeth
(100-9430)..... ea



MAXIMA® Tissue Forceps

German Stainless Steel, Standard Pattern, 1 x 2 Teeth

- 4½"
(104-6618)..... ea
- 5"
(104-8243)..... ea



MAXIMA® Allis Tissue Forceps

German Stainless Steel

- 6", 4 x 5 Teeth
(101-1848)..... ea
- 7", 5 x 6 Teeth
(100-0228)..... ea



Tissue Forceps

Standard Stainless Steel, Standard Pattern

- 5½", 1 x 2 Teeth
(100-4877)..... ea
- 5½", 2 x 3 Teeth
(100-0389)..... ea
- 6", Fenestrated Handle, 1 x 2 Teeth
(100-9871)..... ea



Iris Tissue Forceps

Standard Pattern, 1 x 2 Teeth, 0.8-mm Wide Tips

- 4", Straight (100-0364)..... ea
- 4", Full Curve (100-8011)..... ea
- 4½", Half Curve (100-7103)..... ea



MAXIMA® Iris Tissue Forceps 4"

German Stainless Steel, Straight, Standard Pattern, 1 x 2 Teeth, 0.8-mm Wide Tips

- (101-6898)..... ea



MAXIMA® Bergh Cilia Forceps 3½"

German Stainless Steel, Horizontal Serrations, 5-mm Wide Jaws

- (101-8226)..... ea



MAXIMA® Eye Dressing Forceps 4"

German Stainless Steel, Standard Pattern, 0.8-mm Wide Tips

- Half Curved (101-7739)..... ea



- Straight (101-6897)..... ea



Eye Dressing Forceps 4"

Straight, Standard Pattern, 0.8-mm Wide Tips

- (100-5515)..... ea



Eye Dressing Forceps 4"

Half Curved, Standard Pattern, 0.8-mm Wide Tips

- (100-0673)..... ea



MAXIMA® Swiss Cilia and Suture/Splinter Forceps 3¾"

German Stainless Steel, Precision Fitted, Diagonal Jaws
Great for hair and suture removal.

- (102-3436)..... ea



Lucae Dressing Forceps 5 1/2"

Standard Stainless Steel, Serrated Tips, Bayonet Shape

(100-3781)..... ea



MAXIMA® Wilde Dressing Forceps 5"

German Stainless Steel, Serrated Tips, Angular

(104-5539)..... ea

Hartman-Noyes Alligator Forceps

Standard Steel, Serrated Jaws



3 1/2"
(100-7328)..... ea



5 1/2"
(100-0886)..... ea



6 1/2"
(100-6937)..... ea



MAXIMA® Hartman-Noyes Ear Forceps 3 1/4"

German Stainless Steel, Alligator Type, Serrated Jaws

(104-5508)..... ea



FOR COMPLETE OFFERING OF STERILIZATION POUCHES SEE PAGE 30

←

HENRY SCHEIN®
BRAND PRODUCTS



MAXIMA® Splinter Forceps

German Stainless Steel, Fine, Serrated Pointed Tips

- 3½" (104-5576).....ea
- 4" (104-6190).....ea



Plain Splinter Forceps 3½"

Standard Stainless Steel

- (100-3799).....ea



Plain Splinter Forceps 4½"

- (100-2111).....ea



MAXIMA® Plain Splinter Forceps

German Stainless Steel

- 3½" (101-7834).....ea



- 4½" (101-6895).....ea



MAXIMA® Walter Splinter Forceps

German Stainless Steel

- 4¼", Straight (102-3562).....ea
- Curved (100-0637).....ea



MAXIMA® Feilchenfeld Splinter Forceps

German Stainless Steel

- 3" (104-9311).....ea



- 4½" (104-7650).....ea



Carmalt Splinter Forceps 4¼"

German Stainless Steel, Curved


- (104-7583).....ea



Carmalt Splinter Forceps 4¼"

German Stainless Steel, Straight

- (104-8498).....ea



Virtus Splinter Forceps 6"
Standard Stainless Steel, Straight
(100-1538)..... ea



Kelly Forceps 5½"
Standard Stainless Steel
Curved
(100-7878)..... ea




Straight
(100-8745)..... ea



MAXIMA® Kelly Forceps 5½"
German Stainless Steel
Curved
(100-8381)..... ea



Halstead Mosquito Hemostat 5"
German Stainless Steel..... ea
Specify:
Straight (100-0846) Curved (100-2376)



Straight
(100-3357)..... ea



Standard Stainless Steel..... ea
Specify:
Straight (100-4662) Curved (100-4862)



MAXIMA® Hartman Mosquito Forceps 3½"
 German Stainless Steel, Straight
 (102-7483).....ea



MAXIMA® Hartman Mosquito Forceps 3½"
 German Stainless Steel, Curved
 (102-4346).....ea




Hartman Mosquito Forceps 3½"
 Standard Stainless Steel
 Straight
 (100-5499).....ea



Curved
 (100-7530).....ea



MAXIMA® Halstead Mosquito Forceps 5"
 German Stainless Steel
 Curved
 (100-2376).....ea




Curved, Extra Delicate
 (104-5506).....ea

Straight
 (100-0846).....ea

Straight, Extra Delicate
 (104-5500).....ea

www.HenrySchein.com/HSCares



Rochester-Ochsner Forceps 8"
 Straight, 1 x 2 Teeth
 (100-3431)..... ea




Rochester-Carmalt Forceps 8"
 Curved, Longitudinal Serrations
 (100-5793)..... ea



Rochester-Carmalt Forceps 8"
 Curved, Longitudinal Serrations
 (100-5793)..... ea



MAXIMA® Utility Sterilizer Forceps 8"
 German Stainless Steel, Straight
 (102-4137)..... ea



Rochester-Pean Hemostat 6 1/4"
 Economy, Straight
 (100-0629)..... ea



MAXIMA® Rochester-Pean Hemostat 6 1/4"
 German Stainless Steel, Curved
 (100-6870)..... ea



Rochester-Pean Hemostat 6 1/4"
 Economy, Curved
 (100-2005)..... ea



MAXIMA® Rochester-Pean Hemostat 6 1/4"
 German Stainless Steel, Straight
 (100-5738)..... ea



MAXIMA® Crile Forceps 5½"
German Stainless Steel, Curved
(100-1117)..... ea



MAXIMA® Crile Forceps 6¼"
German Stainless Steel, Curved
(104-7370)..... ea



MAXIMA® Crile Hemostat 5½"
German Stainless Steel, Straight
(100-8285)..... ea



MAXIMA® Crile Forceps 6¼"
German Stainless Steel, Straight
(104-7413)..... ea




Crile Hemostat Forceps 5½"
Standard Steel
Curved
(100-6966)..... ea



Rochester-Ochsner Forceps 6¼"
Straight, 1 x 2 Teeth
(100-0665)..... ea



Straight
(100-5322)..... ea



Foerster Sponge Forceps 7"
 Standard Stainless Steel, Straight, Serrated Jaws
 (100-6729).....ea



MAXIMA® Foerster Sponge Forceps 7"
 German Stainless Steel, Curved, Serrated Jaws
 (104-9341).....ea




MAXIMA® Foerster Sponge Forceps 9¾"
 German Stainless Steel, Serrated Jaws
 Straight
 (101-1737).....ea



MAXIMA® Fletcher Sponge Forceps 9½"
 German Stainless Steel, Straight
 (104-8249).....ea



Curved
 (104-8859).....ea



MAXIMA® Bozeman Sponge Forceps 10¼"
 German Stainless Steel, Curved
 (104-7369).....ea



MAXIMA® Foerster Sponge Forceps 7"
 German Stainless Steel, Straight, Serrated Jaws
 (104-9217).....ea



MAXIMA® Backhaus Towel Clamp 5¼"
 German Stainless Steel
 (101-3568)..... ea



MAXIMA® Lorna Nonperforating Towel Clamp 3½"
 German Stainless Steel
 (104-8244)..... ea



MAXIMA® Backhaus Towel Clamp
 German Stainless Steel
 3½"
 (101-4934)..... ea
 5¼"
 (101-3568)..... ea



MAXIMA® Lorna Nonperforating Towel Clamp 5¼"
 German Stainless Steel
 (104-8245)..... ea



Backhaus Towel Clamp 5¼"
 Standard Stainless Steel
 (100-0341)..... ea



Jones Towel Clamp 3½"
 Economy
 (100-8003)..... ea



MAXIMA® Castroviejo Needle Holder 5½"

Straight, with Lock

German Stainless Steel, Smooth Jaws

(100-2146)..... ea

Standard Stainless Steel

(953-9789)..... ea



MAXIMA® Castroviejo Needle Holder 5½"

Standard Stainless Steel, Straight, with Lock

(953-9789)..... ea



Collier Needle Holder 5"

Standard Stainless Steel, Fenestrated Jaws

(100-9278)..... ea



MAXIMA® Ryder Needle Holder

German Stainless Steel, Tungsten Carbide, Serrated Jaw

5"
(104-8251)..... ea

10"
(104-9323)..... ea



MAXIMA® Derf Needle Holder 4¾"

German Stainless Steel, Tungsten Carbide

(104-8597)..... ea

German Stainless Steel

(104-7419)..... ea



Derf Needle Holder 4¾"

Standard Stainless Steel

(100-7035)..... ea




Mathieu Needle Holder 5½"
Standard Stainless Steel
(100-0347).....ea



MAXIMA® Mathieu Needle Holder with Spring Handle 5½"
German Stainless Steel
(104-7424).....ea



Webster Needle Holder 5"
Standard Stainless Steel
Delicate Pattern, Smooth Jaws
(100-6371).....ea



Olsen-Hegar Needle Holder/Scissors Combination 5½"
Suture/Scissors Combo, Serrated Jaws
(100-1125).....ea



Olsen-Hegar Needle Holder/Scissors Combination 6½"
Serrated Jaws
(100-3586).....ea



Olsen-Hegar Needle Holder/Scissors Combination 7½"
(100-4564).....ea



MAXIMA® Webster Needle Holder 5"
German Stainless Steel, Delicate Pattern, Smooth Jaws
(101-9003).....ea



Tungsten Carbide
(104-6385).....ea



MAXIMA® Olsen-Hegar Needle Holder 5½"
German Stainless Steel
Suture/Scissor Combo, Serrated Jaws
(100-1395).....ea



MAXIMA® Olsen-Hegar Needle Holder 7¼"
German Stainless Steel, Suture/Scissors Combo
(104-7428).....ea



MAXIMA® Olsen-Hegar Needle Holder 6½"
German Stainless Steel, Suture/Scissor Combo, Serrated Jaws
(102-9182).....ea



MAXIMA® Olsen-Hegar Needle Holder/Scissors Combination 5³/₄"
 German Stainless Steel
 Tungsten Carbide, Serrated Jaws
(100-1432)..... ea



MAXIMA® Crile-Wood Needle Holder
 German Stainless Steel
 Serrated Jaws
 6"
(101-7128)..... ea
 6", Tungsten Carbide
(100-7943)..... ea
 7"
(104-8255)..... ea




Smooth Jaw
 6" Tungsten Carbide, Delicate
(104-5541)..... ea




MAXIMA® Baumgartner Needle Holder 5¹/₂"
 German Stainless Steel
(101-7415)..... ea




Tungsten Carbide, Serrated Jaw
(100-2001)..... ea



MAXIMA® Halsey Needle Holder 5"
 German Stainless Steel
 Serrated Jaws
(104-7420)..... ea



MAXIMA® Boynton Needle Holder 5"
(106-7948)..... ea




German Stainless Steel, Tungsten Carbide
(100-6378)..... ea



MAXIMA® Mayo-Hegar Needle Holder
German Stainless Steel, Tungsten Carbide



- 5¾", Serrated Jaw
(104-5517)..... ea
- 6"
(100-6116)..... ea
- 7¼"
(104-8253)..... ea



Mayo-Hegar Needle Holder
Standard Stainless Steel

- 5"
(100-0959)..... ea
- 6"
(100-2570)..... ea
- 8"
(100-6577)..... ea



Crile Wood Needle Holder 6"
Standard Stainless Steel

- (100-2492)**..... ea



MAXIMA® Halsey Needle Holder
Smooth Jaws

- 5", German Stainless Steel and Tungsten Carbide
(104-5514)..... ea
- 5¾", German Stainless Steel
(102-7437)..... ea



MAXIMA® Mayo-Hegar Needle Holder
German Stainless Steel

- 5¾", Delicate
(100-8670)..... ea
- 6"
(104-5531)..... ea
- 7¼"
(104-7548)..... ea



MAXIMA® Senn Retractor 6 3/8"

German Stainless Steel

Retractor.....ea

Specify:

Sharp (100-1499) 6 3/8", Blunt (100-1810)



Maxima Self-retaining Retractor with 4" Prongs

German Stainless Steel, Blunt Points

(104-6440)..... ea



Senn Retractor 6 3/8"

Standard Stainless Steel, Blunt Points, Double Ended

(100-1694)..... ea



Self-Retaining Retractor 4"

Sharp, Blunt Points


(104-5534)..... ea



MAXIMA® Ragnell Retractor 6"

German Stainless Steel, Delicate, Double Ended, 3mm x 8mm and 5mm x 15mm Blades

(104-8261)..... ea



MAXIMA® Alm Retractor 2 3/4"

German Stainless Steel, 2 1/4 Spread, 1/8-Deep Prongs

(104-7432).....ea



Buck Neurological Hammer 7³/₄"

German Stainless Steel, Chrome, with Brush and Needle (Screw Into Handle)

(104-5498)..... ea



MAXIMA® Weitlaner Retractor 5¹/₂"

German Stainless Steel, Sharp Points, 3 x 4 Teeth

(104-7436)..... ea



MAXIMA® Taylor Percussion Hammer 7"

German Stainless Steel, Solid Handle, Pediatric Size

(104-8462)..... ea



MAXIMA® Frazier Dura Hooks

German Stainless Steel

Sharp Point

(104-8994)..... ea

Blunt Point

(104-8263)..... ea



MAXIMA® Weitlaner Retractor 4"

German Stainless Steel, Blunt Points, 2 x 3 Teeth

(104-9635)..... ea



Surgical Mallet

7¹/₂" mallet with 2" stainless steel head. Used with a chisel to split teeth and/or remove bone.


(100-6094)..... ea



Laryngeal Mirrors

German Stainless Steel

- 20 mm
(104-5767)..... ea
- 24 mm
(104-6638)..... ea



Ear Syringes

With Shield

- 2 oz
(100-1128)..... ea
- 3 oz
(100-3595)..... ea
- 4 oz
(100-3751)..... ea



Vienna Nasal Speculum

Standard Steel, Standard Pattern

- Small
(100-0773)..... ea
- Medium
(100-1290)..... ea
- Large
(100-9171)..... ea



MAXIMA® Double-Action Nail Nipper 6¼"

German Stainless Steel, Sharp, Concave Jaws, Single Barrel
(104-5530)..... ea



Maxima Nail Nipper 4¾"

German Stainless Steel, Concave Jaws, Spring Barrel, Chrome Smooth barrel-spring action and precise sharp blades to cut the thickest toenails. Heavy duty.
(102-6341)..... ea



MAXIMA® Nail Nipper 5½"

German Stainless Steel, Concave Jaws, Double Spring
(102-5405)..... ea



MAXIMA® Nail Splitters

German Stainless Steel, Straight Jaws, Double Spring

- 4"
(104-5512)..... ea
- 5", Delicate
(104-7539)..... ea
- 5"
(104-5526)..... ea



MAXIMA® Nail Splitters 5"

Standard Stainless Steel, English Anvil, Straight Jaws, Double Spring Economy
(104-8507)..... ea

Maxima® Schamberg Comedone Extractor 3¾"

German Stainless Steel



Standard Pattern, Small Crimped Loop
(104-8113)..... ea



Standard Pattern, Crimped Loop, Extra Fine
(104-5520)..... ea



Square Loop Ends
(104-9155)..... ea



MAXIMA® Unna Comedone Extractor

German Stainless Steel, Spoon End

5¾"
(104-7841).....ea



Maxima® Fox Dermal Curettes

German Stainless Steel

2 mm
(104-5584)..... ea
 3 mm
(104-5731).....ea
 4 mm
(104-5548)..... ea



MAXIMA® Excavator Curette 5½"

German Stainless Steel, Double End, Small 1.5 mm and Medium 2 mm
(104-5536)..... ea



Billeau Flexible Ear Loop 6½"

Medium #2
(104-6676).....ea



Disposable Vaginal Specula

Recommended for routine gynecological exams. Made of hard, clear, nontoxic, high strength plastic resin.

- Allows examination of vaginal mucosa at a 360° angle without rotation
- Provides wide range of adjustability (angular and elevation adjustments)
- One-handed manipulation for ease of operation
- Self-holding in various opening positions
- Decreases risk of contamination
- If needed, apply lubrication for better comfort during examination
- Made of hard, clear, nontoxic, high-strength plastic resin

Small (102-5396)	10/pkg
Medium (102-5397)	10/pkg
Large (102-9422)	10/pkg



MAXIMA® Graves Vaginal Specula

German Stainless Steel

Small (104-5540)	ea
Medium (104-5522)	ea
Large (104-5535)	ea



Graves Vaginal Speculum

Original Model

Standard Steel.....	ea
Specify:	
Small (100-1888) Large (100-9184)	
Medium (100-6018) X-Large (953-0600)	



Disposable Vaginal Specula—Light Source Compatible

- Disposable specula made of high-strength, clear, nontoxic plastic
- Allow examination of the vaginal mucosa at 360° angle without rotation
- Provide wide range of adjustability (angular and elevation adjustments), ease of operation and manipulation with one hand
- Self-holding in various opening positions
- Single use of speculum decreases the risk of cross contamination
- Compatible with corded light source of leading brands, eliminating the need for positioning a light source independently

Specula.....	25/case
Specify:	
Small (112-6154) Large (112-6156)	
Medium (112-6155)	



Pederson Vaginal Speculum

Economy

Standard Steel.....	ea
Specify:	
Small (100-1170) Medium (100-3729)	
Large (100-9000) X-Large (953-0074)	



MAXIMA® Pederson Vaginal Speculum

German Stainless Steel, Narrow Blades

Small (104-8971)	ea
Medium (104-5601)	ea
Large (104-5543)	ea
X-Narrow (104-7925)	ea



MAXIMA® Uterine Sound SIMS 13"
 German Stainless Steel, Graduated in Centimeters
 (104-5527).....ea



MAXIMA® Uterine Sound SIMS
 German Stainless Steel, Graduated in Inches
 13"
 (104-7461).....ea



MAXIMA® Baby Tischler Biopsy Punch
 German Stainless Steel
 7¾"
 (104-6825).....ea



Maxima® Tischler Cervical Biopsy Punch Forceps
 German Stainless Steel, 6-mm x 3-mm x 1.5-mm Bite
 (104-5523).....ea



MAXIMA® Schroeder-Braun Uterine Tenaculum Forceps 9¾"
 German Stainless Steel, Straight
 (102-5362).....ea



MAXIMA® Kevorkian-Younge Endocervical Biopsy Curette 12"
 German Stainless Steel
 3-mm x 12-mm Loop with Basket
 (104-8110).....ea



Maxima® Kevorkian-Younge Uterine Biopsy Punch
 German Stainless Steel
 9¾"
 (104-9098).....ea



ER Laceration Tray

Sterile, Single-Use
Not made with natural rubber latex.
Hard plastic tray includes:

- 1-9¾" x 12" large center strip bag
- 1 precut 24" x 24" CSR wrap
- 2-19½" x 23" blue absorbent towels, scrim
- 1-18" x 26" fenestrated drape (black & white)
- 1-10-cc L/F Luer-lock syringe
- 1-18 ga x 1½" needle
- 1-27 ga x 1½" needle
- 1-25 ga x ¾" needle
- 1-5" Webster needle holder with large finger loops and beveled hinges
- 1-4½" pair of straight Iris scissors with large loops and tip protector (sharp/sharp)
- 1-5" curved mosquito hemostat with large finger loops
- 1-4¾" pair of Adson tissue forceps, 1" x 2"
- 2-2-oz medicine cups
- 1 sheet of universal medication labels (10/sheet)
- 5-2" x 2", 8-ply gauze dressings
- 10-4" x 4", 8-ply gauze pads

(112-5823)..... ea



**FOR COMPLETE OFFERING
OF WOUND CARE**
— SEE PAGES 87-89 —>



Suture Removal Kit with Plastic Handle Instruments

Sterile, Single-Use
Kit features plastic-handled Littauer scissors with plastic thumb forceps. Flatter, thinner-blade scissors provide easier suture removal by allowing blade to slip under suture. Not made with natural rubber latex.

- Hard plastic tray includes:
- 1 pair 4½" Littauer plastic handle scissors, hook-tipped
 - 1 pair 4" PosiGrip™ plastic thumb forceps
 - 1-3" x 3" gauze sponge

(900-4687)..... ea



Suture Removal Kit with Littauer Scissors

Sterile, Single-Use
Not made with natural rubber latex.
Hard plastic tray includes:

- 1-4½" metal Littauer scissors, hook-tipped
- 1-4" metal forceps
- 1-3" x 3" gauze sponge

(570-0619)..... ea



Staple Removal Kit

Used to easily remove surgical staples from patients with minimal discomfort. Not made with natural rubber latex.

- Hard plastic tray includes:
- 1 stainless-steel skin staple remover
 - 1-3" x 3" gauze sponge

(570-0620)..... ea



Suture Removal Kit with Iris Scissors

Sterile, Single-Use
Not made with natural rubber latex.
Hard plastic tray includes:

- 1 pair of small-loop, 4½" straight metal iris scissors (sharp/sharp)
- 1 serrated, 4¾" metal forceps
- 1-3" x 3" gauze pad
- 1 large alcohol prep pad

(112-5822)..... ea



Non-Fenestrated Drapes

Sterile, Polylined, Tissue/Poly/Tissue
• Heat sealed along 2 outer edges to ensure a pinhole-free plastic layer that guarantees maximum protection against bacterial strike-through from below
• Features polyethylene lining between 2 layers of absorbent-strength tissue, virtually eliminating risk of bacterial migration and strike-through
• Using blue-tissue side up prevents glare
• All sterile fields are available in convenient dispenser boxes of 50 each


White, 18" x 26"
(900-4686)..... 50/box



Fenestrated Drapes

Sterile, with 3" Round Fenestration
Tissue/Poly/Tissue
• Allows access to the operative site while it protects surrounding area
• 3" round fenestration option is die cut and heat sealed, leaving a smooth, even edge that prevents linting
• Not made with natural rubber latex

White, 18" x 26"
(900-4685)..... 50/box




General Drapes


Sterile

- Trilayered SMS fabric is soft, drapable, and naturally repellent to fluids, and its construction provides an excellent barrier to strikethrough
- Low lint-generating fabric reduces particulates and minimizes airborne particulate contamination
- Excellent flammability properties provide additional protection for the patient and staff when using lasers and electrocautery devices
- Not made with natural rubber latex
- Single use

Half, 44" x 60" (570-1656).....22/case




Medium, 44" x 76" (570-1617).....20/case



Large, 60" x 76" (570-1657).....20/case

Large, Reverse Fold, 60" x 76" (570-1663).....20/case




Chest Drapes

Sterile

- Trilayered SMS fabric is soft, drapable, and naturally repellent to fluids, and its construction provides an excellent barrier to strikethrough
- Low lint-generating fabric reduces particulates and minimizes airborne particulate contamination
- Excellent flammability properties provide additional protection for the patient and staff when using lasers and electrocautery devices
- Not made with natural rubber latex
- Single use

100" x 126" (570-1626).....12/case



Split-Sheet Drapes


Sterile

- Trilayered SMS fabric is soft, drapable, and naturally repellent to fluids, and its construction provides an excellent barrier to strikethrough
- Low lint-generating fabric reduces particulates and minimizes airborne particulate contamination
- Excellent flammability properties provide additional protection for the patient and staff when using lasers and electrocautery devices
- Not made with natural rubber latex
- Absorbent reinforcement
- Single use

Drape size: 53" x 76", split size: 2½" x 30" with hook-and-loop tube/cord holders, 2" adhesive tape (570-1647).....26/case

Drape size: 76" x 110", split size: 2½" x 37" with fabric tube/cord holders, 1" adhesive tape (570-1658).....24/case

Drape size: 76" x 124", split size: 3" x 25" with fabric tube/cord holders and 1" adhesive tape (570-1648).....22/case




Utility Drapes with Tape

Sterile

- Trilayered SMS fabric is soft, drapable, and naturally repellent to fluids, and its construction provides an excellent barrier to strikethrough
- Low lint-generating fabric reduces particulates and minimizes airborne particulate contamination
- Excellent flammability properties provide additional protection for the patient and staff when using lasers and electrocautery devices
- Not made with natural rubber latex
- Single use

15" x 26" (570-1684).....100/case



Back Table Covers, Zone Reinforced

Sterile

- Trilayered SMS fabric reinforcement
- For single use
- Not made with natural rubber latex

44" x 75" (570-1596).....20/case

44" x 90" (570-1597).....20/case

50" x 90" (570-1598).....20/case



Ophthalmic Drapes

- Made of trilayered SMS fabric, making it soft, drapable, and naturally repellent to fluids, and its construction an excellent barrier to strike-through
- Low lint-generating fabric reduces particulates and minimizes airborne particulate contamination
- Excellent antifiammability properties provide additional protection for the patient and staff when using laser and electrocautery devices
- Not made with natural rubber latex
- Single use

With Full Incise, 48" x 51" (570-1637).....30/case

With Full Incise and Pouch (570-1639).....30/case




Mayo Stand Cover

Sterile

- SMS barrier fabric reinforcement
- Not made with natural rubber latex
- Single use

Regular, 23" x 54" (570-1616).....30/pkg

X-large, 29½" x 57" (570-1618).....30/pkg




EENT Drape Pack

Sterile

Contents include:

- 1 EENT drape, 76" x 124", with 3" x 25" Split
- 1 head drape, 50" x 36"
- 1 surgical gown, non-reinforced, set-in sleeve, Large
- 1 surgical gown, non-reinforced, set-in sleeve, Large, overwrap with towel
- 1 back table cover, zone reinforced, 44" x 90"
- 1 mayo stand cover, 23" x 54"
- 1 white hand towel, 15½" x 22"
- 1 disposable bag

(570-1652).....4/case




Basic Drape Pack 2

Sterile

Contents include:

- 4 utility drapes with tape, 15" x 26"
- 1 large drape, 60" x 76"
- 1 surgical gown, non-reinforced, set-in sleeve, Large, overwrap with towel
- 1 surgical gown, non-reinforced, set-in sleeve, X-large
- 1 back table cover, zone-reinforced, 44" x 90"
- 1 mayo stand cover, 23" x 54"

(570-1599).....4/case



Basic Drape Pack 3

Sterile

Contents include:

- 4 utility drapes with tape, 15" x 26"
- 1 half drape, 44" x 60"
- 1 large drape, 60" x 76"
- 1 surgical gown, non-reinforced, set-in sleeve, Large, overwrap with towel
- 1 surgical gown, non-reinforced, set-in sleeve, X-large
- 1 back table cover, zone-reinforced, 44" x 90"
- 1 mayo stand cover, 23" x 54"

(570-1600).....6/case

FIND THE PERFECT PAIR OF GLOVES

SEE PAGES 14-18 for our complete offering of Criterion Gloves ←

CRITERION™

by HENRY SCHEIN®



Single Basin Set

(570-1570)..... 10/case
 Contains: 1-6000-cc ring basin & 1-54" x 54" CSR wrap.



Double Basin Set

(570-1571)..... 10/case
 Contains: 1-6000-cc ring basin & 2-54" x 54" CSR wraps.



Major Set-up Procedure Pack

(570-1573)..... 6/case

Contains:	Qty	Contains:	Qty
Disposal bag	1	1200-cc pitcher	1
700-cc emesis basin	1	Yankauer suction handle	1
6000-cc ring basin	1	60-mL irrigating syringe bulb	1
#15 carbon blade	1	¼" x 12" suction tubing	1
32-oz/1000-mL bowl	1	54" x 54" CSR wrap	1
Cautery pencil	1	#10 carbon blades	2
Cautery-tip cleaner	1	Flexible light handle cover	1
60-cc medicine cup	1	18" x 18" 4-ply X-ray laparotomy sponges	5
Standard-tip skin marker	1	Blue cotton towels	6
Needle counter	1	4" x 4" 16-ply X-ray gauze sponges	10



Minor Set-up Procedure Pack

(570-1572)..... 8/case
 Contains: 1 of each: 700-cc emesis basin, 6000-cc ring basin, 32-oz/1000-mL bowl, 60-cc medicine cup; & 54" x 54" CSR wraps.



Cardiac Catheterization Procedure Pack

(570-1588)..... 6/case

Contains:	Qty	Contains:	Qty
36" x 40" band bag	1	25 ga x 1½" needle	1
#11 carbon blade	1	8 ga x 7 cm needle introducer	1
2500-mL guidewire bowl	1	10-mL Luer-lock 3-ring control syringe	1
32-oz/1000-mL bowl	1	XL surgical gowns	2
Table back cover	1	20-mL Luer-lock syringes	2
Set-up cover	1	White paper towels	2
60-cc medicine cup	1	10-mL Luer-lock syringes	4
Femoral angiography drape	1	5½" forcep towels	4
5" curved Halstead Mosquito	1	Blue cotton towels	4
forcep		4" x 4" 4-ply nonwoven sponges	40
Sharps holder	1		
Waterproof label set	1		
Dual-tip skin marker	1		
18 ga x 1½" needle	1		



Basic Laparotomy Procedure Pack

(570-1589).....3/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	6" ruler	1
700-cc emesis basin	1	Yankauer suction handle	1
Mayo stand cover	1	10-mL Luer-lock syringe	1
32-oz/1000-mL bowl	1	60-mL irrigating bulb syringe	1
Specimen cup with screw lid	1	¼" x 12" suction tubing	1
Laparotomy drape	1	54" x 54" CSR wrap	1
X-large fabric-reinforced surgical gown with towel	1	#10 carbon blades	2
Waterproof label set	1	Larg e film-reinforced surgical gowns	2
Dual-tip skin marker	1	60-cc medicine cups	2
Needle counter	1	Utility drapes	4
22 ga x 1½" needle	1	Blue low-lint towels	4
25 ga x 1½" Eclipse™ needle	1	Blue cotton towels	6
1200-cc pitcher	1	4" x 4" 16-ply X-ray gauze sponges	10
		18" x 18" 4-ply X-ray laparotomy sponges	10



Breast Procedure Pack

(570-1690).....3/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	6" ruler	1
700-cc emesis basin	1	Yankauer suction handle	1
6000-cc ring basin	1	10-mL Luer-lock syringe	1
#15 carbon blade	1	60-mL irrigating bulb syringe	1
32-oz/1000-mL bowl	1	¼" x 12" suction tubing	1
Table back cover	1	54" x 54" CSR wrap	1
Mayo stand cover	1	#10 carbon blades	2
Specimen cup with screw lid	1	60-cc medicine cups	2
Chest drape	1	Blue cotton towels	4
XL surgical gown	1	Blue low-lint towels	4
XL surgical gown with towel	1	Utility drapes	4
Waterproof label set	1	4" x 4" 16-ply X-ray gauze sponges	10
Dual-tip skin marker	1	18" x 18" 4-ply laparoscopic X-ray sponges	10
Needle counter	1		
22 ga x 1½" needle	1		
25 ga x 1½" Eclipse™ needle	1		
1200-cc pitcher	1		



Laparoscopy Procedure Pack

(570-1590).....3/case

Contains:	Qty:	Contains:	Qty:
Antifog solution kit	1	22 ga x 1½" needle	1
Disposal bag	1	25 ga x 1½" Eclipse™ needle	1
700-cc emesis basin	1	6" ruler	1
6000-cc ring basin	1	#11 scalpel	1
#15 carbon blade	1	10-mL Luer-lock syringe	1
Back table cover	1	Insufflation tubing with filter	1
Mayo stand cover	1	¼" x 12" suction tubing	1
32-oz (1000-mL) bowl	1	54" x 54" CSR wraps	2
Specimen cup with screw lid	1	60-cc medicine cups	2
Laparotomy drape	1	X-large surgical gowns	2
X-large surgical gown with towel	1	Blue low-lint towels	2
Waterproof label set	1	Blue cotton towels	4
Dual-tip skin marker	1	4" x 4" 16-ply X-ray gauze sponges	10
Needle counter	1	18" x 18" 4-ply X-ray lap sponges	10



Basic Eye Procedure Pack

(570-1583).....8/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	Eye pad	1
Self-closure bag	1	Eye shield	1
32-oz (1000-mL) bowl	1	3-mL Luer-lock syringe	1
Back table cover	1	1-cc Luer-slip syringe	1
Mayo stand cover	1	2-oz bulb syringe	1
Ophthalmic drape	1	Instrument wipe	1
X-large surgical gown	1	60-cc medicine cups	2
X-large surgical gown with towel	1	Blue low-lint towels	4
Waterproof label set	1	3" applicators	10
Standard-tip skin marker	1	4" x 4" 12-ply gauze sponges	10
Needle counter	1	Eye spears	20
25 ga x ⅝" needle	1		



Plastics Procedure Pack

(570-1584).....6/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	Frazier suction handle	1
700-cc emesis basin	1	10-mL Luer-lock syringe	1
#15 stainless blade	1	3-cc Luer-lock syringe	1
32-oz/1000-mL bowl	1	60-mL irrigating bulb syringe	1
Cautery pencil	1	0-mL Luer-lock control syringe	1
Table back cover	1	¼" x 12" suction tubing	1
Mayo stand cover	1	60-cc medicine cups	2
Head drape	1	XL surgical gowns	2
Split drape	1	Utility drapes	4
XL surgical gown with towel	1	Blue cotton towels	6
Dual-tip skin marker	1	6" applicators	10
Needle counter	1	Blank labels	10
1-pocket diathermy pouch	1	Eye spears	10
6" ruler	1	4" x 4" 16-ply gauze X-ray sponges	10
#11 safety scalpel	1		



Hand Procedure Pack

(570-1581).....6/case

Contains:	Qty:	Contains:	Qty:
6½" x 12" suture bag	1	X-large surgical gown	1
6½" x 10" ziplock bag	1	X-large surgical gown with towel	1
4" x 108" Esmark bandage	1	Standard-tip skin marker	1
700-cc emesis basin	1	Needle counter	1
#15 carbon blade	1	25 ga x 1½" needle	1
32-oz (1000-mL) bowl	1	6" ruler	1
Back table cover	1	2-ply stockinette	1
Mayo stand cover	1	10-mL Luer-lock syringe	1
60-cc medicine cup	1	30-mL Luer-lock syringe	1
Hand drape	1	60-mL irrigating bulb syringe	1
Impervious split drape	1	4" x 4" 16-ply gauze	1
1" x 8" Xeroform dressing	1	Blue cotton towels	6



Head and Neck Procedure Pack

(570-1586).....2/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	8" x 3" Telfa pad	1
700-cc emesis basin	1	1200-cc pitcher	1
6000-cc ring basin	1	6" ruler	1
#10 carbon blade	1	Yankauer suction handle	1
#15 carbon blade	1	10-mL Luer-lock syringe	1
Back table cover	1	60-mL irrigating bulb syringe	1
Mayo stand cover	1	¼" x 12" suction tubing	1
Specimen cup with screw lid	1	54" x 54" CSR wrap	1
Half-drape	1	32-oz (1000-mL) bowls	2
Head drape	1	60-cc medicine cups	2
Split drape	1	X-large surgical gowns	2
X-large surgical gown with towel	1	Blue low-lint towels	4
Waterproof label set	1	8" x 18" 4-ply X-ray lap sponges	5
Dual-tip skin marker	1	Blue cotton towels	6
Needle counter	1	4" x 4" 16-ply X-ray gauze sponges	20
18 ga x 1½" needle	1		
25 ga x 1½" needle	1		



Knee Arthroscopy Procedure Pack

(570-1575).....4/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	18 ga x 3½" spinal needle	1
6" x 5-yd elastic bandage	1	6" x 4-yd cast padding	1
6" x 12" Esmark bandage	1	6" ruler	1
4" x 5-yd tan cohesive bandage	1	Large impervious stockinette	1
#11 carbon blade	1	10-mL Luer-lock syringe	1
32-oz (1000-mL) bowl	1	30-mL Luer-lock syringe	1
Back table cover	1	9" x 5" x 2" half-deep tray	1
Mayo stand cover	1	¼" x 12" suction tubing	1
60-cc medicine cup	1	30" x 30" poly wrap	1
Arthroscopy drape	1	Self-closure bags	2
Large drape	1	Breathable X-large surgical gowns	2
Impervious split drape	1	Blue cotton towels	8
1" x 8" Xeroform dressing	1	4" x 4" 16-ply X-ray gauze sponges	10
Waterproof label set	1	4" x 4" 12-ply gauze sponges	10
Standard-tip skin marker	1		
Needle counter	1		
18 ga x 1½" needle	1		



Nasal & Sinus Procedure Pack

(570-1587).....4/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	Needle counter	1
Self-closure bag	1	25 ga x 1½" needle	1
700-cc emesis basin	1	6" ruler	1
32-oz (1000-mL) bowl	1	Yankauer suction handle	1
Back table cover	1	3-oz bulb syringe	1
Mayo stand cover	1	9" x 5" x 2" half-deep tray	1
Specimen cup with screw lid	1	¼" x 12" suction tubing	1
6" tongue depressor	1	60-cc medicine cups	2
Head drape	1	Blue low-lint towels	2
Split drape	1	Blue cotton towels	4
X-large surgical gown	1	6" applicators	6
X-large surgical gown with towel	1	4" x 4" 16-ply X-ray gauze sponges	10
Waterproof label set	1	4" x 4" 12-ply gauze sponges	10
Dual-tip skin marker	1		



Shoulder Arthroscopy Procedure Pack

(570-1578).....6/case

Contains:	Qty:	Contains:	Qty:
6½" x 12" suture bag	1	18 ga x 3½" spinal needle	1
4" x 5-yd tan cohesive bandage	1	6" ruler, medium impervious	1
#11 carbon blade	1	Stockinette	
#15 carbon blade	1	10-mL Luer-lock syringe	1
Table back cover	1	30-mL Luer-lock syringe	1
Mayo stand cover	1	9" x 5" x 2" half-deep tray	1
Large drape	1	¼" x 12" suction tubing	1
Split shoulder drape	1	U-drape	1
Impervious split drape	1	4" x 4" 16-ply gauze X-ray sponge	10
XL surgical gown with towel	1	60-cc medicine cups	2
Standard-tip skin marker	1	XL X-long breathable surgical gowns	2
Needle counter	1	Blue cotton towels	8
18 ga x 1½" needle	1		



ENT Basic Procedure Pack

(570-1585).....6/case

Contains:	Qty:	Contains:	Qty:
Disposal bag	1	Needle counter	1
#15 carbon blade	1	27 ga x 1¼" aluminum-hub A-bevel needle	1
16-oz (500-mL) bowl	1	6" ruler	1
32-oz (1000-mL) bowl	1	10-mL Luer-lock control syringe	1
Back table cover	1	¼" x 12" suction tubing	1
Mayo stand cover	1	Blue cotton towels	3
Half-drape	1	60-cc medicine cups	4
Split drape	1	Blank labels	10
X-large surgical gown	1	4" x 4" 16-ply X-ray gauze sponges	10
X-large surgical gown with towel	1		
Dual-tip skin marker	1		



Podiatry Procedure Pack

(570-1582).....4/case

Contains:	Qty:	Contains:	Qty:
6½" x 12" suture bag	1	Large impervious stockinette	1
4" x 108" Esmark bandage	1	10-mL Luer-lock syringe	1
4" x 5-yd tan cohesive bandage	1	30-mL Luer-lock syringe	1
Back table cover	1	60-mL irrigating bulb syringe	1
Mayo stand cover	1	¼" x 12" suction tubing	1
Lower-extremity drape	1	18" x 18" 4-ply X-ray lap sponge	1
X-large surgical gown	1	4" x 4" 16-ply X-ray gauze sponge	10
X-large surgical gown with towel	1	4" x 4" 12-ply gauze sponges	10
Waterproof label set	1	#15 carbon blades	2
Dual-tip skin marker	1	60-cc medicine cups	2
Needle counter	1	Blue low-lint towels	2
18 ga x 1½" needle	1	6½" ziplock bags	3
25 ga x 1½" needle	1	Blue cotton towels	4
6" ruler	1		



Extremity Procedure Pack

(570-1576).....6/case

Contains:	Qty	Contains:	Qty
6½" x 12" suture bag	1	Standard-tip skin marker	1
4" x 5-yd elastic bandage	1	Needle counter	1
4" x 108" Esmark bandage	1	23 ga x 1½" needle	1
700-cc emesis basin	1	6" ruler	1
#10 carbon blade	1	Large impervious stockinette	1
32-oz (1000-mL) bowl	1	10-mL Luer-lock syringe	1
Cautery pencil with holster	1	2-oz bulb syringe	1
Back table cover	1	¼" x 12-ft suction tubing	1
Mayo stand cover	1	4" x 4" 16-ply gauze	10
60-cc medicine cup	1	4" x 4" 12-ply gauze	10
Large drape	1	6½" x 10" ziplock bags	2
Lower-extremity drape	1	#15 carbon blades	2
Impervious split drape	1	Flexible light-handle covers	2
Cautery electrode	1	Waterproof label sets	3
X-large surgical gown	1	24" x 17" blue cotton towels	7
X-large surgical gown with towel	1		



Upper Extremity Procedure Pack

(570-1580).....4/case

Contains:	Qty	Contains:	Qty
Disposal bag	1	6" ruler, large impervious stockinette	1
Back table cover	1	10-mL Luer-lock syringe	1
Mayo stand cover	1	60-mL irrigating bulb syringe	1
60-cc medicine cup	1	9" x 5" x 2" half-deep tray & ¼" x 12"	1
Lower-extremity drape	1	Suction tubing	1
Three-quarter drape	1	1" x 8" Xeroform dressing	1
1" x 8" Xeroform dressing	1	#15 carbon blades	2
X-large surgical gown	1	Blue cotton towels	6
X-large surgical gown with towel	1	4" x 4" 16-ply X-ray gauze sponges	10
Standard-tip skin marker	1	18" x 18" 4-ply X-ray lap sponges	10
Needle counter	1		
25 ga x 1½" needle	1		



Shoulder Beach Chair Procedure Pack

(570-1579).....4/case

Contains:	Qty	Contains:	Qty
Disposal bag	1	18 ga x 1½" needle	1
4" x 5 yd tan cohesive bandage	1	18 ga x 3½" spinal needle	1
#11 carbon blade	1	6" ruler, medium Impervious stockinette	1
#15 carbon blade	1	10-mL Luer-lock syringe	1
Table back cover	1	30-mL Luer-lock syringe	1
Mayo stand cover	1	9" x 5" x 2" half-deep tray	1
Large drape	1	¼" x 12" suction tubing	1
Beach chair shoulder drape	1	60-cc medicine cups	2
Impervious split drape	1	XL X-long breathable surgical gowns	2
XL surgical gown with towel	1	Blue cotton towels	8
Standard-tip skin marker	1	4" x 4" 16-ply gauze	10
Needle counter	1	X-ray sponges	

HENRY SCHEIN
PRACTICE GREEN

Care for Your Patients and the Planet

www.henryschein.com/practicegreen

WEAR IT WITH
CONFIDENCE.

CRITERION™

A full range of PPE options just right for healthcare professionals. CRITERION™ provides a smart selection of PPE products that let you outfit your practice from **head-to-toe** with quality you can rely on.

FROM HERE



TO HERE



For complete offering
of PPE see pages 14–26.

CRITERION™
COMFORT IN PROTECTION™

by HENRY SCHEIN®

Standard Series Surgical Gowns, AAMI 3

- AAMI 3 performance in the critical zones
- Premium SMS barrier fabric
- Nonreinforced
- Set-in sleeves
- Adjustable hook-and-loop closures
- Includes scrim hand towel



Medium (570-1630).....	32/case
Large (570-1688).....	30/case
X-large (570-1629).....	28/case
2X-large (570-1631).....	26/case
3X-large (570-1662).....	24/case

SOFT Series Surgical Gowns, AAMI 3

- Non-Reinforced, Set-In Sleeves
- AAMI 3 performance in the critical zones
- SOFT SMS barrier fabric
- Adjustable hook-and-loop closures
- Scrim hand towel



Medium (570-1664).....	32/case
Large (570-1665).....	30/case
X-large (570-1666).....	28/case
2X-large (570-1667).....	26/case
3X-large (570-1668).....	24/case

Surgical Gowns, AAMI 4

- AAMI 4 performance in the critical zones
- Breathable viral-barrier fabric
- Raglan sleeves
- Sleeve-n-sleeve construction
- Adjustable hook-and-loop closures
- Includes scrim hand towel



Regular Length

Small/Medium (570-1633).....	30/case
Large (570-1634).....	28/case
X-Large (570-1636).....	26/case
2X-large (570-1638).....	24/case

Long Length

Large (570-1640).....	26/case
X-large (570-1641).....	24/case
2X-large (570-1645).....	22/case



O.R. Towels

- Prewashed, 100% high-quality cotton
- Highly absorbent and virtually lint-free
- All sterile towels are CSR-wrapped and have a special corner fold that allows for aseptic presentation into a sterile field
- Not made with natural rubber latex
- 17" x 27"
- Blue

Non-Sterile (570-1113).....	100/case
Sterile 2s (570-1117).....	2/pkg
4s (570-1116).....	4/pkg
6s (570-1114).....	6/pkg



Laparotomy Sponges

- Prewashed, 100% high-quality cotton with tight weave
- Cross-stitched with super-locked edges to minimize linting and loose threads
- X-ray detectable
- Banded 5s
- Not made with natural rubber latex

18" x 18" (570-1115).....	5/pkg
------------------------------	-------



Double-Sided Shave Prep Razors

- Disposable
- Sharp stainless steel blades for smooth shaving
- Nonslip grip handle

(900-7665).....	100/box
-----------------	---------



Surgical Skin Markers

These surgical marking pens provide an accurate and reliable way to mark the incision site prior to surgery.

- Prep-resistant
- Nonabrasive tip
- Large ink reservoir

Non-Sterile, Fine Tip, Ruler, and Up to 9 Labels (900-7018).....	100/box
---	---------

Sterile, Fine Tip, Ruler, and Up to 6 Labels (900-7019).....	50/box
---	--------

Mini (WriteSite Jr.), Dark Blue (900-7020).....	100/box
--	---------



Alcohol Prep Pads

Sterile

Soft, absorbent, and non-woven pads impregnated with 70% isopropyl alcohol USP. Primarily used as a topical antiseptic; disinfects prior to administering injections.

Medium, 30 mm x 30 mm, Folded (570-2219).....	200/box
Large, 70 mm x 62.5 mm (570-1307).....	200/box

Povidone-Iodine Swabsticks and Prep Pads

Sterile

Features and benefits:

- Topical antiseptic used for skin preparation that helps to reduce bacteria that can cause infection
- Saturated with 10% USP povidone-iodine prep solution (equivalent to 1% titratable iodine)
- Not made with natural rubber latex

Single Pack, 4"L (570-1549).....	25 pkg/box
Contains: 1 swabstick per pkg.	

Triple Pack, 4"L (570-1548).....	25 pkg/box
Contains: 3 swabsticks per pkg.	

Single Pouches, 1½" x 1½", Folded (570-1543).....	1000/case
Contains: 1 per pouch.	



WOUND CLOSURE



TOPICAL SKIN ADHESIVE - OCTYL

High-Strength 2-Octyl Cyanoacrylate Formulation

Advantages

- Rapid set time
- Flexibility without a plasticizer
- Ergonomic ampule designed for ease of handling and targeted application
- Provides an effective microbial barrier and water-resistant protection
- High tensile strength
- Unique formulation will not solidify in ampule, allowing all .5 ml to be used
- Gravity fed applicator delivery system - no need to squeeze ampule



Glass-free reservoir containing 0.5mL of Topical Skin Adhesive

Rigid sleeve helps prevent bending of applicator during adhesive application

Angled applicator tip changes color in the presence of adhesive

Perforation in the middle, allowing ease of adhesive flow

Additional Information

Polymerization starts when adhesive comes in contact with tissue or fluid.

All 0.5 ml of adhesive can be used on the patient, as adhesive will not solidify in the ampule.

Topical Skin Adhesive will remain attached to the skin surface until natural cellular regeneration occurs (approx. 5-7 days). If needed, any standard adhesive remover product may be used to release the adhesive from the skin. Peel off the film, do not pull the skin apart.

Henry Schein® Topical Skin Adhesive - Octyl
0.5mL ampules, Clear, 12/box.....(570-3554)

Order today at HENRYSCHEIN.COM/MEDICAL



EXCLUSIVE BRANDS & PRODUCTS



QUALITY YOU CAN RELY ON



SMART SAVINGS

Wound Closure Strips



- Adhesive wound closures with an exceptionally porous, spun-bond fiber backing include bilateral nylon filaments for high tensile strength
- To be used in addition to or in place of sutures
- Each size is in a color-coded package for quick and easy size identification
- Medical-grade adhesive properties: Affords excellent adhesion, resistance to moisture, superior longevity, and is gentle to the skin
- Hypoallergenic: Designed for patient comfort and longer-term exposure
- Not made with natural rubber latex

- 3/8" x 3" (900-4317).....250/box
Contains: 5 strips per envelope, 50 envelopes per box.
- 1/4" x 3" (900-4315).....150/box
Contains: 3 strips per envelope, 50 envelopes per box.
- 1/4" x 4" (900-4318).....500/box
Contains: 10 strips per envelope, 50 envelopes per box.
- 3/4" x 1 1/2" (900-4316).....300/box
Contains: 6 strips per envelope, 50 envelopes per box.
- 1/2" x 4" (900-4314).....300/box
Contains: 6 strips per envelope, 50 envelopes per box.
- 1" x 4" (900-7506).....100/box
Contains: 4 strips per envelope, 25 envelopes per box.



Topical Skin Adhesive—Butyl

Sterile, 0.5-mL Single-Patient-Use Plastic Ampules

The proven choice of wound closure with over 40 years of clinical experience, this FDA-approved liquid adhesive consists of n-butyl-2 cyanoacrylate.

- Highest bond strength of all topical tissue adhesives on the market
- Will not clog or dry out when opened
- Provides a precise, controlled application
- Inhibits growth of bacteria
- Available in dyed and clear versions
- Flexible version provides more flexible coverage for longer incisions

Nonflexible, without Tip.....ea
Specify:
Clear (112-6106) Blue (112-6111)

Flexible Ampule with Tip, Blue (570-0351).....ea

Flexible, with Tip, Blue (112-7158).....10/box

Topical Skin Adhesive-Octyl



This unique formulation will not solidify in ampule, allowing all 0.5 mL to be used. The ergonomic ampule is designed for ease of handling and targeted application. Sterile, single-patient use.

Made with 2-octyl cyanoacrylate, this formulation provides:

- Rapid set time
- High tensile strength
- Flexibility without a plasticizer
- An effective microbial barrier



0.5-mL Plastic Ampule, Clear (570-3554).....12/box

Disposable Skin Stapler

- Enables user to work with higher efficiency and assists in avoiding cross infection effectively
- Specially designed for surgical skin suturing
- Novel structure and application security allow quick speed when sewing up during operation, modest histological reaction, and good wound matching and healing
- No pain during removal
- Scarring is minimal



Skin Stapler, 35 Wide (570-2956).....ea

Skin Stapler, 35 Regular (570-2957).....ea



Skin Staple Remover (900-7011).....ea

Maxima™ PTFE SURGICAL SUTURES

- Maxima PTFE (polytetrafluoroethylene) sutures are an excellent choice for dental bone grafting and implant procedures.
- Proprietary monofilament construction reduces the potential for bacterial migration into the surgical site.
- Ultrasoft, high-density material is designed with a low surface porosity.
- Black premium needles are made from 300-series stainless steel.



ITEM CODE	OSTEOGENICS CYTOPLAST**	NEEDLE	SIZE	LENGTH	RECOMMENDED USAGE	12/BOX
(112-7253)	CS0418	FS-2 (19mm, 3/8 reverse cut)	2-0	18	Larger suture size with a popular needle for procedures where a longer needle is desired	
(112-7254)	CS051819	FS-2 (19mm, 3/8 reverse cut)	3-0	18	A popular suture size and standard sized needle for general dental implant and bone grafting	
(112-7255)	CS0518	FS-3 (16mm, 3/8 reverse cut)	3-0	18	A popular suture size and smaller needle for dental implant and bone grafting	
(112-7256)	CS0618 RC	FS-3 (16mm, 3/8 reverse cut)	4-0	18	Most popular suture size for dental implant and bone grafting	
(112-7257)	CS0618 PERIO	TF (13mm, 1/2 taper point)	4-0	18	Most popular suture size used for soft tissue grafts and delicate tissues requiring taper point needle	
(112-7258)	-	PS-2 (19mm 3/8 premium reverse cut)	4-0	18	Most popular suture size with the most popular sized premium needle	
(112-7259)	CS0618 PREM	P-3 (13mm, 3/8 premium reverse cut)	4-0	18	Most popular suture size used for soft tissue grafts and delicate tissues requiring cutting point needle	
(112-7260)	-	P-3 (13mm, 3/8 premium reverse cut)	5-0	18	Small suture size and premium needle option for delicate tissues to reduce tissue trauma	
(112-7261)	-	P-3 (13mm, 3/8 premium 6-0 reverse cut)	6-0	18	Smaller suture size and premium needle option for delicate tissues and the least amount of tissue trauma	

**Cytoplast is a registered trademark of Osteogenics.

FS-2 Reverse Cutting 19mm 3/8 Circle	FS-3 Reverse Cutting 16mm 3/8 Circle	P-3 Premium Reverse Cutting 13mm 3/8 Circle	PS-2 Premium Reverse Cutting 19mm 3/8 Circle	TF Taper Point 13mm 1/2 Circle

Made in USA



BRAND PRODUCTS

Surgical Sutures

From A Name You Can Trust.

Maxima® Sutures by Henry Schein®

POLYLITH CARD
Helps reduce moisture content, an advantage with absorbable suture materials that degrade by hydrolysis.

NEEDLE PARK
Two needle parks to accommodate a wide range of needle types, sizes, and double-arm sutures. Needle park also protects needles from damage.

RACETRACK FEATURES AND BENEFITS

LARGER DIAMETER
Larger overall diameter to ensure less suture memory and minimize the number of suture rotations around track.

PRESS-PUSH NEEDLE LIFT
The tab, when pressed from behind, lifts for easy needle removal.



- Silicone-coated 300 series stainless needles give finer point geometry for superior penetration
- Precision needles are strong enough to resist bending, yet ductile enough to resist breaking
- Racetrack packaging helps reduce material memory as the suture is released from the track
- Large breadth of suture and needle combination for a range of surgeries
- Significant cost savings for your clinic compared to the leading brand

Precision Reverse Cutting Needles			PGA ABSORBABLE SUTURE					
			Material	Color	Length	5-0	4-0	3-0
3/8 circle 13mm	NP-3		Braided	Undyed	18" 45 cm	570-1048 12/box	570-1049 12/box	
3/8 circle 19mm	NPS-2		Braided	Undyed	18" 45 cm		570-1045 12/box	570-1046 12/box
		Violet		570-1399 12/box				
3/8 circle 19mm	NFS-2		Braided	Undyed	30" 75 cm	570-1397 12/box	570-1041 12/box	570-1047 12/box
		Violet						

Precision Reverse Cutting Needles			NYLON NON-ABSORBABLE SUTURE					
			Material	Color	Length	5-0	4-0	3-0
3/8 circle 13mm	NP-3		Monofilament	Black	18" 45 cm	570-1039 12/box	570-1040 12/box	
3/8 circle 19mm	NPS-2		Monofilament	Black	18" 45 cm		570-1042 12/box	570-1044 36/box
		Black		570-1043 12/box				
3/8 circle 24mm	NFS-1		Monofilament	Black	18" 45 cm			570-1050 12/box

Precision Reverse Cutting Needles			POLYPROPYLENE NON-ABSORBABLE SUTURE					
			Material	Color	Length	5-0	4-0	3-0
3/8 circle 13mm	NP-3		Monofilament	Blue	18" 45 cm	570-1591 12/box	570-1593 12/box	
3/8 circle 19mm	NPS-2		Monofilament	Blue	18" 45 cm		570-1592 12/box	570-1594 36/box
		Blue		570-1595 12/box				
3/8 circle 19mm	NFS-2		Monofilament	Blue	18" 45 cm			

**SUPERIOR POINT**

Needles are honed with a unique multiedged point to provide greater ease of penetration.

SUPERIOR ATTACHMENT

Needles are drilled for a seamless needle-suture juncture, providing a smoother transition with less tissue disruption.

SUPERIOR PACKAGING

Unique double-sterile barrier packaging ends costly waste.

**NON-ABSORBABLE SUTURES** 12/Box

MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Silk Black Braided	C-9 24mm 1/2 Circle	18"/45cm				632G* (100-6320)	
	C-43 15mm 1/2 Circle	18"/45cm				(100-7762)	
	C-3 13mm 3/8 Circle	18"/45cm		640G* (101-4148)	641G* (101-2590)		
	C-6 18mm 3/8 Circle	18"/45cm		682G* (101-3207)	683G* (100-6830)	1679H* (100-7772)	
	C-7 24mm 3/8 Circle	18"/45cm				684G* (101-2049)	
	C-22 16mm 3/8 Circle	18"/45cm	(101-2636)				
Nylon-Black Monofilament	C-1 11mm 3/8 Circle	18"/45cm	697G* (900-7482)				
	C-3 13mm 3/8 Circle	18"/45cm	1698G* (900-7473)	698G* (900-7470)	699G* (900-7466)		
	C-6 18mm 3/8 Circle	18"/45cm 10"/25cm		661G* (101-7137)	662G* (101-1209) (900-7464)		
	C-7 24mm 3/8 Circle	18"/45cm				663G* (101-3728)	
	C-8 26mm 3/8 Circle	18"/45cm					664G* (900-7476)
	C-17 12mm 3/8 Circle	10"/25cm	(900-7460)	(900-7459)	(900-7461)		
	C-22 16mm 3/8 Circle	18"/45cm	1665G* (101-6053)				
	PC-31 18mm 3/8 Circle	18"/45cm		1666G* (900-7472)	1667G* (900-7487)		
	DSCM-16 16mm 3/8	18"/45cm		(570-0458)	(570-0459)		
	DSM-19 18mm 3/8 Circle	18"/45cm				(570-0479)	
	DSM-24 24mm 3/8 Circle	18"/45cm				(570-0480)	
Polypropylene-Blue Monofilament	C-1 11mm 3/8 Circle	18"/45cm	8697G* (900-7483)				
	C-3 13mm 3/8 Circle	18"/45cm	8695G* (900-7479)	8698G* (900-7484)	8699G* (900-7486)		
	C-6 18mm 3/8 Circle	18"/45cm		8661G* (101-6409)	8683G* (101-7477)		
	C-7 24mm 3/8 Circle	18"/45cm				8684G* (101-7700)	
	C-22 16mm 3/8 Circle	18"/45cm	(101-4378)				
	DS-18 18mm 3/8 Circle	18"/45cm				8665G* (900-7481)	
	PC-31 18mm 3/8 Circle	18"/45cm			8682G* (900-7478)		
	DSM-16 16mm 3/8 Circle	18"/45cm		(570-0471)			
	DSM-19 18mm 3/8 Circle	18"/45cm		(570-0468)		(570-0467)	
	DSM-24 24mm 3/8 Circle	18"/45cm				(570-0470)	
	DSM-26 26mm 3/8 Circle	18"/45cm					(570-0473)
	DSCM-13 13mm 3/8	18"/45cm	(570-0472)				
	DSCM-16 16mm 3/8	18"/45cm		(570-0469)			
Polyester-Green Braided	D-14 17mm 1/2 Circle	18"/45cm		X951H* (100-7773)			

*Comparable to Ethicon™ material.




ABSORBABLE SUTURES 12/Box

MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Plain Gut	C-9 24mm 1/2 Circle	18"/45cm				612G* (100-6120)	
	C-3 13mm 3/8 Circle	18"/45cm		686G* (101-3192)			
	C-6 18mm 3/8 Circle	18"/45cm			(100-7771)		
	C-6 18mm 3/8 Circle	27"/70cm			H821H* (101-0668)	H822H* (101-2976)	
	DSCM-13 13mm 3/8 circle	18"/45cm	1916G* (900-7485)	1915G* (900-7471)			
Chromic Gut	C-3 13mm 3/8 Circle	18"/45cm		687G* (101-8824)			
	C-6 18mm 3/8 Circle	18"/45cm			1637G* (900-7637)	1638H* (900-7636)	
	C-6 18mm 3/8 Circle	27"/70cm			635H* (101-2787)	636H* (101-1582)	
	HR-26 26mm 1/2 Circle	30"/76cm				G122H* (570-0466)	
	DSM-13 13mm 3/8 Circle	18"/45cm			1654G* (570-0475)		
	DSM-19 19mm 3/8 Circle	18"/45cm			1637G* (570-0476)		
Polyglycolic Acid (PGA) Undyed-Braided	C-3 13mm 3/8 Circle	18"/45cm	J492G* (900-7468)	J493G* (900-4417)	J494G* (900-4418)		
	C-6 18mm 3/8 Circle	27"/70cm			J422H* (900-7463)	J423H* (900-7462)	
	PC-30 18mm 1/2 Circle	18"/45cm			J507G* (900-7477)		
	PC-31 18mm 3/8 Circle	18"/45cm			J496G* (900-4419)	J497G* (900-4421)	
	DSM-24 24mm 3/8 Circle	27"/70cm					(570-0482)
	HR-26 26mm 1/2 Circle	27"/70cm				(570-0481)	
Polyglycolic Acid (PGA) Undyed-Braided-Fast Absorbing	DSM-13 13mm 3/8 Circle	18"/45cm		VR493* (570-0477)	VR494* (570-0478)		
Polyglycolic Acid (PGA) Violet	C-6 18mm 3/8 Circle	18"/45cm			J392H* (900-4420)		
	DSM-19 19mm 3/8 Circle	18"/45cm				J393H* (570-0474)	
PGA PCL Undyed-Mono	DSM-13 13mm 3/8 Circle	18"/45cm		Y493G* (900-7467)	Y494G* (900-7465)		
	DSM-19 19mm 3/8 Circle	18"/45cm			Y496G* (900-7488)	Y497G* (900-7474)	
	DSM-13 13mm 3/8 Circle	18"/45cm	Y492G* (570-0465)				
	DSM-19 19mm 3/8 Circle	27"/70cm			(570-0460)	(570-0461)	
	DSM-24 24mm 3/8 Circle	27"/70cm				(570-0463)	
PGA PCL Violet-Mono	DSM-13 13mm 3/8 Circle	18"/45cm		Y463G* (570-0462)	(570-0464)		


* Comparable Ethicon™ material

NEEDLE CHART

REVERSE CUTTING NEEDLES



PRECISION REVERSE CUTTING NEEDLES



PRECISION CONVENTIONAL CUTTING NEEDLES



TAPER CUTTING NEEDLES





Self-Adherent Bandages

Non-Sterile, 5-yd Rolls

- Self-adherent: sticks to itself, not the skin
- Comfortable porous material: pain-free removal that will not damage skin

Latex

Tan, 1" ea
(570-0113)

Tan, 2" ea
(900-4322)

Tan, 3" ea
(900-4323)

Tan, 4" ea
(900-4324)

Tan, 6" 12 rolls/pkg
(900-4328)

Multi, 1" 30 rolls/pkg
(112-7366)

Multi, 2" 36 rolls/pkg
(900-4331)

Multi, 3" 24 rolls/pkg
(900-4332)

Non-Latex

Tan, 1" ea
(112-7355)

Tan, 2" ea
(112-7357)

Tan, 3" ea
(112-7358)

Tan, 4" ea
(112-7359)

Tan, 6" ea
(112-7362)

Multicolored, 1" 30 rolls/pkg
(112-7363)

Multicolored, 2" 30 rolls/pkg
(112-7364)

Multicolored, 3" 24 rolls/pkg
(112-7365)



Elastic Bandages with Clips

Non-Sterile, 4½-yd Rolls

- High cotton content to keep patients cool and comfortable
- Durable bandages can be laundered, sterilized, and reused
- Ideal for support and compression applications, such as treatment of varicosities, sprains, strains, and torn or injured ligaments
- Not made with natural rubber latex

2" 10 rolls/pkg
(900-4677)

3" 10 rolls/box
(900-4678)

4" 10 rolls/box
(900-4679)

6" 10 rolls/box
(900-4680)



Elastic Bandages with VELCRO®

Non-Sterile, 4½-yd Rolls

- Elastic bandage with Velcro® brand fasteners
- Provides even compression with strength, comfort, and breathability
- Premium honeycomb-knitted material
- Not made with natural rubber latex

2" 10/box
(900-4681)

3" 10/box
(900-4682)

4" 10/box
(900-4683)

6" 10/box
(900-4684)



Flexible Fabric Elastic Strips

These strips stretch, conform, and permit free movement. Will not stick to wound. Not made with natural rubber latex.

¾" x 3" 100/box
(112-6132)

1" x 3" 100/box
(112-6133)

Knuckle, 1½" x 3" 100/box
(112-6144)



Sheer Plastic Bandage Strips

Sterile

Individually wrapped adhesive bandages for the protection of minor cuts, punctures, and abrasions. Highly absorbent pad will not stick to wound. Air vents let air in and keep moisture out.

¾" x 3" 100/box
(112-6138)

1" x 3" 100/box
(112-6142)

¾" Spots 100/box
(112-6143)

Professional Tray Pack, ¾" x 3" 1500/case
(112-6134)

Professional Tray Pack, 1" x 3" 1500/case
(112-6135)



Crayon Adhesive Bandages

Sterile

- Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
- Pad will not stick to wounds
- Permits wound to drain
- Not made with natural rubber latex

¾" x 3" 100/box
(112-6995)



Animal Adhesive Bandages

Sterile

- Exclusive design for Henry Schein
- Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
- Nonstick, highly absorbent pad will not stick to wound
- Air vents let air in and keep moisture out
- Not made with natural rubber latex

Spots, ¾" 100/box
(900-4506)

¾" x 3" 100/box
(900-4507)



All-Gauze Sponges

Non-Sterile
Soft, pliable, and absorbent. Available in 3 sizes in both 8-ply and 12-ply styles. Non-sterile.

- Precision sizes and folding: consistency within functions
- Soft and pliable: nonabrasive touch
- Highly absorbent: soaks up unwanted fluid
- Can be used for general cleaning, debriding, and prepping

8 Ply, 2" x 2" (570-2657)	200/pkg
12 Ply, 2" x 2" (570-2658)	200/pkg
8 Ply, 3" x 3" (570-2650)	200/pkg
12 Ply, 3" x 3" (570-2648)	200/pkg
8 Ply, 4" x 4" (570-2656)	200/pkg
12 Ply, 4" x 4" (570-2630)	200/pkg



Premium All-Gauze Sponges

Manufactured from 100% USP type-VII cotton gauze that is bleached and processed to established standards

- Soft and pliable
- Protects as a wound covering as it provides maximum absorbency

Non-Sterile 8 Ply, 2" x 2" (570-2647)	200/pkg
12 Ply, 2" x 2" (570-2628)	200/pkg
12 Ply, 3" x 3" (570-2640)	200/pkg
8 Ply, 4" x 4" (570-2646)	200/pkg
12 Ply, 4" x 4" (570-2645)	200/pkg
Sterile 8 Ply, 2" x 2" (570-2642)	100/pkg
12 Ply, 3" x 3" (570-2631)	80/pkg
8 Ply, 4" x 4" (570-2644)	50/pkg
12 Ply, 4" x 4" (570-2637)	50/pkg

**Essentials™
Non-Woven Sponges**

Non-sterile, poly/rayon blend, non-woven sponges for general use. Soft and absorbent, low-linting 4-ply material.

2" x 2" (112-6757)	5000/case Contains: 25 sleeves per case.
4" x 4" (112-6758)	2000/case



Exodontia Sponges

100% cotton gauze sponges with maximal absorbency cotton inner liner

- Manufactured under strictest standards: superior and consistent product quality
- Available in 3 sizes: convenient for different uses
- 8-ply
- Not made with natural rubber latex

Non-Sterile—Filled 2" x 2" (570-2636)	5000/case
3" x 3" (570-2602)	4000/case
4" x 4" (570-2651)	2000/case
Sterile—Filled 2" x 2" (570-2653)	200/pkg



Premium Exodontia Sponges

Manufactured from all USP 100% cotton. Soft, absorbent, and pure white. Highest quality. 2" x 2", 8-ply.

Non-Sterile—Filled (570-2635)	5000/case
Sterile—Filled (570-2652)	200/pkg



**CRITERION™
Non-Woven Sponges**

Non-Sterile, 4 Ply

- High-quality, finely constructed, rayon/poly-blend that is soft and gentle to the touch
- Multiuse sponge that can be used for all applications of wound-care treatments
- Thicker and more absorbent than standard sponges
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2" (104-2930)	200/pkg
3" x 3" (104-2410)	200/pkg
4" x 4" (104-2411)	200/pkg



**CRITERION™
Non-Woven Sponges**

Sterile 2s, 4 Ply

- High-quality, finely constructed, rayon/poly-blend that is soft and gentle to the touch
- Multiuse sponge that can be used for all applications of wound-care treatments
- Thicker and more absorbent than standard sponges
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2" (104-7897)	50/pkg
3" x 3" (104-4132)	50/pkg
4" x 4" (104-8665)	50/pkg



**Premium
Non-Woven Sponges**

Non-Sterile, 4 Ply

- General-use, non-woven rayon/polyester-blend material
- Quarter-fold design
- Less linting: minimizes lint contamination
- Soft and pliable: prevents irritation to mucosa
- Greater absorbency: reduces number of sponges required
- Closed weave
- Not made with natural rubber latex

2" x 2" (101-9278)	200/pkg
3" x 3" (101-8748)	200/pkg
4" x 4" (101-4336)	200/pkg



**Non-Woven
Cotton Sponges**

Non-sterile, 4 Ply

- Highly absorbent 100% cotton
- 4-ply gauze sponges: thick, soft gauze
- Water bleach: no discoloration after sterilization
- See-through hygienic package imprinted with a sterilization indicator: easy identification
- Tear notch: convenient bag opening

2" x 2" (100-1090)	100/pkg
3" x 3" (100-3835)	100/pkg
4" x 4" (100-4420)	100/pkg



**Syngauze 50™
Non-sterile, 4 Ply**

- General-use nonwoven material with a 70/30% blend of rayon and polyester. 47–50% more absorbent than 12-ply gauze
- Absorbs 10 times its weight: greater absorbency reduces the number of sponges required
- Soft and pliable: prevents irritation to mucosa; allows shaping for a better fit
- Sterilizable: lint-free mesh minimizes lint contamination
- Nonraveling: less mess
- Not made with natural rubber latex

2" x 2" (100-2125)	200/pkg
3" x 3" (101-9668)	200/pkg
4" x 4" (101-4583)	200/pkg
(107-9502)	2000/case



Type VII Gauze

- Sterile
- Not made with natural rubber latex
- 4" x 4"

8 Ply (570-2633)	10/pkg
12 Ply (570-2655)	10/pkg
16 Ply (570-2638)	10/pkg
12 Ply, X-ray Detectable (570-2632)	10/pkg
16 Ply, X-ray Detectable (570-2634)	10/pkg



Compare to Steri-Pad Gauze Pads

Sterile
Can be used for general cleaning, debriding, and prepping. Sterile 1's in peel-back packages. 12-ply.

- Highly absorbent
- 100% woven cotton
- Box is perforated at bottom for easier gauze removal

2" x 2" (570-2659)	100/pkg
3" x 3" (570-2629)	100/pkg
4" x 4" (570-2649)	100/pkg



Combine ABD Pads

- Highly absorbent multilayer ABD (abdominal) dressing
- Soft, nonwoven outer facing wicks moisture to the cellulose core, which quickly absorbs and disperses fluid
- Moisture-resistant barrier prevents fluid strikethrough
- All edges are sealed to prevent linting and leakage
- Not made with natural rubber latex
- Sterile

5" x 9" (570-1470)	25/box
8" x 7½" (570-1472)	20/box
8" x 10" (570-1471)	20/box



Conforming Stretch Gauze Rolls

- 4.1-yd Rolls
- Economical knitted bandage in various widths
 - Made from high-quality materials: nylon, polyester, and cotton
 - Provides light compression as it cushions and protects
 - Offers softness, breathability, and alleviation of irritation
 - Low lint content reduces contamination to wound site
 - Improved packaging is resealable, durable, easy to store, and remains clean
 - Not made with natural rubber latex

Non-Sterile

1" (104-2704)	12 rolls/bag
2" (104-8825)	12 rolls/bag
3" (104-8997)	12 rolls/bag
4" (104-2728)	12 rolls/bag
6" (104-9424)	6 rolls/bag

Sterile

2" (104-2739)	12 rolls/bag
3" (104-9724)	12 rolls/bag
4" (104-6708)	12 rolls/bag
6" (104-4293)	6 rolls/bag



Woven Gauze Unna's Boot

10-yd Roll

- Cotton gauze formulated with a zinc-oxide emulsion
- Used for treating venous leg ulcers, edema, and dermatitis
- Soothes inflammation and rashes
- Not made with natural rubber latex

3" (100-3444)	ea
4" (100-3256)	ea
3" with Calamine (900-7455)	ea
4" with Calamine (900-7456)	ea



Premium Hypoallergenic Surgical Cloth Tape

10-yd Rolls

- Reliable adhesion
- Highly porous
- Thin, flexible backing
- Gentle removal
- Not made with natural rubber latex
- Non-sterile

½" (101-3025)	24/box
1" (101-2323)	12/box
2" (101-5979)	6/box



Premium Hypoallergenic Transparent Surgical Tape

10-yd Rolls

- Excellent adhesion
- Porous, with flexible backing
- Tears easily in both directions
- Not made with natural rubber latex
- Non-sterile

½" (101-6655)	24/box
1" (101-4645)	12/box
2" (101-9226)	6/box



Premium Hypoallergenic Surgical Paper Tape

10-yd Rolls

- Excellent adhesion
- Pinked, nonfraying edges
- Strong woven backing
- Gentle removal
- Not made with natural rubber latex
- Non-sterile

½" (101-4305)	24/box
1" (101-9673)	12/box
2" (101-1332)	6/box



Pneumatic Walking Boot

- Intuitive sole to promote natural gait
- Quick-press release valve
- High-grade polymer material
- Shock-absorption feature in sole
- Air bladder that covers larger surface area
- Rivets on struts for better fitting
- Manufactured with medical-grade materials
- Comfort pads to minimize pressure points
- Available in Tall and Ankle heights

Tall Walking Boot..... ea

Specify:

X-small, Men's Size 2-4; Women's Size 3½-5½ (570-1007)

Small, Men's Size 4½-7; Women's Size 6-8½ (570-1008)

Medium, Men's Size 7½-10; Women's Size 9-11½ (570-1009)

Large, Men's Size 10½-13; Women's Size 12-14½ (570-1010)

X-large, Men's Size 13½+; Women's Size 15+ (570-1006)

Ankle Walking Boot..... ea

Specify:

X-small, Men's Size 2-4; Women's Size 3½-5½ (570-1011)

Small, Men's Size 4½-7; Women's Size 6-8½ (570-1012)

Medium, Men's Size 7½-10; Women's Size 9-11½ (570-1013)

Large, Men's Size 10½-13; Women's Size 12-14½ (570-1014)

X-large, Men's Size 13½+; Women's Size 15+ (570-1015)



Ankle Brace, Stirrup

- Provides medial and lateral stability while allowing plantar and dorsiflexion
- Air chambers decrease edema and help break the pain cycle
- Anatomically contoured semirigid shells conform closely to the ankle with minimal bulk in the shoe
- Adjustable heel strap ensures comfortable and functional fit
- Supports injured and unstable ankles

Universal Fit

(570-1016)..... ea



Knee Brace, Hinged

- Removable, dual-axis polycentric hinges mimic the true anatomy of the knee and provide both medial and lateral support
- Wraparound configuration and strap provide warmth and equal compression
- Front closure provides easy application
- Condyle pads for additional comfort
- Brushed neoprene and seamed edges increase patient comfort and ease of washing

Knee Brace.....ea

Specify:

X-Small, Thigh: 16½"W, Calf: 14"W (570-1029)

Small, Thigh: 17½"W, Calf: 15"W (570-1030)

Medium, Thigh: 18½"W, Calf: 16"W (570-1031)

Large, Thigh: 16½"W, Calf: 14"W (570-1032)

X-Large, Thigh: 21"W, Calf: 18"W (570-1033)

2X-Large, Thigh: 23"W, Calf: 20"W (570-1034)

3X-Large, Thigh: 25"W, Calf: 22"W (570-1035)



Surgical Shoe, Round Toe

- Upper section is made from a soft, breathable black nylon mesh with hook-and-loop closures for uniform compression
- Reinforced sole is flared to provide better medical and lateral stability
- Reinforced, padded heel and soft tricot lining add comfort and support
- Semirigid with padded collar
- Primarily used postsurgically, post-trauma, for foot wounds including burns, and for diabetic and arthritic patients

Surgical Shoe..... ea

Specify:

Women's

X-Small, Size 3½ (570-1025)

Small, Size 4-6 (570-1026)

Medium, Size 6½-8½ (570-1027)

Large, Size 9+ (570-1028)

Men's

Small, Size 6-8 (570-1021)

Medium, Size 8½-10 (570-1022)

Large, Size 10½-12 (570-1023)

X-large, Size 12½+ (570-1024)



Surgical Shoe, Square Toe

- A medical shoe used to protect the foot and toes after an injury or surgery
- Upper section is made from a soft, breathable black nylon mesh with hook-and-loop closures for uniform compression
- Reinforced, padded heel and soft tricot foam lining for additional comfort and support
- Black

Surgical Shoe..... ea

Specify:

Women's

Small (570-1535)

Medium (570-1536)

Large (570-1537)

Men's

Small (570-1538)

Medium (570-1539)

Large (570-1540)

X-Large (570-1541)



Night Splint, Plantar

- 3 padded straps with buckles to ensure immobilization and patient comfort
- Bilateral flexion straps provide 10°-90° dorsiflexion for gentle stretching
- Crafted, 3-dimensional Lycra®-lined soft, cool interior for enhanced patient comfort
- Optional flex wedge included

Night Splint.....ea

Specify:

Small, Men's Size Up to 6½;

Women's Size Up to 8 (570-1017)

Medium, Men's Size 7-9½; Women's Size 8½-11 (570-1018)

Large, Men's Size 10-13;

Women's Size 11½-14½ (570-1019)

X-Large, Men's Size 13½+; Women's Size 15+ (570-1020)



Aluminum Crutches

- Push-pin adjustment makes it easy to adjust crutch height in 1" increments
- Nonskid, jumbo-size, vinyl-contoured tips provide excellent traction
- Easy wing-nut adjustments for hand grips
- 300-lb weight capacity on aluminum crutch

Crutches..... pair

Specify:

Youth, Patient Height: 4'6"-5'2", Crutch Height: 37"-46" (570-1466)

Adult, Patient Height: 5'2"-5'10", Crutch Height: 45"-53" (570-1464)

Tall Adult, Patient Height: 5'10", Crutch Height: 53"-61" (570-1465)



Synthetic Cast Padding

Non-Sterile, 4-yd Rolls

- Made of synthetic fibers
- Provides a safe environment under a hard cast
- Comfortable feel for the patient
- 100% polyester and water-resistant; not waterproof
- Retains consistency under wet or dry conditions
- Conforms to all body contours
- Works with synthetic or plaster casting material
- Not made with natural rubber latex

2"		
(900-4997)	24/bag
3"		
(900-4998)	12/bag
4"		
(900-4999)	12/bag
6"		
(900-7001)6/bag



Synthetic Stockinette

Non-Sterile, 25-yd Roll

- Essential component of a high-quality cast process
- Underlayment for all standard casting applications
- Suitable for use with all casting materials, including plaster and synthetic
- Not made with natural rubber latex

2"		
(112-5513)	ea
3"		
(112-5514)	ea
4"		
(112-5515)	ea
6"		
(112-5516)	ea



Orthopedic Splint, Fiberglass

Non-Sterile, 15-ft Rolls

A device used for supporting or immobilizing a limb. For use in postsurgical procedures and acute trauma treatment of fractures, sprains, soft tissue injuries, and more.

- Water repellent and hydrophilic
- Soft feeling
- Antibiotic and deodorized
- Clean and neat external appearance
- Good strength/molding
- Not made with natural rubber latex

2"		
(112-7054)	ea
3"		
(112-7060)	ea
4"		
(112-7061)	ea
5"		
(112-7062)	ea



Orthopedic Casting Tape, Polyester

Non-Sterile, 4-yd Rolls

- Composed of a specially knitted, multidirectional stretch polyester fabric impregnated with polyurethane resin
- Designed for use in the production of definitive orthopedic casts and splints
- Suitable for use in specialized prosthetic and orthotic devices
- Not made with natural rubber latex

2"			
			10/box
Specify:			
Pink	(112-7037)	Black	(112-7038)
3"			10/box
Specify:			
White	(112-7039)	Red	(112-7044)
Green	(112-7041)	Pink	(112-7045)
Blue	(112-7042)	Black	(112-7046)
Purple	(112-7043)		
4"			10/box
Specify:			
White	(112-7047)	Purple	(112-7050)
Green	(112-7048)	Red	(112-7051)
Blue	(112-7049)	Pink	(112-7052)



Fiberglass Casting Tape

Non-Sterile, 4-yd Rolls

- Fiberglass polyurethane formulation
- Radiotranslucent
- No special glove requirement
- Highly conformable
- Impact resistant and moisture resistant
- Easy to unroll
- Not made with natural rubber latex

2"			
			10 rolls/box
Specify:			
White	(900-4635)	Red	(900-4646)
Blue	(900-4638)	Purple	(900-4648)
Black	(900-4641)	Green	(900-4650)
Pink	(900-4644)		
3"			10 rolls/box
Specify:			
White	(900-4636)	Red	(900-4647)
Blue	(900-4639)	Purple	(900-4649)
Black	(900-4642)	Green	(900-4651)
Pink	(900-4645)		
4"			10 rolls/box
Specify:			
White	(900-4637)	Pink	(900-7002)
Blue	(900-4640)	Purple	(900-7003)
Black	(900-4643)		



Instant Cold Compresses

Disposable, single-use cold pack that provides immediate cold therapy relief without refrigeration. Ready to activate for quick use in emergencies or first aid situations. Compact, lightweight, and flexible for comfortable application. Ideal for temporary relief of pain and swelling caused by sprains, contusions, minor burns, toothaches, and insect bites. Non-sterile.

5" x 5 1/2"		
(112-6147)	50/case
5" x 7"		
(112-6146)	48/case
6" x 9" (15.24 cm x 22.86 cm)		
(112-6145)	24/case



Instant Heat Compresses

Single-use, disposable hot pack that delivers soothing warmth for instant relief—no heating equipment required. Flexible and conforming for comfort, with a ready-to-activate design that heats rapidly. Ideal for temporary relief of frostbite, arthritis discomfort, or aches in cold or wet weather. Portable, compact, and lightweight for on-the-go use. Non-sterile.

6" x 8 1/4" (15.24 cm x 22.86 cm)		
(499-0813)	24/case



Reusable Cold or Hot Gel Packs

Versatile gel pack designed for both cold and hot therapy. Simply freeze for cold therapy or microwave for heat. Remains soft and flexible—even when frozen—for added comfort. Reusable and easy to clean for long-lasting value. Provides temporary relief of arthritis, joint pain, muscle strain, and stiffness. Non-toxic and not made with natural rubber latex. Non-sterile.

4" x 6 1/2" (10.16 cm x 16.51 cm)		
(570-2140)	ea
5 1/4" x 10" (13.34 cm x 25.4 cm)		
(570-2139)	24/case



OneStep+ Pro Hb Analyzer

This point-of-care hemoglobin analyzer is designed to provide quantitative, lab-quality hemoglobin results in only 3 seconds.

- Reagent-free microcuvettes
- Microcuvette shelf life of up to 2 years, even after the canister has been opened
- Cleared for use on patients 6 months old or older
- Measurement range 5 to 25.6 g/dL
- Small, 15-µL sample volume
- 3¼" color LCD display
- CLIA waived

Starter Kit (570-1314)..... ea
 Contains: 1 analyzer, 1 box of 100 microcuvettes, 1 optical system check, 4 AA batteries & 1 USB cable.

Analyzer (570-1311)..... ea
 Contains: 1 analyzer, 4 AA batteries & 1 USB cable.

Microcuvettes (570-1312)..... 200/box
 Contains: 4 canisters of 50 cuvettes.

Controls Set (570-1313)..... ea
 Contains: 3 levels of controls: Low, Mid & High.

Optical System Check (570-1322)..... ea
 Contains: ID chip & microcuvette.



True Metrix Pro Meter

- No coding
- Tiny, 0.5-µL sample size
- Stores 500 results
- As fast as 4 seconds
- BG range: 20-600
- Data management
- Download capabilities
- Glucose control detection
- Strip-release button

#RE4099P-43 (570-3418)..... ea
 Contains: 1 True Metrix Pro meter, installed 3-volt battery, instructions for use, logbook & compact carrying case.



CardioChek Plus Analyzer

Fast, Accurate, Efficient Lipid and Glucose Results at the Point of Care

Recognizing the strong connection between heart attack, stroke, and diabetes, the CardioChek Plus Analyzer provides rapid, on-site blood-chemistry values and measurements to help identify at-risk individuals.

- Fast: On-site lipid and glucose results in as little as 90 seconds; simple, quality-control testing
- Portable: compact and lightweight, easy to transport and store; handheld and battery-operated; and room-temperature test strips
- Patient-friendly: fingerstick with sample size of 1.1 µL to 40 µL
- Wireless communication: easy data transfer
- Measures total and HDL cholesterol, triglycerides, and glucose
- Calculates LDL and non-HDL cholesterol as well as TC/HDL and LDL/HDL ratios
- Certified accuracy
- CLIA waived
- Meets NCEP guidelines for accuracy and precision
- FCC- and ISO-certified wireless communications
- CRMLN-certified for CHOL and HDL tests
- FDA cleared

#PTS-2700-HS (570-0411)..... ea



True Metrix Pro Strips

Advanced technologies—the meter, a complex algorithm, chemistry, and electrodes—on the test strip work together as part of the True Metrix Pro system to produce accurate results. Featuring Triple Sense Technology, the system provides proven accuracy and confidence in results.

- No-coding technology eliminates the need for coding of meter
- Code is assigned during manufacturing process
- Code is printed onto contacts on end of test strip
- Meter reads code on insertion

#R3099P-450, Strips (570-0327)..... 50/box

Controls #R5H01-1, Level 1 (125-4034)..... ea

#R5H01-2, Level 2 (125-4035)..... ea

#R5H01-3, Level 3 (125-4036)..... ea



PTS Panels Starter Kit

#PTS-4060-HS (570-1036)..... ea
 Contains: 1 analyzer, 60 lipid/glucose tests, multichemistry controls & 40-µL cap tubes.



PTS Panels eGlucose Test Strips

#PTS-2713-HS (570-0415)..... 50/box
 Contains: 50 test strips.



PTS Panels Lipid/Glucose Smart Bundle 2-Pack

#PTS-2729-HS (570-0412)..... 2/pkg
 Contains: 1 of each: 15-count lipid test-strip vials & 15-count vial of glucose test strips.



PTS Panels Lipid Test Strips

#PTS-1710-HS (570-0414)..... 15/box
 Contains: 15 test strips.



PTS Panels 40-µL Capillary Tubes

#PTS-2866-HS (570-0413)..... 16/bag



Henry Schein A1CNow[®]+ makes it convenient for diabetes patients to have A1C tested in physicians offices.



(570-3428)

A1CNow[®]+

- ✓ Improve patients' in-office experience
- ✓ Increase efficiency
- ✓ Help drive better health outcomes for diabetes patients



(137-8628)

Ensure Accurate A1C Results with

A1CNow[®]+ Controls[™]

ONESTEP+ DRUG OF ABUSE RAPID URINE TEST CUP

Features & Benefits

- CLIA waived
- Designed to detect a wide range of substances: AMP, BAR, BUP, BZO, COC, MDMA, EDDP, MTD, MET, MOP, OPI, OXY, PCP, PPX, TCA, and THC*
- Simple to use, no need to tilt cup to activate test
- Self-contained cup is ideal for sending preliminary positive specimens to the lab for confirmation
- Lid click feature ensures the cup is sealed properly, maintaining integrity of the test
- Results available in 5 minutes
- Temperature strip helps to monitor specimen tampering
- Adulteration strip offers initial screening and assists in monitoring any tampering of the specimen

*AMP-Amphetamine, BAR-Barbiturate, BUP-Buprenorphine, BZO-Benzodiazepines, COC-Cocaine, MDMA-Ecstasy, EDDP-Methadone Metabolite, MET-Methamphetamine, MTD-Methadone, MOP-Morphine, OPI-Opiates, OXY-Oxycodone, PCP-Phencyclidine, PPX-Propoxyphene, TCA-Tricyclic Antidepressants, THC-Marijuana



Henry Schein[®] OneStep+ Drug of Abuse Rapid Urine Test Cup 25/Box

Kit Contents: (25) OneStep+ Drug of Abuse Rapid Urine Test Cups, (2) Specimen Validity/Adulteration Strips (if applicable), (30) Disposable Gloves, (25) Security Seal Labels, (1) Instructions For use, (1) Procedure Card

Test Cup, 12 Panel—With Adulteration Strip..... (570-3927)
AMP 1000, BAR300, BZO300, COC300, MDMA500, MET1000, MOP300, MTD300, OXY100, PCP25, TCA1000, THC50 (CR, SG, pH)

Test Cup, 14 Panel..... (570-3928)
AMP 1000, BAR300, BUP10, BZO300, COC300, EDDP300, MDMA500, MET1000, MOP300, MTD300, OXY100, PCP25, TCA1000, THC50

Order today at HENRYSCHIEIN.COM/MEDICAL

Urispec™ Plus Urinalysis Reader & Urispec™ 11-Way Urine Strip Promotion



Provides dependable urine status results to detect early stages of many diseases, such as diabetes, kidney disease, and urinary tract infections. Instrument-read results have long proven to be advantageous for both busy health care professionals and patients. Comprehensive interface options and the optimized printouts minimize risk associated with manual transcriptions. Reliable results can be obtained virtually immediately at the point of care.

- CLIA-waived
- Readings eliminate subjectivity of visual color interpretation, which makes urinalysis simpler and more reliable
- Features no-touch operation: After strip placement, measurement starts automatically, which makes urinalysis faster and more hygienic
- Easy switch between different test strips: Urispec Plus is for use with Urispec 11-Way, 5-Way Plus, 4-Way, and GP+A urine test strips

(570-1309)..... ea
Contains: 1 Urispec™ Plus Urinalysis Reader & 12 bottles of Urispec™ 11-Way Urine Strips.

Small Reader Paper Rolls



Can be used on small analyzers like the Urispec™ Plus and other manufacturers' units, including the NanoEnTek FRENDS™. 57 mm x 32 ft

(570-5469)..... 6 rolls/pkg

Urispec™ 4-Way Urine Reagent Strips



- Protein, nitrite, glucose, and leukocytes
- Ideal for OB/GYN
- Infection, diabetes, and kidney function
- Dip and read in 60 seconds
- Strips for visual read or for use with the Urispec Plus Urinalysis Reader

(107-7157).....100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*

Urispec™ 5-Way Plus Urine Reagent Test Strips



- Tests for glucose, blood, protein, nitrite, leukocytes, and urobilinogen
- Strips for visual read or for use with the Henry Schein Urispec Plus Urinalysis Reader only

(900-7508)..... 100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*

Urispec™ GP+A Urine Reagent Test Strips



- Test for glucose and protein, with an ascorbic acid test for accuracy.
- Ideal for metabolic, kidney, and diabetic monitoring
- Dip and read in 60 seconds
- Strips for visual read or for use with Urispec Plus Urinalysis Reader

(100-6567)..... 100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*

Compare to Hydratrend Urispec™ 11-Way Urine Reagent Test Strips



- Unique "read all" feature at 60 seconds
- Similar to Multistix® 10 SG Plus Ascorbic Acid Test and all brands of 10 SG products
- Glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes
- Reader available
- Ascorbic acid pad to ensure uncompromised results
- Room-temperature storage
- Strips for visual read or for use with Urispec Plus Urinalysis Reader

(102-2285)..... 100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*

Urispec™ Pro 2 GP Urine Strips



- Strips for Visual Read
- Tests for glucose and protein.
- Made in the USA

(900-4996)..... 100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*

Urispec™ 5 OB Urine Strips



- Strips for Visual Read
- Tests for glucose, protein, blood, nitrites, and leukocytes.
- Made in the USA

(900-4995)..... 100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*

Urispec™ 10SG Urine Strips



- Strips for Visual Read
- Tests for glucose, bilirubin, ketone, urobilinogen, specific gravity, blood, pH, protein, nitrites, and leukocytes.
- Made in the USA

(900-4994)..... 100/btl

CPT™ Codes: 81000*/81001*/81002*/81003*



Urispec™ Urinalysis Controls

For Most Brands of Routine Urine Dip Sticks and Pregnancy Tests

A bilevel product packaged in plastic tubes and dropper bottles. Controls are liquid, ready-to-use, and are designed to validate urinalysis reagent strips. Packaged in plastic dip tubes and dropper bottles, they are convenient and easy to use. Long open- and closed-vial stability at 2°C–8°C for 2 years unopened, or 90 days once the Dip Control is opened, and 180 days once the Drop Control is opened. This allows the user ample time to utilize the product before expiration, resulting in cost savings.

- Values for all pregnancy test brands
- Assays for common brands of urine dip for routine urinalysis (2 parameter-10SG)
- Includes values for microalbumin, creatinine, and hCG
- Two distinct levels (normal and abnormal)
- Closed-vial dating up to 2 years with refrigerator storage
- 30-day stability at room temperature

3 Bilevel Sets, 15 mL/Tube (6 Vials)
Dip stick in tube, with open-vial stability of 90 days at 2°C–8°C.

(112-6841)..... ea

2 Bilevel Sets, 25 mL/Dropper Bottle (4 Bottles)
Open-vial stability of 180 days at 2°C–8°C.

(112-6842)..... ea

NEW!

Onestep+ Drug of Abuse Rapid Urine Test Cup Kit



- CLIA waived
- Offers a fast and convenient way to detect multiple drugs and their metabolites in urine
- Lid features a click mechanism that ensures the cup is sealed properly, maintaining the integrity of the test
- Temperature strip is included to help detect specimen tampering
- Adulteration strip provides initial screening and assists in detecting any tampering attempts
- Urine specimen
- Fill cup is 130 mL
- Results available in 5 minutes
- Non-sterile

12-Panel Kit, with Adulteration Strip
AMP1000, BAR300, BZO300, COC300, MDMA500, MET1000, MOP300, MTD300, OXY100, PCP2 5, TCA1000, and THC50+(CR/SG/pH)
(570-3927)..... ea

Contains: 25 OneStep+ Drug of Abuse Rapid Urine Test Cups, 2 specimen validity/adulteration strips, 30 disposable gloves, 25 security-seal labels, 1 set of instructions for use & 1 procedure card.

14-Panel Kit
AMP1000, BAR300, BUP10, BZO300, COC300, EDDP300, MDMA500, MET1000, MOP300, MTD3 00, OXY100, PCP25, TCA1000, and THC50
(570-3928).....

Contains: 25 OneStep+ Drug of Abuse Rapid Urine Test Cups, 30 disposable gloves, 25 security-seal labels, 1 set of instructions for use & 1 procedure card.



OneStep+ iFOBT Occult Blood Monoclonal Antibody Kits

Innovative dry-base collection method that utilizes 1 sampling card for the collection of 2 fecal specimens, and requires only 1 immunochemical specimen preparation tube/test cassette system for analysis.

- Meets the American Cancer Society and the US Task Force for Colorectal Cancer (CRC) Foundation recommendation of multiple-day testing
- The sampling tab technology replaces unpleasant liquid buffer tubes and eliminates human error due to patient oversampling
- No dietary or drug restrictions
- 30-day stability/24-month shelf life from date of manufacture
- CLIA waived
- Made in the USA

25-Test Complete Kit

(900-4413).....25 tests/kit
 Contains: 25 test cassettes, 25 buffer tubes & 25 patient return mailers.

Cassettes

(900-7458).....30/box
 Contains: 30 test cassettes & 30 buffer tubes.

20-Patient Mailer Kits

(900-4416).....20/pkg
 Contains: 20 postal-approved envelopes, each containing 1 sampling card slide, 2 applicators & patient instructions.

25-Test Complete Kit

CPT[®] Codes: 82274QW[†]/G0328QW[†]



OneStep Pro+ FIT Tests

This fecal immunochemical test (FIT) is a simple strip procedure that provides qualitative results in 5 minutes. It utilizes colloidal gold to enhance test-line intensity and clarity and requires no dietary or medicinal restrictions, while ensuring quick and accurate results.

Personal Use Kit

(570-0607).....20/box
 Contains: 20 sampling bottles & 20 mailing envelopes.

Complete Kit

(570-0609).....50/box
 Contains: 50 test strips, 50 sampling bottles & 50 collection papers.

Strips Kit

(570-0610).....50/box
 Contains: 50 test strips.



OneStep+ Pregnancy Cassette Tests

- Urine test
- CLIA waived
- 25 mIU/mL detection level
- Sensitivity: 99%
- Room-temperature storage

25-Patient Kit

(900-4073).....ea
 Contains: 25 pouched tests per kit & 1 package insert.

100-Patient Kit

(900-4074).....ea
 Contains: 100 pouched tests per kit & 1 package insert.

CPT[®] Codes: Urine 81025[†]/Serum 84703[†]



OneStep+ Pregnancy Combo Cassette Tests

- Urine or serum test
- CLIA waived for urine only
- 25 mIU/mL detection level
- Sensitivity: 99%
- Room-temperature storage

25-Patient Kit

(900-4075).....25/box
 Contains: 25 pouched tests per kit & 1 package insert.

CPT[®] Codes: 81025QW[†] (Urine); 84703[†] (Serum)



OneStep+ Pregnancy Strip Tests

Rapid, chromatographic immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine to aid in the early detection of pregnancy.

- Urine test
- CLIA waived
- 25 mIU/mL detection level
- Sensitivity: 99%
- Room-temperature storage
- Read time 3 minutes or less

(900-4072).....25/box
 Contains: 25 pouched tests per kit & 1 package insert.

CPT[®] Code: 81025QW[†] (Urine)



One Step+ ER Tests

A rapid, convenient, and noninvasive qualitative method for detecting occult blood in the stool. It is intended for professional use as an aid in the diagnosis of asymptomatic gastrointestinal conditions that may manifest themselves by the presence of occult blood in the stool.

- Patented grid design minimizes the number of unreadable slides due to patient oversampling
- Easier to read borderline specimens
- Provides an intense blue reaction, resulting in better readability than traditional guaiac test kits
- 21-day on-slide stability
- 3-year shelf life with room temperature storage
- Familiar patient collection
- Made in the USA

OneStep+ ER Triple Slides Kit

(112-7141).....50/box
 Contains: 3 x 10-mL ER Developer and 50 mailing envelopes containing triple slide, 3 applicators & patient instructions.

OneStep+ ER Single Slides Kit

(112-7140).....100/box

OneStep+ ER Developer

(112-7142).....ea



Specimen Containers

These 4-oz specimen containers have blue, tamper-evident, leak-resistant tops and are suitable for routine urine or specimen collection and transport. Container includes patient information label.

Sterile, with Tab, Unwrapped, 4 oz


(900-7433).....100/case

Sterile, Individually Wrapped, 4 oz

(900-7431).....100/case

Non-Sterile, Bulk, 4 oz

(900-7432).....500/case



SEASONAL RAPID TESTING SOLUTIONS

HENRY SCHEIN®
BRAND PRODUCTS



OneStep™+ Ultra Flu Test

- CLIA waived
- Results in 10 minutes
- Sensitivity A/B 90.3%/88%
- Specificity A/B 96.7%/99.2%
- For use with nasal and nasopharyngeal swab specimens

25/box.....(570-1886)

Contains: 25 test sticks, 25 sterile nasal swabs, 27 extraction buffer vials, 1 influenza A control swab & 1 influenza B control swab.

CPT® Codes: Influenza A, 87804QW†/Influenza B, 87804QW-59



OneStep™+ Ultra Mono Test

- CLIA-waived for whole blood
- 3-step process, dipstick test
- Results in 5 minutes
- Whole blood, serum, or plasma sample type
- Moderate complexity for serum and plasma
- Performance [reconciled with EBV serology]: 100% sensitivity; 95.9% specificity

25/box.....(570-1728)

Contains: 25 each of: test sticks in a container, test tubes, transfer pipettes, 50-µL capillary pipettes; and 1 each of: diluent, mono positive control, mono negative control, workstation & package insert.

CPT® Code: 86308QW† (whole blood); 86308† (serum and plasma)



OneStep™ Pro+ Strep A Cassette Kit

- CLIA waived
- Easy-to-use flip design; all-in-one extraction, incubation, and test
- Results in 5–10 minutes
- Sensitivity of 96.2% and specificity of 98.7%
- External negative and positive controls included

25/box.....(900-7570)

CPT® Code: 87880QW†



OneStep™ Strep A Test Strips

- CLIA waived
- Sensitivity: 97%; specificity: 95%
- Room-temperature storage
- 25 tests including external negative and positive controls
- Canister packaging for ease of multiple testing

25/box.....(107-4050)

50/box.....(570-0122)

CPT® Code: 87880QW†

Order today at HENRYSCHEIN.COM/MEDICAL



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS





OneStep+ Ultra Flu Test

Intended for the qualitative detection of influenza type A and type B nucleoprotein antigens directly from nasal and nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection.

- CLIA waived
- Results in 10 minutes
- Sensitivity (A/B): 90.3%/88%
- Specificity (A/B): 96.7%/99.2%
- For use with nasal and nasopharyngeal swab specimens
- Made in the USA

(570-1886).....25/box

Contains: 25 test sticks, 25 sterile nasal swabs, 27 extraction buffer vials, 1 influenza A control swab & 1 influenza B control swab.

CPT® Codes: Influenza A, 87804QW[†]/Influenza B, 87804QW-59[†]



OneStep+ Strep A Test Strips

Convenient, compact packaging with everything needed to run 25 tests including positive and negative controls.

- CLIA waived
- Sensitivity: 97%; specificity: 95%
- Room-temperature storage
- Canister packaging for ease of multiple testing
- Walk-away procedure
- 25 tests including external negative and positive controls

(107-4050)..... 25/box

(570-0122)..... 50/box

CPT® Code: 87880QW[†]



OneStep Pro+ Strep A Cassette Kit

- CLIA-waived
- Easy-to-use flip design; all-in-one extraction, incubation, and test
- Results in 5–10 minutes
- Sensitivity of 96.2% and specificity of 98.7%
- External negative and positive controls included
- Made in the USA

(900-7570)..... 25 tests/kit

CPT® Code: 87880QW[†]



OneStep+ Ultra Mono Test

- 3-step process
- 5-minute time to result
- Dipstick test
- Whole blood, serum, or plasma sample type
- CLIA-waived for whole blood
- Moderate complexity for serum and plasma
- Performance (reconciled with EBV serology): 100% sensitivity; 95.9% specificity

(570-1728)..... 25/box

Contains: 25 each of: test sticks in a container, test tubes, transfer pipettes, 50-µL capillary pipettes; and 1 each of: diluent, mono positive control, mono negative control, workstation & package insert.

CPT® Code: 86308QW[†] (whole blood); 86308[†] (serum and plasma)



OneStep+ H. pylori Test Kit

Rapid immunoassay for the qualitative detection of anti-*Helicobacter pylori* antibodies in human whole blood, serum, or plasma.

- May be used with whole blood, serum, or plasma
- CLIA waived for whole blood
- Sensitivity of 95.9% with 91.6% overall agreement to biopsy histology
- Easy-to-use single reagent procedure
- Results in 10 minutes
- Requires only 25 µl of whole blood or 10 µl serum/plasma
- Convenient fingerstick sample collection

(570-2778)..... 20/box

Contains: 20 test devices (test device contains a membrane strip coated with *H. pylori* antigen and a pad with indicator conjugates in a protein matrix), 30 capillary tubes, 1 dropper bottle of developer solution (contains 0.1% sodium azide), instructions for use & procedure card.

CPT® codes: 86318QW[†] (whole blood), 86318[†] (serum/plasma)

HENRY SCHEIN

Cares

www.HenrySchein.com/HSCares





PowerSpin DXB Centrifuge

Featuring a variable-speed range of 500 rpm–5500 rpm to ensure a clean and complete separation. Quiet and dependable, with a maintenance-free motor and stainless-steel tube shields, it can be used to process blood, urine, fecal, and other specimens.

- The all-digital front panel has 3 preprogrammed options: 3 minutes at 4000 g (for BD Barricor™ tubes); 5 minutes at 3200 g; and 10 minutes at 2000 g
- Custom settings are easy and intuitive with an all-digital panel and oversized display
- Zero rpm safety latching system prevents lid from being opened while rotor is spinning
- Quiet running—only 50 dB at 3500 rpm
- 99-minute digital timer with 1-second increments
- 500 rpm–5500 rpm in 100 rpm increments
- Remembers last setting for quick repeat runs
- Adjustable acceleration and deceleration
- Stainless-steel rotor chamber for easy cleaning

#C8763
(570-0483).....ea



PowerSpin™ Centrifuge

Features a 6-place angled rotor designed for 2–10 mL tubes (maximum tube size up to 16 mm x 110 mm) and a fixed speed of 3400 rpm, generating 1350 RCF (g) for superior specimen separation. Set the 0–30-minute timer to desired run length.

- Maintenance-free, brushless motors run quietly at a fixed speed of 3400 rpm, emitting only 58dB, all with a low temperature rise
- Suction-cup rubber feet hold unit firmly in place while virtually eliminating vibration
- Lid features a gasket to reduce aerosol emissions, a window for strobe speed checks, and a safety interlock switch that cuts power to motor if lid is raised
- Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb
- 115 V/60 Hz

Fixed Speed, 6 Place
(112-6905).....ea



PowerSpin™ LX Centrifuge

Spin blood, urine, fecal materials, or semen at the best speed for each type of specimen with this powerful, variable-speed, 6-place or 8-place centrifuge. Maintenance-free, brushless motor runs quietly at a fixed speed up to 4000 rpm, emitting only 58 dB, all with low-temperature rise. Suction-cup rubber feet hold unit firmly in place and virtually eliminate vibration. Lid features a gasket to reduced aerosol emission, a window for strobe speed check, and a safety interlock switch that cuts power to the motor if lid is raised. Angled rotor accepts 2–10-mL tubes (maximum tube size 16 mm x 110 mm).

- Speed range: 300–4000 rpm
- 30-minute timer
- Stainless steel rotor chamber
- Lid latch safety switch
- Dimensions: 10.2" x 11.8" x 7.9"
- Weight: 15 lb
- 115 V/60 Hz

Variable Speed, 6 Place
(112-6907).....ea

Variable Speed, 8 Place
(112-6909).....ea



PortaSpin™ Mobile Centrifuge

- Plugs into a 12-volt car power outlet
- Fixed speed at 3500 rpm
- 6-place fixed angle rotor
- 30-minute timer
- Stainless steel carrying handle
- Optional AC power adapter

#C826
(112-6906).....ea



Micro-Hematocrit Centrifuge

This powerful centrifuge, with all-metal housing, spins up to 24 microhematocrit tubes at 12,000 rpm. Tubes are held in place with a screw-down metal lid and replaceable rubber gasket/ring. The 5-minute mechanical timer has a bell tone to signal a completed run, while the electronic brake can stop the unit on demand. Metal lid features a safety-lid latch that cuts engine power if lid is lifted. Unit comes complete with a metal hematocrit reader tray.

#C-MH30
(112-6910).....ea



Rock-IT™ Blood Tube Rocker

A versatile rocker for hematology (blood) tubes prior to counting; also excellent for vials and gentle mixing of injectables. Nearly silent while gently completing 25 mixing cycles per min. Platform is designed for 15 tubes and has a reversible pad that can be autoclaved. Dimensions: 13.7" x 4.7" x 5"

#L-TTR200, 15 Place
(112-6911).....ea



Vortex Mini Mixer

This economical vortex mixer is rugged and sturdy for the most demanding applications. Dual-mode switch allow user to vortex continuously or intermittently using the touch mode: simply press tube to the soft rubber mixing cup. Features speed control to adjust the vortexing velocity to the user's requirements.

- Powerful and quiet
- Variable speed up to 3000 rpm
- Switch from Auto to Touch mode
- Power-on LED indicator

#L-VM1000
(112-6913).....ea



G380 Microscope

Designed for viewing comfort, this microscope provides crisp, clear imaging on a versatile, ergonomic platform. Standard 10x high-eyepoint eyepieces are engineered to reduce eyestrain for users who wear glasses. Each eye tube features independent focusing diopter adjustment to compensate for the natural vision differential between each eye. Base design adds stability to the body and allows extra room for hand access to the oversized coarse-focus and fine-focus controls and easy access to the coaxial mechanical stage controls.

- Siedentopf binocular head
- Dual diopter adjustments
- 10x high eyepoint widefield eyepieces
- 4X, 10X, 40X & 100X (oil) objectives
- Variable LED illumination
- NA 1.25 Abbe condenser

Microscope, Binocular
(112-6914).....ea
 Contains: 1 dust cover & 2 extra fuses.



Advanced Clinical Microscope, G700 Series

- Tilting binocular head
- Extra-wide-field PL10X eyepieces
- Infinity corrected plan optics: 2X, 4X, 10X, 40X, 100X Oil
- 5-place reversed nosepiece
- Double layer mechanical stage
- True Kohler illumination
- Swing-out Abbe condenser for 2X objective
- Phase contrast available
- Dark field, polarization options

Microscope, Binocular
(570-0437).....ea
 Microscope, Trinocular
(570-0438).....ea



IP730 Series Microscope

Designed to provide crisp, clear imaging on a versatile and ergonomic platform.

- **HEADS:** Binocular or trinocular heads. Featuring 55mm–75mm interpupillary distance adjustment. Both heads rotate 360°, one eye tube with diopter adjustment
- **EYEPIECES:** With 10X, 20mm FN extra-wide field
- **STAGE:** Graduated mechanical stage with low-position coaxial X-Y controls. Ball-bearing movement allows for smooth, precise adjustments. Stage size is 180mm x 150mm and the range of X-Y stage movement is 80mm x 50mm
- **OBJECTIVES:** All models are equipped with either achromatic or plan 4X, 10X, 40XR and 100XR (Oil) objectives
- **CONDENSER:** Focusable and centerable N.A. 1.25 Abbe condenser with graduated iris diaphragm and filter holder
- **ILLUMINATION:** Koehler illumination with field iris diaphragm featuring a 6-V/20-W halogen light source. Variable intensity rheostat wheel is located on the left side of the base
- **ELECTRICAL:** 110 V-60 Hz to 240 V-50 Hz
- 5-year limited warranty, 1 year on electrical system; excludes bulbs and fuses

(112-6915).....ea



Glass Microscope Slides

22 mm x 22 mm; thickness: 0.17 mm–0.25 mm

Plain Slide
(102-7277)..... 72/box
 Frosted End Slide
(102-5605)..... 72/box
 Cover Glass #2, 1-oz Package
(102-6803)..... ea



Top Hat Specimen Collector

800 cc
(112-5810)..... 25/box



Lens Paper

For fine lenses and optical instruments.

- Special high wet strength
- Soft and dust-free
- Will not scratch glass
- Ideal if liquid cleaner is used

4" x 6"
(100-8362)..... 50/pkg



Phlebotomy Tray

- Large carrying capacity
- Includes racks (sharps and alcohol dispenser not included)
- Handy pull-out drawer

#48700-HSI, 9" x 6" x 14"
(112-6916).....ea



Specimen Containers

These 4-oz specimen containers have blue, tamper-evident, leak-resistant tops and are suitable for routine urine or specimen collection and transport. Container includes patient information label.

Sterile, with Tab, Unwrapped, 4 oz
(900-7432)..... 100/case
 Sterile, Individually Wrapped, 4 oz
(900-7431)..... 100/case
 Non-Sterile, Bulk, 4 oz
(900-7432)..... 500/case



A&D Ointment, Unit Dose

Helps treat chafed skin associated with diaper rash and helps seal out wetness. Temporarily protects minor cuts, scrapes, and burns, and temporarily protects and relieves chapped and cracked skin and lips. Also helps prevent and protect from the drying effects of wind and cold weather.


5-gm Packets
(570-3126)..... 144/box
 Contains: 144–5-gm packets per box.



Compare to Active Ingredient in Cortaid®
Hydrocortisone Cream 1%, USP

Itch relief for minor skin irritations, inflammation, and rashes. Active ingredient: hydrocortisone, UPS, 10 mg.

1-oz Tube
(570-3124)..... ea



Compare to Active Ingredient in Vaseline® Petroleum Jelly
White Petrolatum, Unit Dose


Temporarily protects minor cuts, scrapes, burns, and cracked lips. Active ingredient: white petrolatum.

5 gm
(570-3127).....144/box
 Contains: 144 individual 5-gm packets in a dispenser box.



Triple Antibiotic Ointment, USP

1-oz Tube
(570-3130)..... ea



0.9-gm Unit-Dose Packets
(570-3128)..... 144/box



Bacitracin Zinc Ointment, USP

Helps prevent infection in minor cuts and burns. Active ingredient: bacitracin zinc, 500 units.

1-oz Tube
(570-3129)..... ea



0.9-gm Unit-Dose Packs
(570-3125)..... 144/box



Lubricating Jelly

Sterile
 3-gm Unit-Dose Packs
(570-3105)..... 144/box
 4-oz Flip-Top Tube
(570-3106)..... ea



Compare to Tylenol®
APAP Regular Strength Tablets, Unit Dose

Nonaspirin pain reliever and fever reducer.

- Original pain reliever
- 2 tablets per pouch
- Made in the USA

325 mg
(900-4689).....500/box
 Contains: 250 pouches of 2 tablets.



Compare to Advil® or Motrin®
Ibuprofen Tablets, Unit Dose

Each 200-mg tablet provides powerful pain relief from pain and inflammation.

200 mg
(570-0177)..... 500/box
 Contains: 250 pouches of 2 tablets.



Silver Nitrate Applicators, 6" Wood

- Applications: skin or mucous membrane cauterization; removing warts and granulated tissue
- 75% silver nitrate; 25% potassium nitrate
- 50 mg active ingredients/applicator
- Silver nitrate tips are soluble in all secretions
- Control strength by dipping in sterile, distilled water
- Rigid wooden applicator for precision placement
- Convenient tube vial

(112-6994)..... 100/vial

**Hydrogen Peroxide Solution,
USP 3%**

8-oz Bottle
(112-7012)..... ea



16-oz Bottle
(112-7069)..... ea



1-gal Bottle
(112-7070)..... ea



Isopropyl Alcohol ♦

70% Isopropyl
16-oz Bottle
(112-7067)..... ea



32-oz Bottle
(112-7011)..... ea

1-gal Bottle
(112-7068)..... ea



Denatured
For burning. Not for hard-surface disinfection.
16-oz Bottle
(100-0034)..... ea

1-gal Bottle
(100-0108)..... ea

♦ Hazardous material

Distilled Water

- Vapor distilled: clean, natural, and healthy water
- Sturdy PET bottle: excellent for shipping individual cases
- Convenient 3-liter size with molded grip for easy handling
- Pure water—great for use in the medical sector
- Satisfy the demand for distilled water that many medical devices/procedures require
- 2-year shelf life



3-Liter Bottles
(570-2874)..... 4/case

Vaccine Temperature Data Logger

- High accuracy
 - 24 million readings
 - Maximum of 170 files (.csv)
 - Internal (ambient) and external (fridge/freezer) sensors reading display
 - Minimum and maximum reading memories
 - Low and high alarm limits setting
 - Flashing light alarm indicator
 - °C or °F scale
 - Real-time display
 - Waterproof external sensor
 - Low battery indicator
 - Removable SD memory card or USB interface for data retrieval
 - Optional external power
- (570-2722)..... ea



Ear Wash System

A safe and effective option for removing earwax (cerumen).

- Flexible hose and rigid tube options include the bottle of ear wash fluid, trigger mechanism, and 3 tips
- Accessories, including a basin, additional tips, and splash shields for the rigid tube unit, are sold separately*
- Basin is ergonomically molded to fit under the ear
- Simply fill the bottle with warm water, secure the trigger mechanism onto the bottle, twist a disposable tip into the end of the nozzle, and press the trigger to pump the water through the nozzle into the ear
- User should place the basin under the ear to catch the water being flushed out

Ear Wash Kit with Flexible Hose and 3 Tips
(570-1884)..... ea

Ear Wash Kit with Rigid Tube and 3 Tips
(570-1885)..... ea

Ear Wash Tips
(570-1883)..... 20/bag

Plastic Ear Wash Basin
(570-1903)..... ea

Ear Wash Splash Shields for Rigid Tube
(570-1917)..... 10/bag

*There is no option for Ear Wash Splash Shields for the Flexible Hose unit.

Contact our pharmaceutical and vaccine team and let us do the research.

Henry Schein's DxRxSolutions is a single point of contact for Henry Schein customers for answers to various types of questions related to pharmaceuticals, and vaccines. With a mission of helping practices diagnose, prevent, and treat disease, the DxRx team is available as a trustworthy and accessible source of information at your convenience. The DxRx team is available via e-mail, 8:30am–5:00pm (et).

Features:

- DxRx can answer questions related to:
- Latest updates on supply issues
 - Indications and usage as provided on package inserts
 - Brand-to-generic comparisons
 - CPT® and J codes*
 - Product storage requirements
 - ACIP recommendations for vaccinations
 - Flu vaccine and flu-related topics

Benefits:

- Complimentary resource for Henry Schein customers
- Easy-to-access and use resource
- E-mail for answers to pharmaceutical and vaccination questions

*For informational purposes only. Source: www.reimbursementcodes.com or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstance.

Because So Many Rely on You...You Can Rely on Us.



Email: DxRxSolutions@henryschein.com
Visit: www.henryschein.com/drxrsolutions




Lidocaine HCl 2% with Epinephrine Injection, USP

1.7-mL Anesthetic Cartridges

- Validated siliconization process: ensures smooth injection
- Meets all FDA requirements: allows total flexibility in brand selection
- Color-coded stoppers and labeling: easy identification of each anesthetic formula
- Mylar® wrap: prevents breakage
- Blister packaging: no glass/glass contact

1:100,000, Red
(465-1205)..... 50/box
 Contains: 5 blister packs of 10 cartridges per package.

1:50,000, Green
(465-0015)..... 50/box
 Contains: 5 blister packs of 10 cartridges per package.



Articaine 4% with Epinephrine

Indicated for local, infiltrative, or conductive anesthesia in both simple and complex dental procedures. It is available in 2 formulations for you to decide which epinephrine concentration is right for your patient and the procedure.

Articaine 4%..... 50/box
 Specify:

1:100,000 **(112-5510)** 1:200,000 **(112-5511)**



Mepivacaine HCl 3% Plain Injection, USP

1.7-mL Anesthetic Cartridges

- Validated siliconization process: ensures a smooth injection
- Meets all FDA requirements: allows total flexibility in brand selection
- Color-coded stoppers and labeling: easy identification of each anesthetic formula
- Mylar® wrap: prevents breakage
- Blister packaging: no glass/glass contact

(467-3005)..... 50/box
 Contains: 5 blister packs of 10 cartridges per package.

Approved for sale only to licensed practitioners in the USA.

ACCLEAN® 5% Sodium Fluoride Varnish

ACCLEAN®



5% sodium-fluoride varnish is simple to apply. Premeasured unit-dose cups fit easily into enclosed finger ring, and varnish applies effortlessly with included Microbrush® applicator. This easy-to-use, pearlized white varnish sets on contact with saliva and dries to a natural tooth color. Formulated with only natural sweeteners and xylitol. Delivers the highest fluoride concentration available (compared to nonvarnish fluorides).

Unit Dose, 0.5 mL..... 50/box
 Specify:

Bubblegum **(112-4869)** Mint **(112-7169)**
 Cherry **(112-4868)** Raspberry **(112-5963)**

Contains: 50 unit-dose cups, 50 Microbrush® applicators & 1 finger ring.

Bulk Package, 0.5 mL..... 200/box
 Specify:

Bubblegum **(112-5999)** Mint **(112-7170)**
 Cherry **(112-5998)** Raspberry **(112-6001)**



ACCLEAN® 5% Sodium Fluoride Varnish

ACCLEAN®

Brush & Tray System

5% sodium-fluoride varnish in convenient 0.4 mL unit-dose packaging with an application brush. This easy-to-use, translucent varnish sets on contact with saliva and dries to a natural tooth color. Delivers the highest fluoride concentration available with a smooth and even application. Available in 5 great-tasting flavors.

Unit Dose—0.4 mL..... 50/box
 Specify:

Bubblegum **(570-0528)** Mint **(570-0531)**
 Caramel **(570-0529)** Raspberry **(570-0532)**
 Cherry **(570-0530)**

Bulk Package—0.4 mL..... 200/box
 Specify:

Bubblegum **(570-0523)** Mint **(570-0526)**
 Caramel **(570-0524)** Raspberry **(570-0527)**
 Cherry **(570-0525)**

ACCLEAN® Chlorhexidine Gluconate 0.12% Oral Rinse

ACCLEAN®



Antimicrobial oral rinse with 0.12% chlorhexidine gluconate, used for treating gingivitis. Helps reduce tissue inflammation, bleeding, and plaque accumulation. Mint flavor.

4 oz. Bottle—Mint
(570-2250)..... ea

16 oz. Bottle—Mint
(570-2249)..... ea

ACCLEAN® Disposable Prophy Angles

ACCLEAN®



Grip designed for overall control and tactile sensitivity. Not made with rubber latex.

- With cups for high torque and vibration-free operation
- Available with outer ridges for improved interproximal stain removal
- Tapered brush angles for maneuvering around hard-to-reach areas

Fits E-Type
 Angles..... 100/box

Specify:
 Soft, Short, 4 Web, Ribbed **(570-3190)**
 Firm, Short, 4 Web, Ribbed **(570-3189)**
 Soft, Short, 4 Web, Ribbed, with Outer Ridges **(570-3192)**
 Tapered Brush **(570-3191)**

Fits E-Type
 Angles..... 500/box

Specify:
 Soft, Short, 4 Web, Ribbed, with Outer Ridges **(570-3170)**
 Soft, Short, 4 Web, Ribbed **(570-3169)**
 Firm, Short, 4 Web, Ribbed **(570-3168)**



ACCLEAN® Prophy Paste with TCP

ACCLEAN®

1.23% fluoride-ion prophy paste with tri-calcium phosphate for excellent stain removal and polishing performance. Also features enhanced fluoride activity with the addition of TCP, which may help enhance remineralization of the tooth surface and prevent demineralization. Rinses easily and completely and delivers consistent, splatter-free performance formula. Large selection of grits and flavors. Gluten-free. Made in the USA.

Unit-Dose Cups..... 200/box
 Specify:

Extra-Coarse Mint **(570-3188)**
 Coarse Bubblegum **(570-3180)**
 Coarse Cherry **(570-3179)**
 Coarse Mint **(570-3177)**
 Medium Bubblegum **(570-3184)**
 Medium Cherry **(570-3183)**
 Medium Mint **(570-3181)**
 Fine Bubblegum **(570-3185)**
 Fine Cherry **(570-3186)**
 Fine Mint **(570-3187)**

Assorted Pack—Coarse
 Contains: Bubblegum & Cherry **(570-3178)**
 Assorted Pack—Medium
 Contains: Bubblegum & Cherry **(570-3182)**

A

A & D Ointment.....	100
Acclean	102
Adhesive–Skin, Bandages.....	87
Adson–Brown Forceps.....	50
Adson–Dressing Forceps, Tissue Forceps, Micro.....	50–51
Air Spirometers	5
Alcohol Prep Pads.....	12, 81
Alcohol, Isopropyl.....	101
All-Gauze Sponges.....	88
Allis Tissue Forceps.....	52
Alm Retractor.....	66
Analyzers–HB, Lipid, Glucose.....	92
Aneroids	5
Ankle Brace	90
Antibacterial Soap.....	27
Antimicrobial Soap.....	27
APAP Tablets	100
Applicators–Cotton Tipped, Wooden.....	12
Apron	26
Articaine.....	102
Autoclave Bags.....	30
Autoclave Cleaners.....	32
Autoclave–Bags, Pouches.....	30

B

Bacitracin Zinc Ointment.....	100
Back Table Covers	73
Backhaus Towel Clamp	61
Bags–Infectious Waste, Emesis, Trash	39
Bandage and Utility Scissors.....	43
Bandage Shears.....	44
Bandages–Adhesive, Elastic, Self-Adherent, Plastic, Fabric.....	87
Basic Drape Packs	73
Basin, Ring	74
Bergh Cilia Forceps	53
Bibs	10
Billeau Ear Loops.....	69
Biological Indicator Test Pack	32
Biological Monitoring	32
Biopsy–Forceps, Curette, Punch	40
Blades–Surgical, Handles.....	41
Bone-Cutting Forceps.....	50
Boot–Walking, Pneumatic.....	90
Bouffant Caps.....	24
Bowie-Dick Test.....	32
Boynton Needle Holder.....	64
Bozeman Sponge Forceps.....	60
Brace–Ankle, Knee.....	90
Buck Neurological Hammer	67

C

Capes, Exam.....	8
Capillary Tubes	92
CardioChek.....	92
Carmalt Splinter Forceps.....	55
Casting.....	91
Casting Scissors	43

Cast–Padding, Tape	91
Castroviejo Needle Holder.....	62
Cauteries	7
Centrifuges	98
Cervical Biopsy Forceps.....	71
C-Fold Towels.....	13
Chest Drapes.....	73
Chloroprene Gloves	17
Clamps, Towel–Backhaus, Jones, Lorna.....	61
Cleaning Tablets.....	32
Cleanser, Skin.....	27
Coats, Lab	25
Code-A-Color	42
Cold Compresses.....	91
Cold Gel Packs.....	91
Collier Needle Holder	62
Colposcopes.....	6
Combine ABD Pads.....	88
Compresses–Cold, Heat.....	91
Contour Bibs.....	10
Controls–Urinalysis.....	94
Cotton Balls.....	12
Cotton Rolls.....	12
Cotton-Tipped Applicators	12
Cover Gowns.....	26
Cover, Mayo Stand.....	72
Covers, Shoe	24
Cream–Hydrocortisone	100
Crile Forceps	59, 64, 65
Crile Hemostat Forceps.....	51
Crile-Wood Needle Holder	64
Crutches	90
Cryosurgical Kits.....	7
CSR Wraps	30
Cups–Paper, Plastic, Waxed, Medicine	11
Curette–Dermal, Excavator, Biopsy, Ear.....	69

D

Data Logger.....	101
Dean Scissors.....	46
Derf Needle Holder	62
Desiccator	7
Detergents	33
Diagnostic Equipment.....	3–5
Digital Scales.....	6
Digital Thermometers	3
Dispensers, Glove	19
Dispensers, Sanitizer.....	27
Dispensers, Soap	27
Disposable Tray Inserts	42
Disposable Vaginal Speculum	70
Distilled Water	101
Drape Pack, EENT.....	73
Drape Packs.....	73–74
Drape Packs, Basic	73
Drapes, Exam	8
Drapes–Non-Fenestrated, Fenestrated, Split-Sheet, Utility, Half, Chest	72, 73
Dressing Forceps	50–54
Dri-Gard Towels.....	10
Drugs of Abuse Cups.....	93–94
DxRx.....	101

E

Ear Loop Face Masks with Eye Shields.....	24
Ear Loops.....	69
Ear Syringes.....	68
Ear Wash System.....	101
ECG Electrodes.....	7
ECG Paper.....	7
EENT Drape Pack.....	73
EKG Paper	7
Elastic Bandages	87
Electrosurgical Generators	7
Emesis Bags	39
Emesis Basins.....	12
ENT Forceps.....	54
ER Laceration Tray	72
Esmarch Heavy-Duty Bandage and Cast Shears.....	43
Essentials Bibs	10
Essentials Headrest Covers	9
Exam Capes.....	8
Exam Drapes	8
Exam Gloves.....	14–20
Exam Gown.....	8
Exam Shorts.....	8
Exam Stool	3
Excavator Curette.....	69
Exodontia Sponges.....	88
Extractors.....	69
Eye Dressing Forceps	53
Eyewear, Protective	24

F

Face Masks	23–24
Face Shields	24
Fecal Occult.....	95
Feilchenfeld Splinter Forceps.....	55
Fenestrated Drapes.....	72
Fiberglass–Casting Tape, Splint.....	91
FIT Tests.....	95
Fletcher Sponge Forceps	60
Flu Tests.....	97
Fluoride Varnish.....	102
Foerster Sponge Forceps.....	60
Forceps.....	55–60, 71
Fox Dermal Curette.....	69
Frazier Dura Hooks.....	67

G

Gauze Shears.....	44
Gauze–Rolls, Pads, Sponges, Syngauze, Type VII.....	88–89
Gel Packs–Cold, Hot	91
General Drapes.....	73
Germicidal Cloths.....	28–29
Glasses, Safety.....	24
Glove Dispensers.....	19
Gloves, Chloroprene.....	17
Gloves, Exam.....	14–20
Gloves, Latex.....	18
Gloves, Nitrile	16–17

Gloves, Polychloroprene	17
Gloves, Polyisoprene	17
Gloves, Surgical	17-18
Gloves, Vinyl	18
Glucose Analyzer	92
Glucose Tests	92
Glutaraldehyde Solution	29
Goldman-Fox Scissors	48
Gown, Exam	8
Gowns, Cover	26
Gowns, Isolation	26
Gowns, Protective	26
Gradle Eye Suture Scissors	45

H

H. Pylori Tests	97
Half Drapes	73
Halsey Needle Holder	64-65
Halstead Mosquito Forceps	56
Halstead Mosquito Hemostat Forceps	56
Hammer	67
Hand Lotion	27
Hand Sanitizer	27
Hand Sanitizer Dispensers	27
Hand Soap	27
Hand Soap Dispensers	27
Hand Wipes	27
Hartmann Mosquito Forceps	56
Hartmann-Noyes Alligator Forceps	54
Hartmann-Noyes Ear Forceps	54
HB Analyzer	92
HBA1C Tests	93
Headrest Covers	9
Heat Compresses	91
Hemoglobin Analyzer	92
Hemostat Forceps	56-59
Hercules Heavy-Duty Bandage and Plaster Shears	43
Hooks	67
Hot Gel Packs	91
Hudson Tissue Forceps	51
Hydrocortisone Cream	100
Hydrogen Peroxide	101

I

Ibuprofen Tablets	100
Ice Packs	91
iFob	95
Indicator Strips	31
Infectious Waste Bags	39
Instrument Cassettes	31
Instrument Identification—Rings, Tape	42
Instrument Tray	42
Iris—Scissors, Tissue Forceps	49, 53
Isolation Gowns	26
Isopropyl Alcohol	101
IV Administration Sets	37
IV Catheters, Safety	37-38
IV Poles	3

J

Jackets, Lab	25
Jelly—White Petrolatum, Lubricating	100
Jones Towel Clamp	61

K

Kelly Angular Scissors	49
Kelly—Scissors, Forceps	49
Kits—Minor Procedure	72, 74
Knee Brace	90
Knowles Bandage Scissors	43

L

Lab Coats	25
Lab Jackets	25
Lancets, Safety	37
Laparotomy Sponges	81
Laryngeal Mirrors	68
Latex Gloves	18
Lens Paper	99
Lidocaine Injection	102
Lipid Analyzer	92
Lipid Tests	92
Lister Bandage Scissors	43
Littauer Jr. Stitch Scissors	44
Littauer Stitch Scissors	44
Lorna Towel Clamp	61
Lotion, Hand	27
Lubricating Jelly	100
Lucae Forceps	54

M

Mallet	67
Markers—Skin	81
Masks, Face	23-24
Mathieu Needle Holder	63
MAXICIDE	29
Maxiclens	27
MaxiSpray	29
MAXITABS	33
MaxiTest	31-32
MAXIWIPE	28-29
Maxizyme Tabs	33
Mayo Dissecting Scissors	46-47
Mayo Stand	3
Mayo Stand Cover	72
Mayo-Hegar Needle Holder	65
Mechanical Scales	6
Medicine Cups	11
Mepivacaine Injection	102
Metzenbaum Scissors	46-47
Micro Adson Forceps	50
Microbore Extension Set	37
Microscope Slides	99
Microscopes	99
Minor Procedure Kits—ER Laceration, Suture Removal, Staple Removal	72

Mirrors—Laryngeal	68
Mixers	98
Mono Tests	97
Multifold Towels	13

N

Nail Scissors	49
Nail—Nipper, Splitter	68
Nasal Speculum	68
Nebulizer—Compressor System, Kit	5
Needle & Syringe Combination, Conventional	34
Needle & Syringe Combination, Safety	34
Needle Holders	63-65
Needle, Safety Insulin Syringes	35
Needle-Free Connector	37
Needles with Tuberculin Syringes	36
Needles, Blood Collection	37
Needles, Blunt Fill	35
Needles, Conventional	35-36
Needles, Multiple Sample Safety Blood Collection	37
Needles, Safety	34
Nipper, Nail	68
Nitrile Gloves	16-17
Non-Fenestrated Drapes	72
Non-Woven Sponges	87
Northbent Scissors	44

O

O.R. Towels	81
O'Brien Scissors	45
Ointment—A & D, Triple Antibiotic, Bacitracin Zinc,	100
Olsen-Hegar Needle Holder	63-64
OPA Solution	29
Operating Room Scissors	47-48
Ophthalmic Drapes	73
Ophthalmoscopes	3
Oral Rinses	102
Otoscope Speculum Tips	3
Otoscopes	3
Overbed Table	3

P

Packing Strips	89
Padding—Cast	91
Palm Aneroids	5
Paper Cups	11
Paper Towels	13
Patient Belonging Bags	8
Patient Footwear	8
Percussion Hammer	67
Peroxide, Hydrogen	101
Phlebotomy Tray	99
Pillowcases	9
Plain Splinter Forceps	55
Plastic Cups	11
Plastic Surgery Scissors	49
Pocket Aneroids	5
Polychloroprene Gloves	17

Polyester Casting Tape.....	91
Polysoprene Gloves.....	17
Portaspin.....	98
Povidone-Iodine Swabsticks and Prep Pads.....	81
Powerspin.....	98
Pregnancy Tests.....	95
Prep Pads, Alcohol.....	12, 81
Prophy Angles.....	102
Prophy Paste.....	102
Protective Eyewear.....	24
Protective Gowns.....	26
PTS Panels.....	92
Pulse Oximeters.....	3-4
Punch, Biopsy.....	41

R

Ragnell Retractor.....	66
Razors.....	81
Reader—Urinalysis, Paper.....	94
Reagents—Urine.....	94
Retractors.....	66-67
Ring Basins.....	74
Rochester-Carmalt Forceps.....	58
Rochester-Ochsner Forceps.....	58-59
Rochester-Pean Forceps.....	58
Rockers.....	98
Rock-It.....	98
Ruskin Bone-Cutting Forceps.....	50
Ryder Needle Holder.....	62

S

Safety Glasses.....	24
Safety Lancets.....	37
Sanitizer, Hand.....	27
Scales—Mechanical, Digital, Wheelchair, Pediatric.....	6
Scalpels—Safety, Disposable.....	41
Schamberg Comedone Extractors.....	69
Schroeder-Braun Uterine Tenaculum Forceps.....	71
Scissors.....	43-49
Self-Adherent Bandages.....	87
Self-Retaining Retractor.....	66
Self-Seal Sterilization Pouches.....	30
Senn Retractor.....	66
Sharps Containers.....	39
Shields, Face.....	24
Shoe Covers.....	24
Shoe—Surgical.....	90
Shortbent Scissors.....	44
Shorts, Exam.....	8
Silver Nitrate Applicators.....	100
Skin Adhesive.....	82-83
Skin Cleanser.....	27
Skin Markers.....	81
Slides, Microscope.....	99
Soap—Hand, Antimicrobial, Antibacterial.....	27
Sound, Uterine.....	71
Specimen Collectors.....	99
Specimen Containers.....	95, 99
Specula—Nasal, Vaginal.....	70

Speculum—Nasal, Vaginal, Disposable, Lighted, Steel.....	70
Spencer Scissors.....	45
Spencer Stitch Scissors.....	45
Splinter Forceps.....	55-56
Splint—Plantar, Fiberglass.....	90
Split-Sheet Drapes.....	73
Splitter, Nail.....	68
Sponge Forceps.....	60
Sponges—Laparotomy, All-Gauze, Non-Woven, Exodontia, Syngauze.....	88-89
Stand, Mayo.....	3
Standard Procedure Packs.....	75-79
Standard Sets.....	75-79
Staple Removal Kits.....	72
Staple Remover.....	82
Staple—Skin, Remover.....	82
Steam Indicator Labels.....	31
Steam Indicator Tape.....	31
Sterilization Integrators.....	31
Sterilization Solutions.....	29
Sterilization Tubing.....	31
Sterilizer Test Pack.....	32
Stethoscopes.....	5
Stevens Tenotomy Scissors.....	45
Stitch Scissors.....	44-45
Stockinette.....	91
Stool, Exam.....	3
Strabismus Scissors.....	47
Strep A Tests.....	97
Stretcher Sheets.....	9
Strips—Instrument Identification, Wound Closure, Test.....	42, 83, 92, 94-95
Surface Disinfectants.....	28-29
Surgical Blades.....	41
Surgical Gloves.....	17-18
Surgical Gowns.....	81
Surgical Handles.....	41
Surgical Instruments.....	40
Surgical Mallet.....	67
Surgical Masks.....	23
Surgical Milk.....	29
Surgical Shoe.....	90
Suture Removal Kits.....	72
Suture Scissors.....	44
Sutures—Absorbable, Non-Absorbable, Spools.....	83-86
Swiss Cilia Suture/Splinter Forceps.....	53
Syngauze.....	87
Synthetic—Cast Padding, Stockinette.....	91
Syringe—Ear.....	68
Syringes.....	35
Syringes with Exchange Tray.....	36

T

Table Paper.....	9
Table, Overbed.....	3
Tablets—Aspirin, Ibuprofen, APAP.....	100
Tape—Cloth, Adhesive, Surgical, Hypoallergenic, Transparent, Paper, Casting, Fiberglass, Polyester.....	89
Temperature Monitoring.....	101
Test Packs.....	32

Test Strips—Lipid, Glucose, Urine, Reagent, Pregnancy, Strep A.....	92-97
Tests—Lipid, Glucose, HBA1C, Urine, Pregnancy, FIT, iFob, Fecal Occult, Flu, Mono, Strep A, H. Pylori.....	92-97
Thermometers.....	3
Thumb Forceps.....	51
Tissue Forceps.....	50-53
Tissues.....	13
Tongue Depressors.....	12
Topical Skin Adhesive.....	82-83
Towel Clamps—Backhaus, Jones, Lorna.....	61
Towels.....	10
Towels, Dri-Gard.....	10
Towels—O.R.....	81
Toxicology.....	93-94
Trash Liners.....	39
Tray Covers.....	10
Triple Antibiotic Ointment.....	100
True Metrix.....	92
Tubes—Capillary.....	92
Type VII Gauze.....	88

U

Ultrasonic Cleaning Solutions.....	33
Ultrasound Paper.....	7
Underpads.....	13
Universal Scissors.....	43
Unna Comedone Extractors.....	69
Unna's Boot.....	89
Urinalysis.....	94
Urinalysis Controls.....	94
Urine Strips.....	94
Urine Tests.....	93-94
Urispec.....	94
Uterine Sound.....	71
Utility Drapes.....	73
Utility Sterilizer Forceps.....	58

V

Vaginal Speculum—Steel Disposable, Lighted.....	70
Varnish.....	102
Vinyl Gloves.....	18
Virtus Splinter Forceps.....	56
Vortex.....	98

W

Walking Boot.....	90
Wall Systems.....	3
Walter Splinter Forceps.....	53
Washcloths.....	13
Waxed Cups.....	11
Webster Needle Holder.....	63
Weitlaner Retractor.....	67
White Petrolatum.....	100
Wilde Forceps.....	54
Wipes, Hand.....	27
Wire-Cutting Scissors.....	43
Wooden Applicators.....	12
Wound Closure Strips.....	82

Terms of Sale

THE HENRY SCHEIN PRICE POLICY:

We endeavor to maintain prices for the duration of a catalog, but we reserve the right to make price adjustments in response to manufacturers' price increases or extraordinary circumstances. Prices are subject to change without notice.

Henry Schein, Inc. ("Henry Schein") and customer agree that the terms and conditions hereinafter set forth shall govern the relationship between Henry Schein and the customer to the extent that the parties do not have a written agreement in effect that conflicts with such terms and conditions. Customer acknowledges and accepts all such terms and conditions by placing an order for goods with Henry Schein, and upon Henry Schein's delivery of the order to the customer.

Choose Your Payment Method

Reduce the cost and administration of paying Henry Schein—

Pay electronically (ACH Debit) or set up AutoPay. Please call Customer Service for details.

For your convenience, we provide several payment alternatives. Orders billed to your account may be paid by ACH Debit, Check by Phone, or Check. If you prefer, you may use your Henry Schein Credit Card, American Express, Visa, MasterCard or Discover Card when placing your order. All sales are subject to our normal terms and conditions.

Unless otherwise instructed, check payments must be mailed to: Henry Schein, Inc. • Dept. CH 10241 • Palatine, IL 60055-0241

All sales are subject to credit approval. Invoices are payable within agreed terms of sale. Henry Schein and its affiliates reserve the right to offer a dual pricing system for payment by check, Henry Schein-approved ACH and/or Henry Schein-approved EFT vs. any other payment method and, to the extent permitted by law, reserves the right to surcharge for payment by credit card and/or other methods incurring a fee to Henry Schein.

Open Accounts Receivable:

All unpaid accounts receivable past due are subject to a 1.5% finance charge.

DELIVERY TERMS: Delivery dates are not guaranteed. Purchase orders are subject to product availability. Customers' sole and exclusive remedy for failure to deliver shall be refund of monies paid for the products at issue. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Except as noted below, title passes at the time the shipment is loaded at the shipper's dock.

California: For all shipments of goods to customers located within California, title will pass upon receipt of goods by California customers.

Continental U.S.: All orders will be subject to a handling charge. This charge includes freight, except for additional carrier charges related to special delivery services and hazardous material shipments. Special orders are subject to additional freight charges.

Alaska, Hawaii & Pacific Protectorates: Standard shipping methods provide direct, reduced cost, expedited air delivery service to all accounts in Alaska and Hawaii. Customers in the Pacific Protectorates are offered direct surface transport, or postal services for reliable delivery. No additional surcharges apply, except when special services are requested. Low-level hazardous items (dangerous goods in accepted quantities and Consumer Commodity ID 8000) are now available via UPS 2nd-day air.

Guam, Puerto Rico, U.S. Trust Territories & Virgin Islands: All orders will be subject to a handling charge. This charge includes freight through the United States Postal Service (USPS). Special delivery orders and hazardous material shipments can be shipped via United Parcel Service (UPS) for an additional charge. No minimum order amount or weight applies. Speak to your International Representative for details.

Outside U.S. (50 States): If your order is being shipped outside the U.S. (50 states), please refer to the International Terms & Conditions at henryschein.com. Unless otherwise agreed, freight terms are FOB Shipper's Dock ("Ex Works" outside North America). Title passes at the time the shipment is loaded at the shipper's dock. Customer is responsible for compliance with any applicable import requirements.

RX PRODUCTS & CONTROLLED SUBSTANCES:

Regulations require us to limit the sale of Rx and controlled substances only to registered, licensed healthcare professionals. If you are a new customer or have recently moved, please furnish us with a copy of your updated state and federal DEA registration (including verifying your shipping address). Please note that all orders for controlled substances are subject to a due diligence review process. Schedule II controlled substances can be ordered electronically or by mail. For more information on our Controlled Substance Ordering System please visit www.henryschein.com/e222; if you prefer to continue using Federal 222 Forms to order Schedule II controlled substances, mail the form to: Henry Schein, Inc. • 5315 West 74th Street • Indianapolis, IN 46268

Henry Schein restricts the sale or other transfer of medications to prisons/correctional facilities for use in lethal injections, based on our manufacturer agreements. The goods Henry Schein sells are intended to be used for their label-approved purposes or applicable standards of care, which do not include human lethal injection.

THE DRUG SUPPLY CHAIN SECURITY ACT (DSCSA):

(MN, DM, WH, M2) The Drug Supply Chain Security Act (DSCSA) information related to prescription drug products is available on our website www.henryschein.com/pedigree. If you have any problems accessing our website or would like to receive a copy of DSCSA documentation via fax, mail or email, please contact our customer service department.

REGULATORY REQUIREMENTS:

Local regulatory requirements may apply to use or installation of certain products. Be sure to understand and comply with any such requirements prior to purchase, use, or installation of products.

RETURNS:

WE CANNOT ACCEPT ANY RETURNS WITHOUT PRIOR AUTHORIZATION.

To arrange for a return, simply call our Customer Service department or contact your Sales Consultant. The following conditions must be complied with:

• All returns must be accompanied by a copy of your invoice and a reason for the return • Merchandise

DISCOUNTS, REBATES AND DISCLOSURES: Invoice or statement prices may reflect or be subjected to a bundled discount or rebate pursuant to purchase offer, promotion or discount program. You must fully and accurately report to Medicare, Medicaid, Tricare and/or any other federal or State program, upon request by such program, the discounted price(s) or net price(s) for each invoiced item, after giving effect to any applicable discounts or rebates, which price(s) may differ from the extended prices set forth on your invoice. Accordingly, you should retain your invoice and all relevant information for your records. It is your responsibility to review any agreements or other documents, including offers or promotions, applicable to the invoiced products/prices to determine if your purchase(s) are subject to a bundled discount or rebate. Any such discounts must be calculated pursuant to the terms of the applicable purchase offer, promotion or discount program. Participation in a promotional discount program is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the discount program requirements.

must be returned in its original unopened container, unmarked, and properly packaged • Returned products must have been purchased within the previous thirty (30) days. Any returns past thirty (30) days are subject to a restocking fee • Shortages or errors in shipments must be reported within seven (7) days of invoice date to issue credit (if applicable) • Shipping charges will apply on all returns.

Dental equipment returns are subject to Henry Schein's equipment return policy, available upon request.

Exceptions: The following special, customized, or government-regulated items are not returnable:

• Opened handpieces, small equipment, and custom-ordered equipment • Immune globulin products • Special order items (products that we do not ordinarily stock) • Personalized and imprinted items • Opened computer hardware and software • Hazardous/flammable materials • Expired products • Items that cannot be returned to the manufacturer • Any item marked non-returnable or non-refundable • Items required to be shipped and stored frozen • Any drop-shipped products • Invitro diagnostics •

Prescription Drug Returns: Please note that, in order to comply with Federal and State traceability requirements, prescription drugs may be returned providing that the following key elements are met:

- 1) Returns of prescription drugs will only be accepted if Henry Schein is notified within 30 calendar days of shipment date and valid return authorization is issued by Henry Schein.
- 2) The Prescription Drug Marketing Act requires any customer returning prescription drugs to complete and return a Prescription Drug Return Authorization form. Federal law requires that the healthcare entity returning prescription drugs document that the product was kept under proper storage and handling conditions while in their possession and during the return of the product. To get a copy of the form and proper return authorization, please contact Customer Service.
- 3) In addition, traceability regulations require that the healthcare entity returning prescription drugs certifies that the product being returned is the same exact product purchased from Henry Schein.
- 4) Henry Schein will not issue credit for any returned prescription drugs which return was not authorized as provided herein, have been tampered with or where the labeling has been altered in any way.

INSTITUTIONAL, GOVERNMENT & CORPORATE ACCOUNTS:

Terms of Sale follow the same guidelines unless denoted differently in a contract. When applicable, freight is calculated by the weight. There is no minimum order; however, there may be a nominal handling charge on small orders. Some offers and promotions outlined in catalog may not apply. Requests for bids and proposals may be sent to: Henry Schein, Inc. • Bid Department (Dental-E255) (Medical-E270) • 135 Duryea Rd. • Melville, NY 11747.

INTERNATIONAL ORDERS:

We proudly serve healthcare professionals, governments, and dealers throughout the world. To place orders or for inquiries on export terms and conditions please contact the International Department (USA) by phone 1-631-843-5325, fax: 1-631-843-5676, or send us an e-mail at: export@henryschein.com.

ARBITRATION: All Claims related to or arising under or relating to this Agreement except for Claims for non-payment of amounts due to Henry Schein for the sale of products) are to be exclusively and finally determined by binding arbitration in the state of New York, or another location mutually agreeable to the parties. Any and all Claims must be arbitrated on an individual basis, and there shall be no right or authority for any Claims or disputes to be arbitrated on a class action or collective basis. For avoidance of doubt, each party irrevocably waives any right to: (i) have any Claim resolved in connection with any class action or collective action, or (ii) recover any damages or relief directly or indirectly as part of any class action or collective action. The arbitration shall be conducted on a confidential basis pursuant to the Commercial Arbitration Rules of the American Arbitration Association, or if applicable, under its Procedures for Large, Complex Commercial Disputes. Any decision or award as a result of any such arbitration proceeding shall be in writing and shall provide an explanation for all conclusions of law and fact and shall include the assessment of costs, expenses and reasonable attorneys' fees. Any such arbitration shall be conducted by an arbitrator experienced in the disputed subject matter and shall include a written record of the arbitration hearing. The parties reserve the right to object to any individual who shall be or has been at any time employed by or affiliated with a competing organization or entity. An award of arbitration may be confirmed in a court of competent jurisdiction. To the extent that any Claim or dispute is determined to not be subject to arbitration, all other Claims or disputes that would otherwise be subject to arbitration must be arbitrated. As used in this Agreement, "Claims" shall mean any and all liabilities, disputes and expenses whatsoever including, without limitation, claims, adversary proceedings (whether before a court, administrative agency or any other tribunal), damages (whether compensatory, multiple, exemplary or punitive), judgments, awards, penalties, settlements, investigations, costs, responses to subpoenas or other governmental directives and reasonable attorneys' fees and disbursements with respect to any claims that may be sustained, suffered or incurred by a Party hereto.

WARRANTIES: Henry Schein will pass through to the customer, at the time of sale, any transferable product warranties, indemnities and remedies provided to Henry Schein by the applicable manufacturer. EXCEPT AS OTHERWISE PROVIDED HEREIN, TO THE EXTENT PERMITTED BY LAW, HENRY SCHEIN PROVIDES NO WARRANTIES OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION ANY WARRANTY OF MERCHANTABILITY, FITNESS FOR PARTICULAR PURPOSE OR NON-INFRINGEMENT, AND THE CUSTOMER SHALL LOOK TO THE MANUFACTURER OF THE PRODUCT FOR ANY WARRANTY THEREON.

New dental equipment installations include a 90 day labor warranty on such installation.

LIMITATION OF LIABILITY: THE CUSTOMER AGREES TO LOOK SOLELY TO THE MANUFACTURER OF THE PRODUCT FOR ANY CLAIM ARISING DUE TO LOSS, INJURY, DAMAGE OR DEATH RELATED TO THE USE OR SALE OF PRODUCTS. HENRY SCHEIN SHALL NOT BE LIABLE FOR INDIRECT, INCIDENTAL, PUNITIVE, SPECIAL OR CONSEQUENTIAL DAMAGES INCLUDING, BUT NOT LIMITED TO, LOST PROFITS AND LOSS OF GOODWILL, ARISING FROM OR RELATING TO ANY BREACH TO THIS AGREEMENT (OR ANY DUTY OF COMMON LAW, AND WHETHER OR NOT OCCASIONED BY THE NEGLIGENCE OF HENRY SCHEIN OR ITS AFFILIATES), REGARDLESS OF ANY NOTICE OF THE POSSIBILITY OF SUCH DAMAGES.



Quality Products, Affordable Prices, and Outstanding Service.

Why choose Henry Schein® Products?

The products you use in your health care facility have to be high quality, whether it's reliable tools you use for a procedure or supplies that help prevent the spread of infection. But because you're running a business, the price of those products is always a factor when you're making a purchasing decision.

That's why we're proud to offer products that support optimal patient care and help protect your bottom line. Additionally, our products are backed by Henry Schein®, a company customers have known and trusted for more than 93 years.



EXCLUSIVE BRANDS
AND PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS

SOLUTIONS FOR HEALTH CARE PROFESSIONALS



Business Solutions

to help customers successfully manage and grow their business



Clinical Solutions

to stay on the leading edge of patient care



Technology Solutions

to help harness technology for the benefit of the clinical setting and patients served



Supply Chain Solutions

for a complete product selection and inventory to deliver the right products at the right time



To Order: **1-800-772-4346** 8am–8pm (et)

henryschein.com/medical

California customers please use website when ordering for Prop 65 information.

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2026 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products. *For informational purposes only. Source: ReimbursementCodes.com and/or www.cms.gov. Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.